




**ADC® TUNING FORKS**

Non magnetic, corrosion resistant aluminum alloy construction weighs 1/3 of comparable steel tuning forks. Produced from 3/8" x 1" bar stock for superior performance and consistent frequency accuracy. Extra long 2" handle of turned smooth aluminum to facilitate bone conduction tests.

- 500128 Tuning Fork w/fixed weight, 128cps Frequency
- 500256 Tuning Fork w/fixed weight, 256cps Frequency
- 500512 Tuning Fork w/o weight, 512cps Frequency
- 501024 Tuning Fork w/o weight, 1024cps Frequency
- 502048 Tuning Fork w/o weight, 2048cps Frequency
- 504096 Tuning Fork w/o weight, 4096cps Frequency



**MILTEX ALLOY TUNING FORKS**


Aluminum alloy.

- 19-100 C-64 Vibrations Tuning Fork
- 19-102 C-128 Vibrations Tuning Fork
- 19-104 C-256 Vibrations Tuning Fork
- 19-106 C-512 Vibrations Tuning Fork
- 19-108 C-1024 Vibrations Tuning Fork
- 19-110 C-2048 Vibrations Tuning Fork
- 19-112 C-4096 Vibrations Tuning Fork
- 19-120 Set of 5 Tuning Forks, C-128, C-256, C-512, C-1024 & C-2048 in Fitted Soft Case
- 19-120A Fitted Zippered Case Only For Tuning Fork Set

**ADC® NEUROLOGICAL HAMMERS**

Four of the most popular hammers for diagnosis of neurological function.

- 3691 Buck Hammer, 7 1/4", Chrome Plated Handle w/2 sided rubber head, Handle Conceals "screw-in" Brush, Needle Contained Within The Head
- 3693 Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Orange
- 3693BK Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Black
- 3693DG Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Dark Green
- 3693RB Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Royal Blue
- 3695 Wartenberg Pinwheel, 7 1/2", Stainless Steel Handle w/textured grip, Rotating Spur
- 3697 Babinski Hammer, 8 1/2", Octagonal Stainless Steel Handle w/concealed needle, Rubber Head



1-200



1-220

**MILTEX HAMMERS**

- 1-200 Taylor Percussion Hammer, 7", Solid Handle, Stainless Steel
- 1-202 Taylor Percussion Hammer, 7 1/2", Solid Handle, Chrome
- 1-208 Taylor Percussion Hammer, 9", Hoop Handle, Chrome
- 1-210 Dejerine Percussion Hammer, 9 1/2", Chrome
- 1-211 Tromner Percussion Hammer, 9 3/4" w/large & small rubber tips, Hollow Tapered Handle w/sharp point, Chrome
- 1-212 Combo-5 Neuro Diagnostic Set Includes Buck Hammer, C-128 Tuning Fork, Pin Wheel, Brush & Needle
- 1-213 Combo-5 Neuro Diagnostic Set Includes Buck Hammer, C-256 Tuning Fork, Pin Wheel, Brush & Needle
- 1-214 Combo-4 Neuro Diagnostic Set Includes Buck Hammer, Pin Wheel, Brush & Needle
- 1-220 Buck Neurological Hammer, 7 3/4" w/brush & needle (screw into handle), Chrome
- 1-221-A Needle Only For Buck Hammer
- 1-221-SM Rubber End For Buck Hammer, Small
- 1-222 Rabiner Neurological Hammer, 9" w/brush & needle (screw into handle), Chrome



**MILTEX WARTENBERG PINWHEEL**

For neurological examination. Stainless steel.

- 1-226 Pinwheel

**MILTEX HEAD MIRROR**

- 1-240 Frankel Head Band & Mirror Set, Includes 3 1/2" Mirror w/1/2" aperture & #1-244 Head Band
- 1-242 Head Mirror, 3 1/2" Dia w/1/2" aperture, Boilable
- 1-243 Head Mirror, 3" Dia w/5/8" aperture, Boilable
- 1-244 Frankel Head Band Only, Black Fibre, Secure Thumb Screw Size Adjustment, Includes Swivel Mirror Connector

**MILTEX MAUNDER ORAL SCREW**

- 2-20 Oral Screw, White Plastic

## MILTEX BINOCULAR MAGNIFYING LOUPE

Lightweight. Coated optics 1.8X magnification. 30 cm focal distance. 10 cm viewing case. Pupillary distance adjustment.


- 1-4000 Magnifying Loupe w/spectacle frame & carrying case
- 1-4006 Magnifying Loupe Only
- 1-4010 Accessory Clip-On Lenses For 1-4000 & 1-4006, Increasing Magnification to 2.5X w/17 cm focal distance & 4.5 cm viewing area

## MILTEX EKG CALIPER

- 1-3015 Economy Line Caliper w/protective black plastic cap, Graduated in mm, Nickel Plated

## MILTEX FINGER NAIL DRILL

- 33-230 4" Finger Nail Drill, Standard Pattern, Complete w/3 drill points stored in handle
- 33-232 4 1/2" Finger Nail Drill w/chuck head, Complete w/3 drill points stored in handle
- 33-234 Set of 3 Drill Points: Small, Medium & Large, Chrome



## MILTEX MICRO INSTRUMENT STERILIZING CASES

With handles and removable cover. Double rows of pins welded to bottom hold instruments in place. Strips of silicone rubber cushions on bottom and cover give added protection to instruments. Perforated cover and bottom assure complete heat penetration with no build-up of condensation.

- 3-550 8" x 8 1/2" x 1 1/4" Instrument Case
- 3-552 8" x 12 1/4" x 1 1/4" Instrument Case
- 3-554 10" x 10" x 1 1/2" Instrument Case
- 3-570 8 1/2" x 8" x 1 1/4" Instrument Case, Double Rows Of Silicone Rubber On Bottom
- 3-572 8 1/2" x 12" x 1 1/4" Instrument Case, Double Rows Of Silicone Rubber On Bottom
- 3-574 10 1/2" x 10" x 1 1/2" Instrument Case, Double Rows Of Silicone Rubber On Bottom
- 3-576 10 1/2" x 15" x 1 1/4" Instrument Case, Double Rows Of Silicone Rubber On Bottom

### Accessories

- 3-550A Silicone Replacements For 3-550
- 3-552A Silicone Replacements For 3-552
- 3-554A Silicone Replacements For 3-554



## MILTEX RING INSTRUMENT HOLDERS

Latch has snap-lock feature to facilitate loading and removing of instruments.

- 3-70 8" Ring Instrument Holder
- 3-71 10" Ring Instrument Holder
- 3-72 12" Ring Instrument Holder
- 3-74 16" Ring Instrument Holder

## MILTEX STAINLESS WIRE BRUSH

For cleaning rasps and files. Autoclavable.


- 21-379 6 1/4" Wire Brush



## MILTEX STERILIZER FORCEPS

Made of stainless steel. 3-prong.

- 3-10 8" Sterilizer Forceps
- 3-12 12" Sterilizer Forceps



## MILTEX TIP-IT™ INSTRUMENT GUARDS

Used to protect instrument tips and edges from damage during instrument handling and sterilization. Non-sterile. Disposable. Steam and gas sterilizable. Radiopaque. Medical grade plastic. Latex free. Packaging: 50/package.

- 3-2501 Size 1 Instrument Guard, 1/16" x 3/4", White
- 3-2502 Size 2 Instrument Guard, 5/64" x 3/4", Blue
- 3-2503 Size 3 Instrument Guard, 7/64" x 3/4", Green
- 3-2504 Size 4 Instrument Guard, 1/8" x 1", Red
- 3-2505 Size 5 Instrument Guard, 3/16" x 1", Yellow
- 3-2506 Size 6 Instrument Guard, 13/32" x 3/4", Maroon
- 3-2507 Size 7 Instrument Guard, 1/16" x 3/8" x 1", Orange
- 3-2508 Size 8 Instrument Guard, 5/64" x 5/8" x 1", Brown
- 3-2509 Size 9 Instrument Guard, 1/8" x 1" x 1", Black
- 3-2510 Assorted Instrument Guard

**MILTEX TIP IT® INSTRUMENT GUARDS**

Protect instrument tips and edges from damage during handling and sterilization. PVC medical grade plastic for use in medical applications. Latex-free material that is non-allergenic. Large holes (vented styles) to reduce vacuum and moisture build up for easier removal from instrument. Steam and gas permeable to allow for effective sterilization of instruments. Four styles in nine colors and sizes for a variety of choices.

3-2501V	Size 1, 1/16" x 3/4" Solid Vented Instrument Guard, White, Non-Sterile, 50/pkg
3-2502V	Size 2, 5/64" x 3/4" Solid Vented Instrument Guard, Blue, Non-Sterile, 50/pkg
3-2503V	Size 3, 7/64" x 3/4" Solid Vented Instrument Guard, Green, Non-Sterile, 50/pkg
3-2504V	Size 4, 1/8" x 1" Solid Vented Instrument Guard, Red, Non-Sterile, 50/pkg
3-2505V	Size 5, 3/16" x 1" Solid Vented Instrument Guard, Yellow, Non-Sterile, 50/pkg
3-2506V	Size 6, 13/32" x 3/4" Solid Vented Instrument Guard, Maroon, Non-Sterile, 50/pkg
3-2507V	Size 7, 1/16" x 3/8" x 1" Solid Vented Instrument Guard, Orange, Non-Sterile, 50/pkg
3-2508V	Size 8, 5/64" x 5/8" x 1" Solid Vented Instrument Guard, Brown, Non-Sterile, 50/pkg
3-2509V	Size 9, 1/8" x 1" x 1" Solid Vented Instrument Guard, Black, Non-Sterile, 50/pkg
3-2510V	Size 10, Assorted Solid Vented Instrument Guards (equal number of sizes plus 5), Non-Sterile, 50/pkg

**MILTEX CUTICLE SCISSORS**

Curved blades. Chrome.

40-435	3 1/2" Cuticle Scissors, Standard
40-445	3 1/2" Cuticle Scissors, Extra Delicate

**MILTEX INGROWN NAIL SPLITTING SCISSORS**

Straight with one serrated blade.

40-260	6" Nail Splitting Scissors
--------	----------------------------

**MILTEX NAIL SCISSORS**

Curved blades. Chrome.

40-425	3 1/2" Nail Scissors
--------	----------------------

**MILTEX NAIL SPLITTER**

40-230	5" Nail Splitter, Anvil Type
--------	------------------------------


**MILTEX PLATYPUS NAIL PULLING FORCEPS**

Stainless.

40-234	5 1/2" Nail Pulling Forceps, Standard Wide Jaws
40-235	5 1/2" Nail Pulling Forceps, Narrow Jaws

**MILTEX TISSUE & CUTICLE NIPPERS**

40-220	4" Cuticle Nipper, Ring Handle, Chrome, Straight Jaws
40-245CH	4" Cuticle Nipper, Chrome, 5 mm, Convex Jaws
40-245SS	4" Cuticle Nipper, Stainless, 5 mm, Convex Jaws
40-250CH	4 1/2" Cuticle Nipper, Chrome, 10 mm, Convex Jaws
40-250SS	4 1/2" Cuticle Nipper, Stainless, 10 mm, Convex Jaws
40-251CH	5" Cuticle Nipper, Chrome, 10 mm, Convex Jaws
40-251SS	5" Cuticle Nipper, Stainless, 10 mm, Convex Jaws
40-252	5" Cuticle Nipper, Heavy Pattern, Stainless, Convex Jaws

**TECH-MED FINGER RING CUTTER**


Emergency tool designed to quickly remove ring from finger. Chrome plated handle with stainless steel blade.

8120	Finger Ring Cutter, Standard
8122	Finger Ring Cutter, Heavy Duty

**COURIER® INSTRUMENTS BANDAGE & UTILITY SCISSORS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407035	7 1/2" Bandage & Utility Scissors
---------	-----------------------------------

**MILTEX ALL PURPOSE HOSPITAL SHEARS**

Plastic handles. Autoclavable.

5-1000	7 1/2" Hospital Shears, Serrated
5-998	6 1/2" Hospital Shears, Serrated



**TECH-MED ALL PURPOSE SHEARS**

Easily cuts through the toughest materials. 7 1/2" with stainless steel blade. Autoclavable.

- 8100 All Purpose Shears, Black
- 8100BL All Purpose Shears, Blue
- 8100G All Purpose Shears, Green
- 8100P All Purpose Shears, Pink
- 8100R All Purpose Shears, Red
- 8100W All Purpose Shears, White



**MILTEX EDGE™ SCISSORS SHARPENER**


Quickly puts a sharp edge on dull scissors (not intended for Supercut or serrated edge scissors). Autoclavable.

- 5-5000 Hand Held Scissors Sharpener

**MILTEX UTILITY SCISSORS**

Curved blunt blades. 15mm long. Ribbon type.

- 18-1419 3 3/4" Utility Scissors



**ADC® SCISSORS**

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 300 Lister Bandage Scissors, 4 1/2"
- 3000G Lister Bandage Scissors, Gold, 4 1/2"
- 3006 Lister Bandage Scissors w/clip, 4 1/2"
- 3007 Lister Bandage Scissors w/clip, 5 1/2"
- 301 Lister Bandage Scissors, 5 1/2"
- 3010G Lister Bandage Scissors, Gold, 5 1/2"
- 302 Lister Bandage Scissors, 7 1/2"
- 308 Littauer Stitch Scissors, 5 1/2"
- 3404SB Operating Scissors, Straight, Sharp/Blunt, 5 1/2"
- 3404SS Operating Scissors, Straight, Sharp/Sharp, 5 1/2"
- 3410 Mayo Dissecting Scissors, Straight, 5 1/2"
- 3424 Iris Scissors, Straight, 4 1/2"
- 3425 Iris Scissors, Curved, 4 1/2"
- 3656 Umbilical Cord Scissors, 4"
- 323BD Listerette Scissors, 5 1/2", Burgundy
- 323BK Listerette Scissors, 5 1/2", Black
- 323DG Listerette Scissors, 5 1/2", Dark Green
- 323M Listerette Scissors, 5 1/2", Magenta
- 323N Listerette Scissors, 5 1/2", Navy
- 323RB Listerette Scissors, 5 1/2", Royal Blue
- 323TL Listerette Scissors, 5 1/2", Teal
- 323V Listerette Scissors, 5 1/2", Purple



N407025




N407030

**COURIER® INSTRUMENTS LISTER BANDAGE SCISSORS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407025 5 1/2" SS Lister Bandage Scissors
- N407030 7 1/4" SS Lister Bandage Scissors


Personal Attention  
Local Service



### MILTEX LISTER BANDAGE SCISSORS

Extra fine quality.

5-510	3 1/2" Bandage Scissors
5-512	4 1/2" Bandage Scissors
5-514	5 1/2" Bandage Scissors
5-516	7 1/4" Bandage Scissors
5-550	8" Bandage Scissors w/one large finger ring



### MILTEX NURSE'S SCISSORS

Feature fluoride coated blades to prevent adhesive from sticking to cutting surface. 5 1/2" long. Autoclavable.

5-700	Nurses Scissors, Pink Handles & Guard w/carry case
5-702	Nurses Scissors, Blue Handles & Guard w/carry case

### MILTEX UNIVERSAL SCISSORS

Different colors with black handles for easy identification. Designed to cut through even the toughest material. Stainless steel with one serrated blade. 7 1/2" length. The heavy plastic handles are contoured to provide maximum comfort. Autoclavable to 143°C/290°F. Will cut bandage and gauze, cast material, metal sheets and tough clothing material. Angled design allow for placement of the lower blade under the bandage for safe and easy removal. Available in high quality grade German stainless steel or in Vantage economy grade.

5-1025	High Quality Red Universal Scissors
5-1026	High Quality Orange Universal Scissors
5-1027	High Quality Blue Universal Scissors
5-1028	High Quality Purple Universal Scissors
5-1029	High Quality Green Universal Scissors
5-1030	High Quality Yellow Universal Scissors
V95-1000	Economy Black Universal Scissors
V95-1025	Economy Red Universal Scissors
V95-1026	Economy Orange Universal Scissors
V95-1027	Economy Blue Universal Scissors
V95-1028	Economy Purple Universal Scissors
V95-1029	Economy Green Universal Scissors
V95-1030	Economy Yellow Universal Scissors



### MILTEX KNOWLES BANDAGE SCISSORS


5-560	5 1/2" Bandage Scissors
5-561	Angled, 5 1/2" Bandage Scissors



### COURIER® INSTRUMENTS O.R. SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407000	5 1/2" Straight Sharp/Blunt O.R. Scissors
---------	---



5-12


5-32

**MILTEX STANDARD PATTERN OPERATING SCISSORS**

- 5-12 4 1/2" Operating Scissors, Sharp-Blunt Points, Straight
- 5-14 5" Operating Scissors, Sharp-Blunt Points, Straight
- 5-16 5 1/2" Operating Scissors, Sharp-Blunt Points, Straight
- 5-16TC 5 1/2" Operating Scissors, Sharp-Blunt Points, Straight, Carb-N-Sert
- 5-17 6" Operating Scissors, Sharp-Blunt Points, Straight
- 5-18 6 1/2" Operating Scissors, Sharp-Blunt Points, Straight
- 5-19 7 1/2" Operating Scissors, Sharp-Blunt Points, Straight
- 5-2 4 1/2" Operating Scissors, Sharp-Sharp Points, Straight
- 5-22 4 1/2" Operating Scissors, Blunt-Blunt Points, Straight
- 5-24 5" Operating Scissors, Blunt-Blunt Points, Straight
- 5-26 5 1/2" Operating Scissors, Blunt-Blunt Points, Straight
- 5-26TC 5 1/2" Operating Scissors, Blunt-Blunt Points, Straight, Carb-N-Sert
- 5-27 6" Operating Scissors, Blunt-Blunt Points, Straight
- 5-28 6 1/2" Operating Scissors, Blunt-Blunt Points, Straight
- 5-32 4 1/2" Operating Scissors, Sharp-Sharp Points, Curved
- 5-34 5" Operating Scissors, Sharp-Sharp Points, Curved
- 5-36 5 1/2" Operating Scissors, Sharp-Sharp Points, Curved
- 5-36TC 5 1/2" Operating Scissors, Sharp-Sharp Points, Curved, Carb-N-Sert
- 5-37 6" Operating Scissors, Sharp-Sharp Points, Curved
- 5-38 6 1/2" Operating Scissors, Sharp-Sharp Points, Curved
- 5-4 5" Operating Scissors, Sharp-Sharp Points, Straight
- 5-42 4 1/2" Operating Scissors, Sharp-Blunt Points, Curved
- 5-44 5" Operating Scissors, Sharp-Blunt Points, Curved
- 5-46 5 1/2" Operating Scissors, Sharp-Blunt Points, Curved
- 5-46TC 5 1/2" Operating Scissors, Sharp-Blunt Points, Curved, Carb-N-Sert
- 5-48 6 1/2" Operating Scissors, Sharp-Blunt Points, Curved
- 5-49 7 1/2" Operating Scissors, Sharp-Blunt Points, Curved
- 5-52 4 1/2" Operating Scissors, Blunt-Blunt Points, Curved
- 5-54 5" Operating Scissors, Blunt-Blunt Points, Curved
- 5-56 5 1/2" Operating Scissors, Blunt-Blunt Points, Curved
- 5-56TC 5 1/2" Operating Scissors, Blunt-Blunt Points, Curved, Carb-N-Sert
- 5-58 6 1/2" Operating Scissors, Blunt-Blunt Points, Curved
- 5-6 5 1/2" Operating Scissors, Sharp-Sharp Points, Straight
- 5-6TC 5 1/2" Operating Scissors, Sharp-Sharp Points, Straight, Carb-N-Sert
- 5-7 6" Operating Scissors, Sharp-Sharp Points, Straight
- 5-8 6 1/2" Operating Scissors, Sharp-Sharp Points, Straight
- 5-SC-16 5 1/2" Operating Scissors, Sharp-Blunt Points, Straight, SuperCut Blade
- 5-SC-46 5 1/2" Operating Scissors, Sharp-Blunt Points, Curved, SuperCut Blade

**MILTEX MAYO DISSECTING SCISSORS**

- 5-120 5 1/2" Straight Dissecting Scissors, Standard Beveled Blades
- 5-120TC 5 1/2" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-122 5 1/2" Curved Dissecting Scissors, Standard Beveled Blades
- 5-122TC 5 1/2" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-124 6 3/4" Straight Dissecting Scissors, Standard Beveled Blades
- 5-124TC 6 3/4" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-126 6 3/4" Curved Dissecting Scissors, Standard Beveled Blades
- 5-126TC 6 3/4" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-128 9" Straight Dissecting Scissors, Standard Beveled Blades
- 5-128TC 9" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-130 9" Curved Dissecting Scissors, Standard Beveled Blades
- 5-130TC 9" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
- 5-136 5 1/2" Straight Dissecting Scissors, Rounded Blades
- 5-136TC 5 1/2" Straight Dissecting Scissors, Carb-N-Sert, Rounded Blades
- 5-138 5 1/2" Curved Dissecting Scissors, Rounded Blades
- 5-138TC 5 1/2" Curved Dissecting Scissors, Carb-N-Sert, Rounded Blades
- 5-140 6 3/4" Straight Dissecting Scissors, Rounded Blades
- 5-140TC 6 3/4" Straight Dissecting Scissors, Carb-N-Sert, Rounded Blades
- 5-142 6 3/4" Curved Dissecting Scissors, Rounded Blades
- 5-142TC 6 3/4" Curved Dissecting Scissors, Carb-N-Sert, Rounded Blades
- 5-SC-120 5 1/2" Straight Dissecting Scissors, SuperCut Blade
- 5-SC-122 5 1/2" Curved Dissecting Scissors, SuperCut Blade
- 5-SC-124 6 3/4" Straight Dissecting Scissors, SuperCut Blade
- 5-SC-126 6 3/4" Curved Dissecting Scissors, SuperCut Blade
- 5-SC-128 9" Straight Dissecting Scissors, SuperCut Blade
- 5-SC-130 9" Curved Dissecting Scissors, SuperCut Blade



**MILTEX MAYO OPERATING SCISSORS**

Heavy pattern. Sharp points.

- 5-115 5 1/2" Operating Scissors, Straight
- 5-116 5 1/2" Operating Scissors, Curved

**MILTEX THOREK SCISSORS**

Full curve.

- 5-244 7 1/4" Thorek Scissors
- 5-245 10" Thorek Scissors


**MILTEX QUIMBY SCISSORS**

15mm blades. Strongly curved and pointed.

- 21-558 5" Quimby Scissors
- 21D-559 5" Quimby Scissors w/one serrated blade

**MILTEX SIMS SCISSORS**

- 5-220 8" Sims Scissors, Sharp-Sharp Points, Straight
- 5-222 8" Sims Scissors, Sharp-Blunt Points, Straight
- 5-224 8" Sims Scissors, Blunt-Blunt Points, Straight
- 5-226 8" Sims Scissors, Sharp-Sharp Points, Curved
- 5-228 8" Sims Scissors, Sharp-Blunt Points, Curved
- 5-230 8" Sims Scissors, Blunt-Blunt Points, Curved




**MILTEX WIRE CUTTING SCISSORS WITH NOTCH**

- 9-120 4 3/4" Wire Cutting Scissors w/one serrated blade, Straight
- 9-122 4 3/4" Wire Cutting Scissors w/one serrated blade, Curved
- 9-124 4 3/4" Wire Cutting Scissors w/one serrated blade, Angled to Side
- 9-124TC 4 3/4" Wire Cutting Scissors w/one serrated blade, Carb-N-Sert, Angled To Side
- 9D-121 4 3/4" Wire Cutting Scissors, Straight
- 9D-123 4 3/4" Wire Cutting Scissors, Curved
- 9D-132 4 1/4" Wire Cutting Scissors w/one serrated blade, Straight
- 9D-133 4 1/4" Wire Cutting Scissors, Straight
- 9D-134 4 1/4" Wire Cutting Scissors w/one serrated blade, Curved
- 9D-135 4 1/4" Wire Cutting Scissors, Curved

**MILTEX WIRE TWISTER**

Carb-N-Sert with tungsten carbide inserts.

- 22-500TC 6" Wire Twister
- 22-502TC 8" Wire Twister



5-316

5-318

**MILTEX BROPHY SCISSORS**

- 5-316 5 1/2" Brophy Scissors, Straight, Sharp Points
- 5-318 5 1/2" Brophy Scissors, Curved, Sharp Points
- 5D-326 5 1/2" Brophy Scissors, Curved, One Serrated Blade


**MILTEX DEAVER OPERATING SCISSORS**

- 5-100 5 1/2" Straight Operating Scissors, Sharp-Sharp Points
- 5-102 5 1/2" Straight Operating Scissors, Straight, Sharp-Blunt Points
- 5-104 5 1/2" Straight Operating Scissors, Straight, Blunt-Blunt Points
- 5-106 5 1/2" Curved Operating Scissors, Sharp-Sharp Points
- 5-108 5 1/2" Curved Operating Scissors, Sharp-Blunt Points

**MILTEX STRULLY DISSECTING SCISSORS**

Slightly curved blades with probe tips.


- 5-262 8" Dissecting Scissors




5-180




5-185




5-188



5-284



5-SC-180



5-SC-284

**MILTEX METZENBAUM SCISSORS**

- 5-179 5 1/2" Straight Metzenbaum Scissors, Standard
- 5-179A-TC 5 1/2" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-179B-TC 5 1/2" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
- 5-179TC 5 1/2" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-180 5 1/2" Curved Metzenbaum Scissors, Standard
- 5-180A-TC 5 1/2" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-180B-TC 5 1/2" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
- 5-180TC 5 1/2" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-181 7" Straight Metzenbaum Scissors, Standard
- 5-181A-TC 7" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-181B-TC 7" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
- 5-181TC 7" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-182 7" Curved Metzenbaum Scissors, Standard
- 5-182A-TC 7" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-182B-TC 7" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
- 5-182TC 7" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-183 8" Straight Metzenbaum Scissors, Delicate
- 5-184 8" Curved Metzenbaum Scissors, Standard
- 5-184TC 8" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-185 8" Curved Metzenbaum Scissors, Delicate
- 5-186 9" Straight Metzenbaum Scissors, Standard
- 5-186TC 9" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-187 9" Straight Metzenbaum Scissors, Delicate
- 5-187TC 9" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-188 9" Curved Metzenbaum Scissors, Standard
- 5-188TC 9" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-189 9" Curved Metzenbaum Scissors, Delicate
- 5-189A-TC 10" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-189TC 9" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
- 5-190 11" Straight Metzenbaum Scissors, Standard
- 5-192 11" Curved Metzenbaum Scissors, Standard
- 5-192TC 11" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-193 11" Curved Metzenbaum Scissors, Delicate
- 5-193A-TC 12" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
- 5-194 14 1/2" Straight Metzenbaum Scissors, Standard
- 5-196 14 1/2" Curved Metzenbaum Scissors, Standard
- 5-283 5" Straight Metzenbaum Scissors, Delicate
- 5-284 5" Curved Metzenbaum Scissors, Delicate
- 5-286 6" Curved Metzenbaum Scissors, Delicate
- 5-288 7" Curved Metzenbaum Scissors, Delicate
- 5-SC-180 5 1/2" Curved Metzenbaum Scissors, SuperCut Blade
- 5-SC-182 7" Curved Metzenbaum Scissors, SuperCut Blade
- 5-SC-182TC 7" Curved Metzenbaum Scissors, SuperCut, Carb-N-Sert Blades
- 5-SC-188 9" Curved Metzenbaum Scissors, SuperCut Blade
- 5-SC-283 4 1/2" Straight Metzenbaum Scissors, Delicate, SuperCut Blade
- 5-SC-284 4 1/2" Curved Metzenbaum Scissors, Delicate, SuperCut Blade

**MILTEX POTTS-SMITH SCISSORS**

- 5-246 7 1/2" Potts-Smith Scissors, Angled On Side 25°, 20mm Blades
- 5-247 7 1/2" Potts-Smith Scissors, Angled On Side 45°, 20mm Blades
- 5-248 7 1/2" Potts-Smith Scissors, Angled On Side 60°, 20mm Blades
- 5-249A 7" Potts-Smith Scissors, Delicate, Angled On Side 25°, 13mm Blades
- 5-249B 7" Potts-Smith Scissors, Delicate, Angled On Side 45°, 13mm Blades
- 5-249C 7" Potts-Smith Scissors, Delicate, Angled On Side 60°, 13mm Blades

**MILTEX REYNOLDS DISSECTING SCISSORS**

Curved. Tenotomy type dissecting tips.

- 5-172 6" Dissecting Scissors
- 5-174 6" Dissecting Scissors, One Serrated Blade
- 5-176 7" Dissecting Scissors
- 5-178 7" Dissecting Scissors, One Serrated Blade


**MILTEX RICHTER DISSECTING SCISSORS**

Angular blades with one sharp and one blunt point.

- 34-92 5" Dissecting Scissors

**MILTEX PLASTIC SURGERY SCISSORS**


- 5-270 4 3/4" Straight Surgery Scissors, Sharp-Sharp Points
- 5-271 4 3/4" Straight Surgery Scissors, Sharp-Sharp Points, Serrated Blade
- 5-272 4 3/4" Straight Surgery Scissors, Sharp-Blunt Points
- 5-273 4 3/4" Straight Surgery Scissors, Sharp-Blunt Points, Serrated Blade
- 5-274 4 3/4" Straight Surgery Scissors, Blunt-Blunt Points
- 5-275 4 3/4" Straight Surgery Scissors, Blunt-Blunt Points, Serrated Blade
- 5-276 4 3/4" Curved Surgery Scissors, Sharp-Sharp Points
- 5-277 4 3/4" Curved Surgery Scissors, Sharp-Sharp Points, Serrated Blade
- 5-278 4 3/4" Curved Surgery Scissors, Sharp-Blunt Points
- 5-279 4 3/4" Curved Surgery Scissors, Sharp-Blunt Points, Serrated Blade
- 5-280 4 3/4" Curved Surgery Scissors, Blunt-Blunt Points
- 5-281 4 3/4" Curved Surgery Scissors, Blunt-Blunt Points, Serrated Blade
- 5-282 4 3/4" Surgery Scissors, Angled On Side, Sharp-Sharp Points
- 5-282A 4 3/4" Surgery Scissors, Angled On Side, Sharp-Sharp Points, Serrated Blade



**COURIER® INSTRUMENTS LITTAUER STITCH SCISSORS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407045 4 1/2" Littauer Stitch Scissors



**MILTEX LITTAUER JUNIOR STITCH SCISSORS**

- 9-102 4 1/2" Stitch Scissors, Slender Pattern




**MILTEX LITTAUER STITCH SCISSORS**

- 9-104 5 1/2" Stitch Scissors

**MILTEX NORTHBENT STITCH SCISSORS**

Curved. Light pattern.

- 9-103 4 3/4" Stitch Scissors



**MILTEX O'BRIEN STITCH SCISSORS**

- 9-110 3 3/4" Stitch Scissors, Angled, Sharp Points

**MILTEX RAGNELL DISSECTING SCISSORS**

Flat tip and curved.

- 5-290 5" Dissecting Scissors
- 5-291 7" Dissecting Scissors

**MILTEX SHORTBENT STITCH SCISSORS**


9-101 3 1/2" Stitch Scissors, Curved, Delicate



**COURIER® INSTRUMENTS SPENCER STITCH SCISSORS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407040 3 1/2" Spencer Stitch Scissors



**MILTEX SPENCER STITCH SCISSORS**

9-100 3 1/2" Stitch Scissors



**MILTEX STITCH SCISSORS**

9-108 4 1/2" Stitch Scissors, Delicate, Angular Blade



18-1486

**MILTEX WESTCOTT STITCH SCISSORS**

Curved with sharp points.

18-1486 4 1/2" Stitch Scissors



N407100



N407105

**COURIER® INSTRUMENTS DRESSING FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407100 5" Dressing Forceps
- N407105 5 1/2" Dressing Forceps



**MILTEX DRESSING FORCEPS**

Serrated tips.

- 6-10 6" Dressing Forceps, Standard w/serrated handles
- 6-100 5" Dressing Forceps, Curved Narrow Tips
- 6-12 7" Dressing Forceps, Standard w/serrated handles
- 6-14 8" Dressing Forceps, Standard w/serrated handles
- 6-16 10" Dressing Forceps, Standard w/serrated handles
- 6-18 12" Dressing Forceps, Standard w/serrated handles
- 6-2 4" Dressing Forceps, Standard w/serrated handles
- 6-24 4 1/2" Dressing Forceps, Delicate w/fluted handles
- 6-26 5" Dressing Forceps, Delicate w/fluted handles
- 6-28 5 1/2" Dressing Forceps, Delicate w/fluted handles
- 6-28XL 5 1/2" Dressing Forceps, Delicate Lightweight w/fenestrated handles
- 6-30 6" Dressing Forceps, Delicate w/fluted handles
- 6-30TC 6" Dressing Forceps, Delicate w/fluted handles, Carb-N-Sert
- 6-4 4 1/2" Dressing Forceps, Standard w/serrated handles
- 6-6 5" Dressing Forceps, Standard w/serrated handles
- 6-8 5 1/2" Dressing Forceps, Standard w/serrated handles
- 6-8XL 5 1/2" Dressing Forceps, Lightweight w/fenestrated handles



**COURIER® INSTRUMENTS ADSON DRESSING FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407120 4 3/4" Adson Dressing Forceps



**MILTEX ADSON DRESSING FORCEPS**

Serrated.

- 17-2510 4 3/4" Dressing Forceps, Micro, Delicate, 0.5mm Wide Tips
- 6-118 4 3/4" Dressing Forceps, Delicate
- 6-118A 4 3/4" Dressing Forceps, Delicate, Angled
- 6-118TC 4 3/4" Dressing Forceps, Carb-N-Sert
- 6-118XL 4 3/4" Dressing Forceps, Fenestrated Handle



**COURIER® INSTRUMENTS BROWN-ADSON FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407130 4 3/4" Brown-Adson Forceps, 7x7 Teeth

**MILTEX ADSON BAYONET DRESSING FORCEPS**

- 6-198 8 1/4" Dressing Forceps, Serrated
- 6-198TC 8 1/2" Dressing Forceps, Carb-N-Sert

**MILTEX POTTS-SMITH DRESSING FORCEPS**

Serrated.

- 6-154 7" Dressing Forceps
- 6-154TC 7" Dressing Forceps, Carb-N-Sert
- 6-156 8 1/4" Dressing Forceps
- 6-156TC 8 1/4" Dressing Forceps, Carb-N-Sert
- 6-158 9 1/2" Dressing Forceps
- 6-158TC 9 3/4" Dressing Forceps, Carb-N-Sert
- 6-159 12" Dressing Forceps



**MILTEX HUDSON (EWALD) FORCEPS**

- 6-112 4 3/4" Hudson Forceps, Dressing, Serrated
- 6-114 4 3/4" Hudson Forceps, Tissue, 1x2 Teeth



**MILTEX CUSHING DRESSING FORCEPS**

Serrated with scraper end.

- 6-166 7" Dressing Forceps w/Gutch Handle
- 6-166TC 7" Dressing Forceps, Carb-N-Sert
- 6-188 7" Dressing Forceps
- 6-190 7 1/4" Dressing Forceps, Bayonet



**MILTEX SEMKIN FORCEPS**

- 6-104 5" Semkin Forceps, Dressing, Serrated Tips
- 6-106 5" Semkin Forceps, Tissue, 1x2 Teeth
- 6-107 5" Semkin Forceps, Tissue, Curved, 1x2 Teeth
- 6-108 5" Semkin Forceps, Tissue, 2x3 Teeth
- 6-109 6" Semkin Forceps, Dressing, Serrated Tips
- 6-109A 6" Semkin Forceps, Tissue, 1x2 Teeth



6-308



6-310

**MILTEX CARMALT SPLINTER FORCEPS**

- 6-308 4 1/4" Splinter Forceps, Straight
- 6-310 4 1/4" Splinter Forceps, Curved



**MILTEX FEILCHENFELD SPLINTER FORCEPS**

- 6-316 3" Splinter Forceps
- 6-318 4 1/2" Splinter Forceps

**MILTEX PEET SPLINTER FORCEPS**

Double curve.

- 6-350 4 3/4" Splinter Forceps

**MILTEX PHYSICIANS' SPLINTER FORCEPS**

- 6-330 5" Splinter Forceps

**MILTEX SPLINTER LIBERATOR**

Trocar point in reversible handle.

- 6-400 5" Splinter Liberator




N407145


**COURIER® INSTRUMENTS SPLINTER FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407135 3 1/2" Splinter Forceps
- N407140 4" Splinter Forceps
- N407145 4 1/2" Splinter Forceps



6-300



6-302

**MILTEX PLAIN SPLINTER FORCEPS**


Fine points.

- 6-300 3 1/2" Splinter Forceps
- 6-302 4" Splinter Forceps
- 6-304 4 1/2" Splinter Forceps

**MILTEX SWISS JEWELER STYLE FORCEPS**

Non-magnetic stainless steel.

- 17-301 Style 1 Jeweler Style Forceps, Fine, 4 3/4"
- 17-302 Style 2 Jeweler Style Forceps, Strong, Fine, 4 3/4"
- 17-303 Style 3 Jeweler Style Forceps, Narrow, Fine, 4 3/4"
- 17-303C Style 3C Jeweler Style Forceps, Narrow, Fine, 4 3/8"
- 17-303X Style 3 Jeweler Style Forceps, Micro-Fine, 4 3/4"
- 17-304 Style 4 Jeweler Style Forceps, Fine, 4 3/8"
- 17-305 Style 5 Jeweler Style Forceps, Super Fine, 4 3/8"
- 17-305X Style 5 Jeweler Style Forceps, Micro-Fine, 4 3/8"
- 17-306 Style 6 Jeweler Style Forceps, Pickup, 4 1/2"
- 17-307 Style 7 Jeweler Style Forceps, Fine Curved, 4 1/2"
- 17-307X Style 7 Jeweler Style Forceps, Micro-Fine, Curved, 4 1/2"




**MILTEX VIRTUS SPLINTER FORCEPS**

- 6-336 6" Splinter Forceps, Straight
- 6-337 6" Splinter Forceps, 45° Angle

**MILTEX WALTHER SPLINTER FORCEPS**

- 6-312 4 1/4" Splinter Forceps, Straight
- 6-314 4 1/4" Splinter Forceps, Curved



**MILTEX WILLIAMS SPLINTER FORCEPS**


- 6-338 5 1/2" Splinter Forceps, Straight
- 6-340 5 1/2" Splinter Forceps, 45° Angle
- 6-342 5 1/2" Splinter Forceps, Right Angle



**COURIER® INSTRUMENTS ADSON TISSUE FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407125 4 3/4" Adson Tissue Forceps, 1x2 Teeth



17-2500 6-120A 6-121

**MILTEX ADSON TISSUE FORCEPS**

- 17-2500 4 3/4" Tissue Forceps, 1x2 Teeth, 0.5 mm Wide Tips, Delicate
- 6-120 4 3/4" Tissue Forceps, 1x2 Teeth, Delicate
- 6-120A 4 3/4" Tissue Forceps, 1x2 Teeth, Angled, Delicate
- 6-120XL 4 3/4" Tissue Forceps, 1x2 Teeth, Fenestrated Handles, Lightweight
- 6-121 4 3/4" Tissue Forceps, 1x2 Teeth, Cross Serrated Jaws
- 6-122 4 3/4" Tissue Forceps, 2x3 Teeth, Delicate
- 6-123 4 3/4" Tissue Forceps, 1x2 Teeth w/suturing platform
- 6-123TC 4 3/4" Tissue Forceps, 1x2 Teeth, Cross Serrated Typing Platform, Carb-N-Sert, Lightweight
- 6-123XL 4 3/4" Tissue Forceps, 1x2 Teeth, Fenestrated Handles w/platform, Lightweight




N407110 N407115

**COURIER® INSTRUMENTS TISSUE FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407110 5" Tissue Forceps 1x2 Teeth
- N407115 5 1/2" Tissue Forceps 1x2 Teeth



6-124 6-124XL

**MILTEX BROWN-ADSON TISSUE FORCEPS**

7 x 7 side grasping teeth.

- 6-124 4 3/4" Tissue Forceps
- 6-124TC 4 3/4" Tissue Forceps, Cross Serrated Typing Platform, Carb-N-Sert
- 6-124XL 4 3/4" Tissue Forceps, Lightweight, Fenestrated Handles
- 6-126 4 3/4" Tissue Forceps, Angular



**MILTEX TISSUE FORCEPS**

- 6-40 4 1/2" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-42 5" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-44 5 1/2" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-44XL 5 1/2" Tissue Forceps, 1x2 Teeth, Standard Lightweight w/fenestrated handles
- 6-46 6" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-47 7" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-48 8" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-50 10" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-52 12" Tissue Forceps, 1x2 Teeth, Standard w/serrated handles
- 6-60 4 1/2" Tissue Forceps, 2x3 Teeth, Standard w/serrated handles
- 6-62 5" Tissue Forceps, 2x3 Teeth, Standard w/serrated handles
- 6-64 5 1/2" Tissue Forceps, 2x3 Teeth, Standard w/serrated handles
- 6-66 6" Tissue Forceps, 2x3 Teeth, Standard w/serrated handles
- 6-72 5" Tissue Forceps, 3x4 Teeth, Standard w/serrated handles
- 6-74 5 1/2" Tissue Forceps, 3x4 Teeth, Standard w/serrated handles
- 6-76 6" Tissue Forceps, 3x4 Teeth, Standard w/serrated handles
- 6-82 5" Tissue Forceps, 4x5 Teeth, Standard w/serrated handles
- 6-84 5 1/2" Tissue Forceps, 4x5 Teeth, Standard w/serrated handles
- 6-86 6" Tissue Forceps, 4x5 Teeth, Standard w/serrated handles
- 6-90 4 1/2" Tissue Forceps, 1x2 Teeth, Delicate w/fluted handles
- 6-92 5" Tissue Forceps, 1x2 Teeth, Delicate w/fluted handles
- 6-94 5 1/2" Tissue Forceps, 1x2 Teeth, Delicate w/fluted handles
- 6-96 6" Tissue Forceps, 1x2 Teeth, Delicate w/fluted handles



**MILTEX POTTS-SMITH TISSUE FORCEPS**

1 x 2 teeth and serrated tips.

- 6-160 7" Tissue Forceps
- 6-162 8 1/4" Tissue Forceps
- 6-164 9 1/2" Tissue Forceps
- 6-164A 12" Tissue Forceps



**MILTEX BONNEY TISSUE FORCEPS**

- 6-148 7" Tissue Forceps, 1 x 2 Teeth, Serrated Tips
- 6-150 7" Tissue Forceps, 2 x 3 Teeth

**MILTEX MAYO RUSSIAN TISSUE FORCEPS**

Fenestrated jaws and slotted grooved handles.

- 6-210 9" Mayo Russian Tissue Forceps

**MILTEX RUSSIAN TISSUE FORCEPS**

- 6-142 6" Tissue Forceps
- 6-144 8" Tissue Forceps
- 6-145 10" Tissue Forceps

**MILTEX WANGENSTEEN TISSUE FORCEPS**


Rounded jaws with fine cross serrations.

- 16-80 9" Wangensteen Tissue Forceps

**MILTEX STONE TISSUE FORCEPS**


4 x 5 teeth.

- 6-132 6" Tissue Forceps



**MILTEX ALLIS TISSUE FORCEPS**


- 16-10 7 1/4" Allis Tissue Forceps, 5 x 6 Teeth
- 16-12 10" Allis Tissue Forceps, 5 x 6 Teeth
- 16-20 6" Allis Tissue Forceps, 4 x 5 Teeth, Angular Jaws
- 16-22 6" Allis Tissue Forceps, 5 x 6 Teeth, Angular Jaws
- 16-24 7 1/4" Allis Tissue Forceps, 5 x 6 Teeth, Angular Jaws
- 16-25 Allis-NC Tissue Forceps, 6" w/double row of non-traumatic teeth
- 16-26 Allis-NC Tissue Forceps, 7 1/4" w/double row of non-traumatic teeth
- 16-3 Baby Allis Tissue Forceps, 5", 4 x 5 Teeth, Extra Delicate Jaws 2.5mm Wide
- 16-30 Allis-Adair Tissue Forceps, 6 1/2", 9 x 10 Teeth
- 16-31 Boys-Allis Tissue Forceps, 6", 5 x 6 Teeth, Heavy Pattern
- 16-33 Judd-Allis Tissue Forceps, 6", 3 x 4 Teeth, Delicate
- 16-34 Judd-Allis Tissue Forceps, 7 1/2", 3 x 4 Teeth, Delicate
- 16-38 Thoms-Allis Tissue Forceps, 8", 6 x 7 Teeth
- 16-4 Baby Allis Tissue Forceps, 5 1/2", 4 x 5 Teeth, Delicate Jaws, 4mm Wide
- 16-6 6" Allis Tissue Forceps, 4 x 5 Teeth
- 16-8 6" Allis Tissue Forceps, 5 x 6 Teeth



## ADC® SPECIALTY FORCEPS


Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 315 Magill Catheter Forceps, Child, 8"
- 316 Magill Catheter Forceps, Adult, 9 3/4"
- 317 Young Tongue Seizing Forceps, 6 1/4"



## MILTEX MAGILL ENDOTRACHEAL CATHETER INTRODUCING FORCEPS

- 2-70 9" Catheter Introducing Forceps, Adult
- 2-72 7" Catheter Introducing Forceps, Child, Extra Delicate



## MILTEX UTILITY STERILIZER FORCEPS


Made of stainless steel.

- 3-20 8" Straight Sterilizer Forceps
- 3-22 8" Curved Sterilizer Forceps
- 3-24 11" Straight Sterilizer Forceps
- 3-26 11" Curved Sterilizer Forceps

## MILTEX WILL HOSPITAL UTILITY FORCEPS

Cross serrated tying platform. 1.2mm wide.


- 18-970 4" Utility Forceps



## RUSCH MAGILL INTUBATING FORCEPS

Made of lightweight, polished stainless steel. Non-sterile.


- 008500050 Closed Jaw Forceps, Child
- 008501300 Open Jaw Forceps, Adult
- 008502100 Open Jaw Forceps, Child
- 008505300 Closed Jaw Forceps, Adult



## ADC® HEMOSTATIC FORCEPS

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 310 Kelly Forceps, Straight, 5 1/2"
- 3101 Crile Forceps, Straight, 5 1/2"
- 3102 Crile Forceps, Curved, 5 1/2"
- 311 Kelly Forceps, Curved, 5 1/2"
- 312 Kelly Forceps, Straight, 6 1/4"
- 313 Kelly Forceps, Curved, 6 1/4"
- 314 Halstead Mosquito Forceps, Straight, 5"
- 3141 Halstead Mosquito Forceps, Curved, 5"



N407210




N407215

## COURIER® INSTRUMENTS KELLY FORCEPS


Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407210 5 1/2" Straight Kelly Forceps
- N407215 5 1/2" Curved Kelly Forceps




**MILTEX KELLY FORCEPS**

- 7-36 5 1/2" Kelly Forceps, Straight
- 7-38 5 1/2" Kelly Forceps, Curved



**MILTEX CRILE FORCEPS**

- 7-42 5 1/2" Crile Forceps, Straight
- 7-44 5 1/2" Crile Forceps, Curved
- 7-46 6 1/4" Crile Forceps, Straight
- 7-48 6 1/4" Crile Forceps, Curved
- 7-50 5 1/2", Baby Crile Forceps, Straight, Extra Delicate
- 7-52 5 1/2", Baby Crile Forceps, Curved, Extra Delicate



7-11  
7-12

**MILTEX PACKER MOSQUITO FORCEPS**

Flat jaws.

- 7-11 5" Mosquito Forceps, Straight
- 7-12 5" Mosquito Forceps, Curved



N407200




N407205

**COURIER® INSTRUMENTS HALSTED MOSQUITO FORCEPS**


Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407200 5" Straight Halsted Mosquito Forceps
- N407205 5" Curved Halsted Mosquito Forceps




**MILTEX HALSTEAD MOSQUITO FORCEPS**

- 18-1934 5" Halstead Mosquito Forceps, Straight, Non-Magnetic
- 18-1936 5" Halstead Mosquito Forceps, Curved, Non-Magnetic
- 7-10 5" Halstead Mosquito Forceps, Curved, Extra Delicate
- 7-14 5" Halstead Mosquito Forceps, Straight, 1 x 2 Teeth
- 7-16 5" Halstead Mosquito Forceps, Curved, 1 x 2 Teeth
- 7-2 5" Halstead Mosquito Forceps, Straight
- 7-20 5" Halstead Mosquito Forceps, Extra Delicate, Jaws 45° Angled On Flat
- 7-4 5" Halstead Mosquito Forceps, Curved
- 7-8 5" Halstead Mosquito Forceps, Straight, Extra Delicate




## MILTEX HARTMAN & HARTMAN-NOYES ALLIGATOR FORCEPS

- 19-434 Hartman Alligator Forceps, 3 1/4" Shaft, Serrated Jaws
- 20-140 Hartman-Noyes Alligator Forceps, 5 1/2" Shaft, Serrated Jaws
- 20-142 6 1/2" Shaft Alligator Forceps, Serrated Jaws
- 20-144 5 1/2" Shaft Alligator Forceps, 1x2 Teeth



## MILTEX HARTMAN MOSQUITO FORCEPS

- 7-24 3 1/2" Hartman Mosquito Forceps, Straight
- 7-24A 3 1/2" Hartman Mosquito Forceps, Straight Jaws, Shanks Angled On Flat
- 7-26 3 1/2" Hartman Mosquito Forceps, Curved



## MILTEX KOCHER HEMOSTATIC FORCEPS

- 16-120 10 1/2" Hemostatic Forceps, Straight, Longitudinal Serrations
- 16-122 10 1/2" Hemostatic Forceps, Curved, Longitudinal Serrations
- 7-66 5 1/2" Hemostatic Forceps, Straight, 1x2 Teeth
- 7-68 5 1/2" Hemostatic Forceps, Curved, 1x2 Teeth

## MILTEX BABY OCHSNER FORCEPS

1 x 2 teeth. Extra delicate.

- 7-70 5 1/2" Baby Ochsner Forceps, Straight
- 7-72 5 1/2" Baby Ochsner Forceps, Curved

## MILTEX BABY PEAN FORCEPS

Extra delicate.

- 7-110 5 1/2" Baby Pean Forceps, Straight
- 7-112 5 1/2" Baby Pean Forceps, Curved


## MILTEX PEAN FORCEPS

Curved with longitudinal serrations.

- 30-1805 9" Pean Forceps
- 30-1815 10 1/4" Pean Forceps

## MILTEX PENNINGTON HEMOSTATIC FORCEPS

- 7-220 6" Hemostatic Forceps




## MILTEX RANKIN FORCEPS

- 7-92 6 1/4" Rankin Forceps, Straight
- 7-94 6 1/4" Rankin Forceps, Curved

## MILTEX ROCHESTER-CARMALT FORCEPS

- 7-170 6 1/4" Rochester-Carmalt Forceps, Straight
- 7-172 6 1/4" Rochester-Carmalt Forceps, Curved
- 7-180 8" Rochester-Carmalt Forceps, Straight
- 7-182 8" Rochester-Carmalt Forceps, Curved



## MILTEX ROCHESTER-OSHSNER FORCEPS

- 7-150 6 1/4" Straight Rochester-Ochsner Forceps
- 7-160 6 1/4" Curved Rochester-Ochsner Forceps
- 7-166 9" Curved Rochester-Ochsner Forceps
- 7-152 7 1/4" Straight Rochester-Ochsner Forceps
- 7-154 8" Straight Rochester-Ochsner Forceps
- 7-156 9" Straight Rochester-Ochsner Forceps
- 7-157 10" Straight Rochester-Ochsner Forceps
- 7-162 7 1/4" Curved Rochester-Ochsner Forceps
- 7-164 8" Curved Rochester-Ochsner Forceps
- 7-167 10" Curved Rochester-Ochsner Forceps




## MILTEX ROCHESTER-PEAN FORCEPS

- 7-118 5 1/2" Straight Rochester-Pean Forceps
- 7-120 6 1/4" Straight Rochester-Pean Forceps
- 7-122 7 1/4" Straight Rochester-Pean Forceps
- 7-124 8" Straight Rochester-Pean Forceps
- 7-126 9" Straight Rochester-Pean Forceps
- 7-128 10 1/4" Straight Rochester-Pean Forceps
- 7-130 12" Straight Rochester-Pean Forceps
- 7-136 5 1/2" Curved Rochester-Pean Forceps
- 7-138 6 1/4" Curved Rochester-Pean Forceps
- 7-140 7 1/4" Curved Rochester-Pean Forceps
- 7-142 8" Curved Rochester-Pean Forceps
- 7-144 9" Curved Rochester-Pean Forceps
- 7-146 10 1/4" Curved Rochester-Pean Forceps
- 7-148 12" Curved Rochester-Pean Forceps


## MILTEX DANDY FORCEPS

- 7-102 5 1/2" Dandy Forceps, Curved Sideways




## MILTEX HALSEY NEEDLE HOLDERS


- 8-8 5" Halsey Needle Holder, Smooth, Regular Jaws
- 8-8A-TC 5" Halsey Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-8CG 5" Halsey Needle Holder, Ceram-A-Grip, Smooth
- 8-8TC 5" Halsey Needle Holder, Smooth, Carb-N-Sert T.C. Jaws
- 8D-36 5" Halsey Needle Holder, Serrated Jaws



N407300



N407305



N407310



N407315




N407320

## COURIER® INSTRUMENTS NEEDLE HOLDERS


Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407300 Webster 5" Needle Holder
- N407305 Halsey 5" Needle Holder
- N407310 Baumgartner 5 1/2" Needle Holder
- N407315 Mayo-Hegar 5" Needle Holder
- N407320 Mayo-Hegar 6" Needle Holder



## MILTEX CRILE-WOOD NEEDLE HOLDERS

- 8-50 6" Crile-Wood Needle Holder, Serrated Jaws
- 8-50A-TC 6" Crile-Wood Needle Holder, Smooth Jaws, Carb-N-Sert T.C. Jaws
- 8-50CG 6" Crile-Wood Needle Holder, Ceram-A-Grip
- 8-50TC 6" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-51 7" Crile-Wood Needle Holder, Serrated Jaws
- 8-51CG 7" Crile-Wood Needle Holder, Ceram-A-Grip
- 8-51TC 7" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-52 8" Crile-Wood Needle Holder, Serrated Jaws
- 8-52CG 8" Crile-Wood Needle Holder, Ceram-A-Grip
- 8-52TC 8" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws, Carb-N-Sert T.C. Jaws
- 8-53 9" Crile-Wood Needle Holder, Serrated Jaws
- 8-53TC 9" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-55TC 10 1/2" Crile-Wood Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-57TC 12" Crile-Wood Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-750TC 6" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
- 8-751TC 7" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
- 8-752TC 8" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
- 8-753TC 9" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws



## MILTEX BAUMGARTNER NEEDLE HOLDERS

- 8-40 5 1/2" Baumgartner Needle Holder, Regular Jaws
- 8-40CG 5 1/2" Baumgartner Needle Holder, Ceram-A-Grip
- 8-40TC 5 1/2" Baumgartner Needle Holder, Carb-N-Sert T.C. Jaws

## MILTEX METZENBAUM NEEDLE HOLDER

Fenestrated jaws.

- 8-74 7 1/4" Metzenbaum Needle Holder, Straight
- 8-76 7 1/4" Metzenbaum Needle Holder, Curved

## MILTEX WANGENSTEEN NEEDLE HOLDER

- 8-134 10 3/4" Wangensteen Needle Holder, Long Jaws
- 8-134TC 10 3/4" Wangensteen Needle Holder, Carb-N-Sert, Serrated
- 8-136 10 1/2" Wangensteen Needle Holder, Short Jaws




## MILTEX MAYO-HEGAR NEEDLE HOLDERS

- 8-130 10 1/2" Mayo-Hegar Needle Holder
- 8-130TC 10 1/2" Mayo-Hegar Needle Holder, Carb-N-Sert T.C. Jaws
- 8-42 5" Mayo-Hegar Needle Holder, Serrated Jaws
- 8-42CG 5" Mayo-Hegar Needle Holder, Ceram-A-Grip
- 8-42TC 5" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-44 6" Mayo-Hegar Needle Holder, Serrated Jaws
- 8-44CG 6" Mayo-Hegar Needle Holder, Ceram-A-Grip
- 8-44TC 6" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-46 7" Mayo-Hegar Needle Holder, Serrated Jaws
- 8-46CG 7" Mayo-Hegar Needle Holder, Ceram-A-Grip
- 8-46TC 7" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-48 8" Mayo-Hegar Needle Holder, Serrated Jaws
- 8-48CG 8" Mayo-Hegar Needle Holder, Ceram-A-Grip
- 8-48TC 8" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
- 8-49TC 12" Mayo-Hegar Needle Holder, Carb-N-Sert T.C. Jaws
- 8-744TC 6" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
- 8-746TC 7" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
- 8-748TC 8" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws

## MILTEX STRATTE NEEDLE HOLDER

Double curved jaws and bent shanks.

- 8-102 9" Stratte Needle Holder
- 8-102TC 9" Stratte Needle Holder, Carb-N-Sert, Serrated



## MILTEX DERF NEEDLE HOLDERS

- 8-4 4 3/4" Derf Needle Holder, Regular Jaws
- 8-4CG 4 3/4", Ceram-A-Grip
- 8-4TC 4 3/4" Derf Needle Holder, Carb-N-Sert T.C. Jaws

## MILTEX HEANEY NEEDLE HOLDER

Curved jaws.

- 8-100 8 1/2" Heaney Needle Holder
- 8-100TC 8 1/2" Heaney Needle Holder, Carb-N-Sert, Serrated



## MILTEX WEBSTER NEEDLE HOLDERS

- 8-6 5" Webster Needle Holder, Smooth
- 8-6CG 5" Webster Needle Holder, Ceram-A-Grip
- 8-6TC 5" Webster Needle Holder, Smooth, Carb-N-Sert
- 8-7 4 3/4" Webster Needle Holder, Smooth, Very Delicate
- 8-706TC 5" Webster Needle Holder w/quick-release ratchet, Carb-N-Sert, Smooth
- 8-7A-TC 4 3/4" Webster Needle Holder, Carb-N-Sert, Extra Delicate, Serrated
- 8-7CG 4 3/4" Webster Needle Holder, Ceram-A-Grip
- 8-7TC 4 3/4" Webster Needle Holder, Smooth, Very Delicate, Carb-N-Sert



**MILTEX OLSEN-HEGAR NEEDLE HOLDER-SCISSORS**

- 8-15 5 1/2" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-15CG 5 1/2" Olsen-Hegar Needle Holder, Serrated, Ceram-A-Grip
- 8-15TC 5 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-16 6 1/2" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-16CG 6 1/2" Olsen-Hegar Needle Holder, Serrated, Ceram-A-Grip
- 8-16TC 6 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-17 7 1/4" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-17TC 7 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-14 4 3/4" Olsen-Hegar Needle Holder, Serrated, Extra Delicate, Regular Jaws
- 8-14A 4 3/4" Olsen-Hegar Needle Holder, Smooth, Extra Delicate, Regular Jaws
- 8-14A-TC 4 3/4" Olsen-Hegar Needle Holder, Smooth, Extra Delicate, Carb-N-Sert T.C. Jaws
- 8-14TC 4 3/4" Olsen-Hegar Needle Holder, Serrated, Extra Delicate, Carb-N-Sert T.C. Jaws

**MILTEX CASTROVIEJO NEEDLE HOLDER**


- 18-1816TC 4 1/8" Castroviejo Needle Holder, Carb-N-Sert, Round Handles w/lock, Delicate, Smooth, Straight
- 18-1820 5 1/2" Castroviejo Needle Holder, Straight
- 18-1824 5 1/2" Castroviejo Needle Holder, Curved
- 18-1828 5 1/2" Castroviejo Needle Holder w/lock, Straight
- 18-1828A-TC 5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Serrated, Straight
- 18-1828TC 5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Smooth, Straight
- 18-1831 5 1/2" Castroviejo Needle Holder, Extra Delicate w/lock, Straight
- 18-1832 5 1/2" Castroviejo Needle Holder w/lock, Curved
- 18-1832A-TC 5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Serrated, Curved
- 18-1832TC 5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Smooth, Curved
- 18-1833 5 1/2" Castroviejo Needle Holder, Extra Delicate w/lock, Curved



**COURIER® INSTRUMENTS BACKHAUS TOWEL CLAMPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407220 3 1/2" Backhaus Towel Clamp



**MILTEX BACKHAUS TOWEL CLAMP**

- 7-504 3 1/2" Towel Clamp
- 7-506 5 1/4" Towel Clamp

**MILTEX JONES TOWEL CLAMP**

- 7-500 3 1/2" Towel Clamp, Standard Size
- 7-501 2 1/4" Towel Clamp, Lightweight

**MILTEX ROEDER TOWEL CLAMP**

- 7-510 5 1/4" Towel Clamp w/ball stops



**MILTEX ALM RETRACTOR**


- 11-12 2 3/4" Alm Retractor, Spread 2 1/4", Prongs 1/8" Deep
- 11-13 3 7/8" Alm Retractor, Spread 3", Prongs 1/2" Deep

**MILTEX BALFOUR RETRACTORS**

- 11-420 7" Spread Abdominal Retractor w/fenestrated side blades 2 1/2"D & Solid Center Blade 3" x 1 3/4"D
- 11-421 10" Spread Abdominal Retractor w/fenestrated side blades 2 1/2"D & solid center blade 3" x 1 3/4"D
- 11-424 7" Spread Abdominal Retractor w/solid side blades 2 3/4"D & solid center blade 3" x 1 3/4"D
- 11-428 7" Spread Abdominal Retractor w/fenestrated side blades 4"D solid center blade 2 3/4" x 4"D
- 11-430 4" Spread Pediatric Abdominal Retractor w/fenestrated side blades & solid center blade 1"D, Light Pattern
- 11-432 7" Spread Abdominal Retractor w/ratchet bar, 2 Pair Interchangeable Fenestrated Side Blades 2 1/4" & 3 1/2"D, 2 Solid Center Blades 2 5/8" x 1 3/4" & 3 1/4" x 2 3/4"
- 11-440 Center Blade, 3" x 1 3/4"D
- 11-442 Center Blade, 2 3/4" x 4"D
- 11-444 Center Blade, 2 5/8" x 1 3/4"D
- 11-446 Center Blade, 3 1/4" x 2 3/4"D
- 11-448 Center Blade, 2" x 5"D, Malleable

**MILTEX DEAVER RETRACTORS**

11-200	1" x 9" Deaver Retractor
11-202	1" x 12" Deaver Retractor
11-204	1" x 13" Deaver Retractor
11-206	1 1/2" x 12" Deaver Retractor
11-208	2" x 12" Deaver Retractor
11-210	3" x 12" Deaver Retractor
11-218	5/8" x 8 1/2" Baby Deaver Retractor
11-219	3/4" x 7" Baby Deaver Retractor
11-220	7/8" x 8" Baby Deaver Retractor
11-221	1" x 12" Deaver Retractor w/hollow grip handle
11-222	1" x 13" Deaver Retractor w/hollow grip handle
11-223	1 1/2" x 12" Deaver Retractor w/hollow grip handle
11-224	2" x 12" Deaver Retractor w/hollow grip handle
11-224A	3" x 12" Deaver Retractor w/hollow grip handle
11-225	4" x 12" Deaver Retractor w/hollow grip handle

**MILTEX GELPI RETRACTORS**

11-365	5 1/2" Pediatric Gelpi Retractor w/grip lock
11-366	6 3/4" Gelpi Retractor w/grip lock
11-368	6 3/4" Gelpi Retractor w/ball stops & grip lock

**MILTEX MEYERDING RETRACTORS**

11-240	9 1/2" Retractor, Large Blade, 3 1/2" x 2"
11-242	9 1/2" Retractor, Medium Blade, 3" x 1"
11-244	9 1/2" Retractor, Small Blade, 2" x 5/8"
11-650	Self-Retaining, 1" x 2 3/8" Blades
11-652	Self-Retaining, 1" x 2 1/8" Blades
11-654	Self-Retaining, 1" x 1 3/4" Blades
11-656	Self-Retaining, 3/4" x 2 3/8" Blades
11-82	6 1/4" Skin Hook & Retractor, Double End
11-88	7" Finger Retractor, Style A, Blade 1/4" x 11/16"
11-90	7" Finger Retractor, Style B, Blade 3/16" x 5/8"
11-92	7" Finger Retractor, Style C, Blade 1/8" x 1/4"
11-94	7" Finger Retractor, Style D, Blade 3/32" x 3/8"
11-96	7" Finger Retractor, Style E, Blade 1/4" x 7/8"

**MILTEX NECK RAKE RETRACTORS**

11-30	6" Neck Rake Retractor, 1 Sharp Prong, Flexible
11-32	6" Neck Rake Retractor, 2 Sharp Prongs, Flexible
11-34	6" Neck Rake Retractor, 3 Sharp Prongs, Flexible
11-36	6" Neck Rake Retractor, 1 Plunt Prong, Flexible
11-38	6" Neck Rake Retractor, 2 Blunt Prongs, Flexible
11-40	6" Neck Rake Retractor, 3 Blunt Prongs, Flexible
11-50	6" Neck Rake Retractor, 1 Sharp Prong, Rigid
11-52	6" Neck Rake Retractor, 2 Sharp Prongs, Rigid
11-54	6" Neck Rake Retractor, 3 Sharp Prongs, Rigid
11-56	6" Neck Rake Retractor, 1 Blunt Prong, Rigid
11-58	6" Neck Rake Retractor, 2 Blunt Prongs, Rigid
11-60	6" Neck Rake Retractor, 3 Blunt Prongs, Rigid
11D-27	5" Neck Rake Retractor, 3 Sharp Prongs, Flexible

**MILTEX RAGNELL RETRACTOR**

Double end 3 x 8mm and 5 x 15mm blades. Delicate.

11-73	6" Ragnell Retractor
-------	----------------------

**MILTEX RIBBON RETRACTORS**

Malleable.

11-130	3/4" x 13" Ribbon Retractor
11-132	1" x 13" Ribbon Retractor
11-134	1 1/4" x 13" Ribbon Retractor
11-136	1 1/2" x 13" Ribbon Retractor
11-140	2" x 13" Ribbon Retractor
11-144	5/8" x 8" Baby Ribbon Retractor
11-146	1" x 8" Baby Ribbon Retractor

**MILTEX RICHARDSON RETRACTORS**

11-246	9 1/2" Richardson Retractor w/loop handle, 1" x 3/4"
11-247	9 1/2" Richardson Retractor w/loop handle, 1 1/4" x 1"
11-248	9 1/2" Richardson Retractor w/loop handle, 1 1/2" x 1 1/2"
11-249	9 1/2" Richardson Retractor w/loop handle, 2" x 3/4"
11-250	9 1/2" Richardson Retractor w/hollow grip handle, 1" x 3/4"
11-252	9 1/2" Richardson Retractor w/hollow grip handle, 1 1/4" x 1"
11-254	9 1/2" Richardson Retractor w/hollow grip handle, 1 1/2" x 1"
11-256	9 1/2" Richardson Retractor w/hollow grip handle, 2" x 3/4"
11-260	Richardson-Eastman Retractor, Double End, 2/set (10" & 11")
11-260A	Richardson-Eastman Retractor, Double End, 11"
11-260B	Richardson-Eastman Retractor, Double End, 10"

**MILTEX RIGBY APPENDECTOMY RETRACTOR**

11-374 6 3/4" Appendectomy Retractor w/grip lock

**MILTEX SENN RETRACTORS**

11-74 6 3/8" Senn Retractor, 3 Sharp Prongs

11-76 6 3/8" Senn Retractor, 3 Blunt Prongs

**MILTEX VOLKMAN FINGER RETRACTORS**

11-42 4 1/2" Finger Retractor, 2 Prongs, Sharp

11-43 4 1/2" Finger Retractor, 3 Prongs, Sharp

11-44 4 1/2" Finger Retractor, 4 Prongs, Sharp

11-45 4 1/2" Finger Retractor, 6 Prongs, Sharp

11-46 4 1/2" Finger Retractor, 2 Prongs, Blunt

11-47 4 1/2" Finger Retractor, 3 Prongs, Blunt

11-48 4 1/2" Finger Retractor, 4 Prongs, Blunt

11-49 4 1/2" Finger Retractor, 6 Prongs, Blunt

**MILTEX VOLKMAN RETRACTORS**

11-500 8 1/2" Volkman Retractor, Hook, Sharp, Hollow Handle

11-502 8 1/2" Volkman Retractor, Hook, Sharp, Ring Handle

11-504 8 1/2" Volkman Retractor, 2 Prongs, Sharp, Hollow Handle

11-506 8 1/2" Volkman Retractor, 3 Prongs, Sharp, Hollow Handle

11-508 8 1/2" Volkman Retractor, 4 Prongs, Sharp, Hollow Handle

11-510 8 1/2" Volkman Retractor, 6 Prongs, Sharp, Hollow Handle

11-512 8 1/2" Volkman Retractor, 2 Prongs, Blunt, Hollow Handle

11-514 8 1/2" Volkman Retractor, 3 Prongs, Blunt, Hollow Handle

11-516 8 1/2" Volkman Retractor, 4 Prongs, Blunt

11-530 8" Volkman Retractor, 2 Prongs, Sharp, Ring Handle

11-532 8" Volkman Retractor, 3 Prongs, Sharp, Ring Handle

11-534 8" Volkman Retractor, 4 Prongs, Sharp, Ring Handle

11-536 8" Volkman Retractor, 6 Prongs, Sharp, Ring Handle

11-538 8" Volkman Retractor, 2 Prongs, Blunt, Ring Handle

11-540 8" Volkman Retractor, 3 Prongs, Blunt, Ring Handle

11-542 8" Volkman Retractor, 4 Prongs, Blunt, Ring Handle

11-544 8" Volkman Retractor, 6 Prongs, Blunt, Ring Handle

11-546 8 3/4" Volkman Retractor, 2 Prongs, Sharp, USA Pattern Ring Handle

11-547 8 3/4" Volkman Retractor, 3 Prongs, Sharp, USA Pattern Ring Handle

11-548 8 3/4" Volkman Retractor, 4 Prongs, Sharp, USA Pattern Ring Handle

11-549 8 3/4" Volkman Retractor, 6 Prongs, Sharp, USA Pattern Ring Handle

11-551 8 3/4" Volkman Retractor, 2 Prongs, Blunt, USA Pattern Ring Handle

11-552 8 3/4" Volkman Retractor, 3 Prongs, Blunt, USA Pattern Ring Handle

11-553 8 3/4" Volkman Retractor, 4 Prongs, Blunt, USA Pattern Ring Handle

11-554 8 3/4" Volkman Retractor, 6 Prongs, Blunt, USA Pattern Ring Handle

27-230 9" Volkman Retractor, Hook, Sharp, 15mm Deep, Ring Handle

27-232 9" Volkman Retractor, Hook, Blunt, 15mm Deep, Ring Handle

**MILTEX WEITLANER RETRACTORS**

11-600 4" Weitlaner Retractor, 2 x 3 Teeth, Sharp

11-602 4" Weitlaner Retractor, 2 x 3 Teeth, Blunt

11-608 5 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp

11-610 5 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt

11-611 5 1/2" Weitlaner Retractor, Solid Blades, Blunt

11-614 6 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp

11-616 6 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt

11-617 8" Weitlaner Retractor, 3 x 4 Teeth, Sharp

11-618 8" Weitlaner Retractor, 3 x 4 Teeth, Blunt

11-620 9 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp

11-622 9 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt

**MILTEX SPRING WIRE RETRACTOR**

11-19 3 1/4" Wire Retractor w/blades 10mm W x 13mm D

**MILTEX WIENER HOOKS**

18-418 5" Wiener Hook, Sharp

18-420 5" Wiener Hook, Blunt

**MILTEX PAYR PYLORUS CLAMPS**

16-229 6" Baby Payr Pylorus Clamp

16-230 8" Payr Pylorus Clamp

16-232 11" Payr Pylorus Clamp

16-234 13 3/4" Payr Pylorus Clamp

**MILTEX PERCY INTESTINAL FORCEPS**

Serrated jaws.

16-62 6" Percy Intestinal Forceps

**MILTEX WILLIAMS INTESTINAL FORCEPS**

Cross serrated jaws.

16-56 6 1/2" Williams Intestinal Forceps

**MILTEX BABCOCK INTESTINAL FORCEPS**

16-42 Baby, 5 1/2" Intestinal Forceps, Delicate Jaws 6mm Wide

16-44 6 1/4" Intestinal Forceps, Jaws 9mm Wide

16-46 8 1/4" Intestinal Forceps, Jaws 10mm Wide

16-48 9 1/2" Intestinal Forceps, Jaws 15mm Wide

**MILTEX BAINBRIDGE FORCEPS**

Longitudinal serrations with cross tip serrations.

7-230 6" Bainbridge Forceps, Straight

7-232 6" Bainbridge Forceps, Curved

7-234 7 1/4" Bainbridge Forceps, Straight

7-236 7 1/4" Bainbridge Forceps, Curved

**MILTEX BALLENGER SWIVEL KNIVES**

20-216 7 3/4" Swivel Knife, Straight, 3mm Width

20-218 7 3/4" Swivel Knife, Straight, 4mm Width

20-220 7 3/4" Swivel Knife, Straight, 5mm Width

20-222 7 3/4" Swivel Knife, Bayonet, 3mm Width

20-224 7 3/4" Swivel Knife, Bayonet, 4mm Width

**MILTEX BALLENGER V-SHAPED CHISEL**

20-278	6 1/4" Ballenger V-Shaped Chisel, 2mm Wide
20-280	6 1/4" Ballenger V-Shaped Chisel, 4mm Wide
20-282	6 1/4" Ballenger V-Shaped Chisel, 6mm Wide
20-284	6 1/4" Ballenger V-Shaped Chisel, 8mm Wide

**MILTEX BALLENTINE FORCEPS**

Longitudinal serrations. Single tooth. Straight.

30-1718	8 1/2" Ballentine Forceps
---------	---------------------------

**MILTEX BETHUNE RIB SHEARS**

Double curved handles.

25-350	13 3/4" Rib Shears
--------	--------------------

**MILTEX BLAKE GALL STONE FORCEPS**

14-84	8 1/2" Gall Stone Forceps, Straight
14-86	8 1/2" Gall Stone Forceps, Curved

**MILTEX BRONCHOSCOPIC FORCEPS**

Double articulated jaws.

23-600	Bronchoscopic Forceps, Serrated Grasping Jaws, 30cm Shaft, 4mm Pass Diameter
23-622	Bronchoscopic Forceps, Kahler Alligator Grasping Jaws, 40cm Shaft, 4mm Pass Diameter
23-640	Bronchoscopic Forceps, Serrated Jaws 3 x 12mm, 30cm Shaft, 4mm Pass Diameter
23-642	Bronchoscopic Forceps, Serrated Jaws 3 x 12mm, 40cm Shaft, 4mm Pass Diameter
23-644	Bronchoscopic Forceps, Serrated Jaws 3 x 12mm, 50cm Shaft, 4mm Pass Diameter
23-680	Bronchoscopic Forceps, Jackson Cup Jaws 2mm Diameter Straight, 30cm Shaft
23-690	Bronchoscopic Forceps, Jackson Cup Jaws 4mm Diameter Straight, 30cm Shaft
23-720	Bronchoscopic Forceps, Jackson Cup Jaws 4mm Diameter Angular, 30cm Shaft
23-804	Bronchoscopic Forceps, Serrated Jaws 2 x 8mm, Miniature Size, 20cm Shaft, 3mm Pass Diameter

**MILTEX MENGHINI BIOPSY NEEDLES**

Improved model with built-in stop plug that never requires replacement.

13-130	1.2mm Biopsy Needle, 18G For Liver Biopsy, Adult Size, 70mm Long
13-132	1.4mm Biopsy Needle, 17G For Liver Biopsy, Adult Size, 70mm Long
13-156	1.8mm Biopsy Needle, 16G For Kidney Biopsy, 120mm Long

**MILTEX PATTERSON TROCARS**

8 1/2". Trocar points, cannula and handle are stainless.

13-10	Size 3 Patterson Trocar, 10.5mm Dia.
13-11	Size 4 Patterson Trocar, 12mm Dia.
13-8	Size 1 Patterson Trocar, 6mm Dia.
13-9	Size 2 Patterson Trocar, 7.5mm Dia.

**MILTEX JOHNS HOPKINS BULLDOG CLAMP**

7-280	1 1/2" (3.8 cm) Bulldog Clamp, Straight
7-282	2" (5 cm) Bulldog Clamp, Straight
7-284	2 1/2" (6.3 cm) Bulldog Clamp, Straight
7-286	3" (7.6 cm) Bulldog Clamp, Straight
7-290	1 1/2" (3.8 cm) Bulldog Clamp, Curved
7-292	2" (5 cm) Bulldog Clamp, Curved
7-294	2 1/2" (6.3 cm) Bulldog Clamp, Curved
7-296	3" (7.6 cm) Bulldog Clamp, Curved
7-298	3 1/2" (8.9 cm) Bulldog Clamp, Curved

**MILTEX BULLDOG CLAMP APPLYING FORCEPS**

25-1420	9 1/2" Bulldog Clamp Applying Forceps
---------	---------------------------------------

**MILTEX BUNNELL HAND DRILL**

Lightweight and chrome plated with 5/32" capacity stainless Jacobs chuck. Cannulated for use with Kirschner wires and thin drills. Chuck key stores in handle.

27-180	Hand Drill
27-202	.035" Dia., 12" Kirschner Wires, Trocar Point w/flat chuck end
27-204	.045" Dia., 12" Kirschner Wires, Trocar Point w/flat chuck end
27-206	.0625" Dia., 12" Kirschner Wires, Trocar Point w/flat chuck end
27-208	.035" Dia., 7 1/2" Kirschner Wires, Trocar Point At Both Ends
27-210	.045" Dia., 7 1/2" Kirschner Wires, Trocar Point At Both Ends
27-212	.0625" Dia., 7 1/2" Kirschner Wires, Trocar Point At Both Ends

**MILTEX BUNNELL TENDON STRIPPER**

6"	
27-1020	3mm I.D. Tendon Stripper
27-1022	4mm I.D. Tendon Stripper
27-1024	5mm I.D. Tendon Stripper
27-1026	6mm I.D. Tendon Stripper
27-1028	7mm I.D. Tendon Stripper

**MILTEX SPURLING-KERRISON CERVICAL RONGEURS**

26-330	Up-Bite Cervical Rongeur, 3mm, 6" Shaft
26-331	Up-Bite Cervical Rongeur, 5mm, 6" Shaft
26-332	Down-Bite Cervical Rongeur, 3mm, 6" Shaft
26-335	40° Forward Angle Cervical Rongeur, 5mm, 6" Shaft
26-336	Up-Bite Cervical Rongeur, 3mm, 8" Shaft
26-337	Up-Bite Cervical Rongeur, 5mm, 8" Shaft
26-338	Down-Bite Cervical Rongeur, 3mm, 8" Shaft
26-339A	40° Forward Angle Cervical Rongeur, 3mm, 8" Shaft
26-339B	40° Forward Angle Cervical Rongeur, 5mm, 8" Shaft

**MILTEX-KERRISON CERVICAL RONGEURS**

7" shaft. Improved handle.

26-390	5mm Cervical Rongeur, 40° Forward Angle
26-393A	2mm Down-Bite Cervical Rongeur
26-388	2mm Cervical Rongeur, 40° Forward Angle
26-389	3mm Cervical Rongeur, 40° Forward Angle
26-392	5mm Up-Bite Cervical Rongeur
26-392A	2mm Up-Bite Cervical Rongeur
26-392B	3mm Up-Bite Cervical Rongeur
26-393	5mm Down-Bite Cervical Rongeur
26-393B	3mm Down-Bite Cervical Rongeur

**MILTEX STILLE PLASTER SHEARS**

Serrated blades.

27-3060	11" Plaster Shears, Medium
27-3062	15" Plaster Shears, Standard

**MILTEX STOOKEY RONGEUR**

Angled on side. 3mm bite.

26-200	7 3/4" Stookey Rongeur
--------	------------------------

**MILTEX U.S. PATTERN CAST KNIFE**

Wood handle. 1 1/2" blade.

27-3002	5 1/2" Cast Knife
---------	-------------------

**MILTEX WILLIGER BONE CURETTE**

Double end. 3mm and 4mm oval cups.

21-322	5 1/2" Williger Bone Curette
--------	------------------------------

**MILTEX WOLFE-BOEHLER CAST BREAKER**

Jaws 15mm wide with deep serrations.

27-3080	10" Cast Breaker
---------	------------------

**MILTEX SPURLING PITUITARY RONGEURS**

4 x 10mm cup jaws.

26-430	5" Shank Pituitary Rongeur, Straight
26-435	7" Shank Pituitary Rongeur, Straight
26-436	7" Shank Pituitary Rongeur, Curved-Up
26-437	7" Shank Pituitary Rongeur, Curved-Down

**MILTEX BLOUNT KNEE RETRACTOR**

27-223	7" Knee Retractor, Blade 1 3/4"D, 5/8"W Tapering to 3/16"
--------	---

**MILTEX BRUNS PLASTER SHEARS**

Serrated blade.

27-3070	9 1/4" Plaster Shears
---------	-----------------------

**MILTEX REINER CAST KNIFE**

Metal handle. 1 1/2" blade.

27-3004	7" Cast Knife
---------	---------------

**MILTEX VOLKMAN DOUBLE END CURETTES**

27-920	5 1/2" Double End Curette, Oval Cups 5 x 10mm & 6 x 20mm
27-924	5" Double End Curette, Oval Cup 6 x 10mm & Round Cup 7mm Dia.
27-926	6 1/2" Double End Curette, Oval Cup 8 x 14mm & Round Cup 10mm Dia.
27-928	8" Double End Curette, Oval Cup 9 x 18mm & Round Cup 11mm Dia.

**MILTEX SMITH-PETERSEN GOUGES**

Overall length 8".

27-549	1/4" Gouge, Straight
27-550	3/8" Gouge, Straight
27-551	1/2" Gouge, Straight
27-553	5/8" Gouge, Straight
27-554	3/4" Gouge, Straight
27-556	1" Gouge, Straight
27-558	1 1/4" Gouge, Straight

**MILTEX SMITH-PETERSEN LAMINECTOMY RONGEURS**

3mm bite.

26-255	7 1/2" Laminectomy Rongeur, Straight Jaws
26-256	7 1/2" Laminectomy Rongeur, Slightly Curved Jaws
26-258	9 1/2" Laminectomy Rongeur, Straight Jaws

**MILTEX SMITH-PETERSEN OSTEOTOMES**

Overall length 8".

27-510	1/4" Osteotome, Straight
27-512	3/8" Osteotome, Straight
27-514	1/2" Osteotome, Straight
27-516	5/8" Osteotome, Straight
27-518	3/4" Osteotome, Straight
27-520	1" Osteotome, Straight
27-522	1 1/4" Osteotome, Straight
27-530	1/4" Osteotome, Curved
27-532	3/8" Osteotome, Curved
27-534	1/2" Osteotome, Curved
27-536	5/8" Osteotome, Curved
27-538	3/4" Osteotome, Curved
27-540	1" Osteotome, Curved
27-542	1 1/4" Osteotome, Curved

**MILTEX RANEY SCALP CLIP APPLYING FORCEPS**

26-590	6 1/4" Scalp Clip Applying Forceps
26-591	Reusable Scale Clips, Black Finish, 1 dz/pkg

**MILTEX SCHWARTZ CLIP APPLYING FORCEPS & CLIPS**

26-1040	8" Forceps
26-1041	1" Vessel Clip, Straight, Jaws 1.7 x 8mm
26-1042	1" Vessel Clip, Slightly Angled, Jaws 1.7 x 8mm
26-1043	1" Vessel Clip, Strong Angle, Jaws 1.7 x 8mm

**MILTEX PENFIELD DISSECTORS**

26-1450	7 1/4" Dissector, Style No. 1, Double End, 6mm Sharp Cup & 6mm Blunt Dissector
26-1451	7 3/4" Dissector, Style No. 2, Double End, Wax Packer & 6mm Blunt Dissector, Slightly Curved
26-1452	7 1/2" Dissector, Style No. 3, Double End, Wax Packer & 6mm Blunt Dissector, Fully Curved
26-1453	8 3/4" Dissector, Style No. 4, Single End, Light Pattern, 3mm Blunt Dissector, Slightly Curved
26-1454	11 1/2" Dissector, Style No. 5, Double End, Slightly Curved, Blunt Dissector Blades 7mm & 8mm Wide

**MILTEX SMITHWICK HOOKS**

26-1170	7" Smithwick Hook, Blunt, 10mm Deep
26-1174	12" Smithwick Hook, Dissector

**MILTEX MEAD MALLET**

19-796 7 1/2" Mead Mallet, Head 8 oz w/2 replaceable nylon caps

**Accessories**

19-797 Nylon Cap Replacement For 19-796

**MILTEX SWISS PATTERN OSTEOTOMES**

Overall length 5".

- 27-599 2mm Osteotome
- 27-600 4mm Osteotome
- 27-601 6mm Osteotome
- 27-601A 8mm Osteotome
- 27-602 10mm Osteotome
- 27-602A 12mm Osteotome
- 27-603 14mm Osteotome
- 27-604 18mm Osteotome
- 27-605 20mm Osteotome
- 27-610 Set of 9 Osteotomes 2mm To 20mm w/case

**Accessories**

27-609 Case Only, 5 1/2" x 6 1/2" x 1/2", Fitted To Hold 9 Osteotomes 2mm To 20mm

**MILTEX GERZOG MALLET**

8 oz head lead filled with stainless jacket.

19-790 7 1/2" Gerzog Mallet

**MILTEX POOLE SUCTION TUBE**

- 10-310 23 Fr Poole Suction Tube, Curved
- 10-312 30 Fr Poole Suction Tube, Straight

**MILTEX SCHWEIGER EXTRA-CAPSULAR FORCEPS**

4 x 5 teeth.

18-1058 3 3/4" Schweiger Extra-Capsular Forceps

**MILTEX VIRCHOW BRAIN SECTIONING KNIFE**

Double edge.

- 34-50 7 1/2" Brain Sectioning Knife, Blade 180mm x 20mm
- 34-52 10" Brain Sectioning Knife, Blade 250mm x 40mm



**MILTEX 4" KEYES DERMAL PUNCHES**

- 33-20 2 mm Dermal Punch
- 33-22 3 mm Dermal Punch
- 33-24 4 mm Dermal Punch
- 33-25 5 mm Dermal Punch
- 33-26 6 mm Dermal Punch
- 33-28 8 mm Dermal Punch
- 33-30 Keyes Dermal Punch Set w/3 1/2" Handle, Screw-In Keyes Knives 2,3,4,5,6 & 7 mm Mounted on Storage & Sterilization Rack Handle For Dermal Punch
- 33-30-1 2 mm Replacement Head For Dermal Punch
- 33-30-2 3 mm Replacement Head For Dermal Punch
- 33-30-3 4 mm Replacement Head For Dermal Punch
- 33-30-4 5 mm Replacement Head For Dermal Punch
- 33-30-5 6 mm Replacement Head For Dermal Punch
- 33-30-6 7 mm Replacement Head For Dermal Punch
- 33-30-7 7 mm Replacement Head For Dermal Punch



**MILTEX STERILE DISPOSABLE BIOPSY PUNCHES**

Seamless, stainless steel cutting edge with ribbed handle for maximum control.


- 33-31 2 mm Biopsy Punch, 50/bx
- 33-31A 1.5 mm Biopsy Punch, 50/bx
- 33-32 3 mm Biopsy Punch, 50/bx
- 33-33 3.5 mm Biopsy Punch, 50/bx
- 33-34 4 mm Biopsy Punch, 50/bx
- 33-35 5 mm Biopsy Punch, 50/bx
- 33-36 6 mm Biopsy Punch, 50/bx
- 33-37 8 mm Biopsy Punch, 50/bx
- 33-31-P/25 2 mm Biopsy Punch w/plunger, 25/bx
- 33-31A-P/25 1.5 mm Biopsy Punch w/plunger, 25/bx
- 33-31AA 1mm Biopsy Punch, 50/bx
- 33-31AA-P/25 1 mm Biopsy Punch w/plunger, 25/bx
- 33-31B 2.5mm Biopsy Punch, 50/bx
- 33-38 Assorted Biopsy Punch, 50/bx

**MILTEX DERMAL CURETTES**

6 1/4". Oval spoon.

- 33-12 Size 0 Dermal Curette
- 33-13 Size 1 Dermal Curette
- 33-14 Size 2 Dermal Curette
- 33-15 Size 3 Dermal Curette


Ask your representative  
to show you illustrations  
of instruments not  
shown here.



### MILTEX DISPOSABLE DERMAL CURETTES

Features a unique symmetrically ground beveled edge. Manufactured of surgical grade stainless steel providing a clean, precise edge every time guaranteed. Individually packaged in sterile pouches.

33-52	2mm Dermal Curette, 50/bx
33-53	3mm Dermal Curette, 50/bx
33-54	4mm Dermal Curette, 50/bx
33-55	5mm Dermal Curette, 50/bx
33-57	7mm Dermal Curette, 50/bx



### MILTEX FOX DERMAL CURETTES

5 1/2", round.

33-2	6 mm Dermal Curette
33-3	5 mm Dermal Curette
33-4	4 mm Dermal Curette
33-5	3 mm Dermal Curette
33-5A	2 mm Dermal Curette
33-5B	1 mm Dermal Curette




### MILTEX PIFFARD DERMAL CURETTES

33-10	Size 3 Dermal Curette, 6 1/4" Oval
33-11	Size 4 Dermal Curette, 6 1/4" Oval
33-11A	Size 5 Dermal Curette, 6 1/4" Oval
33-16-0	Size 0 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-00	Size 00 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-1	Size 1 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-2	Size 2 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-3	Size 3 Dermal Curette, 5 1/2" Oval Narrow Handles
33-6	Size 1 Dermal Curette, 6 1/4" Oval
33-8	Size 2 Dermal Curette, 6 1/4" Oval

### MILTEX SAALFIELD COMEDONE EXTRACTOR

With acne lancet. Professional pattern.

33-204	5 3/4" Comedone Extractor
33-206	4 1/4" Comedone Extractor



### MILTEX SCHAMBERG COMEDONE EXTRACTOR

3 3/4".

33-200	Square Loop Ends Comedone Extractor
33-201	Standard Pattern Comedone Extractor w/crimped small loop
33-202	Standard Pattern Comedone Extractor w/extra fine loops, Small Loop Crimped

### MILTEX COMEDONE EXTRACTORS

33-203	4" Comedone Extractor, Double End, Flat & Wire End Expressors
33-214	6" Comedone Extractor, One Fenestrated Cup & One Delicate Curved Lancet w/protective cap

### MILTEX GREG COMEDONE EXTRACTOR

One fenestrated cup and chuck with interchangeable No. 65 lancet blade, protected by cap.

33-215	5" Comedone Extractor
<b>Accessories</b>	
33-216	Replacement Blade Only



### MILTEX GRUBER EAR SPECULA

19-10	Set of 4 Ear Specula Sizes 1, 2, 3 & 4 w/screw top case, Light Short Pattern, Chrome
19-12	Set of 4 Ear Specula Sizes 1, 2, 3 & 4 w/o case, Light Short Pattern, Chrome
19-14	Single Size 1 Ear Specula, 2x3mm, Light Short Pattern, Chrome
19-16	Single Size 2 Ear Specula, 3x4mm, Light Short Pattern, Chrome
19-18	Single Size 3 Ear Specula, 4x5mm, Light Short Pattern, Chrome
19-20	Single Size 4 Ear Specula, 5x6mm, Light Short Pattern, Chrome
19-21	Set of 4 Ear Specula Sizes 1, 2, 3 & 4 Oval, Standard Adult Sizes, Long Pattern, Chrome
19-22	Single Size 1 Ear Specula, 4x5mm Oval, Standard Adult Sizes, Long Pattern, Chrome
19-23	Single Size 2 Ear Specula, 5x6mm Oval, Standard Adult Sizes, Long Pattern, Chrome
19-24	Single Size 3 Ear Specula, 6x7mm Oval, Standard Adult Sizes, Long Pattern, Chrome
19-25	Single Size 4 Ear Specula, 7x8mm Oval, Standard Adult Sizes, Long Pattern, Chrome

#### Accessories

19-10A	Case Only For Set of 4 Specula (Child & Infant Sizes), Light Short Pattern, Chrome
--------	--



**MILTEX TOYNBEE EAR SPECULA**

Round. Chrome.

- 19-30 Set of 3 Ear Specula Sizes 1, 2 & 3
- 19-32 Single Size 1 Ear Specula, 4mm
- 19-34 Single Size 2 Ear Specula, 5mm
- 19-36 Single Size 3 Ear Specula, 6mm

**MILTEX BOUCHERON EAR SPECULA**

Chrome. Round.

- 19-2 Set of 4 Ear Specula: Sizes 1, 2, 3 & 4
- 19-2-1 Single Size 1 Ear Specula, 4.5mm
- 19-2-2 Single Size 2 Ear Specula, 5.5mm
- 19-2-3 Single Size 3 Ear Specula, 6.5mm
- 19-2-4 Single Size 4, 7.5mm


**MILTEX FARRIOR EAR SPECULA**

Oval and angled end. Nonglare stainless.

- 19-50-1 Single Size 1 Ear Specula, 3 x 4mm
- 19-50-2 Single Size 2 Ear Specula, 4 x 5mm
- 19-50-3 Single Size 3 Ear Specula, 4.5 x 5.5mm
- 19-50-4 Single Size 4 Ear Specula, 5 x 6mm
- 19-50-5 Single Size 5 Ear Specula, 5.5 x 6.5mm
- 19-50-6 Single Size 6 Ear Specula, 6 x 7mm
- 19-50-7 Single Size 7 Ear Specula, 6.5 x 7.5mm
- 19-50-8 Single Size 8 Ear Specula, 7 x 8mm
- 19-50-9 Single Size 9 Ear Specula, 7.5 x 8.5mm

**MILTEX SCHUHKNECHT POSTAURICULAR RETRACTOR**

- 19-655 5" Postauricular Retractor, 3 x 3 Teeth, 3" Opening



**MILTEX BUCK EAR CURETTES**

6 1/2"

- 19-260 00 Ear Curette, Sharp, Straight Size
- 19-262 0 Ear Curette, Sharp, Straight Size
- 19-264 1 Ear Curette, Sharp, Straight Size
- 19-266 2 Ear Curette, Sharp, Straight Size
- 19-268 3 Ear Curette, Sharp, Straight Size
- 19-270 00 Ear Curette, Sharp, Angled Size
- 19-272 0 Ear Curette, Sharp, Angled Size
- 19-274 1 Ear Curette, Sharp, Angled Size
- 19-276 2 Ear Curette, Sharp, Angled Size
- 19-278 3 Ear Curette, Sharp, Angled Size
- 19-280 00 Ear Curette, Blunt, Straight Size
- 19-282 0 Ear Curette, Blunt, Straight Size
- 19-284 1 Ear Curette, Blunt, Straight Size
- 19-286 2 Ear Curette, Blunt, Straight Size
- 19-288 3 Ear Curette, Blunt, Straight Size
- 19-290 00 Ear Curette, Blunt, Angled Size
- 19-292 0 Ear Curette, Blunt, Angled Size
- 19-294 1 Ear Curette, Blunt, Angled Size
- 19-296 2 Ear Curette, Blunt, Angled Size
- 19-298 3 Ear Curette, Blunt, Angled Size

**MILTEX DISPOSABLE EAR CURETTES**

Each size is color coded for easy identification. The handle is designed with a scored area for finger placement. Single use packaged in an individual clear sleeve.

- 19-320 Ring Tip Ear Curette, Pink, 50/bx
- 19-321 Loop Tip Ear Curette, White, 50/bx
- 19-322 Small Spoon Tip Ear Curette, Gray, 50/bx
- 19-323 Large Spoon Tip Ear Curette, Peach, 50/bx



**MILTEX BILLEAU FLEXIBLE EAR LOOPS**

Overall length 6 1/2".

- 19-314 Small Ear Loop, Size No. 1
- 19-316 Medium Ear Loop, Size No. 2
- 19-318 Large Ear Loop, Size No. 3

**MILTEX SPRATT (BRUN) MASTOID CURETTES**

Oval cups. 6 1/2".

- 19-700 Size 5/0 Mastoid Curette
- 19-702 Size 4/0 Mastoid Curette
- 19-704 Size 3/0 Mastoid Curette
- 19-706 Size 2/0 Mastoid Curette
- 19-708 Size 0 Mastoid Curette
- 19-710 Size 1 Mastoid Curette
- 19-712 Size 2 Mastoid Curette
- 19-714 Size 3 Mastoid Curette
- 19-716 Size 4 Mastoid Curette
- 19-718 Size 5 Mastoid Curette
- 19-720 Size 6 Mastoid Curette


**MILTEX BUCK (SPENCER) EAR PROBE**

- 19-210 5 1/2" Ear Probe, Sterling Shaft

**MILTEX QUIRE MECHANICAL FINGER**


For removal of foreign bodies. Spring handle.

- 19-400 4" Quire Mechanical Finger



**MILTEX LATHBURY APPLICATORS**

- 19-190 5 1/2" Lathbury Applicator, Triangular Tip
- 19-192 5 1/2" Lathbury Applicator, Cross Serrated Tip



**MILTEX EAR SYRINGES**

Chrome-plated, all brass syringes with precision ground plunger, and shaft graduated in 1/2 ounce. Complete with 1 tapered and 1 short bulbous tip. Easily disassembled for sterilization.

- 19-380 2 oz Ear Syringe w/o shield
- 19-381 2 oz Ear Syringe w/shield
- 19-382 3 oz Ear Syringe w/o shield
- 19-383 3 oz Ear Syringe w/shield
- 19-384 4 oz Ear Syringe w/o shield
- 19-385 4 oz Ear Syringe w/shield
- 19-390 Ear Syringe Short Bulbous Tip, Chrome
- 19-394 Ear Syringe Taper Tip, Chrome
- 19-396 Ear Syringe Shield, Chrome

**MILTEX EAR SUCTION TUBE**

Two interchangeable tips 20 ga and 22 ga.

- 19-1820 Ear Suction Tube



**MILTEX BARON SUCTION TUBES**

7.5 cm tube with finger cut-off.

- 19-580 3 Fr Suction Tubes
- 19-582 5 Fr Suction Tubes
- 19-584 7 Fr Suction Tubes



**MILTEX FRAZIER FERGUSON SUCTION TUBES**

With finger cut-off.

- 19-570 6 Fr Suction Tubes, Angled
- 19-572 7 Fr Suction Tubes, Angled
- 19-573 8 Fr Suction Tubes, Angled
- 19-574 9 Fr Suction Tubes, Angled
- 19-575 10 Fr Suction Tubes, Angled
- 19-576 11 Fr Suction Tubes, Angled
- 19-577 12 Fr Suction Tubes, Angled
- 20-460 9 Fr Suction Tubes, Straight

**MILTEX HOUSE SUCTION-IRRIGATION TUBES**


9cm working length.

- 19-562 1.8mm Suction Tube & 1.8mm Irrigating Tube
- 19-563 2.5mm Suction Tube & 2mm Irrigating Tube
- 19-564 3mm Suction Tube & 2mm Irrigating Tube
- 19-565 4.5mm Suction Tube & 2mm Irrigating Tube

**MILTEX HARTMAN EAR FORCEPS**

Delicate serrated jaws 2.5 x 6mm.

- 19-410 5" Hartman Ear Forceps



**MILTEX JANSEN (GRUENWALD) EAR FORCEPS**

6-197 6 1/2" Ear Forceps, Bayonet Shape



**MILTEX LUCAE EAR FORCEPS**

19-370 5 1/2" Serrated Ear Forceps, Bayonet Shape


**MILTEX MICRO ALLIGATOR EAR FORCEPS**

- 19-2081-B Ear Forceps, Serrated Jaws 0.6mm Wide, Ebony Finish
- 19-2090-B Ear Forceps, Oval Cup Jaws 0.5mm Wide, Ebony Finish
- 19-2100 Ear Forceps, 3 1/4" Shaft, Serrated Jaws 0.8mm Wide
- 19-2100-B Ear Forceps, Serrated Jaws 0.8mm Wide, Ebony Finish
- 19-2110 Ear Forceps, Oval Cup Jaws 0.8mm Wide
- 19-2110-B Ear Forceps, Oval Cup Jaws 0.8mm Wide, Ebony Finish

**MILTEX NOYES ALLIGATOR FORCEPS**

Spring handle.

19-404 4" Noyes Alligator Forceps



**MILTEX STRUEMPEL EAR FORCEPS**

Alligator type. Fenestrated oval spoon jaws.

19-438 3" Shaft Struempel Ear Forceps



**MILTEX WEINGARTNER EAR FORCEPS**


Alligator type. Concave serrated jaws.

19-437 3" Shaft Weingartner Ear Forceps



**MILTEX WILDE EAR FORCEPS**

- 19-350 5" Ear Forceps, Angular w/serrated tips
- 19-356 5 1/4" Ear Forceps, Cross Action, Angular w/serrated tips
- 20-120 6 1/2" Ear Forceps, Angular




**MILTEX WULLSTEIN EAR FORCEPS**

- 19-2120 3" Shaft Ear Forceps, Oval Cup, Jaws 1mm Wide
- 19-2130 3" Shaft Ear Forceps, Round Cup, Jaws 1mm Diameter

**MILTEX SEXTON EAR KNIFE**

- 19-500 7" Ear Knife, Straight
- 19-502 7" Ear Knife, Bayonet
- 19-504 7" Ear Knife, Angled
- 19-520 7 1/2" Ear Knife, In Reversible Swivel Handle



19-2060

**MILTEX WULLSTEIN EAR SCISSORS**

Straight and very delicate.

19-2060 3" Shaft Ear Scissors

**MILTEX PERKINS ENDAURAL RETRACTORS**

- 19-668 5 1/4" Perkins Endaural Retractor, Left Serrated Solid Blade, Right Side 3 Teeth
- 19-669 5 1/4" Perkins Endaural Retractor, Right Serrated Solid Blade, Left Side 3 Teeth

**MILTEX ZAUFEL-JANSEN RONGEUR**

Double action. Slightly curved jaws 5mm wide.

19-840 7" Zaufel-Jansen Rongeur

**MILTEX HARTMAN EAR SPECULA**

Round.

- 19-40 Set of 3 Ear Specula Sizes 3, 4 & 5mm
- 19-48 Set of 4 Ear Specula Sizes 4, 5, 6 & 7mm
- 19-48-A Single Size 3mm Ear Specula
- 19-48-B Single Size 4mm Ear Specula
- 19-48-C Single Size 5mm Ear Specula
- 19-48-D Single Size 6mm Ear Specula
- 19-48-E Single Size 7mm Ear Specula

**VANTAGE VIENNA NASAL SPECULUM**

- V920-2 5 1/2" Nasal Specula, Large
- V920-3 5 1/2" Nasal Specula, Medium
- V920-4 5 1/2" Nasal Specula, Small



**MILTEX VIENNA NASAL SPECULA**

- 20-2 5 3/4" Nasal Specula, Large
- 20-3 5 3/4" Nasal Specula, Medium
- 20-4 5 3/4" Nasal Specula Small
- 20-6 Lightweight, 5 3/4" Nasal Specula, Adult Size
- 20-7 Lightweight, 5 3/4" Nasal Specula, Child Size

**MILTEX BECKMAN NASAL SPECULUM**

Angled on side.

20-24 6" Beckman Nasal Speculum



20-28 to 32

**MILTEX KILLIAN SEPTUM SPECULUM**

6" with set screw.

- 20-28 2" Long Blades Septum Speculum
- 20-29 2 1/2" Long Blades Septum Speculum
- 20-30 3" Long Blades Septum Speculum
- 20-32 3 1/2" Long Blades Septum Speculum

**MILTEX SONNENSCHN NASAL SPECULUM**

20-20 6" Nasal Speculum

**MILTEX TIECH (HALLE) INFANT NASAL SPECULUM**

20-12 5 3/4" Infant Nasal Speculum

**MILTEX WIENER ANTRUM RASP**

With trocar tip.

- 20-378 8" Wiener Antrum Rasp, No 1 Small
- 20-380 8" Wiener Antrum Rasp, No 2 Medium
- 20-382 8" Wiener Antrum Rasp, No 3 Large

**MILTEX COTTLE SEPTUM SPECULUM**

Ebony finish. With set screw.

- 20-36-35 6" Septum Speculum, Narrow 35mm Long Blades
- 20-36-50 6" Septum Speculum, Narrow 50mm Long Blades
- 20-36-75 6" Septum Speculum, Narrow 75mm Long Blades
- 20-36-90 6" Septum Speculum, Narrow 90mm Long Blades
- 20-37-45 5 1/2" Septum Speculum, Extra Delicate, Very Thin 45mm Long Blades
- 20-37-75 5 1/2" Septum Speculum, Extra Delicate, Very Thin 75mm Long Blades



20-124

**MILTEX HARTMAN NASAL DRESSING FORCEPS**

Serrated tips.

- 20-124 7" Nasal Dressing Forceps, Improved Delicate Pattern
- 20-126 7 1/4" Nasal Dressing Forceps, Regular Model

**MILTEX PIERCE ANTRUM TROCAR**

5mm diameter cannula and tubing connection.

20-404 7" Pierce Antrum Trocar

**MILTEX PIERCE DOUBLE END ELEVATOR**

20-230 8 1/2" Double End Elevator w/right &amp; left blades

**MILTEX WALSHAM SEPTUM STRAIGHTENER**

21-657 9 1/4" Septum Straightener, Universal w/2 convex blades  
 21-658 9 1/4" Septum Straightener, Right w/1 convex & 1 concave blade  
 21-659 9 1/4" Septum Straightener, Left w/1 convex & 1 concave blade

**MILTEX COTTLE SEPTUM ELEVATOR**

Double end. Graduated in cm.

21-60 9" Cottle Septum Elevator

**MILTEX PENNINGTON SEPTUM ELEVATOR**

Double end for right and left side.

20-246 8 1/4" Pennington Septum Elevator

**MILTEX COUPLAND SUCTION TUBE**

With 2 interchangeable tips 1.5 and 2mm. Chrome.

20-466 10" Coupland Suction Tube

**MILTEX BALLENGER TONSIL SEIZING FORCEPS**

Slight curve. 3 x 3 prongs. Open ring.

22-560 8 1/2" Tonsil Seizing Forceps

**MILTEX RONIS ADENOID PUNCH**

Shaft curved upwards.

20-910 5" Adenoid Punch, Small Triangular Jaws  
 20-912 5" Adenoid Punch, Large Triangular Jaws

**MILTEX BARNHILL ADENOID CURETTES**

Overall length 8 1/2".

20-800 Size 0 Adenoid Curette, 11mm Blade  
 20-802 Size 1 Adenoid Curette, 13mm Blade  
 20-804 Size 2 Adenoid Curette, 15mm Blade  
 20-806 Size 3 Adenoid Curette, 17mm Blade  
 20-808 Size 4 Adenoid Curette, 19mm Blade

**MILTEX VOGEL ADENOID CURETTE**

Infant size. 7mm blade.

20-809 8 1/2" Vogel Adenoid Curette

**MILTEX TISSUE TWIST HOOK**

21-430 5 1/2" Tissue Twist Hook

**MILTEX BLAIR CLEFT PALATE ELEVATOR**

22-300 7 3/4" Cleft Palate Elevator, L-Shape, Blunt  
 22-302 7 3/4" Cleft Palate Elevator, Right Angle, Semi-Sharp Tip  
 22-304 8" Cleft Palate Elevator, Curved Blade, Blunt

**MILTEX DAVIS MOUTH GAG**

2-120 6 1/4" Davis Mouth Gag, Double Bite, Left, w/o blades  
 2-122 6 1/4" Davis Mouth Gag, Double Bite, Right, w/o blades  
 2-130 Size 1 Blade, 3/4" x 1"  
 2-132 Size 2 Blade, 1" x 2 3/4"  
 2-132A Size 2 Blade, 1" x 2 3/4" w/tube  
 2-134 Size 3 Blade, 1" x 3 1/4"  
 2-134A Size 3 Blade, 1" x 3 1/4" w/tube  
 2-136 Size 4 Blade, 1" x 3 3/4"  
 2-136A Size 4 Blade, 1" x 3 3/4" w/tube  
 2-138 Size 5 Blade, 1" x 4"  
 2-138A Size 5 Blade, 1" x 4" w/tube  
 2-140 Size 2 Ring Tongue Retractor Blades, 7/8" x 2 1/4", Left  
 2-141 Size 3 Ring Tongue Retractor Blades, 1 1/8" x 3", Left  
 2-142 Size 4 Ring Tongue Retractor Blades, 1 1/4" x 3 1/2", Left  
 2-143 Size 5 Ring Tongue Retractor Blades, 1 1/2" x 3 3/4", Left  
 2-144 Size 2 Ring Tongue Retractor Blades, 7/8" x 2 1/4", Right  
 2-145 Size 3 Ring Tongue Retractor Blades, 1/18" x 3", Right  
 2-146 Size 4 Ring Tongue Retractor Blades, 1 1/4" x 3 1/2", Right  
 2-147 Size 5 Ring Tongue Retractor Blades, 1 1/2" x 3 3/4", Right

**MILTEX LARYNGEAL MIRRORS**

Boilable, thin model.

23-1-000 Size 000 Mirror w/octagon threaded handle  
 23-10-3 Size 3 Mirror w/octagon threaded handle  
 23-12-4 Size 4 Mirror w/octagon threaded handle  
 23-14-5 Size 5 Mirror w/octagon threaded handle  
 23-16-6 Size 6 Mirror w/octagon threaded handle  
 23-18-7 Size 7 Mirror w/octagon threaded handle  
 23-2-00 Size 00 Mirror w/octagon threaded handle  
 23-29-000 Size 000 Mirror Only w/threaded stem  
 23-30-00 Size 00 Mirror Only w/threaded stem  
 23-32-0 Size 0 Mirror Only w/threaded stem  
 23-34-1 Size 1 Mirror Only w/threaded stem  
 23-36-2 Size 2 Mirror Only w/threaded stem  
 23-38-3 Size 3 Mirror Only w/threaded stem  
 23-4-0 Size 0 Mirror w/octagon threaded handle  
 23-40-4 Size 4 Mirror Only w/threaded stem  
 23-42-5 Size 5 Mirror Only w/threaded stem  
 23-44-6 Size 6 Mirror Only w/threaded stem  
 23-46-7 Size 7 Mirror Only w/threaded stem  
 23-6-1 Size 1 Mirror w/octagon threaded handle  
 23-62 Universal Chuck Handle, Round, Chrome  
 23-8-2 Size 2 Mirror w/octagon threaded handle

**Accessories**

23-60 Octagon Threaded Handle

# Mini-Torr Plus®

## Non-Invasive Blood Pressure Monitor



- Cost effective NIBP
- Quick, accurate readings
- Proven reliable, oscillometric technology
- Easy-to-read LED display with adjustable brightness
- Optional pulse oximetry, integrated printer, temperature, and lightweight roll stand or IV pole

From the people who brought you  
the **BCI® 3301**

We have the  
solutions &  
accessories  
you need



---

smiths

... a leader in  
pulse oximetry  
for more than  
15 years




Endometrial Cannula



Endometrial Sampling Set



Monsel's Solution



Os Locator




Z-Type Forceps

**NEW!**


## MILTEX GYNECOLOGY PRODUCTS

Miltex, Inc., a name you can trust in providing your practice with the highest quality hand-held surgical instruments and disposables, is pleased to offer the following comprehensive collection of Gynecology products.

- **Sterile, Disposable SoftFlex™ and FirmFlex™ Endometrial Biopsy Cannula**, ideal for accurate and economical in-office specimen collection of endometrial tissue.
- **Sterile, Disposable Endometrial Sampling Set**, ideal for sample collection of clearly differentiated endometrial tissue. Set includes curette and syringe.
- **Monsel's and Lugol's Solutions**, conveniently packaged in single-use 8 ml bottles with applicator swabs for ready-to-use, one-step application.
- **Os Cervical Dilator Set**, excellent for cervical dilation prior to curette introduction for endometrial sampling, hysteroscopic and other related gynecological procedures. Three convenient patterns provided within one reusable and autoclavable set.
- **Z-Type Hysterectomy & Parametrium Scissors and Forceps**  
Nine Miltex Z-Type Scissors patterns with precise form and bend to cut tissue lateral to the uterus are configured to precisely match sixteen Miltex Z-Type Forceps, which provide a strong, yet atraumatic grip of the uterine tissue.



FirmFlex™  
SoftFlex™



Endometrial Sampling Set



Monsel's & Lugol's Solutions



Os Cervical Dilator Set



Z-Type Scissors & Forceps

*Miltex*  
Redefining Excellence  
www.miltex.com

**MILTEX BLUMENTHAL ORAL SURGERY RONGEUR**

22-480	6" Oral Surgery Rongeur, Beaks At 30° Angle
22-481	7" Oral Surgery Rongeur, Beaks At 30° Angle
22-482	6" Oral Surgery Rongeur, Beaks At 90° Angle
22D-420	6" Oral Surgery Rongeur, Beaks At 45° Angle
22D-422	4 1/2" Oral Surgery Rongeur, Beaks At 30° Angle
22D-424	4 1/2" Oral Surgery Rongeur, Beaks At 45° Angle

**MILTEX ORAL SURGERY RONGEUR**

22-486	5 1/2" Oral Surgery Rongeur, No. 4 Pattern, Slightly Curved Beaks
22-488	6 1/2" Oral Surgery Rongeur, No. 4A Pattern, Slightly Curved Beaks
22-490	6 3/4" Oral Surgery Rongeur, No. 5 Pattern, Side Cutting, Slightly Curved
22-492	6 1/2" Oral Surgery Rongeur, No. 6A Pattern, End Cutting, Curved Beaks
22D-426	6 1/2" Oral Surgery Rongeur, No. 5A Pattern, Side Cutting, Slightly Curved
22D-428	6" Oral Surgery Rongeur, No. 5S Pattern, Side Cutting, Slightly Curved
22D-432	6 1/2" Oral Surgery Rongeur, No. 1 Pattern, Strong Curved Beaks

**MILTEX TIVNEN TONSIL SEIZING FORCEPS**

Angled 3 x 3 prongs. Open ring.

22-586	8" Tonsil Seizing Forceps
--------	---------------------------

**MILTEX WHITE TONSIL SEIZING FORCEPS**

Medium curve and open ring.

22-550	7 1/2" Tonsil Seizing Forceps
22-552	9 1/2" Tonsil Seizing Forceps

**MILTEX BECK-SCHENCK TONSIL SNARE**

22-1010	10" Tonsil Snare, Complete w/1 straight & 1 vedder tip
---------	--

**Accessories**

22-1016	Medium Vedder Tip Only, 23 x 19mm
22-1020	Straight Tip Only, 2 1/2"

**MILTEX TONGUE SEIZING FORCEPS**

2-50	Collin Tongue Seizing Forceps, 6 3/4" Jaws, 25mm Wide
2-60	Young Tongue Seizing Forceps, 6 1/2" w/soft rubber jaws

**Accessories**

2-62	Rubber Jaws Only For 2-60
------	---------------------------

**MILTEX UNIVERSITY OF MINNESOTA CHEEK RETRACTOR**


Double end 14mm and 20mm wide.

22D-209	5 1/2" Cheek Retractor
---------	------------------------

**RUSCH OBERTO MOUTH PROP**


A specialized product for use in electro-convulsive therapy and epileptic system. Latex with central air vent. Non-sterile.

567500	Oberto Mouth Prop
--------	-------------------

**MILTEX JACKSON TRACHEA TUBES WITH ROTATING LOCK**

Improved model with rotating lock. Stainless steel.

23-1139SS	Size 00 Trachea Tube
23-1140SS	Size 0 Trachea Tube
23-1142SS	Size 1 Trachea Tube
23-1144SS	Size 2 Trachea Tube
23-1146SS	Size 3 Trachea Tube
23-1148SS	Size 4 Trachea Tube
23-1150SS	Size 5 Trachea Tube
23-1152SS	Size 6 Trachea Tube
23-1154SS	Size 7 Trachea Tube
23-1156SS	Size 8 Trachea Tube
23-1158SS	Size 9 Trachea Tube
23-1159SS	Size 10 Trachea Tube

**MILTEX BARRAQUER EYE SPECULUM**

18-35	1 3/4" Eye Speculum, Large Blades 15mm Wide
18-36	1 3/4" Eye Speculum, Small Blades 10mm Wide

**MILTEX WILLIAMS EYE SPECULUM**

18-22	3 1/4" Williams Eye Speculum, Large Blades 15 x 5mm
18-24	2 3/4" Williams Eye Speculum, Small Blades 11 x 5mm

**MILTEX SAUER EYE SPECULUM**

Infant size. 10mm blades.

18-47	1 1/4" Sauer Eye Speculum
-------	---------------------------

**MILTEX PLAIN WIRE EYE SPECULUM**

Blades 15mm wide.

18-34	2 1/2" Wire Eye Speculum
-------	--------------------------

**MILTEX KNAPP EYE SPECULUM**

Blades 15 x 6mm.

18-4	3 1/4" Knapp Eye Speculum
------	---------------------------

**MILTEX MURDOCK EYE SPECULUM**

Blades 15mm wide. Chrome.

18-50 2 1/4" Murdock Eye Speculum

**MILTEX COOK EYE SPECULUM**

18-63 1 7/8" Cook Eye Speculum, Infant Size, 9mm Blades

18-64 1 7/8" Cook Eye Speculum, Child Size, 11mm Blades



**MILTEX DESMARRES CHALAZION FORCEPS**

18-1200 3 1/2" Desmarres Chalazion Forceps, Small, Inside Ring Measures 11 x 17mm

18-1202 3 1/2" Desmarres Chalazion Forceps, Medium, Inside Ring Measures 12 x 23mm

18-1204 3 1/2" Desmarres Chalazion Forceps, Large, Inside Ring Measures 17 x 28mm

**MILTEX SPENCER CHALAZION FORCEPS**

Inside ring diameter 10mm with slide lock.

18-1188 4" Spencer Chalazion Forceps

**MILTEX CASTROVIEJO EYE SPECULUM**

18-10 4" Castroviejo Eye Speculum, Large Blades 18 x 5mm O.D.


18-12 3 3/4" Castroviejo Eye Speculum, Medium Blades 15 x 5mm O.D.



**MILTEX EYE MAGNET & FOREIGN BODY LOOP WITH MAGNIFIER**

Adjustable magnifier attached. Reversible plastic handle, plated metal parts. Plastic case.

18-398M Eye Magnet/Loop w/magnifier



**MILTEX EYE SUTURE SCISSORS**

Slightly curved with sharp pointed short blades.

18-1652 3 3/4" Eye Suture Scissors, Gradle

18-1653 4" Eye Suture Scissors, Ribbon Type

**MILTEX VANNAS CAPSULOTOMY SCISSORS**

Very delicate sharp points.

18-1620 3 1/4" Capsulotomy Scissors, Straight

18-1622 3 1/4" Capsulotomy Scissors, Curved

18-1623 3 1/4" Capsulotomy Scissors, Angled To Side

18-1624 3 1/4" Capsulotomy Scissors, Angled On Flat



N407005



N407010



N407015



N407020

**COURIER® INSTRUMENTS IRIS SCISSORS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407005 4" Straight Iris Scissors

N407010 4" Curved Iris Scissors

N407015 4 1/2" Straight Iris Scissors

N407020 4 1/2" Curved Iris Scissors



### MILTEX IRIS SCISSORS

18-1396	3 1/2" Straight Iris Scissors, Delicate w/20mm blades
18-1398	3 1/2" Curved Iris Scissors, Delicate w/20mm blades
18-1404	4" Straight Iris Scissors, Heavy Pattern
18-1405	4" Curved Iris Scissors, Heavy Pattern
18-1415	4" Straight Iris Scissors, Standard, Ribbon Type
18-1416	4" Curved Iris Scissors, Standard, Ribbon Type
18-1417	3 3/4" Straight Iris Scissors, Miniature Pattern w/14mm blades, Ribbon Type
18-1418	3 3/4" Curved Iris Scissors, Miniature Pattern w/14mm blades, Ribbon Type
5-300	4" Straight Iris Scissors
5-302	4" Curved Iris Scissors
5-304	4 1/2" Straight Iris Scissors
5-304TC	4 1/2" Straight Iris Scissors, Carb-N-Sert
5-306	4 1/2" Curved Iris Scissors
5-306TC	4 1/2" Curved Iris Scissors, Carb-N-Sert
5-308	4 1/2" Iris Scissors, Angled On Side
5-SC-304	4 1/2" Straight Iris Scissors, SuperCut Blade
5-SC-304TC	4 1/2" Straight Iris Scissors, SuperCut Carb-N-Sert Blades
5-SC-306	4 1/2" Curved Iris Scissors, SuperCut Blade
5-SC-306TC	4 1/2" Curved Iris Scissors, SuperCut Carb-N-Sert Blades
5D-303	4 1/2" Curved Iris Scissors (curved sideways)



### MILTEX IRIS TISSUE FORCEPS

1 x 2 teeth.

18-785	4" Tissue Forceps, Straight, Extra Delicate, 0.5mm Wide Tips
18-786	4" Tissue Forceps, Straight, Standard, 0.8mm Wide Tips
18-787	4" Tissue Forceps, Straight, Non Magnetic, Standard, 0.8mm Wide Tips
18-788	4" Tissue Forceps, Half Curved, Standard, 0.8mm Wide Tips
18-789	4" Tissue Forceps, Half Curved, Extra Delicate, 0.5mm Wide Tips
18-790	4" Tissue Forceps, Full Curved, Standard, 0.8mm Wide Tips
18-791	4" Tissue Forceps, Full Curved, Extra Delicate, 0.5mm Wide Tips

### MILTEX WECKER IRIS SPATULA

Sterling blade 30 x 2mm.

18-584	5" Wecker Iris Spatula
--------	------------------------

### MILTEX SMITH-FISHER IRIS SPATULA

One flat end 2mm wide, one spoon end 2mm. Sterling.

18-566	5" Smith-Fisher Iris Spatula
--------	------------------------------

### MILTEX BERGH CILIA FORCEPS

18-1102	3 1/2" Cilia Forceps
18-1108	4 1/4" Cilia Forceps, 10mm Wide Jaws
18-1110	3 1/2" Cilia Forceps, 5mm Wide Jaws

### MILTEX PIERSE COLIBRI TYPE TISSUE FORCEPS

18-880	3 1/4" Tissue Forceps, Tips 0.12mm
18-881	3 1/4" Tissue Forceps, Tips 0.3mm

### MILTEX CASTROVIEJO SUTURING FORCEPS

With tying platform.

18-950	4" Suturing Forceps, 1 x 2 Teeth 0.9mm
18-951	4" Suturing Forceps, 1 x 2 Teeth 0.3mm
18-952	4" Suturing Forceps, 1 x 2 Teeth 0.6mm
18-953	4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.12mm
18-954	4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.6mm
18-955	4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.9mm

### MILTEX RHOTON MICRO FORCEPS

Straight.

17-1100	7" Rhoton Micro Forceps, Suture w/0.7mm tying platform
17-1150	7" Rhoton Micro Forceps, 0.3mm Tips

### MILTEX STRABISMUS TISSUE FORCEPS

1mm wide tips with 1 x 2 teeth.

18-840	4" Tissue Forceps
--------	-------------------



### MILTEX SWISS CILIA & SUTURE

3 3/4" forceps. Excellent for suture removal.

18-1107	Swiss Cilia & Suture
---------	----------------------

### MILTEX WALDEAU FIXATION FORCEPS

Concave jaws 5mm wide with fine teeth.

18-860	4 3/8" Waldeau Fixation Forceps
--------	---------------------------------

### MILTEX ZIEGLER CILIA FORCEPS


2mm wide smooth jaws.

18-1104	3 1/2" Cilia Forceps
---------	----------------------

**MILTEX SNELLEN ENTROPIUM FORCEPS**

35mm wide plate.

- 18-1160 3 1/2" Snellen Entropium Forceps, Left
- 18-1162 3 1/2" Snellen Entropium Forceps, Right



**MILTEX STEVENS TENOTOMY SCISSORS**

- 18-1460 4 1/8" Straight Stevens Tenotomy Scissors, Short Blades, Sharp Points
- 18-1462 4 1/8" Straight Stevens Tenotomy Scissors, Short Blades, Blunt Points
- 18-1464 4 1/8" Curved Stevens Tenotomy Scissors, Short Blades, Sharp Points
- 18-1466 4 1/8" Curved Stevens Tenotomy Scissors, Short Blades, Blunt Points
- 18-1470 4 1/2" Straight Stevens Tenotomy Scissors, Long Blades, Sharp Points
- 18-1472 4 1/2" Straight Stevens Tenotomy Scissors, Long Blades, Blunt Points
- 18-1474 4 1/2" Curved Stevens Tenotomy Scissors, Long Blades, Sharp Points
- 18-1476 4 1/2" Curved Stevens Tenotomy Scissors, Long Blades, Blunt Points
- 18-1477 3 3/4" Straight Stevens Tenotomy Scissors, Blunt Points, Ribbon Type
- 18-1478 3 3/4" Curved Stevens Tenotomy Scissors, Blunt Points Ribbon Type
- 18-SC-1472 4 1/2" Straight Stevens Tenotomy Scissors, SuperCut Blade
- 18-SC-1476 4 1/2" Curved Stevens Tenotomy Scissors, SuperCut Blade

**MILTEX WESTCOTT TENOTOMY SCISSORS**


Wide handles.

- 18-1487 5 1/4" Tenotomy Scissors, Right
- 18-1488 5 1/4" Tenotomy Scissors, Left



**MILTEX NOYES IRIS SCISSORS**

- 18-1510 4 1/2" Iris Scissors, Straight, Sharp-Sharp Points
- 18-1512 4 1/2" Iris Scissors, Straight, Sharp-Blunt Points
- 18-1514 4 1/2" Iris Scissors, Straight, Blunt-Blunt Points
- 18-1516 4 1/2" Iris Scissors, Curved, Sharp-Sharp Points




**MILTEX STRABISMUS SCISSORS**

- 18-1453 4 1/4" Strabismus Scissors, Straight, Ribbon Type
- 18-1454 4 1/4" Strabismus Scissors, Curved, Ribbon Type
- 5-312 4" Strabismus Scissors, Straight
- 5-312TC 4 1/2" Strabismus Scissors, Straight, Carb-N-Sert
- 5-314 4" Strabismus Scissors, Curved
- 5-314TC 4 1/2" Strabismus Scissors, Curved, Carb-N-Sert

**MILTEX WESTCOTT UTILITY SCISSORS**

Curved and blunt tips.

- 18-1480 4 1/2" Utility Scissors



**MILTEX TONOMETERS**

Autoclavable Schiottz tonometers made with jeweled hammers to minimize friction. Scales are included, permitting reading from above. Scale markings etched for accuracy. Supplied in case with weights and test block.

- 18-132 Tonometer w/3 weights (5.5, 7.5 & 10 gr)
- 18-132A Same as No. 18-132, w/Allen plunger retractor
- 18-139 Sterilizing Case For Tonometers, 6 1/2" x 3 1/2" x 1 1/4", Fitted To Firmly Hold Tonometer, Weights & Test Block



**MILTEX VON GRAEFE CATARACT KNIFE**

- 18-140 No. 1 Von Graefe Cataract Knife, 1.5 x 25mm
- 18-142 No. 2 Von Graefe Cataract Knife, 1.7 x 27mm
- 18-144 No. 3 Von Graefe Cataract Knife, 2 x 30mm
- 18-146 No. 4 Von Graefe Cataract Knife, 2.3 x 33mm
- 18-148 No. 5 Von Graefe Cataract Knife, 2.5 x 35mm

**MILTEX ZIEGLER KNIFE-NEEDLE**

- 18-259 4 1/2" Ziegler Knife-Needle, Size 0, 3mm Blade
- 18-260 4 1/2" Ziegler Knife-Needle, Size 1, 5mm Blade
- 18-262 4 1/2" Ziegler Knife-Needle, Size 2, 6mm Blade
- 18-264 4 1/2" Ziegler Knife-Needle, Size 3, 7mm Blade

**MILTEX WHEELER DISCISSION KNIFE**

- 18-152 1 x 18mm Wheeler Discission Knife

**MILTEX WEST LACRIMAL SAC CHISEL & GOUGE**

- 18-1962 Chisel, 6 1/4", Straight Blade 5mm Wide
- 18-1964 Gouge, 6 1/4", Straight Blade 4.5mm Wide

**MILTEX WILDER LACRIMAL DILATORS**

- 18-694 4" Wilder Lacrimal Dilator, Fine Taper
- 18-696 4" Wilder Lacrimal Dilator, Medium Taper
- 18-698 4" Wilder Lacrimal Dilator, Heavy Taper

**MILTEX CASTROVIEJO LACRIMAL DILATOR**

One end needle point and other end medium.

- 18-691 5 1/2" Castroviejo Lacrimal Dilator

**MILTEX ZIEGLER LACRIMAL DILATOR**

Double end.

- 18-680 5 1/2" Ziegler Lacrimal Dilator


**MILTEX WILDER LENS LOOP**

- 18-484 5 1/4" Wilder Lens Loop, Serrated 4 x 15mm
- 18-486 5 1/4" Wilder Lens Loop, Smooth 4 x 15mm


**MILTEX WILLIAMS LACRIMAL PROBES**

5" sterling.

- 18-722 Size 0000-000 Lacrimal Probe
- 18-724 Size 00-0 Lacrimal Probe
- 18-726 Size 1-2 Lacrimal Probe
- 18-728 Size 3-4 Lacrimal Probe
- 18-730 Size 5-6 Lacrimal Probe
- 18-732 Size 7-8 Lacrimal Probe



18-374



18-400

**MILTEX DIX FOREIGN BODY NEEDLE & SPUD**

- 18-374 Spud, 4 3/4", 1.7 x 15mm Blade
- 18-400 Needle & Spud, Protected In Sliding Handle
- 18-402 Needle & Spud, Protected In Reversible Screw Handle
- 18-402A Needle w/screw socket
- 18-402B Spud w/screw socket
- 18-404 Needle & Gouge, Protected In Reversible Screw Handle



**MILTEX CASTROVIEJO SPATULA**

Double end 3.5mm wide angled and 2mm wide slightly curved. Sterling.

- 18-564 5 1/2" Castroviejo Spatula

**MILTEX WALTON FOREIGN BODY GOUGE**

Round end.

- 18-376 4 1/2" Walton Foreign Body Gouge

**MILTEX WALTER CORNEAL SPUD**

- 18-302 4 5/8" Walter Corneal Spud, Straight Blade 4mm
- 18-304 4 5/8" Walter Corneal Spud, Curved Blade 4mm



**MILTEX DESMARRS CORNEAL DISSECTOR**

3 x 8mm blade.


- 18-331 4 3/4" Desmarres Corneal Dissector



**MILTEX VON GRAEFE (BECKER) CYSTOTOME**

1mm blade. Right angle. Malleable shaft.

- 18-314 4 3/4" Von Graefe Cystotome



**MILTEX VON GRAEFE HOOKS**

- 18-450 5 1/2" Von Graefe Hook, Small, 8mm Deep
- 18-452 5 1/2" Von Graefe Hook, Medium, 10mm Deep
- 18-454 5 1/2" Von Graefe Hook, Large, 11mm Deep

**MILTEX PENCIL-TYPE MICRO SCISSORS/FORCEPS**

7 1/4" overall length. Extremely delicate working ends.

- 17-100 1mm Wide Serrated Jaws, Straight Shaft Forceps
- 17-115 4mm Long Blades, Straight Shaft Forceps
- 17-215 4mm Long Blades, Angled Shaft Forceps
- 17-220 8mm Long Blades, Angled Shaft Forceps

**MILTEX DITTEL URETHRAL SOUNDS**

11 1/2". Straight.

- 29-8-10 10 Fr Urethral Sound
- 29-8-12 12 Fr Urethral Sound
- 29-8-14 14 Fr Urethral Sound
- 29-8-16 16 Fr Urethral Sound
- 29-8-18 18 Fr Urethral Sound
- 29-8-20 20 Fr Urethral Sound
- 29-8-22 22 Fr Urethral Sound
- 29-8-24 24 Fr Urethral Sound
- 29-8-26 26 Fr Urethral Sound
- 29-8-28 28 Fr Urethral Sound
- 29-8-30 30 Fr Urethral Sound
- 29-8-32 32 Fr Urethral Sound
- 29-8-34 34 Fr Urethral Sound
- 29-8-8 8 Fr Urethral Sound

**MILTEX LEFORT URETHRAL SOUNDS**

11" with standard 1-56 thread for filiform (not included).

- 29-10-10 10 Fr Urethral Sound
- 29-10-12 12 Fr Urethral Sound
- 29-10-14 14 Fr Urethral Sound
- 29-10-16 16 Fr Urethral Sound
- 29-10-18 18 Fr Urethral Sound
- 29-10-20 20 Fr Urethral Sound
- 29-10-22 22 Fr Urethral Sound
- 29-10-24 24 Fr Urethral Sound
- 29-10-26 26 Fr Urethral Sound
- 29-10-28 28 Fr Urethral Sound
- 29-10-30 30 Fr Urethral Sound
- 29-10-8 8 Fr Urethral Sound

**MILTEX RANDALL KIDNEY STONE FORCEPS**

- 29-284 9 1/4" Kidney Stone Forceps, Quarter Curved
- 29-286 8 1/4" Kidney Stone Forceps, Half Curved
- 29-288 7 3/4" Kidney Stone Forceps, Three Quarter Curved
- 29-290 7 1/2" Kidney Stone Forceps, Full Curved



**MILTEX WALTHER FEMALE DILATOR-CATHETERS**

5 1/4".

- 29-33-12 12 Fr Female Dilator-Catheter
- 29-33-14 14 Fr Female Dilator-Catheter
- 29-33-16 16 Fr Female Dilator-Catheter
- 29-33-18 18 Fr Female Dilator-Catheter
- 29-33-20 20 Fr Female Dilator-Catheter
- 29-33-22 22 Fr Female Dilator-Catheter
- 29-33-24 24 Fr Female Dilator-Catheter
- 29-33-26 26 Fr Female Dilator-Catheter
- 29-33-28 28 Fr Female Dilator-Catheter
- 29-33-30 30 Fr Female Dilator-Catheter
- 29-33-32 32 Fr Female Dilator-Catheter
- 29-33-34 34 Fr Female Dilator-Catheter
- 29-33-36 36 Fr Female Dilator-Catheter
- 29-33-38 38 Fr Female Dilator-Catheter

**MILTEX WALTHER URETHRAL SOUNDS**

11". Shaft tapered downward toward handle.

- 29-18-16 16 Fr Urethral Sound
- 29-18-24 24 Fr Urethral Sound
- 29-18-18 18 Fr Urethral Sound
- 29-18-20 20 Fr Urethral Sound
- 29-18-22 22 Fr Urethral Sound
- 29-18-26 26 Fr Urethral Sound
- 29-18-28 28 Fr Urethral Sound
- 29-18-30 30 Fr Urethral Sound
- 29-18-32 32 Fr Urethral Sound
- 29-18-34 34 Fr Urethral Sound
- 29-18-36 36 Fr Urethral Sound

**MILTEX YOUNG RENAL PEDICLE CLAMP**

Longitudinal serrations with cross serrations at tip.

- 29-300 9" Renal Pedicle Clamp, Curved

**MILTEX NO SCALPEL VASECTOMY INSTRUMENTS**

- 29-805 No Scalpel Vasectomy Ring Clamp, 3.5mm ID, 5 1/2"
- 29-850 No Scalpel Vasectomy Hemostat, Sharp Point, Smooth Jaws, Curved, 5 1/2"
- 29-855 No Scalpel Vasectomy Instrument Set

**RUSCH NO SCALPEL VASECTOMY DISSECTING FORCEPS**


Sharp pointed hemostat is designed to puncture the scrotal sac, pierce the vas, and deliver a portion of vas from the sac where it is clipped and divided. Stainless steel. Smooth jaws. Gold handles.

- 244000 5 1/2" Dissecting Forceps, Curved

**RUSCH NO SCALPEL VASECTOMY FIXATION RING CLAMP**

Used to isolate and immobilize the vas within the scrotum and then secure it prior to clipping. Stainless steel. 3.5mm I.D. Gold handles. Non-sterile.

- 244001 5 1/2" Vasectomy Fixation Ring Clamp



**COOPER SURGICAL SANI-SPEC® SPECULA**

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Provides protection of the doctor, the nurse and the patient from cross contamination. Plastic material avoids unpleasant temperature shock associated with metal. Smooth pre-lubricated blades with 6 position grooves.

- 82329 Small Specula, 1"W x 4 1/2"L, 10/bg, 10 bg/cs
- 82330 Medium Specula, 1 1/4"W x 4 3/4"L, 10/bg, 10 bg/cs
- 82331 Large Specula, 1 1/3"W x 5 1/3"L, 10/bg, 10 bg/cs
- 82340 Small Specula, Individually Wrapped, 50/cs
- 82341 Medium Specula, Individually Wrapped, 50/cs

**COOPER SURGICAL GALENICA GYNE-LITE™ SYSTEM ACCESSORIES**

**Accessories**

- 9220 Adaptor, 10/bg
- 9230 Sheaths, 100/bx, 4 bx/cs
- 9240 Bulb



**STARLINE VAGINAL SPECULA**

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength.

- 16108309 Small Vaginal Specula, 10/bg, 10 bg/cs
- 16108310 Medium Vaginal Specula, 10/bg, 10 bg/cs
- 16108311 Large Vaginal Specula, 10/bg, 10 bg/cs



**COOPER SURGICAL GALENICA VAGINAL SPECULA**

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength. Unit Pack (poly bag) non-sterile. Unit Pack (peeling pouch) sterile.

- 8309 Small Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8310 Medium Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8311 Large Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8615 Small Unit Pack Vaginal Specula, 50/cs
- 8616 Medium Unit Pack Vaginal Specula, 50/cs
- 8617 Large Unit Pack Vaginal Specula, 50/cs

**Accessories**

- 9310 Smoke Evacuator For All Specula, 40/cs



N407400



N407405

**COURIER® INSTRUMENTS GRAVES VAGINAL SPECULA**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407400 Small Vaginal Specula
- N407405 Medium Vaginal Specula
- N407410 Large Vaginal Specula



30-10



30-10B

**MILTEX GRAVES VAGINAL SPECULA**

- 30-10 Small Vaginal Specula, 3/4" x 3", Original Model
- 30-10B Small Vaginal Specula, 3/4" x 3", Ebony Finish
- 30-15 Medium Vaginal Specula, 1 3/8" x 4", Original Model
- 30-15B Medium Vaginal Specula, 1 3/8" x 4", Ebony Finish
- 30-20 Large Vaginal Specula, 1 1/2" x 4 1/2", Original Model
- 30-20B Large Vaginal Specula, 1 1/2" x 4 1/2", Ebony Finish
- 30-25 Medium Vaginal Specula, 1 3/8" x 4", Wide Angle Blades, Improved
- 30-30 Medium Vaginal Specula, 1 3/8" x 4", Wide Angle Blades, Open Side

**VANTAGE GRAVES VAGINAL SPECULUM**

- V930-10 Small Vaginal Speculum
- V930-15 Medium Vaginal Speculum
- V930-15A Screw Nuts For Speculum
- V930-15X Center Nut For Speculum
- V930-20 Large Vaginal Speculum
- V930-22 X-Large Vaginal Speculum



**COURIER® INSTRUMENTS PEDERSON VAGINAL SPECULA**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407415 Small Vaginal Specula
- N407420 Medium Vaginal Specula
- N407425 Large Vaginal Specula



**MILTEX PEDERSON VAGINAL SPECULA**

Narrow blades.

- 30-50 Small Vaginal Specula, 1/2" x 3"
- 30-55 Medium Vaginal Specula, 7/8" x 4"
- 30-60 Large Vaginal Specula, 1" x 4 3/4"
- 30-62 Extra Narrow Vaginal Specula, 5/8" x 4"

**MILTEX LEEP PEDERSON SPECULUM**

Has narrow blades for patients with a smaller vaginal opening. The integrated smoke evacuation tube provides excellent smoke reduction.

- 31-551 Medium Pederson Speculum
- 31-601 Large Pederson Speculum
- 31-611 X-Large Pederson Speculum


**MILTEX LEEP VIEW-MORE SPECULUM**

The wider 1 1/2" yoke area on this speculum provides a better view of the cervix and allows for easier access with other instruments. The integrated smoke evacuation tube provides excellent smoke reduction.

- 31-66 Medium View-More Speculum
- 31-67 Large View-More Speculum

**VANTAGE PEDERSON VAGINAL SPECULUM**

- V930-50 Small Vaginal Speculum
- V930-55 Medium Vaginal Speculum
- V930-60 Large Vaginal Speculum
- V930-61 X-Large Vaginal Speculum



**MILTEX LATERAL VAGINAL RETRACTOR**

Thumb ratchet retracts up to 2 1/4" for extra wide viewing area. Fits inside any vaginal speculum.

- 30-1355 8 1/4" Lateral Vaginal Retractor

**MILTEX SIMS VAGINAL SPECULUM**

- 30-174 Single End Vaginal Specula, Small 1 1/16" x 2 1/2"
- 30-176 Single End Vaginal Specula, Medium 1 1/4" x 3"
- 30-178 Single End Vaginal Specula, Large 1 1/2" x 3 1/2"
- 30-210 Double End Vaginal Specula, Small 1" x 2 1/2" & 1 1/4" x 3"
- 30-215 Double End Vaginal Specula, Medium 1 1/8" x 3" & 1 1/4" x 3 1/2"
- 30-220 Double End Vaginal Specula, Large 1 3/8" x 3 1/2" & 1 1/2" x 4"

**MILTEX COLLIN VAGINAL SPECULUM**

- 30-120 Small Vaginal Specula, 1 1/4" x 3 3/4"
- 30-125 Medium Vaginal Specula, 1 3/8" x 4"
- 30-130 Large Vaginal Specula, 1 5/8" x 4 1/2"

**MILTEX CUSCO VAGINAL SPECULA**

Folding handles.

- 30-100 Large Vaginal Specula, 1 3/8" x 4 1/2" w/side screw
- 30-70 Extra Narrow Vaginal Specula, 3/4" x 3 1/4" w/side screw
- 30-75 Small Vaginal Specula, 1 1/4" x 3 1/2" w/center screw
- 30-80 Medium Vaginal Specula, 1 3/8" x 4" w/center screw
- 30-85 Large Vaginal Specula, 1 3/8" x 4 1/2" w/center screw
- 30-90 Small Vaginal Specula, 1 1/4" x 3 1/2" w/side screw
- 30-95 Medium Vaginal Specula, 1 3/8" x 4" w/side screw

**MILTEX DEVILBISS VAGINAL SPECULA**

Two sets of screws permit parallel or bivalve opening of the blades for examination or operative use.

- 30-140 Small Vaginal Specula, 3 1/2" x 7/8"
- 30-141 Medium Vaginal Specula, 4" x 1 1/4"
- 30-142 Large Vaginal Specula, 5" x 1 3/8"

**MILTEX AUVARD WEIGHTED VAGINAL SPECULA**

9".

- 30-186 2 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
- 30-187 2.5 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
- 30-188 3 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
- 30-193 2.5 lbs Weighted Vaginal Specula, Blade 4" x 1 3/8", Tapering To 1 3/4" At Tip, Blade At Acute Angle 45°
- 30-195 2.5 lbs Weighted Vaginal Specula, Steiner-Auvard, Slightly Curved Blade 5 1/2" x 1 1/4" At 90° Angle

**MILTEX PHANEUF FORCEPS**

Horizontal serrated jaws. 1 x 2 teeth.

- 30-1735 8" Phaneuf Forceps, Straight
- 30-1740 8" Phaneuf Forceps, Angled On Flat

**MILTEX BOZEMAN SPONGE FORCEPS**


10 1/2".

- 7-622 10 1/2" Double Curved Sponge Forceps



**MILTEX BOZEMAN UTERINE DRESSING FORCEPS**


- 7-620 10 1/2" Straight Dressing Forceps
- 7-624 Double Curved Dressing Forceps w/1 large finger ring



**COURIER® INSTRUMENTS FOERSTER FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407225 9 1/2" Straight Serrated Foerster Forceps



**MILTEX SCHROEDER UTERINE VULSELLUM FORCEPS**

30-1000 10" Uterine Vulsellum Forceps, Straight  
 30-1005 10" Uterine Vulsellum Forceps, Curved On Flat  
 30-1010 10" Uterine Vulsellum Forceps, Curved Sideways


**MILTEX WINTER PLACENTA FORCEPS**

30-2530 11" Placenta Forceps, Cup Jaws, Straight  
 30-2535 11" Placenta Forceps, Cup Jaws, Curved

**MILTEX BRAUN EPISIOTOMY SCISSORS**

Angled to side. Guarded lower blade.

30-2190 5 1/2" Episiotomy Scissors  
 30-2195 8 1/2" Episiotomy Scissors



**PRECISION DYNAMICS SECURLINE™ AMNIOTIC MEMBRANE PERFORATOR**

Designed to snap-not poke or prick membrane. Overall length of 10 1/2" provides ample reach for practitioner. Perforator is sterile. 100 individual pouches are placed in a convenient dispenser box that contains a level window that allows one to check product levels at a glance.

3550-00-PDG Amniotic Membrane Perforator, Sterile, 100/bx


**MILTEX BARRETT UTERINE TENACULUM FORCEPS**

30-985 7 1/4" Uterine Tenaculum Forceps, Straight



**MILTEX HEGAR UTERINE DILATORS**

- 30-515 Uterine Dilator, 7", Single End, Set of 14: Sizes 4 to 17 mm In Khaki Roll
- 30-515A Khaki Roll Only For Set of 14 Uterine Dilator
- 30-520 Uterine Dilator, 7", Single End, Set of 26: Sizes 1 to 26 mm In Khaki Roll
- 30-520A Khaki Roll Only For Set of 26 Uterine Dilator
- 30-530-1 1 mm Uterine Dilator, 7", Single End
- 30-530-10 10 mm Uterine Dilator, 7", Single End
- 30-530-105 10.5 mm Uterine Dilator, 7", Single End
- 30-530-11 11 mm Uterine Dilator, 7", Single End
- 30-530-115 11.5 mm Uterine Dilator, 7", Single End
- 30-530-12 12 mm Uterine Dilator, 7", Single End
- 30-530-125 12.5 mm Uterine Dilator, 7", Single End
- 30-530-13 13 mm Uterine Dilator, 7", Single End
- 30-530-135 13.5 mm Uterine Dilator, 7", Single End
- 30-530-14 14 mm Uterine Dilator, 7", Single End
- 30-530-145 14.5 mm Uterine Dilator, 7", Single End
- 30-530-15 1.5 mm Uterine Dilator, 7", Single End
- 30-530-150 15 mm Uterine Dilator, 7", Single End
- 30-530-155 15.5 mm Uterine Dilator, 7", Single End
- 30-530-16 16 mm Uterine Dilator, 7", Single End
- 30-530-165 16.5 mm Uterine Dilator, 7", Single End
- 30-530-17 17 mm Uterine Dilator, 7", Single End
- 30-530-175 17.5 mm Uterine Dilator, 7", Single End
- 30-530-18 18 mm Uterine Dilator, 7", Single End
- 30-530-185 18.5 mm Uterine Dilator, 7", Single End
- 30-530-19 19 mm Uterine Dilator, 7", Single End
- 30-530-2 2 mm Uterine Dilator, 7", Single End
- 30-530-20 20 mm Uterine Dilator, 7", Single End
- 30-530-21 21 mm Uterine Dilator, 7", Single End
- 30-530-22 22 mm Uterine Dilator, 7", Single End
- 30-530-23 23 mm Uterine Dilator, 7", Single End
- 30-530-24 24 mm Uterine Dilator, 7", Single End
- 30-530-25 2.5 mm Uterine Dilator, 7", Single End
- 30-530-250 25 mm Uterine Dilator, 7", Single End
- 30-530-26 26 mm Uterine Dilator, 7", Single End
- 30-530-3 3 mm Uterine Dilator, 7", Single End
- 30-530-35 3.5 mm Uterine Dilator, 7", Single End
- 30-530-4 4 mm Uterine Dilator, 7", Single End
- 30-530-45 4.5 mm Uterine Dilator, 7", Single End
- 30-530-5 5 mm Uterine Dilator, 7", Single End
- 30-530-55 5.5 mm Uterine Dilator, 7", Single End
- 30-530-6 6 mm Uterine Dilator, 7", Single End
- 30-530-65 6.5 mm Uterine Dilator, 7", Single End
- 30-530-7 7 mm Uterine Dilator, 7", Single End
- 30-530-75 7.5 mm Uterine Dilator, 7", Single End
- 30-530-8 8 mm Uterine Dilator, 7", Single End
- 30-530-85 8.5 mm Uterine Dilator, 7", Single End
- 30-530-9 9 mm Uterine Dilator, 7", Single End
- 30-530-95 9.5 mm Uterine Dilator, 7", Single End
- 30-550 Uterine Dilators, 7 1/2", Double End, Set of 8: 3-5 mm & 17-18 mm In Khaki Roll
- 30-550A Khaki Roll Only For Set of 8 Uterine Dilator
- 30-555-1112 11-12 mm Uterine Dilator, 7 1/2", Double End
- 30-555-1314 13-14 mm Uterine Dilator, 7 1/2", Double End
- 30-555-1516 15-16 mm Uterine Dilator, 7 1/2", Double End
- 30-555-1718 17-18 mm Uterine Dilator, 7 1/2", Double End
- 30-555-34 3-4 mm Uterine Dilator, 7 1/2", Double End
- 30-555-56 5-6 mm Uterine Dilator, 7 1/2", Double End
- 30-555-78 7-8 mm Uterine Dilator, 7 1/2", Double End
- 30-555-910 9-10 mm Uterine Dilator, 7 1/2", Double End



**MILTEX PRATT UTERINE DILATORS**

11-1/2". Solid. Double ended.

- 30-560-1315 13-15 Fr Pratt Uterine Dilator
- 30-560-1719 17-19 Fr Pratt Uterine Dilator
- 30-560-2123 21-23 Fr Pratt Uterine Dilator
- 30-560-2527 25-27 Fr Pratt Uterine Dilator
- 30-560-2931 29-31 Fr Pratt Uterine Dilator
- 30-560-3335 33-35 Fr Pratt Uterine Dilator
- 30-560-3739 37-39 Fr Pratt Uterine Dilator
- 30-560-4143 41-43 Fr Pratt Uterine Dilator
- 30-561 Set of 8 Pratt Uterine Dilator, sizes: 13-15 Fr to 41-43 Fr in a Khaki Roll

**Accessories**

- 30-561A Khaki Roll Only For Set of 8



**MILTEX SCHROEDER-BRAUN UTERINE TENACULUM FORCEPS**

30-965 9 3/4" Uterine Tenaculum Forceps, Straight



**MILTEX SIMS UTERINE SOUND**

13". Malleable. Silver plated.

- 30-650 Sims Uterine Sound Graduated in Centimeters
- 30-655 Sims Uterine Sound Graduated in Inches



**MILTEX EMMETT UTERINE TENACULUM HOOKS**

9". Sharp points.

- 30-951 Style 1 Uterine Tenaculum Hook, Half Curved, Small
- 30-952 Style 2 Uterine Tenaculum Hook, Acute Angle
- 30-953 Style 3 Uterine Tenaculum Hook, Right Angle
- 30-954 Style 4 Uterine Tenaculum Hook, Half-Curved, Large
- 30-955 Style 5 Uterine Tenaculum Hook, Double Angle



N407435



N407440

**COURIER® INSTRUMENTS TISCHLER BIOPSY PUNCH**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.


- N407435 Baby-Tischler Biopsy Punch
- N407440 Tischler Biopsy Punch 3x7mm



**COURIER® INSTRUMENTS KEVORKIAN-YOUNGE BIOPSY FORCEPS**

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.


- N407430 Kevorkian-Youngue Biopsy Forceps



**MILTEX KEVORKIAN-YOUNGE BIOPSY FORCEPS**

Shaft interchangeable on pistol grip spring handle. 3.5 mm x 8 mm bite, teeth on lower jaw. Complete with pistol grip handle. Rotating punch.

- 30-1490 9" Shaft Biopsy Forceps
- 30-1491 Biopsy Forceps Shaft Only



30-1482

**MILTEX KEVORKIAN-YOUNGE UTERINE BIOPSY FORCEPS**

3.5 mm x 8 mm bite, teeth on lower jaw.

30-1482 9 1/2" Shaft Uterine Biopsy Forceps

**MILTEX SCHUBERT UTERINE BIOPSY FORCEPS**

30-1420 8 1/2" Uterine Biopsy Forceps, Straight

30-1425 8 1/2" Uterine Biopsy Forceps, Curved

30-1430 11" Uterine Biopsy Forceps, Angled Shanks

30-1435 11" Uterine Biopsy Forceps, S-Shape, Triangular Bite

**MILTEX EPPENDORFER UTERINE BIOPSY FORCEPS**

30-1485 8 1/2" Shaft Uterine Biopsy Forceps

**MILTEX SOMER UTERINE ELEVATING FORCEPS**


Curved jaws.

30-1190 9" Uterine Elevating Forceps

**MILTEX TISCHLER BIOPSY FORCEPS**

30-1494 3 x 7mm Bite Biopsy Forceps, Complete w/pistol grip handle

30-1495 Shaft Only For Biopsy Forceps



30-1492

**MILTEX TISCHLER CERVICAL BIOPSY PUNCH FORCEPS**

30-1442 8 1/4" Shaft (3 x 7 mm bite) Biopsy Punch Forceps


30-1443 Burke Modified Baby Tischler Biopsy Punch Forceps (4.5 x 3 x 1.5 mm bite)

**MILTEX WITTNER UTERINE BIOPSY FORCEPS**

8 1/2" shaft with teeth on lower jaw.

30-1486 Standard Model Uterine Biopsy Forceps w/angled tip

30-1487 Straight Tip Uterine Biopsy Forceps



30-1490  
30-1491

30-1492  
30-1493

30-1490 to 1495  
25-400

**MILTEX EPPENDORFER-KRAUSE BIOPSY FORCEPS**

Shaft interchangeable on pistol grip spring handle. Oval 4.5 mm x 6.5 mm bite. Complete with pistol grip handle. Rotating punch.

30-1492 9" Shaft Biopsy Forceps

30-1493 Biopsy Forceps Shaft Only

**OWEN MUMFORD AMIELLE CARE VAGINAL DILATORS**

Following radiotherapy to the vagina, cervix or pelvis, the woman may experience some internal scarring (adhesions) to the treated area, which can make the vagina narrower, drier and less elastic. Preventing the formation of scar tissue within the vagina helps ensure that it remains supple so that sexual intercourse and vaginal examinations can be performed without discomfort or pain. Designed in conjunction with the UK National Forum of Gynecological Oncology Nurses, Amielle Care vaginal dilators help to reduce the risk of scarring and offer an effective way for women to take an active role in their treatment.

SM2150 Vaginal Dilator

**OWEN MUMFORD AMIELLE COMFORT VAGINAL DILATORS**

Vaginismus is a condition where the vaginal muscles go into spasm, causing the vagina to feel tight and often making sexual intercourse impossible. This may be due to a number of reasons including fear and pain of penetration, gynecological surgery or trauma associated with childbirth. Dyspareunia can also inhibit sexual intercourse as it may cause pain and discomfort upon penetration. This is often associated with vaginal dryness or infection. Designed in conjunction with psychosexual therapists, Amielle Comfort offers an effective way for women to take an active role in the treatment of these and other conditions.

SM2100 Vaginal Dilator



**MILTEX THOMAS UTERINE CURETTES**

11". Blunt blades on malleable shanks.

30-1225-1 Size 1 Thomas Uterine Curette

30-1225-2 Size 2 Thomas Uterine Curette

30-1225-3 Size 3 Thomas Uterine Curette

30-1225-4 Size 4 Thomas Uterine Curette

30-1225-5 Size 5 Thomas Uterine Curette

30-1225-6 Size 6 Thomas Uterine Curette

**MILTEX SIMS UTERINE BIOPSY CURETTE**

Extra small sharp loop. 2.5 x 5mm outside measurement.

30-1380 10" Uterine Biopsy Curette



**MILTEX SIMS UTERINE CURETTES**


11". Sharp blades on malleable shanks.

- 30-1205-0 Size 0 Sims Uterine Curette
- 30-1205-00 Size 00 Sims Uterine Curette
- 30-1205-1 Size 1 Sims Uterine Curette
- 30-1205-2 Size 2 Sims Uterine Curette
- 30-1205-3 Size 3 Sims Uterine Curette
- 30-1205-4 Size 4 Sims Uterine Curette
- 30-1205-5 Size 5 Sims Uterine Curette
- 30-1205-6 Size 6 Sims Uterine Curette

**MILTEX BUMM UTERINE CURETTES**

12". Extra large sizes on malleable shanks.

- 30-1270-1 Size 1 Uterine Curette, Sharp
- 30-1270-2 Size 2 Uterine Curette, Sharp
- 30-1270-3 Size 3 Uterine Curette, Sharp
- 30-1270-4 Size 4 Uterine Curette, Sharp
- 30-1275-1 Size 1 Uterine Curette, Blunt
- 30-1275-4 Size 4 Uterine Curette, Blunt



**MILTEX KOGAN ENDOCERVICAL SPECULUM**

- 30-1350 9 1/2" Endocervical Speculum w/gauge & set screw, Offset Ring Handles, Fenestrated Jaws 7 x 25mm
- 30-1351 11" Endocervical Speculum, Improved Pattern w/ratchet & gauge, Extra Delicate Fenestrated Jaws 3mm Wide At Tip x 25mm Long



**MILTEX NOVAK ENDOCERVICAL SUCTION BIOPSY CURETTES**

9-3/4", Luer hub.

- 30-1390 4 mm Dia. Biopsy Curette, 12 Fr
- 30-1391 2 mm Dia. Biopsy Curette, Extra Delicate
- 30-1392 3 mm Dia. Biopsy Curette, Delicate

**MILTEX ENDOMETRIAL BIOPSY CANNULA**

Offered in either a SoftFlex™ or FirmFlex™ style for easy one-step sampling. The SoftFlex™ is ideal for standard or routine patient sampling applications where a more flexible cannula is desired by the clinician. The FirmFlex™ is a slightly more rigid cannula ideal for more controlled insertion and sampling in patient cases where atrophic endometrium or post-menopausal patients with stenosis of the cervix is a concern. Both the SoftFlex™ and the FirmFlex™ patterns are sterilized via ETO, packaged individually for single-use applications, and provided in convenient quantities of 25 per box. The cannula are made of a clean, latex-free and flexible polypropylene outer sleeve with inner plunger. The cannula sleeve, with plunger in closed position, is inserted through the cervix and into the uterine cavity. Once properly positioned, the plunger is withdrawn and the design generates sufficient vacuum for maximum endometrial sample yield to be captured within the cannula sleeve without the need for external suction devices. Detailed directions for use are packaged with cannula and should be reviewed thoroughly and carefully prior to each use or demonstration.

- 30-3010 SoftFlex™ Endometrial Biopsy Cannula, Sterile, Single Use, 25/bx
- 30-3012 FirmFlex™ Endometrial Biopsy Cannula, Sterile, Single Use, 25/bx

**MILTEX ENDOMETRIAL SAMPLING SET**

A disposable plastic curette and syringe system which is very similar in function and appearance to the stainless steel Randall endometrial suction biopsy curette pattern. Includes a disposable curette and a 10ml Twist-and-Lock syringe. The curette is made of clear, latex-free polypropylene and has a semi-rigid design with a slightly curved and tapered Randall-type cutting edge tip which optimizes access to the uterine contours and efficient sample collection. The semi-rigid design of the curette allows for controlled insertion in female patients with stenosis of the cervix, and the syringe system creates 21" Hg of negative pressure for optimized suction and sample capture in situations where the endometrium is atrophic and more difficult to capture sufficient sample yield.

- 30-3015 Endometrial Sampling Set, Includes: Disposable Curette w/a Twist-and-Lock Syringe, Sterile, Single Use, 10/bx

**MILTEX RANDALL ENDOMETRIAL SUCTION BIOPSY CURETTE**

Luer hub. 4mm diameter.

- 30-1395 9 3/4" Endometrial Suction Biopsy Curette

**COOPERSURGICAL GALENICA ANOSCOPES**

Easy manipulation for septic anal/rectal examinations. Decreases the risks of cross-infections. 56° bevel at the insertion end. Adaptable to fiber optic lighthead. Made from hard non-toxic plastic resin of fine medical grade. Disposable and lightweight.

- 8607 Anoscope, 25/bx, 4 bx/cs

If you do not see it here, ask your representative.

**MILTEX HIRSCHMANN ANOSCOPES**

- 28-20 Small Anoscope, 2 1/4" Long, 9/16" I.D.  
 28-22 Medium Anoscope, 2 5/8" Long, 11/16" I.D.  
 28-24 Large Anoscope, 3 1/2" Long, 7/8" I.D.

**MILTEX KELLY ANOSCOPE**

- 28-40 2" Long x 7/8" O.D. Anoscope

**COOPERSURGICAL SANI-SCOPE®**

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Allows protection for the doctor, the nurse and the patient from cross contamination. Single use obturator.

- 82420 9/10" Dia x 3 3/8" L Sani-Scope®, 25/bg, 4 bg/cs

**MILTEX BARR-SHUFORD RECTAL SPECULUM**

- 28-52 4 1/2" Tapered From 5/8" To 2", Chrome Rectal Speculum

**MILTEX PENNINGTON RECTAL SPECULUM**

- 28-10 3 3/8" Long x 1" O.D. Rectal Speculum w/boilable obturator,  
 4 3/4" Long Handle, Chrome  
 28-11 Obturator Only

**MILTEX PRATT RECTAL SPECULUM**

1" x 3 1/2" blades, with set screw. Stainless steel.

- 28-1 Pratt Rectal Speculum

**MILTEX SAWYER RECTAL SPECULUM**

- 28-70 Sawyer Rectal Speculum, 11" (27.9 cm), Blade 7/8" (2.2 cm) x 2 1/2" (6.4 cm)  
 28-72 Sawyer Rectal Speculum, 11" (27.9 cm), Blade 1 1/8" (2.8 cm) x 3" (7.6 cm)  
 28-74 Sawyer Rectal Speculum, 11" (27.9 cm), Blade 1 1/2" (3.8 cm) x 3 1/2" (8.9 cm)

**MILTEX SIMS RECTAL SPECULUM**


Fenestrated blades. Stainless steel.

- 28-2 Sims Rectal Speculum

**MILTEX BODENHAMMER RECTAL SPECULUM**

Blades 3 1/2" long x 3/4" wide at distal end.

- 28-4 6" Rectal Speculum w/set screw

**MILTEX BRINKERHOFF RECTAL SPECULUM**

- 28-30 Small Rectal Speculum, 3 1/2" L, 1" Tapering to 1/2"  
 28-32 Medium Rectal Speculum, 4" L, 1" Tapering to 5/8"  
 28-34 Large Rectal Speculum, 4 3/4" L, 1 1/4" Tapering to 3/4"

**MILTEX CHELSEA EATON ANAL SPECULUM**

- 28-79 2 5/8" Long x 7/8" O.D. Anal Speculum  
 28-80 2 5/8" Long x 1 1/4" O.D. Anal Speculum

**MILTEX FANSLER (IVES) RECTAL SPECULUM**

3 1/4" long x 1" largest diameter with plastic obturator (boilable). Chrome-plated.

- 28-12 Fansler (Ives) Rectal Speculum  
 28-14 Plastic Obturator Only

**MILTEX HIRSCHMAN PROCTOSCOPE**

- 28-48 5" Long x 15/16" O.D. Proctoscope

**MILTEX PRATT CRYPT HOOKS**

- 28-120 8 1/2" Pratt Crypt Hook, Narrow Style  
 28-122 8 1/2" Pratt Crypt Hook, Short Style  
 28-124 8 1/2" Pratt Crypt Hook, Long Style

## MILTEX BARR CRYPT HOOK

28-140 10" Barr Crypt Hook, Short & Long Hook Ends

## MILTEX BARR DOUBLE ENDED PROBE

28-109 7" Double Ended Probe, Rigid

## MILTEX BARR FISTULA PROBE

28-104 11 1/2" Overall Fistula Probe w/8" sterling shaft, Eye At Tip

## MILTEX ROSSER CRYPT HOOK

28-131 6 3/4" Rosser Crypt Hook


## MILTEX PRATT RECTAL PROBE

28-108 11" Overall Rectal Probe w/7 3/4" sterling shaft

## MILTEX PRATT RECTAL SCISSORS

S-shaped with blunt points.

28-244 8 1/4" Rectal Scissors



## MILTEX MCGIVNEY HEMORRHOIDAL LIGATOR

Offset handle for better vision and comfort. 7" working length. Can be disassembled entirely for cleaning. Supplied with drum loading cone and 100 latex O-rings.


28-150 Hemorrhoidal Ligator, 7", w/drum loading cone & O-rings

28-152 Hemorrhoid Grasping Forceps, 7 1/2", Angled Shanks

28-153 Loading Cone For 28-150

### Accessories

28-151 Latex O-Rings For 28-150, 100/pkg



28-360

## MILTEX RECTAL COTTON CARRIER

28-360 15" Shaft Rectal Cotton Carrier w/spring action




## MILTEX SMITH (BUIE) ANAL RETRACTORS

28-64 Self-Retaining Anal Retractor, Swivel Blades 2 1/2" x 7/8" Open to 4"


28-65 Self-Retaining Anal Retractor, Extra Long Swivel Blades 4" x 7/8" Open to 4"

## MILTEX BARR ANAL RETRACTOR


28-6 8" Anal Retractor, Self-Retaining, Blades 2 3/4" x 7/8"



28-320



28-323



28-341

## MILTEX TURRELL BIOPSY FORCEPS

Rotating shafts. Angled jaws 4 x 8mm bite with basket.

28-320 10" Biopsy Forceps w/ring handle

28-322 14" Biopsy Forceps w/ring handle

28-323 14" Biopsy Forceps w/ring handle & teeth

28-324 16" Biopsy Forceps w/ring handle

28-326 10" Shaft Only

28-328 14" Shaft Only


28-329 14" w/teeth Shaft Only

28-330 16" Shaft Only

### Accessories

28-340 Ring Handle


28-341 Pistol Grip Handle




### MILTEX IMPROVED BIOPSY FORCEPS WITH ROTATING SHAFTS

In pistol-grip handle. Shafts easily disassemble from handle. Constructed without small assembly parts. All shafts and handles are interchangeable. Jaws 3.5 mm x 10 mm.

28-380	25 cm (10") Straight Complete Biopsy Forceps w/pistol grip handle
28-382	35 cm (14") Straight Complete Biopsy Forceps w/pistol grip handle
28-386	35 cm (14") Angled Complete Biopsy Forceps w/pistol grip handle
28-400	Pistol Grip Spring Handle Only
28-410	25 cm (10") Straight Shaft Biopsy Forceps Only
28-412	35 cm (14") Straight Shaft Biopsy Forceps Only
28-416	35 cm (14") Angled Shaft Biopsy Forceps Only



28-300



28-303

### MILTEX YEOMAN BIOPSY FORCEPS

28-294	9" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket
28-296	13" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket
28-298	15" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket
28-300	10" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket
28-302	14" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket
28-303	14" Rotating Shaft Biopsy Forceps w/teeth & ring handle, Straight 4 x 8mm Bite w/basket
28-304	16" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket
28-306	10" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket
28-308	14" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket
28-309	14" Rotating Shaft Only w/teeth, Straight 4 x 8mm Bite w/basket
28-310	16" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket

### MILTEX WATSON-WILLIAMS POLYP FORCEPS

4 1/4" shaft.

20-504	7 3/4" Polyp Forceps, Small Jaws 6mm Wide
20-506	7 3/4" Polyp Forceps, Large Jaws 9mm Wide

### RUSCH O'CONNOR RECTAL FINGER

Latex. Non-sterile.

372000	Rectal Finger
--------	---------------

### MILTEX BUIE SUCTION TUBE

With finger valve.

28-370	16" Buie Suction Tube, 15 Fr
--------	------------------------------

### M-PACT AUTOPSY SAW

0295-270 Autopsy Saw

### M-PACT BONE DUST MOBILE VACUUM

0295-420 Bone Dust Vacuum

#### Accessories

0295-410 Disposable Hepa Filter Cartridge

### MILTEX POST MORTEM HAMMER

34-220 9 1/2" Post Mortem Hammer w/hook handle

### MILTEX VETERINARY SURGICAL NEEDLES

MS1281-1	Size 1 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-2	Size 2 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-3	Size 3 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-4	Size 4 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-5	Size 5 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-6	Size 6 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-7	Size 7 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1281-8	Size 8 Vet Surgical Needle, 3/8 Circle Cutting Edge, 12/pkg
MS1282-1	Size 1 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-3	Size 3 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-4	Size 4 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-5	Size 5 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-6	Size 6 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-7	Size 7 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg
MS1282-8	Size 8 Vet Surgical Needle, 1/2 Circle Cutting Edge, 12/pkg

Can't find what  
you're looking for?  
Call your  
representative.



## Look to ALLEVYN<sup>◇</sup> for a wide range of responsive dressings

With two unique, advanced absorption technologies—hydrocellular polyurethane foam and sodium polyacrylate—ALLEVYN dressings can manage all exudate levels, from none to extra-high. This achieves optimal moisture balance for a rapid wound healing environment.

ALLEVYN offers 3 adhesive levels—secure, gentle, and non-adhesive—to address all moisture, friction, and skin situations. Plus, our custom-designed shapes can handle any terrain.

For dressings that let you respond to any wound environment, look no farther than ALLEVYN.

Go to [www.allevyn.com](http://www.allevyn.com) for more information.



Respond to the wound healing environment with ALLEVYN

©2005 Smith & Nephew  
\*Trademark of Smith & Nephew. Certain marks reg. U.S. Pat. & TM Off.

