




76100

77600

PRO ADVANTAGE® APPLICATORS

- 75900 Applicator Stick, 6" x 1/12", 1000/bx, 20 bx/cs
- 76100 Cotton-Tipped Applicator, 3" x 1/12", Wooden Shaft, Non-Sterile, 100/bg, 10 bg/bx, 10 bx/cs
- 76200 Cotton-Tipped Applicator, 6" x 1/12", Wooden Shaft, Non-Sterile, 100/bg, 10 bg/bx, 10 bx/cs
- 76700 Cotton-Tipped Applicator, 6" x 1/12", Wooden Shaft, Sterile, 1/pk, 100 pk/bx, 10 bx/cs
- 76701 Cotton-Tipped Applicator, 6" x 1/10", Plastic Shaft, Sterile, 1/pk, 100 pk/bx, 10 bx/cs
- 76800 Cotton-Tipped Applicator, 6" x 1/12", Wooden Shaft, Sterile, 2/pk, 100 pk/bx, 10 bx/cs
- 77100 Cotton-Tipped Applicator, 3" x 1/12", Wooden Shaft, Sterile, 2/pk, 100 pk/bx, 10 bx/cs
- 77600 OB/GYN Applicator, 8", Rayon Tipped, Plastic Shaft, 50/bx, 10 bx/cs



SOLON APPLICATORS

- 35800 Plain Applicator, 12" x 1/12", 1000/bx, 12 bx/cs
- 35901 Wood Applicator Stick, Sterile, 6" x 1/12", 2/pk, 100 pk/bx, 10 bx/cs
- 36100 Cotton-Tipped Applicator, 3" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 36108 Cotton-Tipped Applicator, Extra Cotton, 3" x 1/12", Wood Stick, Large Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36121 Cotton-Tipped Applicator, 3" x 1/12", Wood Stick, Small Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36122 Cotton-Tipped Applicator, 3" x 1/12" Wood Stick, Tapered Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36181 Cotton-Tipped Applicator, Extra Cotton, 3" x 1/10", Plastic Stick, Large Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36200 Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 36208 Cotton-Tipped Applicator, Large Tip, 6" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 36217 Non-Sterile Cotton Tipped, Double End, Wood Shaft, 6" x 1/12", 100/bg, 8 bg/bx, 10 bx/cs
- 36224 Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, Small Tip, 100/bg, 10 bg/bx, 10 bx/cs

SOLON APPLICATORS (Continued)

- 36225 Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, Tapered Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36228 Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, X-Small Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36281 Cotton-Tipped Applicator, Extra Cotton, 6" x 1/10", Plastic Stick, Large Tip, 100/bg, 10 bg/bx, 10 bx/cs
- 36600 Sterile Calcium Alginate Applicator, Nasopharyngeal, 5 1/2" x .035, Regular, 1/pk, 100 pk/bx, 10 bx/cs, Type 1
- 36601 Polyester-Tipped Nasopharyngeal Applicator, Sterile, 5 1/2" x .035", Aluminum Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36602 Rayon-Tipped Applicator, Nasopharyngeal, Sterile, 5 1/2" x .035", Aluminum Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36603 Cotton-Tipped Applicator, Nasopharyngeal, Sterile, 5 1/2" x .35", Aluminum Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36604 Polyester-Tipped Applicator, Sterile, 5 1/2" x .025", Aluminum Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36607 Sterile Calcium Alginate Applicator, Urethro-Genital, 5 1/2" x .025, Fine, 1/pk, 100 pk/bx, 10 bx/cs, Type 4
- 36702 Rayon-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36703 Polyester-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36704 Calcium Alginate-Tipped Applicator, Sterile, 6" x 1/10", Plastic Stick, 1/pk, 100 pk/bx, 10 bx/cs, Type 3
- 36705 Calcium Alginate-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 1/pk, 100 pk/bx, 10 bx/cs, Type 2
- 36716 Polyester-Tipped Applicator, Sterile, 6" x 1/10", Plastic Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36721 Sterile, Rayon-Tipped Applicator, 6" x 1/10", Plastic Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 36800 Cotton Tipped Applicator, Sterile 6" x 1 1/2" Wood, 2/pk, 100 pk/bx, 10 bx/cs
- 36806 Sterile, Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, 2/pk, 1000 pk/cs
- 36815 Polyester-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 2/pk, 100 pk/bx, 10 bx/cs
- 36816 Polyester-Tipped Applicator, Sterile, 6" x 1/10", Plastic Stick, 2/pk, 100 pk/bx, 10 bx/cs
- 37100 Sterile Cotton Tipped, Wood Shaft, 3", 2/Sterile Pack, 100 pk/bx, 10 bx/cs
- 37500 Rayon-Tipped Proctoscopic Applicator, Non-Sterile, Polypropylene Stick, 16" x 5/32", 50/bx, 10 bx/cs
- 37502 Rayon-Tipped Proctoscopic Applicator, Non-Sterile, Polypropylene Stick, 16" x 5/32", 100/bx, 10 bx/cs
- 37504 Rayon-Tipped Proctoscopic Applicator, Non-Sterile, Polypropylene Stick, 16" x 5/32", Bulk, 500/cs
- 37600 Rayon-Tipped OB-GYN Applicator, Non-Sterile, Paper Stick, 8" x 5/32", 50/bx, 10 bx/cs
- 37602 Rayon-Tipped OB-GYN Applicator, Non-Sterile, Paper Stick, 8" x 5/32", 100/bx, 10 bx/cs
- 37603 Rayon-Tipped OB-GYN & Proctoscopic Applicator, Sterile, Polypropylene Stick, 8" x 5/32", 2/pk, 50 pk/bx, 10 bx/cs
- 37604 8" Ob-Gyn Paper Shaft, Rayon-Tipped Applicator, Bulk, 500/cs
- 54100 Solon Care™, Non-Sterile, Wood, 3", Flat, 1000/bx, 10 bx/cs
- 56805 Cotton Tipped Applicators, Wood Shaft, 6", Sterile, 5/pk, 40 pk/bx, 10 bx/cs
- 56810 Cotton Tipped Applicators, Wood Shaft, 6", Sterile, 10/pk, 20 pk/bx, 10 bx/cs
- 56820 Cotton Tipped Applicators, Wood Shaft, 6", Sterile, 20/pk, 10 pk/bx, 10 bx/cs

Applicators



BIRCHWOOD SCOPETTES® JR NON-SIZED APPLICATORS

The non-sized swab head is somewhat softer and larger than the sized swab. Many clinicians prefer the softer, fluffier non-sized swab. 5/8" lint free rayon swab head with a 3/16" cushion tip. Rolled paper stem.

34-7027-12	Single Ended, Large Tip 8" Applicator, 100/bx, 12bx/cs
34-7027-8	Single Ended, Large Tip 8" Applicator, Bulk, 500/cs
93103	Large Tip Applicators, 16", Sterile, 50/pk, 2 pk/cs



COVIDIEN/KENDALL Q-TIPS COTTON-TIPPED APPLICATORS

Single-tipped, wood stick applicators. Soft, absorbent tips. Versatile.

8884540400	3" Cotton-Tipped Applicator, 100/bg, 10 bg/bx, 10 bx/cs
8884540500	6" Cotton-Tipped Applicator, 100/bg, 10 bg/bx, 10 bx/cs
8884541300	6" Cotton-Tipped Applicator, Sterile, 2/pk, 200/cs
8884541400	6" Cotton-Tipped Applicator, Sterile, 2/pk, 200/bx, 10 bx/cs

CROSSTEX ADVANTAGE PLUS COTTON TIPPED APPLICATORS

Similar to the Premium brand, they are the perfect alternative for the price conscious practice.

H3AC	3" Applicator, 100/bg, 10 bg/bx, 10 bx/cs
H6AC	6" Applicator, 100/bg, 10 bg/bx, 10 bx/cs

CROSSTEX PREMIUM COTTON TIPPED APPLICATORS

All wood dowels with the tightest bud in the industry. Packaged in polywrapped box, applicators are sealed in autoclavable bags of 100 tips. Made in USA.

H3C	3" Applicator, 100/bg, 10 bg/bx, 10 bx/cs
H6C	6" Applicator, 100/bg, 10 bg/bx, 10 bx/cs

Your Patients Deserve The Best Quality, You Want The Best Value.

Solon has your answer — with a full range of disposable products for your practices every need.

A few examples (left to right) ►

- #36121 — 3" Small Cotton Tip
- #36600 — Sterile Calcium Alginate Tip
- #36225 — 6" Tapered Cotton Tip
- #36228 — 6" Extra Small Cotton Tip
- #36721 — Sterile Rayon-Tipped Applicator
- #32000 — 6" Tongue Depressor

Quality
Your Practice
Can Count On



For more info, please visit:

www.solonme.com





P150110

P150115



P150120

P150125



P150130

PRO ADVANTAGE® FABRIC ADHESIVE BANDAGE

Ideally suited for areas where flexibility is important. Absorbent, breathable island pad cushions and protects wound. Long lasting adhesive stays secure even when wet. Latex free.

- P150110 Adhesive Bandage, Knuckle Bands, 1 1/2" x 3", 100/bx, 12 bx/cs
- P150115 Adhesive Bandage, Finger Tip, 1 3/4" x 2", 100/bx, 12 bx/cs
- P150120 Adhesive Bandage, Strips, 1" x 3", 100/bx, 12 bx/cs
- P150125 Adhesive Bandage, Strips, 2" x 4", 50/bx, 12 bx/cs
- P150130 Adhesive Bandage, Strips, 3/4" x 3", 100/bx, 12 bx/cs



P150135

P150140



P150145

P150150


PRO ADVANTAGE® SHEER ADHESIVE BANDAGE

Absorbent breathable island pad cushions and protects wound. Long-lasting adhesive stays secure even when wet. Latex free.

- P150135 Adhesive Bandage, Strips, 1" x 3", 100/bx, 12 bx/cs
- P150140 Adhesive Bandage, Strips, 2" x 4", 50/bx, 12 bx/cs
- P150145 Adhesive Bandage, Strips, 3/4" x 3", 100/bx, 12 bx/cs
- P150150 Adhesive Bandage, Spots, 7/8" Dia., 100/bx, 12 bx/cs



CBD3018



CBD3019

ASO CAREBAND™ CLEAR ADHESIVE BANDAGES

Latex free with an island pad to protect the wound on all sides. This transparent tape bandage allows natural skin tones to show through.

- CBD3018 3/4" x 3" Clear Adhesive Strips, Latex Free, 100/bx, 12 bx/cs
- CBD3019 1" x 3" Clear Adhesive Strips, Latex Free, 100/bx, 12 bx/cs

ASO CAREBAND™ SHEER ADHESIVE STRIP BANDAGES - TRAY PRESENTATIONS

Made with translucent tape to blend evenly with a variety of skin tones. This is a convenient tray presentation for ease of bandage access. Latex free.

- ASO2020-010 Sheer Adhesive Bandage, 3/4" x 3", Latex Free, 150/tray, 10 tray/cs
- ASO2021-010 Sheer Adhesive Bandage, 1" x 3", Latex Free, 150/tray, 10 tray/cs



CRA5261



GAR5293



GAR5561



WHS5018




WHS5375

ASO CAREBAND™ DECORATED BANDAGES

Latex free. Decorated bandages offer the same high quality protection as our plastic bandages including the non-stick pad to protect the wound. Easy to open and apply. Stays secure - even when wet.

- CBD5747-012-130 Pink Ribbon Health Awareness Bandage, 3/4" x 3" Strips, Latex Free, 100/bx, 12 bx/cs
- CBD5893 Pink Camouflage Bandage, X-Large, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
- CBD5894 Green Camouflage Bandage, 2" x 3", Latex Free, 50/bx, 12 bx/cs
- CBD5895 Pink Camouflage Bandage, 2" x 3", Latex Free, 50/bx, 12 bx/cs
- CBD5896 Pink Camouflage Bandage, X-Large, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
- CRA5261 Crayola Bandages, 5/8" x 3" Strips, Latex Free, Assorted (red, yellow & blue), 100/bx, 12 bx/cs
- GAR5293 Garfield Bandages, 3/4" x 3" Strips, Latex Free, 100/bx, 12 bx/cs
- GAR5561 Garfield Adhesive Bandages, 7/8" Spots, Latex Free, 100/bx, 12 bx/cs
- WHS5018 Wild Hot Bandages, 3/4" x 3" Strips, Latex Free, 100/bx, 12 bx/cs
- WHS5375 Wild Hot Spot Bandages, 7/8", Latex Free, 100/bx, 12 bx/cs

Adhesive Bandages



CBD4018



CBD4019

ASO CAREBAND™ FABRIC ADHESIVE STRIP BANDAGES

Latex free. Ideally suited for areas where flexibility and adhesion are important.

CBD4016	X-Large Fabric Strip Bandage, 2" x 4", Latex Free, 50/bx, 12 bx/cs
CBD4018	Fabric Strip Bandage, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
CBD4019	Fabric Strip Bandage, 1" x 3", Latex Free, 100/bx, 12 bx/cs
CBD4022	Fabric Spot Bandage, 7/8", Latex Free, 100/bx, 12 bx/cs
CBD4024	Fabric Knuckle Strips, Latex Free, 100/bx, 12 bx/cs
CBD4025	Fabric Fingertip Strips, Latex Free, 100/bx, 12 bx/cs
CBD4476	Patch Bandage, 2" x 3", Latex Free, 50/bx, 12 bx/cs

ASO CAREBAND™ FOAM STRIP

CBD8095	3/4" x 3" Foam Strip, 80/bx, 12 bx/cs
CBD8096	1" x 3" Foam Strip, 80/bx, 12 bx/cs



CBD2016



CBD2018



CBD2019



CBD2022

ASO CAREBAND™ SHEER ADHESIVE STRIP BANDAGES

Made with translucent tape to blend evenly with a variety of skin tones. Comes in 6 different sizes. Latex free.

CBD2016	X-Large Sheer Strip, 2" x 4", Latex Free, 50/bx, 12 bx/cs
CBD2018	Sheer Strips, 1" x 3", Latex Free, 100/bx, 12 bx/cs
CBD2019	Sheer Strips, 1" x 3", Latex Free, 100/bx, 12 bx/cs
CBD2022	Sheer Spots, 7/8", Latex Free, 100/bx, 12 bx/cs
CBD2027	Sheer Mini Strips, 3/8" x 1 1/2", Latex Free, 100/bx, 12 bx/cs
CBD2760-024-000	Sheer Junior Strips, 5/8" x 2 1/4", Latex Free, 100/bx, 24 bx/cs

ASO CAREBAND™ HYDROCOLLOID BANDAGES

Latex free. Intended for management of lightly exuding wounds by absorbing exudate and by maintaining a moist wound environment. Perfect for blisters and mild burns.

CBD6406-12	Bandage, Hydrocolloid, Medium, Latex Free, 10/bx, 12 bx/cs
CBD6407-12	Bandage, Hydrocolloid, Large, Latex Free, 6/bx, 12 bx/cs

CAREBAND™
ADHESIVE BANDAGES

QWIKstrip™
ADHESIVE BANDAGES



www.asomedical.com



ASO CAREBAND™ HYPOALLERGENIC ADHESIVE STRIP BANDAGES

Latex free with an island pad to protect the wound on all sides. This hypoallergenic tape bandage is ideal for those with sensitive skin.

CBD4744-012-000 Hypo-Allergenic Strips, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs

ASO CAREBAND™ WATERPROOF BANDAGES

Island pad seals out moisture allowing wound to heal in moist environments. 100% latex free.

CBD1321 Reinforced Waterproof Bandage, Tan, 1" x 3", Latex Free, 100/bx, 12 bx/cs
 CBD1851 Reinforced Waterproof Bandage, Tan, 2" x 3", Latex Free, 50/bx, 12 bx/cs
 CBD3321 Reinforced Waterproof Bandage, Clear, 1" x 3", Latex Free, 100/bx, 12 bx/cs
 CBD3823 Reinforced Waterproof Bandage, Clear, 2" x 3", Latex Free, 50/bx, 12 bx/cs

ASO CAREBAND™ PLASTIC ADHESIVE STRIP BANDAGES

100% latex free with absorbent, breathable, non-stick pad to cushion and protect wounds. Stays secure, even when wet.

CBD1018 Plastic Adhesive Strips, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
 CBD1019 Plastic Adhesive Strips, 1" x 3", Latex Free, 100/bx, 12 bx/cs

ASO QWIKSTRIP™ BANDAGES

Latex free. Simple, one-handed application to avoid direct contact with the bandage itself. Developed by a surgeon to reduce the threat of contaminating the sterile pad. Allows healthcare professionals to apply bandages without needing to remove their gloves. Has a highly absorbent, non-stick sterile pad. Adheres to Universal Safety Precautions.

QWK1018 Plastic Strips, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
 QWK1019 1" x 3" Plastic Strips, 100/bx, 12 bx/cs
 QWK2018 3/4" x 3" Sheer Strips, 100/bx, 12 bx/cs
 QWK2019 Sheer Strips, 1" x 3", Latex Free, 100/bx, 12 bx/cs
 QWK2022 Sheer Spots, 7/8", Latex Free, 100/bx, 12 bx/cs
 QWK4018 Fabric Strips, 3/4" x 3", Latex Free, 100/bx, 12 bx/cs
 QWK4019 Fabric Strips, 1" x 3", Latex Free, 100/bx, 12 bx/cs



3M™ NEXCARE™ ACTIVE™ STRIPS FLEXIBLE FOAM BANDAGES

These flexible foam bandages are designed for patients with an active lifestyle. The unique foam stretches, flexes and conforms better than conventional plastic bandages. Advanced water-resistant adhesive provides extra sticking power on damp or sweating skin.

512-36 Bandages, Foam, 1 1/16" x 3" Strips, 30/bx, 12 bx/dz, 2 dz/cs



7028



7031

DUKAL ADHESIVE BANDAGES

Latex free.

7028 Bandage, Plastic Adhesive Strips, 3/4" x 3", 100/bx, 50 bx/cs
 7031 Bandage, Plastic Adhesive Strips, 1" x 3", 100/bx, 50 bx/cs
 7602 Bandage, Flexible Fabric Adhesive Strips, Knuckle, 100/bx, 24 bx/cs
 7604 Bandage, Sheer Adhesive Strips, 3/4" x 3", 100/bx, 24 bx/cs
 7605 Bandage, Flexible Fabric Adhesive Strips, 3/4" x 3", 100/bx, 24 bx/cs
 7608 Bandage, Sheer Adhesive Strips, 1" x 3", 100/bx, 24 bx/cs
 7609 Bandage, Flexible Fabric Adhesive Strips, 1" x 3", 100/bx, 24 bx/cs
 7611 Bandage, Plastic Adhesive Strips, X-Large", 100/bx, 24 bx/cs
 7612 Bandage, Sheer Adhesive Strips, X-Large", 100/bx, 24 bx/cs
 7613 Bandage, Flexible Fabric Adhesive Strips, X-Large", 100/bx, 24 bx/cs
 7614 Plastic Spot Bandage, 100/bx, 24 bx/cs
 7615 Bandage, Sheer Adhesive Spot, 100/bx, 24 bx/cs
 7616 Bandage, Plastic Adhesive Strips, 3/4" x 3", 100/bx, 24 bx/cs
 7617 Bandage, Plastic Adhesive Strips, 1" x 3", 100/bx, 24 bx/cs
 7618 Bandage, Flexible Fabric Adhesive Strips, Fingertip, 100/bx, 24 bx/cs

COVIDIEN/KENDALL CURITY FLEXIBLE BANDAGES

Made of fabric that stretches and conforms with movement. Ideal for fingers, toes, elbows and other bending places. Strong, hypoallergenic adhesive formula sticks even when wet. "Ouchless" TELFA pad promotes healing by not sticking to wound. Individually wrapped in easy-to-open, pull-tab, sterile sleeves. Latex Free.

5024LF Fabric Adhesive Bandage, 3/4", Latex Free, 50/bx, 24 bx/cs
 6131LF Fabric Adhesive Bandage, 3/4", Latex Free, Assorted Pink, Orange and Yellow Neon, 50/bx, 24 bx/cs
 6133LF Fabric Adhesive Bandage, 1", Latex Free, Assorted Pink, Orange and Yellow Neon, 50/bx, 24 bx/cs
 6135LF Fabric Adhesive Bandage, 1", Latex Free, 50/bx, 24 bx/cs
 6136LF Fabric Adhesive Bandage, 2" x 3 1/2", Latex Free, 50/bx, 12 bx/cs
 6137LF Fabric Adhesive Bandage, 3/4" x 3", for Sensitive Skin, Latex Free, 50/bx, 12 bx/cs
 6138LF Fabric Adhesive Bandage, 1" x 3", for Sensitive Skin, Latex Free, 50/bx, 12 bx/cs
 6409LF Fabric Adhesive Bandage, Fingertip, Latex Free, 40/bx, 30 bx/cs
 6411LF Fabric Adhesive Bandage, Knuckle, Latex Free, 30/bx, 40 bx/cs
 6417LF Fabric Adhesive Bandage, 3/8" x 1 1/2" Junior, Latex Free, 100/bx, 36 bx/cs

COVIDIEN/KENDALL CURITY ADHESIVE BANDAGES

Latex free. Made of perforated vinyl film. All bandages feature "ouchless" TELFA pad; promote healing by not sticking to wound. Hypoallergenic adhesive reduces potential for skin irritation. Porous adhesive allows skin to breathe.

3065LF Adhesive Bandage, 3/4" x 3", Sheer, 50/bx, 72 bx/cs
 3270LF Adhesive Bandage, 1" x 3", Sheer, 50/bx, 72 bx/cs
 6193LF Adhesive Bandage, 3/4" x 3", Plastic, 50/bx, 24 bx/cs
 6403LF Adhesive Bandage, 7/8", Sheer Spots, 100/bx, 36 bx/cs
 6416LF Adhesive Bandage, 1 1/2" x 1 1/2", Patch, 50/bx, 24 bx/cs
 6947LFWC Adhesive Bandage, 3/4" x 3", Plastic, 50/bx, 72 bx/cs
 7028LF Adhesive Bandage, 1" x 3", Plastic, 50/bx, 72 bx/cs
 7893LF Adhesive Bandage, 1" x 3", Plastic, 50/bx, 24 bx/cs

Adhesive Bandages



HARTMANN-CONCO FLEX-BAND® FABRIC ADHESIVE BANDAGES

Ideal for dressing minor injuries and as a bandage for sensitive skin. Made with a hypoallergenic adhesive. Offers widthway strength for maximum conformability and flexibility. The full adhesive coverage around the wound pad firmly secures the bandage and provides a protective seal. Firmly adheres to the skin, yet can be removed without leaving sticky residue. Round-beveled edges that prevent premature lifting or peeling of the bandage. Strong fabric that is flexible and conforms to any contour. Absorbent wound pad that will not stick to wound. Water-repellant layer that protects against moisture penetration. Latex free.

46100000	3" x 3" 4-Wing Bandage, 50/bx, 24 bx/cs
46110000	7/8" Round Spots Bandage, 100/bx, 24 bx/cs
46120000	2 1/2" x 1 3/8" Toe Shield Bandage, 100/bx, 24 bx/cs
46130000	2" x 3" Patch Bandage, 50/bx, 24 bx/cs
46140000	1 1/2" x 3" Knuckle Bandage, 100/bx, 24 bx/cs
46150000	1 1/4" x 1" Oval Bandage, 100/bx, 24 bx/cs
46160000	3/4" x 3" Small Strip Bandage, 100/bx, 24 bx/cs
46180000	1" x 3" Large Strip Bandage, 100/bx, 24 bx/cs
46200000	2" x 1 3/4" Fingertip Small Bandage, 100/bx, 24 bx/cs
46210000	2" x 2 1/2" Fingertip Big Bandage, 50/bx, 24 bx/cs

DERMA SCIENCES ECONOMY ADHESIVE BANDAGES

Latex free.

99982	Bandage Sheer Spot, 100/bx, 36 bx/cs
99990	1" x 3" Flexible Bandage, 100/bx, 36 bx/cs
99991	3/4" x 3" Flexible Bandage, 100/bx, 36 bx/cs

J&J BAND-AID® ADHESIVE BANDAGES

005705	X-Large Adhesive Bandages, 10/bx, 24 bx/cs
--------	--

J&J BAND-AID® ADHESIVE PADS

004764	2 1/4" x 3" Adhesive Pads, 10/bx, 24 bx/cs
004768	2 7/8" x 4" Adhesive Pads, 10/bx, 24 bx/cs

J&J BAND-AID® CLEAR ADHESIVE BANDAGES

004670	All One Size Clear Adhesive Bandages, 30/bx, 24 bx/cs
004708	Spots Clear Adhesive Bandages, 50/bx, 24 bx/cs
005701	Assorted Clear Adhesive Bandages, 45/bx, 24 bx/cs

J&J BAND-AID® DECORATED™ ADHESIVE BANDAGES

Latex free.

004416	Barbie™ Assorted Adhesive Bandages, 25/bx, 24 bx/cs
004458	Scooby-Doo™ Assorted Adhesive Bandages, 25/bx, 24 bx/cs
004473	SpongeBob™ SquarePants Assorted Adhesive Bandages, 20/bx, 24 bx/cs
004484	Dora The Explorer Assorted Adhesive Bandages, 25/bx, 24 bx/cs
004522	SpiderMan™ Assorted Adhesive Bandages, 20/bx, 24 bx/cs
004602	Elmo Assorted Adhesive Bandages, 30/bx, 24 bx/cs



J&J BAND-AID® FLEXIBLE FABRIC ADHESIVE BANDAGES

For areas where freedom of movement is desirable, such as across fingers, elbows, knees and other joints.

004430	Assorted Flexible Fabric Adhesive Bandages, 30/bx, 24 bx/cs
004431	All One Size Flexible Fabric Adhesive Bandages, 30/bx, 24 bx/cs
004434	3/4" x 3" Adhesive Bandage Strip, 100/bx, 12 bx/cs
004436	Fingertip Adhesive Bandage, 100/bx, 12 bx/cs
004438	Knuckle Adhesive Bandage, 100/bx, 12 bx/cs
004444	1" x 3" Adhesive Bandage Strip, 100/bx, 12 bx/cs
004452	Knuckle & Fingertip Flexible Fabric Adhesive Bandages, 20/bx, 24 bx/cs
005685	X-Large Flexible Fabric Adhesive Bandages, 10/bx, 24 bx/cs
005693	Assorted Adhesive Bandages, 45/bx, 24 bx/cs
005753	Bandage, Adhesive, Flexible Fabric, Travel Pack 8 Count, 72 bx/cs



J&J BAND-AID® PLASTIC ADHESIVE BANDAGES

Made from a strong, air-vented plastic backing and coated with a long-lasting adhesive for your protection.

005634	3/4" x 3" Adhesive Bandage Strip, 100/bx, 12 bx/cs
005635	Family Pack Plastic Adhesive Bandages, 60/bx, 24 bx/cs
005644	1" x 3" Adhesive Bandage Strip, 100/bx, 12 bx/cs

J&J BAND-AID® ADHESIVE BANDAGE ACTIV FLEX™

Advanced Hydrocolloid Technology, clinically proven to heal faster, minimize the appearance of scars, Unique design for fingers, Waterproof material, Ultra Thin.

004413	Adhesive Bandage Activ-Flex™, Large, 1.75" x 2.75", 6/bx, 24 bx/cs
004414	Adhesive Bandage Activ-Flex™, Regular, .75" x 3", 10/bx, 24/cs
005722	Adhesive Bandage Activ-Flex™ Finger-Wrap™, Assorted, 10/bx, 24/cs

J&J BAND-AID® ADHESIVE BANDAGE ADVANCED HEALING BLISTER

Gel strip relieves pain while protecting the blister for faster healing, Helps minimize the appearance of scars, Unique design for fingers & toes, Waterproof material.

004488	Advanced Healing Blister Adhesive Bandages, 1-4/25" x 2-4/25", 6/bx, 24/cs
004582	Advanced Healing Blister Adhesive Bandages, 2-3/5" x 5/8", 6/bx, 24/cs

J&J BAND-AID® ADHESIVE BANDAGES - ANTIBIOTIC

Convenient one-step infection protection, Minimizes appearance of scars.

SP255	Suture Spool Polyviolene 1-0, Green Braided, 100 yds
005569	Adhesive Bandages, .75" x 3", 20/bx, 24 bx/cs
005570	Adhesive Bandages, Assorted Sizes, 20/bx, 24 bx/cs



J&J BAND-AID® SHEER STRIPS & SPOTS

Helps provide all-purpose protection that lets skin breathe.

004634	Strips, 3/4" x 3", 100/bx, 12 bx/cs
004644	Strips, 1" x 3", 100/bx, 12 bx/cs
004666	All One Size Sheer Adhesive Bandages, 40/bx, 24 bx/cs
004667	Assorted Sheer Adhesive Bandages, 40/bx, 24 bx/cs
004668	Assorted Sheer Adhesive Bandages, 60/bx, 24 bx/cs
004669	Assorted Sheer Adhesive Bandages, 80/bx, 24 bx/cs
004711	Assorted Sheer Strips, 280/bx, 18 bx/cs
005716	Strips 2" x 4 1/2", 50/bx, 18 bx/cs

J&J BAND-AID® SPORT STRIP® ADHESIVE BANDAGES

004723	X-Wide All One Size Adhesive Bandages, 30/bx, 24 bx/cs
--------	--

J&J BAND-AID® TOUGH STRIPS™ ADHESIVE BANDAGES

004408	All One Size Adhesive Bandages, 20/bx, 24 bx/cs
004421	Finger-Care™ Assorted Adhesive Bandages, 15/bx, 24 bx/cs
004424	X-Large Adhesive Bandages, 10/bx, 24 bx/cs


J&J BAND-AID® WATER BLOCK PLUS® ADHESIVE BANDAGES

004446	Finger-Care™ Assorted Adhesive Bandages, 20/bx, 24 bx/cs
004474	Finger-Wrap™ All One Size Adhesive Bandages, 20/bx, 24 bx/cs
004535	Adhesive Bandages, 1" x 2.5" Clear, 30/bx, 24 bx/cs
005658	Finger-wrap®, Large, 2" x 3", 10/bx, 24 bx/cs
005659	Clear Assorted Adhesive Bandages, 30/bx, 24 bx/cs

NUTRAMAX BLUE METAL X-RAY DETECTABLE ADHESIVE BANDAGES

Latex free.

1114025	1" x 3" Adhesive Bandage, Plastic, 1500/cs
1631025	1" x 3" Adhesive Bandage, Flexible Fabric, 1300/cs
1637025	Knuckle Adhesive Bandage, Flexible Fabric, 1800/cs
1647025	3" Adhesive Bandage, Fingertip, Flexible Fabric, 1400/cs



NUTRAMAX CHILDREN'S CHARACTER ADHESIVE BANDAGES

Adhesive bandages imprinted with either cheerful cartoons or bright neon colors. Vented to help prevent skin maceration. Highly absorbent non-stick pad provides a moist environment for proper wound healing. Each bandage individually wrapped and sterilized to provide maximum protection. Latex free.

1065737	Adhesive Bandage, Scooby Doo™, 3/4" x 3", 100/bx, 12 bx/cs
1073737	Looney Tunes™ Bugs Bunny & Tasmanian Devil Adhesive Bandage, 3/4" x 3". 100/bx, 12 bx/cs
1074737	Tweety Spot, 7/8" Adhesive Bandage, 100/bx, 24 bx/cs
1075413	Glitter® Strips Stat Strip®, 100/bx, 12 bx/cs
1075538	Peanuts™, Charlie Brown® & Lucy™ Stat Strip®, 3/4" x 3", 100/bx, 12 bx/cs
1075737	Looney Tunes™ Bugs Bunny™ & Daffy Duck™ Assorted, Stat Strip®, 3/4" x 3", 100/bx, 12 bx/cs
1076413	Assorted Neon Colors Stat Strip®, Adhesive Bandage, 3/4" x 3", 100/bx, 12 bx/cs
1076538	Peanuts™, Snoopy® Sports Stat Strip®, 3/4" x 3", 100/bx, 12 bx/cs
1076737	Looney Tunes™, Wile E. Coyote & Road Runner Adhesive Bandage, 3/4" x 3", 100/bx, 12 bx/cs
1077538	Peanuts™ Snoopy® Assorted Friends & Pals, Stat Strip®, .75" x 3", 100/bx, 12bx/cs

Adhesive Bandages

NUTRAMAX CHILDREN'S CHARACTER ADHESIVE BANDAGES (Continued)

1079797	Adhesive Bandage, Tasmanian Devil™ & Tweety™, 7/8" Spot, 100/bx, 24 bx/cs
1083737	Looney Tunes™ Tweety™ Flowers Assorted Red/Blue, .75" x 3", Stat Strip®, 100/bx, 12 bx/cs
1084737	Looney Tunes™ Characters, 2" x 3" Patch, Stat Strip®, A fun bandage for large cuts and abrasions as well as knees & elbows, 50/bx, 12bx/cs
1085737	Looney Tunes™ Assortment, Bugs Bunny™ Characters, .75" x 3", Stat Strip®, 100/bx, 12 bx/cs
1086737	Adhesive Bandage, Bugs Bunny™ & Sylvester™, 3/4" x 3", 100/bx, 12 bx/cs
1087737	Adhesive Bandage, Spiderman, Wolverine & Ironman, 3/4" x 3", 100/bx, 12 bx/cs
1087837	Adhesive Bandage, Storm, Captain America & Silver Surfer, 3/4" x 3", 100/bx, 12 bx/cs
1087937	Adhesive Bandage, Hulk, The Thing & Thor, 3/4" x 3", 100/bx, 12 bx/cs
1307538-24	Peanuts™ Joe Cool®, 7/8" Spot, 100/bx, 24 bx/cs
15600	Circus® Characters, 3/4" x 3", Stat Strip®, 100/bx, 12 bx/cs
15601	Herbie® The Dinosaur Stat Strip®, 3/4" x 3", 100/bx, 12 bx/cs
15650	Adhesive Bandage, Planets & Stars, 3/4" x 3", 100/bx, 12 bx/cs
15661	Adhesive Bandage, Planets & Stars, 5/8" x 2 1/4", 100/bx, 12 bx/cs
16700	Camouflage Assorted Pink & Blue, Stat Strip®, .75" x 3", 100/bx, 12 bx/cs
CIRCS34B	Circus® Characters, 3/4" x 3", Bulk 1500/cs
GLIAST100	Glitter® Stars & Strips, Assorted Sizes, Stat Strip®, 100/bx, 12 bx/cs

NUTRAMAX DURABAND™ HEAVYWEIGHT FLEXIBLE BANDAGES

Adhesive bandages. Non-adherent pad. Sterile 1's. Latex free.

1638214	3" Adhesive Bandage, Knuckle, Perforated, Bulk, 3500/cs
1648214	3" Adhesive Bandage, Fingertip, Perforated, Bulk, 3000/cs
1653214	2" x 3" Adhesive Bandage, Patch, Perforated, Bulk, 3000/cs
1695214	1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs

NUTRAMAX BLUE NON-METAL FABRIC BANDAGES

Sterile and latex free this American White Cross is made of comfortable flexible fabric so that it conforms to the body. Blue color is ideal for food handlers.

1641025	Bandage Fingertip, 3", Bulk, 1400/cs
1650025	Bandage, 1" x 3", Bulk, 1500/cs
1664025	Bandage Fingertip, 2", Bulk, 2400/cs
1682025	Bandage, Knuckle, Bulk, 1800/cs

NUTRAMAX ECONOMY ADHESIVE BANDAGES

Latex free.

99949	Flexible Knuckle Bandage, 100/bx, 36 bx/cs
99979	2" Flexible Fingertip Bandage, 100/bx, 36 bx/cs
99992	1" x 3" Sheer Bandage, 100/bx, 36 bx/cs
99993	3/4" x 3" Sheer Bandage, 100/bx, 36 bx/cs



NUTRAMAX FIRST AID™ ADHESIVE BANDAGES

Adhesive bandages with non-adherent pad are sterile and lightweight. Stay secure and comfortable, allowing wound to remain protected even when wet. Assorted sizes. Latex free.

1010033	Plastic, 3/4" x 3" Adhesive Bandage, 10/bx, 288 bx/cs
1045033	Plastic, 3/4" x 3" Adhesive Bandage, 60/bx, 24 bx/cs
1070033	Plastic, 2" x 4" Adhesive Bandage, 50/bx, 24 bx/cs
1075033	Plastic, 3/4" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1090033	Plastic, 1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1104000	Plastic, 3/4" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
1114000	Plastic, 1" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
1130000	Plastic, 2" x 4 1/2" Adhesive Bandage, 1000/cs
1250033	Sheer Assorted Sizes Adhesive Bandage, 60/bx, 24 bx/cs
1260033	Sheer Assorted Sizes Adhesive Bandage, 80/bx, 24 bx/cs
1265033	Sheer, 2" x 4" Adhesive Bandage, 10/bx, 36 bx/cs
1275033	Sheer, 3/4" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1290033	Sheer, 1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1298033-24	Sheer Junior, 3/8" x 1 1/2" Adhesive Bandage, 100/bx, 24 bx/cs
1304000	Sheer, 3/4" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
1307033-24	Sheer Spots, 7/8" Adhesive Bandage, 100/bx, 24 bx/cs
1308033	Clear Strip Patch, 1 1/2" x 1 1/2" Adhesive Bandage, 100/bx, 12 bx/cs
1309033	Sheer Patch, 1 1/2" x 1 1/2" Adhesive Bandage, 100/bx, 12 bx/cs
1314000	Sheer, 1" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
1342000	Sheer Junior, 3/8" x 1 1/2" Adhesive Bandage, Bulk, 30000/cs
1415033	Clear Strip Assorted Sizes Adhesive Bandage, 30/bx, 24 bx/cs
1475033	Clear Strip, 3/4" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1490033	Clear Strip, 1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1510033	Flexible, 3/4" x 3" Adhesive Bandage, 10/bx, 288 bx/cs
1570033	Flexible, 2" x 4" Adhesive Bandage, 50/bx, 24 bx/cs
1580033	Flexible, 3/4" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1595033	Flexible, 1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1602033	Flexible Knuckle, 3" Adhesive Bandage, 100/bx, 12 bx/cs
1607033	Flexible Fingertip, 2" Adhesive Bandage, 100/bx, 12 bx/cs
1612033	Flexible Fingertip, 3" Adhesive Bandage, 50/bx, 12 bx/cs
1617033	Flexible Patch, 2" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1622033	Flexible Fabric 4-Wing, 3" x 3" Adhesive Bandage, 50/bx, 12 bx/cs
1626000	Flexible, 3/4" x 3" Adhesive Bandage, 130/tray, 10 trays/cs
1631000	Flexible, 1" x 3" Adhesive Bandage, 130/tray, 10 trays/cs
1637000	Flexible Knuckle, 3" Adhesive Bandage, Bulk, 1800/cs
1642000	Flexible Fingertip, 2" Adhesive Bandage, Bulk, 2400/cs
1650000	Plastic Patch, 2" x 3" Adhesive Bandage, 2000/cs
1652000	Flexible Patch, 2" x 3" Adhesive Bandage, Bulk, 1500/cs
1657000	Flexible Fabric 4-Wing, 3" x 3" Adhesive Bandage, Bulk, 1200/cs
1775033	Tricot, 3/4" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1790033	Tricot, 1" x 3" Adhesive Bandage, 100/bx, 12 bx/cs
1804000	Tricot, 3/4" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
1814000	Tricot, 1" x 3" Adhesive Bandage, 150/tray, 10 trays/cs

NUTRAMAX FLEXIBLE STRIP PADS

Soft, comfortable rayon/nylon fabric. Stretchable for maximum mobility. Sterile and LATEX FREE. Highly absorbent non-stick pad. Seals around entire wound to keep out dirt and contaminants. Water resistant.

FV24D	Small Fingertip Pad, 2", Latex Free, 100/bx, 12 bx/cs
FV3UR	Large Fingertip Pad, 3", Latex Free, 50/bx, 12 bx/cs
NVJED	Knuckle Pad, Latex Free, 100/bx, 12 bx/cs
PV23R	2" x 3" Patch, Latex Free, 50/bx, 12 bx/cs
SV13D	1" x 3" Pad, Latex Free, 100/bx, 12 bx/cs
SV3ER	4-Wing Pad, Latex Free, 50/bx, 12 bx/cs
SVT3D	3/4" x 3" Pad, Latex Free, 100/bx, 12 bx/cs



NUTRAMAX STAT STRIP™ ADHESIVE BANDAGES

Patented packaging and opening system provide the professional applying the strip and the patient with a level of protection unavailable in any other bandage. Unique 'snap open' wrapper is convenient and time saving as application can be performed with gloved hands. The hand never touches the adhesive or absorptive pad. Latex free.

15200	3/4" x 3" Sheer Adhesive Bandage, 100/bx, 12 bx/cs
15205	1" x 3" Sheer Adhesive Bandage, 100/bx, 12 bx/cs
15210	3/4" x 3" Flexible Fabric Adhesive Bandage, 100/bx, 12 bx/cs
15215	1" x 3" Flexible Fabric Adhesive Bandage, 100/bx, 12 bx/cs
FLEXS34B	3/4" x 3" Flexible Fabric Bandage, 1300/cs



NUTRAMAX SUPERFLEX™ ADHESIVE BANDAGES

This foam strip provides the maximum in comfort and conforming available in adhesive bandages today. Water resistant adhesive keeps the bandage in place even after it gets wet. Absorbent, non-stick pad will not stick to the wound and allows for easy removal. Latex free.

1785033	1" x 3" Adhesive Bandage, 150/tray, 10 trays/cs
---------	---



PRO ADVANTAGE® BUTTERFLY WOUND CLOSURE

For closure of small cuts and wounds. Latex free.

P150105	Medium Butterfly Wound Closure, 100/bx, 12 bx/cs
---------	--



B1550



R1540

3M™ STERI-STRIP™ ADHESIVE SKIN CLOSURES

Porous, nonwoven, rayon-backed material. Hypoallergenic adhesive. Designed to permit skin ventilation and inhibit maceration.

B1550	Nonreinforced Skin Closures, Blendtone, 1/8" x 3", 5 strips/envelope, 50 env/bx, 4 bx/cs
B1551	Nonreinforced Skin Closures, Blendtone, 1/4" x 3", 3 strips/envelope, 50 env/bx, 4 bx/cs
B1553	Nonreinforced Skin Closures, Blendtone, 1/4" x 4", 10 strips/envelope, 50 env/bx, 4 bx/cs
B1557	Nonreinforced Skin Closures, Blendtone, 1/2" x 4", 6 strips/envelope, 50 env/bx, 4 bx/cs
B1559	Nonreinforced Skin Closures, Blendtone, 1/2" x 2", 6 strips/envelope, 50 env/bx, 4 bx/cs
R1540	Reinforced Skin Closures, 1/8" x 3", 5 strips/envelope, 50 env/bx, 4 bx/cs
R1541	Reinforced Skin Closures, 1/4" x 3", 3 strips/envelope, 50 env/bx, 4 bx/cs
R1542	Reinforced Skin Closures, 1/4" x 1 1/2", 6 strips/envelope, 50 env/bx, 4 bx/cs
R1546	Reinforced Skin Closures, 1/4" x 4", 10 strips/envelope, 50 env/bx, 4 bx/cs
R1547	Reinforced Skin Closures, 1/2" x 4", 6 strips/envelope, 50 env/bx, 4 bx/cs
R1548	Reinforced Skin Closures, 1" x 5", 4 strips/envelope, 25 env/bx, 4 bx/cs
R1549	Reinforced Skin Closures, 1/2" x 2", 6 strips/envelope, 50 env/bx, 4 bx/cs

Closure Strips



3M™ STERI-STRIP™ ANTIMICROBIAL SKIN CLOSURES

Made of a porous, non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with filaments for added strength. Iodophor in the adhesive controls skin bacteria for up to 72 hours.

A1840	1/8" x 3" Skin Closures, 5 strips/env, 50 env/bx, 4 bx/cs
A1841	1/4" x 3" Skin Closures, 3 strips/env, 50 env/bx, 4 bx/cs
A1842	1/4" x 1 1/2" Skin Closures, 6 strips/env, 50 env/bx, 4 bx/cs
A1846	1/4" x 4" Skin Closures, 10 strips/env, 50 env/bx, 4 bx/cs
A1847	1/2" x 4" Skin Closures, 6 strips/env, 50 env/bx, 4 bx/cs
A1848	1" x 5" Skin Closures, 4 strips/env, 25 env/bx, 4 bx/cs



3M™ STERI-STRIP™ COMPOUND BENZOIN TINCTURE

May be used to increase adhesion of 3M Steri-Strip™ Skin Closures to the skin.

C1544	2/3 cc Vial Benzoin Tincture, 40/bx, 4 bx/cs
-------	--



3M™ STERI-STRIP™ ELASTIC SKIN CLOSURES

Made of a soft, non-woven elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.

E4540	1/8" x 3" Skin Closures, 5 strips/env, 50 env/bx, 4 bx/cs
E4541	1/4" x 3" Skin Closures, 3 strips/env, 50 env/bx, 4 bx/cs
E4542	1/4" x 1 1/2" Skin Closures, 6 strips/env, 50 env/bx, 4 bx/cs
E4546	1/4" x 4" Skin Closures, 10 strips/env, 50 env/bx, 4 bx/cs
E4547	1/2" x 4" Skin Closures, 6 strips/env, 50 env/bx, 4 bx/cs
E4548	1" x 5" Skin Closures, 4 strips/env, 25 env/bx, 4 bx/cs
E4549	1/2" x 2" Skin Closures, 6 strips/env, 50 env/bx, 4 bx/cs



3M™ STERI-STRIP™ SKIN CLOSURE STACKING RACK

For organizing and storing a variety of 3M Steri-Strip™ skin closure product styles and sizes.

1518	Skin Closure Stacking Rack, Small Tray, 24/cs
1519	Skin Closure Stacking Rack, Medium Tray, 12/cs
1555	Skin Closure Stacking Rack



3M™ STERI-STRIP™ WOUND CLOSURE SYSTEM

Includes (1) 3M Tegaderm™ Dressing: 2 3/8" x 1 7/8" and (1) Steri-Strip reinforced skin closure strip: 1 7/8" x 1/2".

W8512	Wound Closure System, 1/env, 25 env/bx, 4 bx/cs
W8514	Wound Closure System, 3/env, 25 env/bx, 4 bx/cs
W8516	Wound Closure System, 5/env, 25 env/bx, 4 bx/cs

DERMA SCIENCES SHUR-STRIP™ ADHESIVE SKIN CLOSURE STRIPS

81117	Adhesive Skin Closure Strips, 1/8" x 3", 50/bx, 4 bx/cs
81118	Adhesive Skin Closure Strips, 1/4" x 3", 50/bx, 4 bx/cs
81119	Adhesive Skin Closure Strips, 1/4" x 1 1/2", 50/bx, 4 bx/cs
81120	Adhesive Skin Closure Strips, 1/4" x 4", 50/bx, 4 bx/cs
81121	Adhesive Skin Closure Strips, 1/2" x 4", 50/bx, 4 bx/cs
81122	Adhesive Skin Closure Strips, 1" x 4", 25/bx, 4 bx/cs



DERMA SCIENCES SUTURE STRIP® PLUS WOUND CLOSURE STRIPS

Made of a flexible, nonwoven material with a skin-friendly adhesive that adheres aggressively to the skin, yet flexes enough on areas of movement to reduce the chance of skin shear or blistering. Requires no additional adhesive enhancers and is water-resistant.

TP1100	1/8" x 3" Wound Closure Strips, 5/pouch, 50 pch/bx
TP1101	1/4" x 3" Wound Closure Strips, 3/pouch, 50 pch/bx
TP1102	1/4" x 4" Wound Closure Strips, 10/pouch, 50 pch/bx
TP1103	1/2" x 4" Wound Closure Strips, 6/pouch, 50 pch/bx
TP1104	1/4" x 1 1/2" Wound Closure Strips, 6/pouch, 50 pch/bx
TP1105	1" x 5" Wound Closure Strips, 4/pouch, 25 pch/bx



DERMA SCIENCES SUTURE STRIP® WOUND CLOSURE STRIPS

Made of flexible, nonwoven material with skin-friendly adhesive. Allows maximum porosity and vapor transmission.

TS3100	1/8" x 3" Wound Closure Strips, 5/pouch, 50 pch/bx
TS3101	1/4" x 3" Wound Closure Strips, 3/pouch, 50 pch/bx
TS3102	1/4" x 4" Wound Closure Strips, 10/pouch, 50 pch/bx
TS3103	1/2" x 4" Wound Closure Strips, 6/pouch, 50 pch/bx
TS3104	1/4" x 1 1/2" Wound Closure Strips, 6/pouch, 50 pch/bx

ASO CAREBAND™ BUTTERFLY CLOSURE BANDAGES

Latex free. Intended for closure of small cuts and wounds. Center of closure won't stick to wounds.

CBD6418 Medium Butterfly Closures, 1 3/8" x 1/3", Latex Free, 100/bx, 12 bx/cs

DUKAL WOUND CLOSURE STRIPS

Can be used in conjunction with or as an excellent alternative to sutures or staples. Each strip consists of thin non-woven reinforced fibers that are arranged in a bilateral pattern for high tensile strength, to provide incision support and maintain closure. This adhesive backed product is translucent, porous, and breathable to promote good healing. Similar to Steri-Strips.

5150 Wound Closure Strips, Sterile, 1/8" x 3", 5/pkg, 50 pkg/bx, 4 bx/cs
 5151 Wound Closure Strips, Sterile, 1/4" x 3", 3/pkg, 50 pkg/bx, 4 bx/cs
 5152 Wound Closure Strips, Sterile, 1/4" x 1 1/2", 6/pkg, 50 pkg/bx, 4 bx/cs
 5156 Wound Closure Strips, Sterile, 1/4" x 4", 10/pkg, 50 pkg/bx, 4 bx/cs
 5157 Wound Closure Strips, Sterile, 1/2" x 4", 6/pkg, 50 pkg/bx, 4 bx/cs
 5158 Wound Closure Strips, Sterile, 1" x 4", 4/pkg, 25 pkg/bx, 4 bx/cs
 7600 Butterfly Wound Closure Strips, Medium, 100/bx, 24 bx/cs
 C-01 Wound Closure Strips, Non-Sterile, 1/8" x 3", 5/card, 500 cd/bx
 C-02 Wound Closure Strips, Non-Sterile, 1/4" x 3", 3/card, 500 cd/bx
 C-03 Wound Closure Strips, Non-Sterile, 1/4" x 1 1/2", 6/card, 500 cd/bx
 C-04 Wound Closure Strips, Non-Sterile, 1/4" x 4", 5/card, 500 cd/bx
 C-05 Wound Closure Strips, Non-Sterile, 1/2" x 4", 3/card, 500 cd/bx
 C-06 Wound Closure Strips, Non-Sterile, 1" x 4", 2/card, 500 cd/bx



DYNAREX SECURE STRIP WOUND CLOSURE STRIPS

Sterile.

3521 1/8" x 3" Wound Closure Strips, 5/card, 1 card/pch, 50 pch/bx, 4 bx/cs
 3522 1/4" x 1 1/2" Wound Closure Strips, 3/card, 2 cards/pch, 50 pch/bx, 4 bx/cs
 3523 1/4" x 3" Wound Closure Strips, 3/card, 1 card/pch, 50 pch/bx, 4 bx/cs
 3524 1/4" x 4" Wound Closure Strips, 5/card, 2 cards/pch, 50 pch/bx, 4 bx/cs
 3525 1/2" x 4" Wound Closure Strips, 3/card, 2 cards/pch, 50 pch/bx, 4 bx/cs
 3526 1" x 4" Wound Closure Strips, 1/card, 4 cards/pch, 50 pch/bx, 4 bx/cs

ETHICON PROXI-STRIP™ SKIN CLOSURES

101K 1/4" x 3" Skin Closures, 3/pk, 50 pk/bx
 103K 1/4" x 4" Skin Closures, 10/pk, 50 pk/bx
 104K 1/2" x 4" Skin Closures, 6/pk, 50 pk/bx

J&J BAND-AID® BUTTERFLY CLOSURES

Secures small wounds and incisions, and reinforces wound site after removing sutures.

005573 Butterfly Closures, 1 3/4" x 3/8", 10 Count Box, 48/cs



NUTRAMAX BUTTERFLY CLOSURES ADHESIVE BANDAGES

Latex free.

1955000 3/8" x 1 13/16" Adhesive Bandage, Medium, Bulk, 16000/cs
 1975033 3/8" x 1 13/16" Adhesive Bandage, Medium, 100/bx, 12 bx/cs
 1980000 1/2" x 2 3/4" Adhesive Bandage, Large, Bulk, 10800/cs
 1990033-24 1/2" x 2 3/4" Adhesive Bandage, Large, 100/bx, 24 bx/cs

SMITH & NEPHEW LEUKOSTRIP WOUND CLOSURE STRIPS

66002876 38 x 4mm, 1/8" x 1 1/2", 50/bx, 4 bx/cs
 66002878 6.4 x 3.8mm, 1/4" x 3", 50/ctn, 4 ctn/cs
 66002879 6.4 x 102mm, 1/4" x 4", 50/bx, 4 bx/cs
 66002880 102 x 13mm, 1/2" x 4", 50/bx, 4 bx/cs



FIRST AID ONLY SPRAY ON BANDAGE

Quick and convenient, The spray on bandage protects cuts and abrasions while it helps fight infection, Water soluble, and contains no fluorocarbons, so it's environmentally safe, Active Ingredients: Ethyl Alcohol, PVP/VA Copolymer, Benocaine 3.2%, and Benzethonium Chloride 0.2%.

M527 Spray on Bandage, 3 oz. Aerosol Can, 12/cs

J&J BAND-AID® LIQUID BANDAGES

003937 Liquid Bandages, 10/bx, 24 bx/cs



FERNDAL MASTISOL® MEDICAL ADHESIVE

Latex free. Clear, non-irritating adhesive that will secure the most difficult dressings, yet acts as a skin protectant to help reduce allergic reactions. Not water soluble and is excellent for diaphoretic patients, patients with oily skin and heavy beards. Helps the effectiveness of strip skin closures and minimizes the use of sutures.

0523-06 2 oz Dispenser Cap Bottle Medical Adhesive
 0523-15 15ml Unit Dose Medical Adhesive
 0523-16 15ml Spray Medical Adhesive
 0523-48 2/3cc Vials Medical Adhesive, 48/bx

Cotton Rolls/Balls



FERNDAL DETACHOL® ADHESIVE REMOVER

Non-irritating adhesive remover is recommended for removal of Mastisol®. Does not contain acetone, benene, chloroform or ether. Does not sting. Excellent for use around the face as the pleasant fragrance will not irritate the patient's eyes or nose.

0513-04	4 oz Dispenser Cap Adhesive Remover
0513-06	2 oz Dispenser Cap Adhesive Remover
0513-15	15ml Unit Dose Adhesive Remover
0513-48	2/3cc Vials Adhesive Remover, 48/bx

GRAHAM PROFESSIONAL CELLUCOTTON COIL

Cotton beauty coil is ideal for facials, skin care, manicures, and perms.

42251	40 ft Cotton Coil, White, 12/cs
44045	3 Pound Dispenser Box

GRAHAM PROFESSIONAL COTTON ROLL

43360	Cotton Roll 1lb Non-Sterile, 14 r/cs
-------	--------------------------------------

TIDI VENTURE™ COTTON ROLLS

Superior quality Venture™ Cotton Rolls. Highly absorbent cotton construction with an economical price.

919131	Medium, 3/8" x 1 1/2", Cotton Roll, 50/bundle, 40bundle/cs
--------	--

CROSSTEX BRAIDED COTTON ROLLS

Manufactured from pure, short fiber, 100% cotton and held together by thin strands of polyester filament, thus assures a very soft, non-irritating. Made from 100% cotton. Highly absorbent. Starch free for maximum softness. Will not adhere to mucous membranes. Tightly braided to maintain shape when wet. Polywrapped box to protect contents

DNCB	Roll, Non-Sterile, Braided, #2 Medium, 1 1/2" x 3/8", 2000/bx, 10 bx/cs
------	---

DUKAL COTTON ROLL

Made of 100% bleached cotton the roll is easy to unwind, and the short fibers make it simple to tear. Our cotton rolls are packaged in a re-sealable poly drawstring bag for your convenience. They are ideal for cleaning wounds, for application of medications, and where padding is needed.

CR1-12	Cotton Roll, 1 lb., 1 r/bg, 12 bg/cs
--------	--------------------------------------



DYNAREX COTTON ROLLS

Non-sterile.

3250	1 1/2" x 3/8" Cotton Roll, 2000/bx, 10 bx/cs
------	--

COVIDIEN/KENDALL CURITY LAKESIDE COTTON ROLL

Soft and absorbent. Nonsurgical bleached cotton. Non-sterile.

6030	12 1/2" x 56" Cotton Roll, 25/cs
------	----------------------------------

CROSSTEX PREMIUM COTTON ROLLS

Recognized as the industry standard, these cotton rolls are soft, pliable, non-linting, non-irritating - absorbing up to 10 times their weight, thus minimizing the amount needed during a procedure. Made from 100% cotton. No cellulose or rayon fibers. Highly absorbent. Will not adhere to mucous membranes. Maintain shape when wet. Polywrapped box to protect contents

DNC	Roll, Non-Sterile, #2 Medium, 1 1/2" x 3/8", 2000/bx, 12 bx/cs
-----	--

CROSSTEX ADVANTAGE PLUS COTTON ROLLS

Similar to our Premium brand, they are the perfect alternative for the price conscious practice.

DNAC	Roll, Non-Sterile, #2 Medium, 1 1/2" x 3/8", 2000/bx, 12 bx/cs
------	--

CROSSTEX DELUXE COTTON ROLL DISPENSER

This deluxe, heavy duty dispenser provides an aseptic barrier against operatory aerosols and accidental spatter. Top of dispenser opens for easy loading of up to 50 cotton rolls; sliding drawer will dispense up to 4 rolls with a single pull. Non-skid base prevents movement when placed on tabletop or counter. Simplifies the use of cotton rolls while providing maximum protection for you, your staff and your patients. Size: 4" x 3.3" x 2" (10 x 9 x 5cm)

PDCBL	Dispenser, 4" x 3.3" x 2", Blue, 1/ctn
PDCWH	Dispenser, 4" x 3.3" x 2", White, 1/ctn



PRO ADVANTAGE® COTTON BALLS

P159025	Cotton Ball, Medium, 2000/bg, 2 bg/cs
P159045	Cotton Ball, Large, 1000/bg, 2bg/cs

DUKAL COTTON BALLS

Made of 100% bleached cotton. Absorbent and clean. They are available in two sizes, medium and large. Soft and uniform size.

801	Medium Cotton Balls, 2000/bg, 2 bg/cs
802	Large Cotton Balls, 1000/bg, 2 bg/cs



COVIDIEN/KENDALL CURITY COTTON PREPPING BALLS

100% cotton. Soft and absorbent. Designed for a variety of uses throughout the hospital. Packed in convenient poly bags. Non-sterile.

2600	Medium Cotton Balls, 500/bg, 8 bg/cs
2601	Large Cotton Balls, 200/bg, 10 bg/cs
2602	Large Cotton Balls, 50/bg, 80 bg/cs

TIDI COTTON BALLS

969152	Large Cotton Balls, Non-Sterile, 1000/bg, 2000/cs
969153	Medium Cotton Balls, Non-Sterile, 2000/bg, 4000/cs
969164	X-Large Cotton Balls, Non-Sterile, 500/bg, 4 bg/cs

TIDI RAYON BALLS

969162	Non-Sterile Rayon Balls, Large, 1000/bg, 2000/cs
969163	Non-Sterile Rayon Balls, Medium, 2000/bg, 4000/cs



P902015

P902020



P902045

P902050

PRO ADVANTAGE® ALCOHOL PREP PADS

Provides rapid-acting skin antiseptics prior to injection in non-sterile and sterile formulas. Contains 70% isopropyl alcohol.

- P902015 Pro Advantage® Prep Pad, Alcohol, Non Sterile, Medium, 1 3/16" x 2 5/8", 200 pk/bx, 15 bx/cs (Unavailable For Shipment To Canada)
- P902020 Pro Advantage® Prep Pad, Alcohol, Non Sterile, Large, 1.75" x 3.5", 100 pk/bx, 10 bx/cs (Unavailable For Shipment To Canada)
- P902045 Pro Advantage® Prep Pad, Alcohol, Sterile, Medium, 1 3/16" x 2 5/8", 200 pk/bx, 20 bx/cs (Unavailable For Shipment To Canada)
- P902050 Pro Advantage® Prep Pad, Alcohol, Sterile, Large, 1.75" x 3.5", 100 pk/bx, 10 bx/cs (Unavailable For Shipment To Canada)

DUKAL ALCOHOL PADS

Alcohol Prep Pad is a 2-ply isopropyl alcohol saturated pad that is expertly sealed to ensure that the pad remains fully saturated.

- 852 Alcohol Prep Pads, Non-Sterile, 200/bx, 20 bx/cs
- 853 Alcohol Prep Pads, Sterile, 200/bx, 20 bx/cs



DYNAREX PREP PADS

- 1103 Alcohol Prep, Medium, 200/bx, 10 bx/cs
- 1104 Alcohol Prep, Medium, 100/bx, 20 bx/cs
- 1106 Alcohol Prep, Large, 100/bx, 10 bx/cs
- 1108 Iodophor PVP Prep, Medium, 100/bx, 10 bx/cs
- 1113 Prep Pad, Alcohol, Medium, Sterile, 200/bx, 10 bx/cs



COVIDIEN/KENDALL CURITY ALCOHOL PREP PADS

Lighter weight pad material provides lower cost alternative to premium WEBCOL prep pads. 70% isopropyl alcohol. Two ply sponge for cleansing and disinfecting. Sealed in airtight, poly-lined foil pouches.

- 5150 Alcohol Prep Pad, Medium, Non-Sterile, 200/bx, 20 bx/cs
- 5750 Alcohol Prep Pad, Medium, Sterile, 200/bx, 20 bx/cs



COVIDIEN/KENDALL WEBCOL ALCOHOL PREPS

Premium pad material provides maximum absorbency for scrubbing and cleansing. 70% isopropyl alcohol. Sterile, nonwoven sponges designed for moderate cleansing. Economical.

- 5033 WEBCOL Alcohol Preps, Large, 1 ply, 200/bx, 20 bx/cs
- 5110 WEBCOL Alcohol Preps, Large, 2 ply, 200/bx, 20 bx/cs



B35900



B60301



B60307

PDI ALCOHOL PREP PADS

- B35900 Acetone Alcohol Prep Pad, Non-Sterile, Small, 1" x 2 1/2", 1/pk, 100 pk/bx, 10 bx/cs
- B60301 Sterile Alcohol Prep Pad, Medium, 1 3/16" x 2 5/8", 1/pk, 200 pk/bx, 10 bx/cs
- B60307 Sterile Alcohol Prep Pad, Medium, 1 3/16" x 2 5/8", 200/bx, 20 bx/cs
- C69900 Sterile Alcohol Prep Pad, Large, 1 3/4" x 3 1/2", 100/bx, 10 bx/cs



MAJOR PREP PADS

- 700333 Alcohol Prep Pads, 100's

COVIDIEN/KENDALL PREPIES SKIN BARRIER WIPES

Forms protective barrier between skin and adhesive products. Helps protect skin from feces, urine, chemical irritants and gastrointestinal fluids.

- 6560 Skin Barrier Wipe, 2-ply, Large, Non-Sterile, 50/bx, 20 bx/cs

Wipes/Prep Pads/Towelettes - Soaps/Cleaners



GRAHAM PROFESSIONAL ESTHETIC WIPES

4-ply medical grade nonwoven wipes are lint-free. Ideal for a variety of facial treatments.

- 52508 SpaEssentials Wipes 4-Ply Nonwoven, White, 2" x 2", 200/pk, 25 pk/cs
- 52509 SpaEssentials Wipes 4-Ply Nonwoven, White, 4" x 4", 200/pk, 10 pk/cs

TRIAD STING RELIEF PADS

- 10-6400 Sting Relief Pad, 100/bx, 10 bx/cs
- 10-6400B Sting Relief Pad, Bulk, 3000/cs

EMERGING NANOSIL SPRAY GEL™

This cutting edge product is a revolutionary silver antimicrobial spray gel that kills bacteria, yeast and fungi. Nanosil Spray Gel™ forms a protective matrix over the affected area that lasts for hours. Use Nanosil Spray Gel™ to treat infection and promote faster healing for: Burns & Radiation Burns, Decubitus Ulcers & Diabetic Ulcers, Pressure Ulcers & Bed Sores, Cuts, Lacerations & Post-Op Surgical Incisions.

- NSG-04 Silver Antimicrobial Wound Care Spray Gel, 4 oz
- NSG-04CS Silver Antimicrobial Wound Care Spray Gel, 4 oz, 24/cs

3M™ CAVILON™ ANTISEPTIC SKIN CLEANSER

Gently removes of soil without irritating the skin and helps prevent the risk of bacterial contamination associated with exposure to urine, feces or body fluids.

- 3380 8 oz Bottle Antiseptic Skin Cleanser, 12/cs



3M™ REMOVER LOTION

A lotion for removing 3M DuraPrep™ Surgical Solution and other iodophor preps from skin.

- 8610 1/2 oz Packet Remover Lotion, 20/bx, 5 bx/cs
- 8611 4 oz Bottle Remover Lotion, 12/cs



GOJO PROVON® ANTIMICROBIAL SKIN CLEANSER

- 4104-12 w/2% CHG, Antimicrobial Skin Cleanser, 32 fl oz Pump Bottle, 12/cs

SMITH & NEPHEW DERMAL® WOUND CLEANSER

An over-the-counter, pH-balanced, non-toxic, non-irritating, no-rinse skin/wound cleanser. It contains surfactants and moisturizers that gently help remove debris and dressing residue from the wound bed without harsh scrubbing. The cleanser is not soap and does not require rinsing.

- 449000 16 oz. Bottle Wound Cleanser, 12/cs
- 59449200 8 oz. Bottle Wound Cleanser, 12/cs



CERTOL EXCELLE ANTISEPTIC SKIN CLEANSER

4% Chlorhexidine Gluconate Solution that is an effective, rapid antimicrobial skin cleanser. Ideal for surgical handwash.

- EXC128 4" Chlorhexidine Gluconate Skin Cleanser, Gallon, 4/cs



KIMBERLY-CLARK ANTISEPTIC SKIN CLEANSER

Antiseptic Skin Cleanser is a 1.25% PCMX antimicrobial formula that kills germs.

- 93036 Kimcare Baby Powder Scent, 500mL, White, 18/cs



7012

7055

MADA SKIN CLEANSERS

- 7012 MadaFoam, Non-Aerosol Skin Cleaner/Moisturizer, 7 oz, 12/cs
- 7012-6 MadaFoam, Non-Aerosol Skin Cleaner/Moisturizer, 7 oz, 6/cs
- 7055 MadaGel Instant Hand Sanitizer Gel, 4 oz Squeeze Bottle, 24/cs

CROSSTEX SANICLENZ® ANTIMICROBIAL SKIN CLEANSER

Great surgical scrub. Persistent, residual antiseptics. Effective against a broad spectrum of bacteria, fungi & yeast. Fast:CHG acts rapidly. Non-toxic. Aloe Vera, Vitamin E, keratin and Liposome's.

- JSCG Cleanser, Gallon, 4/cs
- JSCP Cleanser, 18 oz., 16/cs



MOLNLYCKE HIBICLENS® ANTISEPTIC ANTIMICROBIAL SKIN CLEANSER

Bonds to skin to sustain its antimicrobial activity for up to six hours. Less-irritating than other antiseptics. Effectiveness increases with use due to a cumulative effect. Non-flammable formula. Not altered in the presence of blood. Antimicrobial skin cleanser helps reduce bacteria that can potentially cause disease. Skin wound and general skin cleansing. Surgical hand scrub. Personnel hand wash. Preoperative skin preparation. Contains 4% CHG.

- 57504 Skin Cleanser, 4 oz Liquid, 48/cs
- 57508 Skin Cleanser, 8 oz Liquid, 24/cs
- 57516 Skin Cleanser, 16 oz Liquid, 12/cs
- 57517 Skin Cleanser, 15mL Packets, Liquid, 50/ctn, 8 ctn/cs
- 57532 Skin Cleanser, 32 oz Liquid, 12/cs
- 57591 Skin Cleanser, Gallons Liquid, 4/cs

Accessories

- 59901 Hand Pump For 16 oz Hibiclens (57516)
- 59902 Hand Pump For 32 oz Hibiclens (57532)
- 59903 Hand Pump For Gallon Hibiclens (57591)

MOLNLYCKE HIBICLENS® DISPENSER & ACCESSORIES

- 59994 Foot Wall Dispenser
- 59995 Empty Globe, 12/cs

J&J BAND-AID® WOUND CLEANSING

- 003939 2.3 oz. One-Step Cleansing + Infection Protection Foam, 6/bx, 4 bx/cs
- 004459 6 oz. Hurt-Free™ Antiseptic Wash, 6/bx, 4 bx/cs

3M™ WOUND CLEANSER

- 91101 Wound Cleanser, 4 oz Squeeze Bottle, 12/cs



HARTMANN-CONCO FLEXI-TUBE® TUBULAR GAUZE

Provides protection to delicate skin. Stretches in width directions. Easy to apply, just cut to exact length needed. Holds primary dressings in place. Also used as a first layer for lymphedema bandaging. Latex free.

- 85010000 2.4" x 21.9 yds Tubular Gauze
- 85020000 3.1" x 21.9 yds Tubular Gauze
- 85030000 3.9" x 21.9 yds Tubular Gauze
- 85040000 4.7" x 21.9 yds Tubular Gauze



HARTMANN-CONCO STAT-WRAP® COHESIVE GAUZE BANDAGE

Combines the best features of a quality gauze bandage with a cohesive bandage to make it easy to use and economical in time and expense. Can be comfortably worn for long periods of time and never adheres to skin or hair. Stays in place with less chance of rolling or constricting - even in wet environments. The cohesive properties minimize linting and loose threads and make even the most difficult-to-bandage areas easy to wrap. Adheres only to itself without tape. Secures splints, I.V.'s, leads and catheters. CAUTION: This product contains natural rubber latex which may cause allergic reactions.

- 60100000 1" x 3 yd Bandage, 16 rl/bx, 5 bx/cs
- 60200000 2" x 3 yd Bandage, 8 rl/bx, 5 bx/cs
- 60300000 3" x 4.5 yd Bandage, 8 rl/bx, 5 bx/cs
- 60400000 4" x 4.5 yd Bandage, 8 rl/bx, 4 bx/cs



PRO ADVANTAGE® GAUZE SPONGES - NON-STERILE

Made from 100% premium bleached cotton. These highly absorbent sponges are available in a broad range of sizes and plus. Applicable for dressing, general cleaning, prepping, packing and debriding wounds. Available in sterile and non-sterile.

- P157112 Gauze Sponge, 2" x 2", 12-Ply, Non-Sterile, 200/sleeve, 40 slv/cs
- P157113 Gauze Sponge, 3" x 3", 12-Ply, Non-Sterile, 200/sleeve, 20 slv/cs
- P157114 Gauze Sponge, 4" x 4", 12-Ply, Non-Sterile, 200/sleeve, 10 slv/cs
- P157115 Gauze Sponge, 2" x 2", 8-Ply, Non-Sterile, 200/sleeve, 25 slv/cs
- P157116 Gauze Sponge, 3" x 3", 8-Ply, Non-Sterile, 200/sleeve, 20 slv/cs
- P157117 Gauze Sponge, 4" x 4", 8-Ply, Non-Sterile, 200/sleeve, 20 slv/cs
- P157118 Gauze Sponge, 4" x 4", 16-Ply, Non-Sterile, 200/sleeve, 10 slv/cs

Dressings/Gauze



P157002

P157004

PRO ADVANTAGE® GAUZE SPONGES - STERILE

Made from 100% premium bleached cotton. These highly absorbent sponges are available in a broad range of sizes and plys. Applicable for dressing, general cleaning, prepping, packing and debriding wounds. Available in sterile and non-sterile. Packaged sterile.

P157002	Gauze Sponge, 1's, 2" x 2", 12-Ply, Sterile, 100/bx, 24 bx/cs
P157003	Gauze Sponge, 1's, 3" x 3", 12-Ply, Sterile, 100/bx, 24 bx/cs
P157004	Gauze Sponge, 1's, 4" x 4", 12-Ply, Sterile, 100/bx, 12 bx/cs
P157022	Gauze Sponge, 2's, 2" x 2", 8-Ply, Sterile, 50/tray, 30 tray/cs
P157023	Gauze Sponge, 2's, 3" x 3", 12-Ply, Sterile, 40/tray, 30 tray/cs
P157024	Gauze Sponge, 2's, 4" x 4", 8-Ply, Sterile, 25/tray, 24 tray/cs
P157025	Gauze Sponge, 2's, 4" x 4", 12-Ply, Sterile, 25/tray, 24 tray/cs

AMD-RITMED X-RAY GAUZE SPONGES

84440	Sponge, Sterile, Non-Strung, X-Ray, 2" x 2", 5/pk, 40 pk/cs
-------	---

AMD-RITMED PREMIUM GAUZE SPONGES

A4120	4" x 4", 12 Ply, Sterile 2's, 50/pk, 24pk/cs
A4121	4" x 4", 12 Ply, Sterile 10's, 10/tray, 12 trays/cs
A4121T	4" x 4", 12 Ply, Sterile 10's, 10/Plastic tray, 128tray/cs
A8120	8" x 4", 12 Ply, Sterile 2's, 50/pk, 20pk/cs



A2337



A7632

AMD-RITMED VITAL-GAUZE MULTI-PURPOSE GAUZE SPONGES

A2337	4" x 4" Gauze Sponges, 4-Ply, Sterile 2s, 50/pk, 24 pk/cs
A2338	2" x 2" Gauze Sponges, 4-Ply, Sterile 2s, 50/pk, 60 pk/cs
A2339	3" x 3" Gauze Sponges, 4-Ply, Sterile 2s, 50/pk, 48 pk/cs
A7632	2" x 2" Gauze Sponges, 4-Ply, Non-Sterile, 200/pk, 40 pk/cs
A7633	3" x 3" Gauze Sponges, 4-Ply, Non-Sterile, 200/pk, 20 pk/cs
A7634	4" x 4" Gauze Sponges, 4-Ply, Non-Sterile, 200/pk, 10 pk/cs



A4320



A4420

AMD-RITMED WOUND CARE GAUZE SPONGES

A3301	3" x 3" Gauze Sponge w/cellulose fill, Non-Sterile, 100/pk, 40 pk/cs
A4301	4" x 3" Gauze Sponge w/cellulose fill, Non-Sterile, 100/pk, 20 pk/cs
A4320	4" x 3" Gauze Sponge w/cellulose fill, Sterile 2s, 50/pk, 24 pk/cs
A4401	4" x 4" Gauze Sponge w/cellulose fill, Non-Sterile, 100/pk, 20 pk/cs
A4420	4" x 4" Gauze Sponge w/cellulose fill, Sterile 2s, 50/pk, 24 pk/cs
A4801	8" x 4" Gauze Sponge w/cellulose fill, Non-Sterile, 100/pk, 10 pk/cs

COVIDIEN/KENDALL CURITY GAUZE PADS

Individually wrapped, 12-ply sterile gauze pads. Perforated carton for easy office dispensing.

3381	Sterile Gauze Pads, Peel Back Package, 2" x 2", 12 Ply, 100/bx, 24 bx/cs
6132	Sterile Gauze Pads, Peel-Back Package, 3" x 3", 12 ply, 100/bx, 24 bx/cs
6309	4" x 4" Gauze Pads, 12-Ply, Sterile 1s, 100/bx, 12 bx/cs

**COVIDIEN/KENDALL CURITY COVER SPONGES**

Filled with layer of cellulose wadding covered with nonwoven fabric. Economical, all-purpose sponge with good absorbent capacity. Nonsterile sponges conveniently indexed in 10's for faster, easier handling.

1700	Cover Sponges, Nonsterile, Bulk, 3" x 3", 100/bg, 40 bg/cs
1713	Cover Sponges, Non-Sterile, Bulk, 3" x 4", 100/bg, 20 bg/cs
1792	Cover Sponges, Non-Sterile, Bulk, 4" x 4", 100/bg, 20 bg/cs
2913	Cover Sponges, Sterile 2s in Peel-Back Package, 4" x 4", 50/bg, 24 bg/cs
3157	Cover Sponges, Sterile 2s in Peel-Back Package, 4" x 4", 50/bg, 24 bg/cs


**COVIDIEN/KENDALL KERLIX 4 X 4 SPONGES**

Prewashed, fluff-dried woven gauze with crinkle-weave pattern for loft and bulk to cushion and protect wound areas. Ideal for all-purpose floor use as well as superior primary dressing. Absorbency and texture ideal for prepping, cleansing, packing and debriding wounds.

4032	Kerlix Sponge, Bulk, 12-Ply, 2000/cs
4588	Kerlix Sponge, 16-Ply, 1280/cs
5042	Kerlix Sponge, Bulk, 16-Ply, 2000/cs
5072	4 x 4 Sponges, Sterile 2's in Peel-Back Package, 12 ply, 50/tray, 124 tray/cs
6120	Kerlix Sponge, 12-Ply, 1280/cs



2252



3033

COVIDIEN/KENDALL CURITY GAUZE SPONGES

Made from high-grade 100% US grown. Exceeds USP Type VII gauze requirements. Well suited for variety of applications; wound dressing, general cleaning, minor prepping, wound packing and debriding. Nonsterile sponges conveniently indexed in 10's for faster, easier handling.

1806	2" x 2" Gauze Sponges, 8-Ply, Sterile 2s, 100/tray, 30 tray/cs
1903	3" x 3" Gauze Sponges, 12-Ply, Sterile 2s in Peel-Back Package, 50/bg, 48 bg/cs
2029	3" x 4" Gauze Sponges, 12-Ply, Sterile 2s in Peel-Back Package, 50/bg, 24 bg/cs
2146	2" x 2" Gauze Sponges, 8-Ply, Non-Sterile, 200/bg, 25 bg/cs
2187	4" x 4" Gauze Sponges, 8-Ply, Sterile 2s in Peel-Back Package, 50/bg, 24 bg/cs
2252	2" x 2" Gauze Sponges, 12-Ply, Non-Sterile, Bulk, 200/bg, 40 bg/cs
2259	4" x 8" Gauze Sponges, 12-Ply, Sterile 2s in Peel-Back Package, 50/bg, 20 bg/cs
2346	3" x 3" Gauze Sponges, 12-Ply, Non-Sterile, 200/bg, 20 bg/cs
2438	3 1/2" x 4" Gauze Sponges, 32-Ply, Non-Sterile, 200/bg, 5 bg/cs
2556	4" x 4" Gauze Sponges, 8-Ply, Non-Sterile, 200/bg, 20 bg/cs
2634	4" x 4" Gauze Sponges, 12-Ply, Non-Sterile, 200/bg, 10 bg/cs
2733	4" x 4" Gauze Sponges, 16-Ply, Non-Sterile, 200/bg, 10 bg/cs
2835	4" x 8" Gauze Sponges, 12-Ply, Non-Sterile, 200/bg, 10 bg/cs
3033	4" x 4" Gauze Sponges, 12-Ply, Sterile 2s, 50/tray, 24 tray/cs
3967	4" x 4" Gauze Sponges, 12-Ply, Sterile 5s in Soft Pouch, 100/bg, 16 bg/cs
3968	4" x 4" Gauze Sponges, 8-Ply, Sterile 10s in Soft Pouch, 100/bg, 16 bg/cs
3971	4" x 4" Gauze Sponges, 12-Ply, Sterile 10s in Soft Pouch, 100/bg, 16 bg/cs
3973	4" x 4" Gauze Sponges, 16-Ply, Sterile 10s in Soft Pouch, 100/bg, 16 bg/cs
6318	4" x 4" Gauze Sponges, 8-Ply, Sterile 10s in Plastic Tray, 1280/cs
6939	4" x 4" Gauze Sponges, 12-Ply, Sterile 10s in Plastic Tray, 1280/cs
7267	3" x 4" Gauze Sponges, 16-Ply, Non-Sterile, 200/bg, 10 bg/cs
7605	4" x 4" Gauze Sponges, 16-Ply, Sterile 10s in Plastic Tray, 1280/cs

**Please ask before
ordering a product based
on the pictured item.**

**COVIDIEN/KENDALL VERSALON NONWOVEN ALL-PURPOSE SPONGES**

Ideal for applying ointments, prepping, wiping needles, cleaning slides and more. Rayon/polyester blend construction. Cost effective. Highly absorbent. Virtually lint-free.

8042	Versalon All Purpose Sponge, Peel Back Pkg, 2" x 2", 3000/cs
8043	Versalon All Purpose Sponge, Peel Back Pkg, 3" x 3", 2400/cs
8044	4" x 4" Sponges, 4-Ply, Sterile 2s, 50/tray, 12 tray/cs
8045	Versalon All Purpose Sponge, Soft Pouch, 4" x 4", 800/cs
8046	Versalon All Purpose Sponge, Peel Back Pkg, 4" x 3", 1200/cs
8047	Nonwoven All-Purpose Sponges, Sterile 10's in Soft Pouch, 4" x 4", 4 ply, 1600/cs
9022	Versalon All Purpose Sponge, 4-Ply, 2" x 2", 200/bg, 40 bg/cs
9023	Nonwoven All-Purpose Sponges, Non-Sterile, Bulk, 3" x 3", 4 ply, 200/bg, 20 bg/cs
9024	4" x 4" Sponges, 4-Ply, Non-Sterile, 200/bg, 10 bg/cs
9026	Versalon All Purpose Sponge, 4-Ply, 4" x 3", 2000/cs
9132	Nonwoven All-Purpose Sponges, Non-Sterile, Bulk, 2" x 2", 3 ply, 200/bg, 25 bg/cs
9134	Versalon All Purpose Sponge, 3-Ply, 4" x 4", 4000/cs

**COVIDIEN/KENDALL VISTEC X-RAY DETECTABLE SPONGES**

USP Type VII gauze. Unique dual-colored element combines blue radiopaque strand with white strand which is visible when sponge is saturated with blood. Element is latex free and runs through both sides of sponge to facilitate counting. Precise, double vacuum delinting process assures low lint. Sponges are paper banded in 10's, utilizing state of the art electronic counting eyes to help ensure a count of 10 sponges time after time.

2979	Vistec X-Ray Detectable Sponge, 8" x 4", 2000/cs
4334	Vistec X-Ray Detectable Sponge, 3/5" x 4", 32-Ply, 1000/cs
7114	Vistec X-Ray Detectable Sponge, 4" x 4", 2000/cs
7148	Vistec X-Ray Detectable Sponge, 4" x 4", 16-Ply, 2000/cs
7317	X-Ray Detectable Sponges, Sterile, Banded 10's - Trays, 4" x 4", 16 ply, 1280/cs
7318	Vistec X-Ray Detectable Sponge, 4" x 8", 12-Ply, 800/cs
7319	Vistec X-Ray Detectable Sponge, 4" x 8", 16-Ply, 640/cs
7320	Vistec X-Ray Detectable Sponge, 3.5" x 4", 32-Ply, 480/cs
7321	Vistec X-Ray Detectable Sponge, 4" x 8", 24-Ply, 800/cs
7462	Vistec X-Ray Detectable Sponge, 4" x 8", 12-Ply, 2000/cs
7463	Vistec X-Ray Detectable Sponge, 4" x 8", 16-Ply, 1000/cs
7466	Vistec X-Ray Detectable Sponge, 4" x 8", 24-Ply, 1000/cs

Dressings/Gauze

CROSSTEX NON-WOVEN 4 PLY SPONGES

Unique rayon/polyester blend. Provides maximum absorbency. Less adhesion to wounds for faster healing. Soft, smooth, strong, virtually lint free. Ideal for sustained patient use, disinfecting, general cleaning.

ENC3NW	3" x 3", Sponge, Premium, 40g, Non-Sterile, 4000/cs
ENC3NWLA	3" x 3", Sponge, Advantage, 35g, Non-Sterile, 4000/cs
ENC4NW	4" x 4", Sponge, Premium, 40g, Non-Sterile, 2000/cs
ENC4NWLA	4" x 4", Sponge, Advantage, 30g, Non-Sterile, 2000/cs
ENCNW	2" x 2", Sponge, Premium, 50g, Non-Sterile, 5000/ctn, 4 ctn/cs
ENCNW4K	2" x 2", Sponge, Premium, 50g, Non-Sterile, 4000/cs
ENCNWL	2" x 2", Sponge, Advantage, 35g, Non-Sterile, 5000/ctn, 4 ctn/cs
ENCNWLA	2" x 2", Sponge, Advantage, 30g, Non-Sterile, 5000/cs
ESNW2	2" x 2", Sponge, 35g, Sterile, 25/bx, 60 bx/cs
ESNW3	3" x 3", Sponge, 35g, Sterile, 25/bx, 48 bx/cs
ESNW4	4" x 4", Sponge, 35g, Sterile, 25/bx, 24 bx/cs

CROSSTEX ULTRA GAUZE™ NON-WOVEN SPONGES

More absorbent than traditional sponges. Releases easily from incision/extraction sites. Will not stick to surface.

ENC3NWU	3" x 3", Sponge, 4 Ply, Non-Sterile, 4000/cs
ENC4NWU	4" x 4", Sponge, 4 Ply, Non-Sterile, 2000/cs
ENCNWU	2" x 2", Sponge, 4 Ply, Non-Sterile, 5000/cs
ENCNWUB	2" x 2", Sponge, 4 Ply, Non-Sterile, 8000/cs

CROSSTEX ALL GAUZE PREMIUM NON-STERILE SPONGES

Made from USP 100% gauze. Available in 8 ply or 16 ply. Highly absorbent.

ENC2	2" x 2", Sponge, 8 Ply, 5000/cs
ENC212	2" x 2", Sponge, 12 Ply, 8000/cs
ENC312	3" x 3", Sponge, 12 Ply, 4000/cs
ENC38	3" x 3", Sponge, 8 Ply, 4000/cs
ENC412	4" x 4", Sponge, 12 Ply, 2000/cs
ENC48	4" x 4", Sponge, 8 Ply, 4000/cs

DERMA SCIENCES STERILE PADS

7060033	2" x 2" Gauze Pad, Sterile, 12-Ply, 100/bx, 12 bx/cs
7070033	3" x 3" Gauze Pad, Sterile, 12-Ply, 100/bx, 12 bx/cs
7080033	4" x 4" Gauze Pad, Sterile, 12-Ply, 100/bx, 12 bx/cs
7205327	2" x 2" Non-Woven Pad, Sterile, 4-Ply, 10/bx, 36 bx/cs
7206327	2" x 2" Non-Woven Pad, Sterile, 4-Ply, 25/bx, 24 bx/cs
7215327	3" x 3" Non-Woven Pad, Sterile, 4-Ply, 25/bx, 24 bx/cs
7225327	4" x 4" Non-Woven Pad, Sterile, 4-Ply, 25/bx, 24 bx/cs
7226327	4" x 4" Non-Woven Pad, Sterile, 4-Ply, 10/bx, 24 bx/cs
7260033	2" x 2" Non-Woven Pad, Sterile, 4-Ply, 100/bx, 12 bx/cs
7270033	3" x 3" Non-Woven Pad, Sterile, 4-Ply, 100/bx, 12 bx/cs
7280033	4" x 4" Non-Woven Pad, Sterile, 4-Ply, 100/bx, 12 bx/cs
7515327	3" x 3" Non-Woven Pad, Sterile, 4-Ply, 10/bx, 36 bx/cs

DERMA SCIENCES DUCARE® WOVEN SPONGES

Made from 100% USP Type VII cotton; the gold standard for general use or for debriding, covering and packing wounds. Comes in a variety of sizes and are available sterile in packages of 2 and non-sterile in sleeves of 200. 4x4's also available in packages of 10, packed in trays. Latex free.

82208	Sponge, Gauze, Woven, 2" x 2", 8-Ply, Sterile, 2's, 50 pkg/bx, 30 bx/cs (3000 ea/cs)
82312	Sponge, Gauze, Woven, 3" x 3", 12-Ply, Sterile, 2's, 25 pkg/bx, 48 bx/cs (2400 ea/cs)
82412	Sponge, Gauze, Woven, 4" x 4", 12-Ply, Sterile, 2's, 25 pkg/bx, 24 bx/cs (1200 ea/cs)
90212	Sponge, Gauze, Woven, 2" x 2", 12-Ply, Non-Sterile, 200/slv, 40 slv/cs (8000 ea/cs)
90412	Sponge, Gauze, Woven, 4" x 4", 12-Ply, Non-Sterile, 200/slv, 10 slv/cs (2000 ea/cs)



1212



1312



1412

DUKAL GAUZE PADS

Made of 100% woven cotton, the pads are individually wrapped for your convenience when only one sponge is needed. This highly absorbent and versatile product will meet all market needs and applications.

1212	Gauze Pad, 2" x 2", 12-Ply, 100/bx, 24 bx/cs
1212-10	Gauze Pad, 2" x 2", 12-Ply, 10/bx, 60 bx/cs
1212-25	Gauze Pad, 2" x 2", 12-Ply, 25/bx, 24 bx/cs
1312	Gauze Pad, 3" x 3", 12-Ply, 100/bx, 24 bx/cs
1313-10	Gauze Pad, 3" x 3", 12-Ply, 10/bx, 60 bx/cs
1313-25	Gauze Pad, 3" x 3", 12-Ply, 25/bx, 24 bx/cs
1412	Gauze Pad, 4" x 4", 12-Ply, 100/bx, 12 bx/cs
1414-10	Gauze Pad, 4" x 4", 12-Ply, 10/bx, 60 bx/cs
1414-25	Gauze Pad, 4" x 4", 12-Ply, 25/bx, 24 bx/cs



6114



6122



6124

DUKAL CLINISORB NON-WOVEN SPONGES

The New Sponge is an excellent example of the DUKAL quality. The effective combination of the weight and open aperture allows each sponge to wick quicker and offer superior absorbency. They are made of a soft rayon/poly blend and are available in both sterile and non-sterile packaging.

6112	Sponge, 2" x 2", Non-Woven, Non-Sterile, 4-Ply, 200/bg, 40 bg/cs
6113	Sponge, 3" x 3", Non-Woven, Non-Sterile, 4-Ply, 200/bg, 20 bg/cs
6114	Sponge, 4" x 4", Non-Woven, Non-Sterile, 4-Ply, 200/bg, 10 bg/cs
6122	Sponge, 2" x 2", Non-Woven, Sterile, 4-Ply, 2/pk, 25 pk/bx, 60 bx/cs
6123	Sponge, 3" x 3", Non-Woven, Sterile, 4-Ply, 2/pk, 25 pk/bx, 48 bx/cs
6124	Sponge, 4" x 4", Non-Woven, Sterile, 4-ply, 2/pk, 25 pk/bx, 24 bx/cs
6125	Sponge, 4" x 4", Non-Woven, Sterile, 4-ply, 5/pk, 20 pk/bx, 8 bx/cs



7224P



7334P

DUKAL PREMIUM NON-WOVEN SPONGES

One of the thickest and softest on the market. This rayon/poly blend is more absorbent than our standard non-woven sponges. The heavier fabric offers quick wicking for maximum performance which results in a reduction of the number of sponges used. Ideal for delicate procedures and sensitive skin.

- 7224P Sponge, 2" x 2" Non-Woven, Premium, Non-Sterile, 4-Ply, 2/pk, 50 pk/bx, 30 bx/cs
- 7334P Sponge, 3" x 3" Non-Woven, Premium, Non-Sterile, 4-Ply, 2/pk, 40 pk/bx, 30 bx/cs
- 7444P Sponge, 4" x 4" Non-Woven, Premium, Non-Sterile, 4-Ply, 2/pk, 50 pk/bx, 12 bx/cs



2102



7224



7334



7444

DUKAL NON-WOVEN GAUZE SPONGES

A rayon polyester blend of non-woven material offers a low lint, absorbable substitute to all cotton sponges.

- 2102 Gauze Sponge, 2" x 2", Non-Sterile, 4-Ply, 200/bg, 20 bg/cs
- 2103 Gauze Sponge, 3" x 3", Non-Sterile, 4-Ply, 200/bg, 20 bg/cs
- 2104 Gauze Sponge, 4" x 4", Non-Sterile, 4-Ply, 200/bg, 10 bg/cs
- 7224 Gauze Sponge, 2" x 2", Sterile, 4-Ply, 2/pk, 50 pk/bx, 30 bx/cs
- 7334 Gauze Sponge, 3" x 3", Sterile, 4-Ply, 2/pk, 40 pk/bx, 30 bx/cs
- 7444 Gauze Sponge, 4" x 4", Sterile, 4-Ply, 2/pk, 50 pk/bx, 12 bx/cs
- 7454SP Gauze Sponge, 4" x 4", Sterile, 4-Ply, 5/pk, 20 pk/bx, 80 bx/cs

TRUST



DUKAL's Quality, Value and Performance can be seen in all our wound care & patient care products for all professional markets!



Caliber™ Island Dressings

Caliber™ Non-Woven Island Dressings provide the perfect combination of protection and comfort. The absorbent pad collects excess exudates and protects the wound from trauma to allow for quick healing. Available in 6 sizes, Sterile and Latex Free.

Transparent Semi-Permeable Dressings

The DUKAL Transparent Dressing with easy 2 step application is an excellent dressing for IV sites, minor burns, and wounds. It keeps contamination out and promotes healing while allowing for moisture vapor transfer. This DUKAL quality, clear, adhesive backed dressing is ideal for IV sites and to observe vaccination sites. This product features a removable patient label for the dressing or chart. Compare to Tegaderm®




Caliber™ Adhesive Bandages

DUKAL's **Caliber™ Adhesive Bandages** are available in many styles and sizes to meet all of your market needs. Durable plastic, sheer, or flexible fabric with a soft cushion pad make **Caliber™ Adhesive Bandages** the best in protection and strength. All bandages are sterile and latex free.

Wound Closure Strips

The **Wound Closure Strips** can be used in conjunction with or as an excellent alternative to sutures or staples. Each DUKAL quality strip consists of thin non-woven reinforced fibers that are arranged in a bilateral pattern for high tensile strength, to provide incision support and maintain closure. This adhesive backed product is translucent, pourous, and breathable to promote good healing.

* Compare to Steri-Strips®




Quality, Value & Performance is our PROMISE

DUKAL CORPORATION • www.dukal.com



Dressings/Gauze

**DUKAL WOVEN COTTON GAUZE SPONGES**

Highly absorbent sponges. Made of 100% woven cotton and ideal for a number of applications. Our high quality sponges can be used for general cleaning, dressings, prepping, packing and debriding wounds. Available in sterile and non-sterile.

2028H	Gauze Sponge, 2" x 2", Non-Sterile, Hospital Fold, 8-Ply, 200/bg, 25 bg/cs
2085	Gauze Sponge, 2" x 2", Non-Sterile, 8-Ply, 200/bg, 25 bg/cs
2085H	Gauze Sponge, 2" x 2", Non-Sterile, 8-Ply, Hospital Fold, 200/bg, 25 bg/cs
2128	Gauze Sponge, 2" x 2", Non-Sterile, 12-Ply, 200/bg, 40 bg/cs
2128-100	Gauze Sponge, 2" x 2", Non-Sterile, 12-Ply, 100/bg, 80 bg/cs
2128H	Gauze Sponge, 2" x 2", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 40 bg/cs
2208S	Gauze Sponge, 2" x 2", Sterile, 8-Ply, Cotton Filled, 200/bg, 25 bg/cs
2212H	Gauze Sponge, 2" x 2", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 40 bg/cs
2280	Gauze Sponge, 2" x 2", Non-Sterile, 8-Ply, Cotton Filled, 200/bg, 25 bg/cs
2283	Gauze Sponge, 2" x 2", Sterile, 8-Ply, Cotton Filled, 200/bg, 25 bg/cs
3084	Gauze Sponge, 3" x 3", Non-Sterile, 8-Ply, 200/bg, 20 bg/cs
3084H	Gauze Sponge, 3" x 3", Non-Sterile, 8-Ply, Hospital Fold, 200/bg, 20 bg/cs
3124	Gauze Sponge, 3" x 3", Non-Sterile, 12-Ply, 200/bg, 20 bg/cs
3124-100	Gauze Sponge, 3" x 3", Non-Sterile, 12-Ply, 100/bg, 40 bg/cs
3124H	Gauze Sponge, 3" x 3", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 20 bg/cs
3308H	Gauze Sponge, 3" x 3", Non-Sterile, 8-Ply, Hospital Fold, 200/bg, 20 bg/cs
3312H	Gauze Sponge, 3" x 3", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 20 bg/cs
408-10	Gauze Sponge, 4" x 4", Sterile Tray, 8-Ply, 10/tray, 128 tray/cs
4084	Gauze Sponge, 4" x 4", Non-Sterile, 8-Ply, 200/bg, 20 bg/cs
4084H	Gauze Sponge, 4" x 4", Non-Sterile, 8-Ply, Hospital Fold, 200/bg, 20 bg/cs
412-10	Gauze Sponge, 4" x 4", Sterile Tray, 12-Ply, 10/tray, 128 tray/cs
4122	Gauze Sponge, 4" x 4", Non-Sterile, 12-Ply, 200/bg, 10 bg/cs
4122-100	Gauze Sponge, 4" x 4", Non-Sterile, 12-ply, 100/bg, 20 bg/cs
4122H	Gauze Sponge, 4" x 4", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 10 bg/cs
416-10	Gauze Sponge, 4" x 4", Sterile Tray, 16-Ply, 10/try, 128 tray/cs
4162	Gauze Sponge, 4" x 4", Non-Sterile, 16-Ply, 200/bg, 10 bg/cs
421-10	Gauze Sponge, 8" x 4", Sterile, 12-Ply, Tray, 10/tray, 80 tray/cs
4408H	Gauze Sponge, 4" x 4", Non-Sterile, 8-Ply, Hospital Fold, 200/bg, 20 bg/cs
4412H	Gauze Sponge, 4" x 4", Non-Sterile, 12-Ply, Hospital Fold, 200/bg, 10 bg/cs
5228	Gauze Sponge, 2" x 2", Sterile, 8-Ply, Cotton Filled, 2/pk, 50 pk/bx, 24 bx/cs
6208	Gauze Sponge, 2" x 2", Sterile, 8-Ply, 2/pk, 50 pk/bx, 30 bx/cs
6312	Gauze Sponge, 3" x 3", Sterile, 12-Ply, 2/pk, 50 pk/bx, 24 bx/cs
6408	Gauze Sponge, 4" x 4", Sterile, 8-Ply, 2/pk, 50 pk/bx, 12 bx/cs
6412	Gauze Sponge, 4" x 4", Sterile, 12-Ply, 2/pk, 25 pk/bx, 24 bx/cs
6842	Gauze Sponge, 8" x 4", Sterile, 12-Ply, 2/pk, 25 pk/bx, 20 bx/cs
8412	Gauze Sponge, 8" x 4", Non-Sterile, 12-Ply, 200/bg, 10 bg/cs

**DUKAL X-RAY DETECTABLE GAUZE SPONGES**

Constructed of only the highest quality gauze. Automated manufacturing ensures correct count.

416-2X	Gauze Sponge, 4" x 4" X-Ray Detectable, Sterile Tray, 16-Ply, 10/tray, 128 tray/cs
432-2X	Gauze Sponge, 4" x 4" X-Ray Detectable, Sterile Tray, 32-Ply, 10/tray, 30 tray/cs
4412-2X	Gauze Sponge, 4" x 4" X-Ray Detectable, Non-Sterile Tray, 12-Ply, Banded 10s, 100/bg, 20 bg/cs
4416-2X	Gauze Sponge, 4" x 4" X-Ray Detectable, Non-Sterile Tray, 16-Ply, Banded 10s, 100/bg, 20 bg/cs
4432-2X	Gauze Sponge, 4" x 4" X-Ray Detectable, Non-Sterile, 32-Ply, Banded, 50/pk, 20 pk/cs (Hermitage 9466T & 9466)
4812-2X	Gauze Sponge, 8" x 4" X-Ray Detectable, Non-Sterile, 12-Ply, Banded, 100/pk, 20 pk/cs
482-2X	Gauze Sponge, 8" x 4" X-Ray Detectable, Sterile Tray, 12-Ply, 10/tray, 48 tray/cs
4824-2X	Gauze Sponge, 8" x 4" X-Ray Detectable, Non-Sterile, 24-Ply, Banded, 50/pk, 20 pk/cs
486-2X	Gauze Sponge, 8" x 4" X-Ray Detectable, Sterile Tray, 16-Ply, 10/tray, 48 tray/cs
488-2X	Gauze Sponge, 8" x 4", X-Ray Detectable, Sterile, Tray, 24-Ply, 10/tray, 48 tray/cs

DUKAL BASIC GAUZE SPONGES

An economical version of DUKAL Gauze Sponges. 100% woven cotton sponges in both sterile and non-sterile packaging.

8500	Gauze Sponge, 2" x 2", Sterile, 8-Ply, 2/pk, 50 pk/bx, 30 bx/cs
8501	Gauze Sponge, 3" x 3", Sterile, 12-Ply, 2/pk, 50 pk/bx, 24 bx/cs
8502	Gauze Sponge, 4" x 4", Sterile, 8-Ply, 2/pk, 50 pk/bx, 12 bx/cs
8503	Gauze Sponge, 4" x 4", Sterile, 12-Ply, 2/pk, 50 pk/bx, 24 bx/cs
8505	Gauze Sponge, 2" x 2", Non-Sterile, 8-Ply, 200/bg, 25 bg/cs
8506	Gauze Sponge, 2" x 2", Non-Sterile, 12-Ply, 200/bg, 40 bg/cs
8507	Gauze Sponge, 3" x 3", Non-Sterile, 8-Ply, 200/bg, 20 bg/cs
8508	Gauze Sponge, 3" x 3", Non-Sterile, 12-Ply, 200/bg, 20 bg/cs
8509	Gauze Sponge, 4" x 4", Non-Sterile, 8-Ply, 200/bg, 20 bg/cs
8510	Gauze Sponge, 4" x 4", Non-Sterile, 12-Ply, 200/bg, 10 bg/cs
8511	Gauze Sponge, 4" x 4", Non-Sterile, 16-Ply, 200/bg, 10 bg/cs

**HARTMANN-CONCO ECONOLUX® GAUZE SPONGES**

Woven with high quality yarn to meet low-linting, high-value needs. Ideal for use in wound cleansing, mechanical and enzyme debriding of wounds, wound packing and swabbing. 100% woven cotton. Folded edges prevent linting and unraveling. Easy opening packages.

416103	2" x 2" Gauze Sponge, 8-Ply, Sterile 2s, 2/pk, 50 pk/bx, 30 bx/cs
416104	4" x 4" Gauze Sponge, 8-Ply, Sterile 2s, 2/pk, 50 pk/bx, 12 bx/cs
416105	4" x 4" Gauze Sponge, 12-Ply, Sterile 2s, 2/pk, 25 pk/bx, 24 bx/cs
416814	2" x 2" Gauze Sponge, 8-Ply, Non-Sterile, 200/slv, 25 slv/cs
416815	2" x 2" Gauze Sponge, 12-Ply, Non-Sterile, 200/slv, 40 slv/cs
416816	4" x 4" Gauze Sponge, 8-Ply, Non-Sterile, 200/slv, 20 slv/cs
416817	4" x 4" Gauze Sponge, 12-Ply, Non-Sterile, 200/slv, 10 slv/cs
416818	4" x 4" Gauze Sponge, 16-Ply, Non-Sterile, 200/slv, 10 slv/cs



416100



56810000



56880000

HARTMANN-CONCO STERILUX® PREMIUM GAUZE SPONGES

Made of premium quality gauze with a high yarn count and non-raveling folded edges to provide a super absorbent, soft and low-linting product. Ideal as wound dressings and for use in wound cleaning and debridement, wound packing and general wound care. A low-linting manufacturing process and folded edges ensure debris-free packaging. An easy-to-open packaging envelope is also coated to prevent dust or other contaminants. Latex free.

416100	4" x 4" Gauze Sponge, 8-Ply, 10/tray, 128 tray/cs
416101	4" x 4" Gauze Sponge, 12-Ply, 10/tray, 128 tray/cs
416102	4" x 4" Gauze Sponge, 16-Ply, 10/tray, 128 tray/cs
416106	4" x 4" Gauze Sponge, X-Ray Detectable, 16-Ply, 10/tray, 135 tray/cs
56810000	2" x 2" Non-Sterile Gauze Sponges, 8-Ply, 200/slv, 25 slv/cs
56820000	2" x 2" Non-Sterile Gauze Sponges, 12-Ply, 200/slv, 40 slv/cs
56840000	3" x 3" Non-Sterile Gauze Sponges, 12-Ply, 200/slv, 20 slv/cs
56850000	4" x 4" Non-Sterile Gauze Sponges, 8-Ply, 200/slv, 20 slv/cs
56860000	4" x 4" Non-Sterile Gauze Sponges, 12-Ply, 200/slv, 10 slv/cs
56870000	4" x 4" Non-Sterile Gauze Sponges, 16-Ply, 200/slv, 10 slv/cs
56880000	2" x 2" Gauze Sponges, Sterile 2s, 8-Ply, 50/bx, 30 bx/cs
56890000	3" x 3" Gauze Sponges, Sterile 2s, 12-Ply, 40/bx, 30 bx/cs
56900000	4" x 4" Gauze Sponges, Sterile 2s, 8-Ply, 50/bx, 12 bx/cs
56910000	4" x 4" Gauze Sponges, Sterile 2s, 12-Ply, 25/bx, 24 bx/cs

J&J FIRST AID COVERS STERILE GAUZE PADS

008520	2" x 2" Sterile Gauze Pads, 10/bx, 24 bx/cs
008521	2" x 2" Sterile Gauze Pads, 25/bx, 24 bx/cs
008522	3" x 3" Sterile Gauze Pads, 10/bx, 24 bx/cs
008523	3" x 3" Sterile Gauze Pads, 25/bx, 24 bx/cs
008524	4" x 4" Sterile Gauze Pads, 10/bx, 24 bx/cs
008525	4" x 4" Sterile Gauze Pads, 25/bx, 24 bx/cs

J&J MIRASORB™ SPONGES

A cotton blend sponge with greater absorbency than woven gauze. Also more resilient than woven gauze, with less lint and fraying than woven gauze.

005531	4" x 4" Sponge, Sterile, 4-Ply, 2s, 50/bx, 18 bx/cs
--------	---

J&J NU GAUZE™ GENERAL-USE SPONGES

Unique rayon and polyester nonwoven provides improved strength, absorbency and scrubbing texture than gauze. Less linting.

005536	4" x 4" Sponge, Sterile, 4-Ply, 2s, 24/bx, 24 bx/cs
--------	---



MEDICOM GENERAL PURPOSE SPONGES

Increased patient comfort with soft viscose/poly blend non-woven material, with no chemicals added to avoid allergic reactions. A great economical alternative to cotton gauze sponges.

2100-CH	4" x 4", 4 Ply, Non-Woven, Non-Sterile, 200/sleeve, 10slv/cs
2101-CH	2" x 2", 4 Ply, Non-Woven, Non-Sterile, 200/sleeve, 25 slv/cs
2103-CH	3" x 3", 4 Ply, Non-Woven, Non-Sterile, 200/sleeve, 20slv/cs

MEDICOM SAFE-GAUZE® SPONGES

Saves time and cost by using fewer non-woven sponges per procedure compared to cotton filled sponges. Easy-to-use non-linting and highly absorbent premium weight sponges ideal for restorative scaling procedures.

4522-5	2" x 2", 4 Ply, Non-Woven, Non-Sterile, 200/slv, 25 slv/cs
4544	4" x 4", 4 Ply Non-Woven, Non-Sterile, 200/sleeve, 10slv/cs

TIDI STERILE GAUZE PADS

Sterile, 12-ply pads are ideal for wound cleaning, dressing changes, skin preps and many other uses. Deep infolded edges protect against loose threads. Meets USP VII absorbent gauze specifications. Individually packed in glassine envelopes.

908110	3" x 3" Gauze Pad, 100/bx, 24 bx/cs
908112	2" x 2" Gauze Pad, 100/bx, 36 bx/cs
908115	4" x 4" Gauze Pad, 100/bx, 12 bx/cs



TIDI COTTON-FILLED SPONGES

908223	3" x 3" Cotton-Filled Sponge, Non-Sterile, 100/bg, 4000/cs
908224	4" x 4" Cotton-Filled Sponge, Non-Sterile, 100/bg, 2000/cs
919000	2" x 2" Cotton-Filled Sponge, Non-Sterile, 200/bg, 5000/cs
919010	2" x 2" Cotton-Filled Sponge, Sterile, Bulk, 200/bg, 5000/cs

TIDI NON-WOVEN SPONGES

4-ply.

260002	2" x 2" Venture Nonwoven Open Weave Sponge, Non-Sterile, 200/slv, 40 slv/cs
260003	3" x 3" Venture Nonwoven Open Weave Sponge, Non-Sterile, 200/slv, 20 slv/cs
260004	4" x 4" Venture Nonwoven Open Weave Sponge, Non-Sterile, 200/slv, 10 slv/cs
912022	2" x 2" Nonwoven Closed Weave Sponge, Non-Sterile, 200/bg, 25 bg/cs
912033	3" x 3" Nonwoven Closed Weave Sponge, Non-Sterile, 200/bg, 20 bg/cs
912044	4" x 4" Nonwoven Closed Weave Sponge, Non-Sterile, 200/bg, 10 bg/cs


Dressings/Gauze



TIDI OMNISORB® II NON-WOVEN SPONGES

High quality sponges feature a unique 4 ply nonwoven closed weave design, rayon/polyester fabric for superior wound management. Low-linting characteristics minimize the incidence of dressing related complications and our quarter-fold design offers minimal cut edge exposure.

- 942022 Non-Sterile Nonwoven Sponge, 2" x 2", 200/bg, 25 bg/cs
 942033 Non-Sterile Nonwoven Sponge, 3" x 3", 200/bg, 20 bg/cs
 942044 Non-Sterile Nonwoven Sponge, 4" x 4", 200/bg, 10 bg/cs



TIDI STERILE & NON-STERILE GAUZE SPONGES

High quality gauze sponges meet USP Type VII absorbent gauze specifications.

- 908200 4" x 8" Gauze Sponge, Non-Sterile, 12-Ply, 100/bg, 2000/cs
 908206 3" x 3" Gauze Sponge, Non-Sterile, 8-Ply, 200/bag, 4000/cs
 908270 2" x 2" Gauze Sponge, Sterile 2s, 8-Ply, 50 pkg/tray, 3000/cs
 908271 3" x 3" Gauze Sponge, Sterile 2s, 12-Ply, 40 pkg/tray, 2400/cs
 908272 4" x 4" Gauze Sponge, Sterile 2s, 12-Ply, 25 pkg/tray, 1200/cs
 908276 4" x 4" Gauze Sponge, Sterile 2s, 8-Ply, 25 pkg/tray, 1200/cs
 908282 2" x 2" Gauze Sponge, Non-Sterile, 8-Ply, 200/bag, 5000/cs
 908291 4" x 4" Gauze Sponge, Non-Sterile, 16-Ply, 200/bag, 2000/cs
 908292 4" x 4" Gauze Sponge, Non-Sterile, 8-Ply, 200/bag, 4000/cs
 908293 3" x 3" Gauze Sponge, Non-Sterile, 12-Ply, 200/bag, 4000/cs
 908294 2" x 2" Gauze Sponge, Non-Sterile, 12-Ply, 200/bag, 8000/cs
 908295 4" x 4" Gauze Sponge, Non-Sterile, 12-Ply, 200/bag, 2000/cs



DYNAREX NON-WOVEN SPONGES

Sterile. 2/envelope.

- 3252 2" x 2" Non-Woven Sponge, Non-Sterile, 200/bg, 20 bg/cs
 3253 3" x 3" Non-Woven Sponge, Non-Sterile, 200/bg, 20 bg/cs
 3254 4" x 4" Non-Woven Sponge, Non-Sterile, 200/bg, 10 bg/cs
 3362 2" x 2" Non-Woven Sponge, 4-Ply, 50 pch/tray, 30 trays/cs
 3363 3" x 3" Non-Woven Sponge, 4-Ply, 50 pch/tray, 24 trays/cs
 3364 4" x 4" Non-Woven Sponge, 4-Ply, 25 pch/tray, 24 trays/cs



PRO ADVANTAGE® NON-ADHERENT STERILE PADS

Ideal for lightly draining wounds, these highly absorbent pads provide cushioning protection without sticking to the healing area. Sterile.

- P150155 Non-Adherent Sterile Pads, 2" x 3", 100/bx, 12 bx/cs
 P150160 Non-Adherent Sterile Pads, 3" x 4", 100/bx, 12 bx/cs



ALBA PERMEABLE NON-ADHERING DRESSINGS

Sterile. Made of a permeable knitted fabric that allows wound drainage to flow freely into the absorbent dressing. This helps keep the wetness from the wound, which minimizes maceration and helps promote healing. White petrolatum lines the dressing, providing a comfortable non-adherence to the wound - even during dressing changes. Ideal for draining surgical wound, open ulcers, abrasions, burns and other uses. The sterile clear film packaging offers a unique stay-flat feature when peeled open.

- 220 3" x 3" Non-Adhering Dressing, Single, 50/bx, 6 bx/cs
 221 3" x 8" Non-Adhering Dressing, Single, 24/bx, 6 bx/cs
 222 3" x 8" Non-Adhering Dressing, Three's, 36/bx, 6 bx/cs
 223 3" x 16" Non-Adhering Dressing, Single, 36/bx, 6 bx/cs
 225 1/2" x 4 yd Non-Adhering Dressing, 12 btl/cs



COVIDIEN/KENDALL TELFA OUCHLESS ADHESIVE DRESSINGS

Made of highly absorbent cotton fabric bonded on both sides with perforated non-adherent film. Won't disrupt healing tissue by sticking to wound. Adhesive strips keep dressing intact. Perforated carton for easy dispensing.

- 6017 2" x 3" Adhesive Dressing, Sterile 1s, 100/bx, 24 bx/cs
 7643 3" x 4" Adhesive Dressing, Sterile 1s, 100/bx, 12 bx/cs

**COVIDIEN/KENDALL CURITY NON-ADHERING DRESSINGS**

Oil emulsion blend impregnated into open mesh, knitted fabric dressing. Non-adhering to wound site. Allows free drainage of exudate away from wound. Ideal for donor sites, delicate suture lines, graft sites, draining wounds.

6111	Non-Adhering Dressing, 1/2" x 4", 1/bottle, 12/cs
6112	Non-Adhering Dressing, 3" x 3", 50/bx, 12 bx/cs
6113	Non-Adhering Dressing, 3" X 8", 24/bx, 6 bx/cs
6114	Non-Adhering Dressing, 3" x 16", 36/bx, 6 bx/cs
6115	Non-Adhering Dressing, 3" x 8", 3/pk, 36 pk/bx, 6 bx/cs
6116	Non-Adhering Dressing, 5" x 9", 12/bx, 6 bx/cs
6121	Non-Adhering Dressing, 9" x 18", 25/bx, 4 bx/cs

**COVIDIEN/KENDALL TELFA ISLAND DRESSING**

Convenient, all-in-one dressing. Soft, nonwoven backing conforms around wound and seals on all four sides. Provides protection from external environment. Non-adherent won't disrupt healing tissue. Ideal dressing for lightly draining wounds. Secure dressing for central and peripheral IV sights.

2562	Telfa Adhesive Island Dressing Plus, 4" x 6", 25/bx, 4 bx/cs
2563	Telfa Adhesive Island Dressing Plus, 6" x 7", 25/bx, 4 bx/cs
7539LF	Telfa M.I.S. DRS 2 x 3 3/4, 400/cs
7540	4" x 5" Island Dressing, Sterile 1s, 25/bx, 8 bx/cs
7541	Telfa Adhesive Island Dressing, 4" x 8", 100/cs
7542	Telfa Adhesive Island Dressing, 4" x 10", 100/cs
7544	Telfa Adhesive Island Dressing, 4" x 14", 50/cs
7550	Island Dressing, Sterile 1's in Peel-Back Package, 4" x 4", 25/bx, 8 bx/cs
7551	Telfa Adhesive Island 6" x 6", 100/cs

**COVIDIEN/KENDALL TELFA OUCHLESS NON-ADHERENT DRESSINGS**

Absorbent cotton pad. Won't disrupt healing tissue by sticking to wound. Ideal as primary dressing for lightly draining wounds. Bonded on both sides with perforated non-adherent film. Can be cut to fit any shape without separating. Easy to use and an ideal take home dressing.

1050	Telfa Dressing 3" x 4", 900/cs
1169	Telfa Dressing 3" x 6", 750/cs
1238	3" x 8" Non-Adherent Dressing, Sterile 1s, 50/bx, 12 bx/cs
12388	Telfa Strips, Bulk, 3" x 8", 400/bg, 10 bg/cs, Made to Order Only
1961	2" x 3" Non-Adherent Dressing, Sterile 1s, 100/bx, 24 bx/cs
2132	Ouchless Non-Adherent Dressing, Sterile 1's in Peel-Back Package, 3" x 4", 100/bx, 24 bx/cs
2891	Telfa Dressing 3" x 8", 1000/cs
3279	Ouchless Non-Adherent Dressing, Non-Sterile, Bulk, Mylar Film On One Side Only, 125/bx, 4 bx/cs



123



134



138

DUKAL NON-ADHERENT PADS

Designed to protect to wound site while promoting the healing process. The non-stick outer facing will allow the dressing to absorb, but not interfere with the granulating tissue by adhering to the wound. The easy removal of the pad will minimize trauma to the wound site. These pads can be cut any size wound without separating.

123	Pad Non-Adherent, 2" x 3", 100/bx, 12 bx/cs
134	Pad Non-Adherent, 3" x 4", 100/bx, 12 bx/cs
138	Pad Non-Adherent, 3" x 8", 1/pk, 75 pk/bx, 8 bx/cs

**J&J ADAPTIC™ NON-ADHERING DRESSING**

Non-adhering cellulose acetate impregnated with specially formulated petrolatum emulsion permits flow of wound exudate to the absorbent dressing. Prevents dressing adherences and helps prevent pooling of fluid.

2012	3" x 3" Non-Adhering Dressing, 50/bx, 12 bx/cs
2013	3" x 8" Non-Adhering Dressing, 3/pkg, 36 pkg/bx, 6 bx/cs
2014	3" x 16" Non-Adhering Dressing, 36/bx, 6 bx/cs
2015	3" x 8" Non-Adhering Dressing, 24/bx, 6 bx/cs
2018	3" x 60" Non-Adhering Dressing, 10/cs
2019	5" x 9" Non-Adhering Dressing, 12/bx, 6 bx/cs
005551	3" x 3" Non-Adhering Sterile Dressing, 12/bx, 24 bx/cs

J&J FIRST AID COVERS ADVANCED CARE

008526	3" x 2.25" Advanced Healing Adhesive Pads, 4/bx, 24 bx/cs
--------	---

J&J FIRST AID COVERS NON-STICK PADS

004720	1 1/2" x 2" Non-Stick Pads, 10/bx, 24 bx/cs
004724	2" x 3" Non-Stick Pads, 10/bx, 24 bx/cs
004978	First Aid Covers Non Stick Pads, Triple Layer, Larger, 4" x 3" 10 Count, 24/cs

J&J RELEASE™ NON-ADHERING DRESSING

Permits absorption of exudate into inner layer, does not trap moisture underneath dressing. Non-adherence permits easy removal. Conforms to body contours. Sterile.

005518	4" x 3" Non-Adhering Dressings, 12/bx, 24 bx/cs
--------	---

NUTRAMAX NON-ADHERENT PAD WITH ADHESIVE

Rayon-filled sponge with non-stick surface prevents adhesion to new tissue in the wound bed. Highly absorbent rayon will wick excess exudate from the wound.

7665033	2" x 3" Non-Adherent Pad, 100/bx, 12 bx/cs
7675033	3" x 4" Non-Adherent Pad, 100/bx, 12 bx/cs
7690000	2" x 3" Non-Adherent Pad, Bulk, 1800/cs

Dressings/Gauze

NUTRAMAX NON-ADHERENT STERILE PADS

Sterile pads will not stick to wounds. Supplied in individual peel-open packets in handy dispenser boxes. Latex free.

7565033	2" x 3" Non-Adherent Pad, 100/bx, 12 bx/cs
7575033	3" x 4" Non-Adherent Pad, 100/bx, 12 bx/cs
7590000	2" x 3" Non-Adherent Pad, Bulk, 3600/cs

COVIDIEN/KENDALL CURITY READY CUT-GAUZE

Non-sterile. Convenient, economical for wide variety of basic dressings. 20 x 10 PSI thread count mesh gauze. 100 square yards per package.

2798	Ready-Cut Gauze, 36" x 36", Non-Sterile, 100 sheets/bg, 10 bg/cs
3195	Ready-Cut Gauze, 18" x 18", Non-Sterile, 400 sheets/bg, 10 bg/cs
6122	Ready-Cut Gauze, 12" x 12", Non-Sterile, 900 sheets/bg, 10 bg/cs



COVIDIEN/KENDALL KERLIX GAUZE ROLLS

Made of prewashed, fluff-dried 100% woven gauze with unique crinkle-weave pattern for loft and bulk to cushion and protect wound areas. Six-ply construction with finished edges reduces loose ends and lint. As primary dressing, open-weave design provides fast wicking action, superior aeration and excellent absorbency. As secondary dressing, it provides bulk, cushioning and ultimate conformability. Ideal for bandaging heads, limbs and difficult-to-dress wounds; i.e., burns, plastic and orthopedic wounds.

1801	Kerlix Roll, Small, 2 1/4" x 3 yds, 6-Ply, Bulk, 96 rl/cs
1892	Kerlix Roll, Large, 4.5" x 4.1 yds, 6-Ply, Bulk, 100 rl/cs
3324	Gauze Rolls, Non-Sterile, 6 ply, Bulk & Packed 12's, 4.5" x 4.1 yds, 12/bx, 4 bx/cs
6715	4.5" x 4.1 yds Gauze Roll, Sterile, 100/cs
6716	Kerlix Roll, Large, 4.5" x 3.1 yds, 8-Ply, 100 rl/cs
6720	Gauze Rolls, Sterile, Soft Pouch, 6 ply, 2.25" x 3.0 yds, 12/bx, 8 bx/cs
6725	Kerlix Roll, Medium, 3.4" x 3.6 yds, 6-Ply, 96 rl/cs
6730	Kerlix Roll, Large, 4.5" x 4.1 yds, 6-Ply, 60 rl/cs
6735	Kerlix Roll, Medium, 3.4" x 3.6 yds, 6-Ply, Bulk, 96 rl/cs
6760	Kerlix Burn Pack, 4.5" x 4.1 yds, 80 rl/cs



COVIDIEN/KENDALL KERLIX AMD ROLLS

AMD Antimicrobial Dressings contain PHMB (Polyhexamethylene Biguanide). Reduce bacterial penetration through the dressing. Resists bacterial colonization within the dressing. Broad-spectrum effectiveness provides protection against gram negative, gram positive and fungi/yeast microorganisms including MRSA and VRE. Limits cross-contamination from patient to the environment. Sterile.

3331	4.5" x 4.1 yds, Sterile in Rigid Plastic Tray, 6 Ply, 1 roll/tray, 60 trays/cs
3332	4.5" x 4.1 yds, Sterile in Soft Pouch, 6-Ply, 100/cs

AMD-RITMED KRIMPED GAUZE BANDAGE ROLLS

A454	4.5" x 4.1 yds Gauze Bandage Roll, Non-Sterile, 6-Ply, 100/cs
A455	4.5" x 4.1 yd Gauze Bandage Roll, Sterile 1s, Soft Pouch, 6-Ply, 100/cs



AMD-RITMED VITAL-ROLL GAUZE BANDAGES

A6902	2" x 131" Gauze Bandage, Non-Sterile, 12/pk, 8 pk/cs
A6903	3" x 131" Gauze Bandage, Non-Sterile, 12/pk, 8 pk/cs
A6904	4" x 131" Gauze Bandage, Non-Sterile, 12/pk, 8 pk/cs
A6906	6" x 131" Gauze Bandage, Non-Sterile, 6/pk, 8 pk/cs
A6922	2" x 131" Gauze Bandage, Sterile, 12/pk, 8 pk/cs
A6923	3" x 131" Gauze Bandage, Sterile, 12/pk, 8 pk/cs
A6924	4" x 131" Gauze Bandage, Sterile, 12/pk, 8 pk/cs
A6926	6" x 131" Gauze Bandage, Sterile, 12/pk, 4 pk/cs



COVIDIEN/KENDALL CONFORM STRETCH BANDAGES

One-ply cotton/polyester blend crocheted bandage. Provides softness, conformability, low lint, high absorbency. Ideal for securing dressings, lvs and splints or for providing mild compression or support. Holds securely to any body contour. Allows for movement and some soft tissue swelling. Stays in place with minimal taping.

2230	Stretch Bandages, Sterile, Soft Pouch, Stretched, 1" x 75", 12/bx, 8 bx/cs
2231	Stretch Bandages, Sterile, Soft Pouch, Stretched, 2" x 75", 12/bx, 8 bx/cs
2232	Stretch Bandages, Sterile, 3" x 75", 12/bx, 8 bx/cs
2236	Stretch Bandages, Sterile, 4" x 75", 12/bx, 8 bx/cs
2238	Stretch Bandages, Sterile, Soft Pouch, Stretched, 6" x 82", 12/bx, 4 bx/cs
2239	Stretch Bandages, Non-Sterile, Bulk, Stretched, 1" x 75", 24/bx, 4 bx/cs
2242	Stretch Bandages, Non-Sterile, 2" x 75", 12/bg, 8 bg/cs
2244	Stretch Bandages, Non-Sterile, Bulk, Stretched, 3" x 75", 12/bx, 8 bx/cs
2247	Stretch Bandages, Non-Sterile, Bulk, Stretched, 4" x 75", 12/bx, 8 bx/cs
2249	Stretch Bandages, Non-Sterile, Bulk, Stretched, 6" x 75", 6/bx, 8 bx/cs

DERMA SCIENCES DUTEX® CONFORMING BANDAGES

This 100% cotton 2-ply woven bandage is ideal as a secondary dressing for wrapping legs and arms and to hold dressings in place. The dressing provides light stretch which enables it to conform to difficult body contours, while also providing light compression. Latex free.

76781	Conforming Bandage, 2" x 4.1 yards, 12/bg, 8bg/cs
76782	Conforming Bandage, 3" x 4.1 yards, 12/bg, 8bg/cs

DUKAL FLUFF BANDAGE ROLL

Made of 100% washed fluff cotton gauze, the roll offers superior performance. It is a highly absorbent dressing, with excellent bulk and loft that will conform as needed and hold dressings in place.

545	Bandage Roll, 4 1/2" x 147", 6-Ply, Fluff Non-Sterile, 100 rl/cs
640	Bandage Roll, 3.5" x 3.6 yds, Fluff, Sterile, 6-Ply, 1 rl/bg, 96 bg/cs
645	Bandage Roll, 4 1/2" x 4.1 yds, 6-Ply, Sterile, Fluff, 1 rl/bg, 100 rl/cs
650	Bandage Roll, 4.5" x 3.1 yds, Fluff, Sterile, 8-Ply, 1 rl/bg, 100 bg/cs

DUKAL BASIC CONFORMING STRETCH GAUZE

The rayon/poly blended product is specifically designed not to slip while in use. There are three packaging options available to meet all needs and applications. "Clean" wrap is non-sterile product individually packaged in a polybag.

- 8513 Conforming Stretch Gauze, 2" Sterile, 12 rl/bg, 8 bg/cs
- 8514 Conforming Stretch Gauze, 3" Sterile, 12 rl/bg, 8 bg/cs
- 8515 Conforming Stretch Gauze, 4" Sterile, 12 rl/bg, 8 bg/cs
- 8516 Conforming Stretch Gauze, 6" Sterile, 6 rl/bg, 8 bg/cs
- 8517 Conforming Stretch Gauze, 2" Non-Sterile, 12 rl/bg, 8 bg/cs
- 8518 Conforming Stretch Gauze, 3" Non-Sterile, 12 rl/bg, 8 bg/cs
- 8519 Conforming Stretch Gauze, 4" Non-Sterile, 12 rl/bg, 8 bg/cs
- 8520 Conforming Stretch Gauze, 6" Non-Sterile, 6 rl/bg, 8 bg/cs



602



702

DUKAL CONFORMING STRETCH GAUZE

This knitted stretch gauze bandage provides slight compression to the site while securely holding the dressing in place. The rayon/poly blended product is specifically designed not to slip while in use. There are three packaging options available to meet all needs and applications. "Clean" wrap is a non-sterile product individually packaged in a polybag.

- 602 Gauze, 2" Stretch, Non-Sterile, 12/bg, 8 bg/cs
- 602PB-96 Gauze, 2" Conforming, Non-Sterile Clean, 12 rl/bg, 8 bg/cs
- 603 Gauze, 3" Stretch, Non-Sterile, 12/bg, 8 bg/cs
- 603PB-96 Gauze, 3" Conforming, Non-Sterile Clean, 12 rl/bg, 8 bg/cs
- 604 Gauze, 4" Stretch, Non-Sterile, 12/bg, 8 bg/cs
- 604PB-96 Gauze, 4" Conforming, Non-Sterile Clean, 12 rl/bg, 8 bg/cs
- 606 Gauze, 6" Stretch, Non-Sterile, 6/bg, 8 bg/cs
- 702 Gauze, 2" Stretch, Sterile, 12 rl/bg, 8 bg/cs
- 703 Gauze, 3" Stretch, Sterile, 12 rl/bg, 8 bg/cs
- 704 Gauze, 4" Stretch, Sterile, 12 rl/bg, 8 bg/cs
- 706 Gauze, 6" Stretch, Sterile, 6 rl/bg, 8 bg/cs



402



403



404



406

DUKAL ROLLED GAUZE

Made of 100% woven gauze the Rolled Gauze can be used to effectively hold dressings securely in place. This 2 ply product is available in sterile and non sterile packaging.

- 402 Rolled Gauze, 2", Non-Sterile, 2-ply, 12 rl/bg, 8 bg/cs
- 403 Rolled Gauze, 3", Non-Sterile, 2-Ply, 12 rl/bg, 8 bg/cs
- 404 Rolled Gauze, 4", Non-Sterile, 2-Ply, 12 rl/bg, 8 bg/cs
- 406 Rolled Gauze, 6", Non-Sterile, 6/bg, 8 bg/cs
- 452 Rolled Gauze, 2", Sterile, 2-ply, 1/pouch, 12 pouch/bg, 8 bg/cs
- 453 Rolled Gauze, 3" Sterile, 2-Ply, 1/pouch, 12/bg, 8 bg/cs
- 454 Rolled Gauze, 4" Sterile, 2-Ply, 1/pouch, 12/bg, 8 bg/cs
- 456 Rolled Gauze, 6" Sterile, 2-Ply, 1/pouch, 6/bg, 8 bg/cs



HARTMANN-CONCO STERILUX® BULKY GAUZE BANDAGE

Designed to pad and protect wound areas. Made of 100% cotton. Provides absorbency, comfort and maximum wound protection. Its six-ply construction and finished edges reduce loose ends and lint. Ideal for difficult-to-dress wounds. Latex free.

- 83400000 4.5 x 4.1 yds Gauze Bandage, Non-Sterile, 100/cs
- 83500000 4.5 x 4.1 yds Gauze Bandage, Sterile, 100/cs

J&J FIRST AID COVERS STERILE ROLLED GAUZE

- 008807 2" x 2.5 yds Sterile Rolled Gauze, 48/cs
- 008808 3" x 2.5 yds Sterile Rolled Gauze, 48/cs
- 008809 4" x 2.5 yds Sterile Rolled Gauze, 6/bx, 2 bx/cs

J&J KLING® CONFORMING GAUZE BANDAGES

- 005521 Sterile Conforming Gauze Bandage, 3" x 2.1 yds, 3/bx, 12 bx/cs
- 005522 Sterile Conforming Gauze Bandage, 4" x 2.1 yds, 5/bx, 12 bx/cs

NUTRAMAX CONFORM BANDAGES

Latex free.

- 6105033 1" x 4 yds Conform Bandages, Sterile, 144/cs
- 6115033 2" x 4 yds Conform Bandages, Sterile, 144/cs
- 6120033 3" x 4 yds Conform Bandages, Sterile, 72/cs
- 6125033 4" x 4 yds Conform Bandages, Sterile, 72/cs



DYNAREX STRETCH GAUZE BANDAGES

Conforming.


- 3102 2" Bandage, Non-Sterile, 12/bg, 8 bg/cs
- 3103 3" Bandage, Non-Sterile, 12/bg, 8 bg/cs
- 3104 4" Bandage, Non-Sterile, 12/bg, 8 bg/cs
- 3106 6" Bandage, Non-Sterile, 6/bg, 8 bg/cs
- 3112 2" Bandage, Sterile, 12/bg, 8 bg/cs
- 3113 3" Bandage, Sterile, 12/bg, 8 bg/cs
- 3114 4" Bandage, Sterile, 12/bg, 8 bg/cs
- 3116 6" Bandage, Sterile, 6/bg, 8 bg/cs



DYNAREX KRINKLE GAUZE ROLLS

- 3161 Sterile, 4 1/2" x 4.1 yds Gauze Roll, 6 Ply, Individual Pouch, 100/cs
- 3162 Non-Sterile, 4 1/2" x 4.1 yds Gauze Roll, 6 Ply, Bulk, 100/cs


Dressings/Gauze



PRO ADVANTAGE® NON-ADHESIVE FOAM DRESSINGS

This dressing is designed with cost-effectiveness in mind, yet will provide a clinically effective dressing for management of exuding wounds. It is made of a non-adherent, highly absorbent, lint free polyurethane foam. Featuring a thick and very soft construction, it is extremely conformable and comfortable, especially for painful wounds where some protection and cushioning are desirable. Conformability makes this dressing ideal for difficult to dress areas of the body. Sterile, Polyurethane.


- P157502 Non-Adhesive Foam Dressing, 2" x 2.125", 10/bx
P157504 Non-Adhesive Foam Dressing, 4" x 4.25", 10/bx



PRO ADVANTAGE® ADHESIVE FOAM ISLAND DRESSINGS

This high-density foam dressing is covered by a layer of comfortable, self-adhesive, water-resistant polyurethane backing helps maintain an optimally moist wound environment, helps prevent strike-through, and aids in the prevention of bacterial contamination of the wound. The moisture vapor transmission rate adjusts depending on the fluid level in the wound. The foam pad has micropores for low adherence to wound surface, yet has a higher fluid-handling capacity than other brands. Foam island gently expands as it absorbs, keeping the dressing in contact with the wound surface for improved fluid handling.

- P157524 Adhesive Foam Island Dressing, 4" x 4", w/ 2" x 2" Pad, 10/bx
P157526 Adhesive Foam Island Dressing, 6" x 6" w/ 4.25" x 4.25" Pad, 10/bx



PRO ADVANTAGE NON-ADHESIVE FOAM DRESSINGS WITH BARRIER

The outer, pink layer used on non-adherent dressing is a breathable film layer. This allows the formation and maintenance of a moist wound healing environment at the wound surface, preventing eschar formation. The pink film is also waterproof and assists in preventing bacterial strike-through. Non-adherent dressings are non-adhesive and are therefore recommended for wounds surrounded by sensitive skin, such as venous ulcers. The lack of adhesive allows the use of the most suitable retention method, for example, tape onto intact skin, appropriate bandaging or a film dressing such as transparent film dressing. Non-adherent dressings are comfortable and can be cut to suit various wound shapes and locations. Non-adherent dressings absorb exudate, they may swell slightly, but do not leave debris in the wound. Dressing changes are therefore easier for the nurse and less traumatic for the patient; wound irrigation can be kept to a minimum. The progression of exudate through the dressing is visible through the pink top film. Sterile, Polyurethane.

- P157512 2" x 2.125", 10/bx
P157514 4" x 4.25", 10/bx



3M™ RESTON™ SELF-ADHERING FOAM PRODUCTS

Distributes body weight, helps prevent skin breakdown. Easily cut and shaped. Protects against abrasions, shearing and pressure. Can be stacked for greater thickness. Adheres efficiently to cushion skin, casts, appliances and other devices.

- 1560M Medium Support Foam Pads, 7 7/8" x 11 3/4" (7/16" thick), 10/pk, 5 pk/cs
1561H High Support Foam Pads, 7 7/8" x 11 3/4" (15/16" thick), 5/pk, 5 pk/cs
1563L Light Support Foam Pads, 4" x 196" (3/16" thick), 1/r, 5 r/cs

3M™ TEGADERM™ FOAM ADHESIVE DRESSING

This highly absorbent, breathable wound dressing is constructed from a conformable polyurethane foam pad and an additional absorbent nonwoven layer with a border of transparent adhesive film impermeable to liquids, bacteria and viruses. Provides total fluid management by a combination of fast wicking, high absorbency and breathability and uses an innovative moist skin adhesive technology. The unique spoke-shaped delivery system allows fast, easy application for wounds over body contours.

- 90610 2" x 2" Foam Adhesive Dressing, Square, 10/bx, 4 bx/cs
90611 2 1/2" x 3" Foam Adhesive Dressing, Oval, 10/bx, 4 bx/cs
90612 4" x 4" Foam Adhesive Dressing, Square, 10/bx, 4 bx/cs
90613 4" x 4 1/2" Foam Adhesive Dressing, Oval, 5/bx, 6 bx/cs
90616 5 1/2" x 6 3/4" Foam Adhesive Dressing, Oval, 5/bx, 3 bx/cs
90619 3" x 3" Foam Adhesive Dressing, Heel Design, 5/bx, 4 bx/cs

3M™ TEGADERM™ FOAM DRESSING - NONADHESIVE

The wound dressing is a fast wicking, non-swelling polyurethane foam with a highly breathable film backing which prevents wound exudate strike-through and is a barrier to outside contamination while the dressing remains intact without leakage.

- 90600 2" x 2" Nonadhesive Foam Dressing, 10/bx, 4 bx/cs
90601 4" x 4" Nonadhesive Foam Dressing, 10/bx, 4 bx/cs
90602 4" x 8" Nonadhesive Foam Dressing, 5/bx, 6 bx/cs
90603 8" x 8" Nonadhesive Foam Dressing, 5/bx, 6 bx/cs
90604 3 1/2" x 3 1/2" Nonadhesive Foam Dressing, 10/bx, 4 bx/cs
90605 4" x 24" Nonadhesive Foam Dressing, 6 r/cs

ETHICON SURGIFOAM

Absorbable gelatin sponge. USP sterile. This item is available on a drop ship to end user only basis. Contact NDC Customer Service for more info.

- 1972 Size 12-7 Surgifoam Absorbable Gelatin Sponge, 12/cs
1973 Size 50 Surgifoam Absorbable Gelatin Sponge, 4/cs
1974 Size 100 Surgifoam Absorbable Gelatin Sponge, 6/cs
1975 Size 100C Surgifoam Absorbable Gelatin Sponge, 6/cs
1977 Hemorrhoidectomy Sponge, Surgifoam Absorbable Gelatin Sponge, 5/cs
1978 1 Gram Surgifoam Absorbable Gelatin Powder, 6/cs

J&J BIOPATCH™ ANTIMICROBIAL DRESSING

Biopatch™ Antimicrobial Dressing is a hydrophilic polyurethane absorptive foam with chlorhexidine gluconate (CHG). The foam material absorbs up to eight times its own weight in fluid, while the CHG incorporated into the dressing inhibits bacterial growth under the dressing. Improved upon the ease of use by adding a visual indicator (blue side up). The first and only percutaneous device dressing proved to reduce CRBSI (Catheter Related Blood Stream Infections).

- 3150 1" Disk w/4mm center hole, Sterile, 10/bx, 4 bx/cs
3151 .75" Disk w/1.5mm center hole, Sterile, 10/bx, 4 bx/cs
3152 1" Disk w/7mm center hole, Sterile, 10/bx, 4 bx/cs

SMITH & NEPHEW ALLEVYN™ PLUS CAVITY DRESSING

Conforming dressing/filler for heavily exuding and deep cavity wounds.

66047571	2" x 2 3/8" Cavity Dressing, 10/bx, 10 bx/cs (475710000004)
66047573	4" x 4" Cavity Dressing, 5/bx, 20 bx/cs (475630000004)
66047574	6" x 8" Cavity Dressing, 3/bx, 16 bx/cs (475740000004)

SMITH & NEPHEW ALLEVYN™ PLUS THIN DRESSING

Flexible dressing for wounds with minimal exudate.


66047576	2" x 2" Thin Dressing, 10/bx, 10 bx/cs
66047578	4" x 4" Thin Dressing, 5/bx, 20 bx/cs (475780000004)
66047579	6" x 8" Thin Dressing, 3/bx, 16 bx/cs (475790000004)



SMITH & NEPHEW ALLEVYN™ ADHESIVE DRESSINGS

Sticks and remains in place, even on difficult-to-dress areas of the body. It may absorb up to four times more exudate than hydrocolloids and is ideal for moderately to heavily exuding wounds. The unique porous wound contact layer allows absorption and gaseous exchange. Will not stick to the wound surface or irritate periwound skin. Waterproof and aids in the prevention of bacterial contamination.

66000805	5" x 5", 10/bx, 5 bx/cs
66000806	7" x 7", 10/bx, 5 bx/cs
66020043	3" x 3" Adhesive Dressing, 10/pkg, 4 pkg/cs
66020044	5" x 5" Adhesive Dressing, 10/pkg, 4 pkg/cs
66020045	7" x 7" Adhesive Dressing, 10/pkg, 4 pkg/cs



SMITH & NEPHEW ALLEVYN™ CAVITY WOUND DRESSINGS

Comprised of highly absorbent hydrocellular foam chips covered with a non-adherent wound contact layer. Designed to be pressed loosely into the wound cavity. Exudate is absorbed, yet no fibers or residue remain in the wound. A moist wound environment is maintained and the dressing is very comfortable in the wound.

66027326	2" Circular (3g) Wound Dressing, 10/pkg, 6 pkg/cs
66027327	4" Circular (17g) Wound Dressing, 5/pkg, 4 pkg/cs
66027328	3 1/2" x 1" Tubular (3g) Wound Dressing, 10/pkg, 6 pkg/cs
66027329	4 3/4" x 1 1/2" Tubular (7g) Wound Dressing, 5/pkg, 4 pkg/cs



SMITH & NEPHEW ALLEVYN™ HEEL DRESSING

Designed to fit snugly around the back of the foot ensuring a better anatomical fit. Can be held in place with a light conforming bandage.

66007630	Heel Dressing, 5/pkg, 6 pkg/cs
----------	--------------------------------



SMITH & NEPHEW ALLEVYN™ HYDROCELLULAR FOAM DRESSINGS

A trilaminate structure comprised of a non-adherent wound contact layer, a highly absorbent central layer and an outer layer that is waterproof and aids in the prevention of bacterial contamination. This structure produces a moist wound environment for those wounds that have traditionally proved difficult to treat because the exudate is hard to manage. Very absorbent and can usually be left in place for up to three to four days.

66020093	6" x 6" Foam Dressing, 10/pkg, 3 pkg/cs
66027643	2" x 2" Foam Dressing, 10/pkg, 6 pkg/cs
66927637	4" x 4" Foam Dressing, 10/pkg, 7 pkg/cs
66927638	8" x 8" Foam Dressing, 10/pkg, 2 pkg/cs



SMITH & NEPHEW ALLEVYN™ HYDROCELLULAR FOAM TRACHEOSTOMY DRESSING

Absorbs large amounts of secretions, which may help prevent maceration and maintain the condition of the skin in the peristomal area. High absorbency requires fewer dressing changes.

66027640	3 1/2" x 3 1/2" Tracheostomy Dressing, 10/pkg, 8 pkg/cs
----------	---



HARTMANN-CONCO PERMAFOAM® NON-ADHESIVE & ADHESIVE FOAM DRESSING

Absorbs wound exudate and wicks fluid away from wound bed. High MVTR. Unsurpassed fluid handling. Conformable, comfortable, breathable and waterproof outer membrane. Conforms to hard-to-dress areas and it's easy to cut and shape for special applications. Latex free.

409401	4" x 4" Foam Membrane Dressing, Non-Adhesive, 10/bx
409405	6" x 6" Foam Membrane Dressing, Non-Adhesive, 5/bx
409408	4" x 4" Foam Membrane Dressing, Comfort Adhesive, 10/bx
409412	6" x 6" Foam Membrane Dressing, Comfort Adhesive, 5/bx
409413	8" x 8" for really large pressure ulcers (generally on the hip), Latex Free, 3/bx
409422	7" x 7" Sacral, Latex Free, 3/bx
409424	6.5" x 7", Concave, Foam Membrane Dressing, Comfort Adhesive, Latex Free, 3/bx
409425	4" x 4", Cavity, Foam Membrane Dressing, Non-Adhesive, Latex Free, 3/bx
409426	3.2" x 3.2", Tracheotomy, Latex Free, 10/bx

COVIDIEN/KENDALL COPA™ ISLAND ULTRA-SOFT FOAM WOUND DRESSING

Ultra-soft hydrophilic foam protects and cushions wounds. Small pores decrease chance of tissue migration into foam dressing. Contains a polyurethane backing designed to help prevent leakage and dressing strike-through. Provides an ideal healing environment for a wide variety of wounds.

55544B	Foam Island Dressing, 4" x 4", 2" x 2" Pad, 10/bx, 5 bx/cs
55566B	Foam Island Dressing, 6" x 6", 4" x 4" Pad, 10/bx, 5 bx/cs
55588B	Foam Island Dressing, 8" x 8", 6" x 6" Pad, 10/bx, 5 bx/cs

COVIDIEN/KENDALL COPA™ PLUS ULTRA-SOFT FOAM WOUND DRESSING

Ultra-soft hydrophilic foam protects and cushions wounds. Small pores decrease chance of tissue migration into foam dressing. Contains a polyurethane backing designed to help prevent leakage and dressing strike-through. Provides an ideal healing environment for a wide variety of wounds.

55522P	Foam Wound Dressing, 2" x 2", 25/bx, 4 bx/cs
55533P	Foam Wound Dressing, 3" x 3", 10/bx, 5 bx/cs
55535P	Foam Wound Dressing, 3 1/2" x 3", Fenestrated, 10/bx, 5 bx/cs
55544P	Foam Wound Dressing, 4" x 4", 10/bx, 5 bx/cs
55548P	Foam Wound Dressing, 4" x 8", 10/bx, 5 bx/cs
55555P	Foam Wound Dressing, 5" x 5", 10/bx, 5 bx/cs
55566P	Foam Wound Dressing, 6" x 6", 10/bx, 5 bx/cs
55588P	Foam Wound Dressing, 8" x 8", 10/bx, 5 bx/cs

COVIDIEN/KENDALL COPA™ ULTRA-SOFT FOAM WOUND DRESSING

Ultra-soft hydrophilic foam protects and cushions wounds. Small pores decrease chance of tissue migration into foam dressing. Provides an ideal healing environment for a wide variety of wounds.

55522	Foam Wound Dressing, 2" x 2", 25/bx, 4 bx/cs
55533	Foam Wound Dressing, 3" x 3", 10/bx, 5 bx/cs
55535	Foam Wound Dressing, 3 1/2" x 3", Fenestrated, 10/bx, 5 bx/cs
55544	Foam Wound Dressing, 4" x 4", 10/bx, 5 bx/cs
55548	Foam Wound Dressing, 4" x 8", 10/bx, 5 bx/cs
55555	Foam Wound Dressing, 5" x 5", 10/bx, 5 bx/cs
55566	Foam Wound Dressing, 6" x 6", 10/bx, 5 bx/cs
55588	Foam Wound Dressing, 8" x 8", 10/bx, 5 bx/cs



MOLNLYCKE TENDRA® MEPILEX® BORDER LITE

Self-adherent soft silicone thin foam dressing. Atraumatic to the wound and the surrounding skin on removal. Minimizes pain to the patient upon dressing removal. Minimizes the risk of maceration. Conforms well to the body contour; low profile. Thin and comfortable during wear. Helps maintain a moist wound environment. Can be cut to size.

281000	1.6" x 2" Self-Adherent Absorbent Foam Dressing, 10/bx, 7 bx/cs
281100	2" x 5" Self-Adherent Absorbent Foam Dressing, 65/cs
281200	3" x 3" Self-Adherent Absorbent Foam Dressing, 10/bx, 7 bx/cs
281300	4" x 4" Self-Adherent Absorbent Foam Dressing, 5/bx, 10 bx/cs
281500	6" x 6" Self-Adherent Absorbent Foam Dressing, 10/bx, 5 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPILEX®

Non-adherent to moist wound beds yet adheres ultra-gently to dry skin. Minimizes trauma and pain during dressing changes. Can remain in place for up to seven days. Helps maintain a moist wound environment. Low risk of peri-wound skin maceration and erosion. Conforms well to body contours. Exudate is visible through the surface of the dressing. Can be used under compression bandages. Can be lifted and adjusted without losing its adherent properties. Can be cut to size.

294190	4" x 4" Self-Adherent Absorbent Foam Dressing, 5/bx, 12 bx/cs
294290	4" x 8" Self-Adherent Absorbent Foam Dressing, 5/bx
294390	6" x 6" Self-Adherent Absorbent Foam Dressing, 5/bx, 6 bx/cs
294490	8" x 8" Self-Adherent Absorbent Foam Dressing, 5/bx, 6 bx/cs
294592	8" x 20" Silicone Soft Transfer Foam Dressing, 2/bx, 6 bx/cs
294890	6" x 8" Silicone Soft Transfer Foam Dressing, 5/bx, 7 bx/cs
295200	3" x 3" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 14bx/cs
295300	4" x 4" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 10 bx/cs
295400	6" x 6" Self-Adherent Absorbent Foam Dressing w/border, 5/bx
295600	6" x 8" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 9 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPILEX® LITE

Absorbent soft silicone thin foam dressing. Atraumatic to the wound and surrounding skin on removal. Minimizes pain to the patient upon dressing removal. Minimizes the risk of maceration. Conforms well to body contour; low profile. Thin and comfortable during wear. Helps maintain a moist wound environment. Can be cut to size.

284090	2.4" x 3.4" Soft Silicone Thin Foam Dressing, 5/bx, 14 bx/cs
284190	4" x 4" Soft Silicone Thin Foam Dressing, 5/bx, 10 bx/cs
284390	6" x 6" Soft Silicone Thin Foam Dressing, 5/bx, 10 bx/cs
284590	8" x 20" Soft Silicone Thin Foam Dressing, 2/bx, 6 bx/cs

**SMITH & NEPHEW ALLEVYN™ ADHESIVE SACRUM**

66000451 9" x 9", 10/bx, 8 bx/cs
 66000700 6" x 5 3/4" Adhesive Sacrum, 10/bx, 8 bx/cs

SMITH & NEPHEW ALLEVYN PLUS FOAM DRESSING

Absorbent dressing for compromised skin with moderate exudate.

66000856 Dressing, Adhesive Sacrum, 6.75" x 6.75" (17 x 17cm), 10/bx
 66047581 2" x 2 3/8" Foam Dressing, 10/bx, 10 bx/cs (475810000004)
 66047583 4" x 4" Foam Dressing, 5/bx, 20 bx/cs (475830000004)
 66047584 6" x 8" Foam Dressing, 3/bx, 16 bx/cs (475840000004)

SMITH & NEPHEW ALLEVYN™ GENTLE BORDER GEL DRESSING

Suitable for use on fragile skin. Advanced triple-layered construction. Absorbent hydrocellular pad. Perforated silicone gel adhesive wound contact layer. Adhesive border for gentle but effective fixation. Highly permeable waterproof outer film.

66800276 Dressing, 3" x 3", Gel Adhesive Hydrocellular Foam w/border, 10/bx, 4 bx/cs
 66800279 Dressing, 5" x 5", Gel Adhesive Hydrocellular Foam w/border, 10/bx, 4 bx/cs
 66800280 Dressing, 7" x 7", Gel Adhesive Hydrocellular Foam w/border, 10/bx, 4 bx/cs

SMITH & NEPHEW ALLEVYN™ GENTLE GEL DRESSING


Suitable for use on fragile skin. Advanced triple-layered construction. Absorbent hydrocellular pad. Perforated soft gel adhesive wound contact layer. Highly permeable waterproof outer film.

66800255 Dressing, 2" x 2", Gel Adhesive Hydrocellular Foam, 10/bx, 6 bx/cs
 66800256 Dressing, 4" x 4", Gel Adhesive Hydrocellular Foam, 10/bx, 7 bx/cs
 66800257 Dressing, 6" x 6", Gel Adhesive Hydrocellular Foam, 10/bx, 3 bx/cs
 66800258 Dressing, 8" x 8", Gel Adhesive Hydrocellular Foam, 10/bx, 2 bx/cs

SMITH & NEPHEW ALLEVYN™ AG ADHESIVE DRESSINGS

Provides the advantages of Allevyn™ hydrocellular dressings in an adhesive format for simple and effective fixation.

66020970 Dressing, Hydrocellular, Adhesive, 3" x 3", 10/bx, 4 bx/cs
 66020973 Dressing, Hydrocellular, Adhesive, 5" x 5", 10/bx, 4 bx/cs
 66800115 Dressing, Hydrocellular, Adhesive, Heel, 5/bx

**PRO ADVANTAGE® HYDROCOLLOID DRESSINGS**

This integrated hydrocolloid formula provides world-class wound healing effectiveness. We offer an ultra-thin hydrocolloid at 15mm. Our standard hydrocolloid is capable of besting the competition wear time while reducing profile and weight by 25%. The outer polyurethane layer provides a barrier to water and bacteria, while maintaining a moist environment. The semi-occlusive properties allow for oxygen and moisture vapor transmission creating the optimal conditions for wound healing. Sterile.

P157602 Hydrocolloid Dressings, Thin, 2" x 2.125", 20/bx
 P157604 Hydrocolloid Dressings, Thin, 4" x 4.25", 10/bx
 P157614 Hydrocolloid Dressings, Regular Thickness, 4" x 4.25", 5/bx

**3M™ TEGADERM™ HYDROCOLLOID DRESSING**

Sterile wound dressings which consist of a hypoallergenic, hydrocolloid adhesive with an outer clear adhesive cover film impermeable to liquids, bacteria and viruses. The inner layer of hydrocolloid adhesive rapidly absorbs exudate faster than the leading competitive hydrocolloid in the first 48 hours. The breathable outer film layer provides a consistently high rate of moisture vapor transmission. Together, these features ensure an optimal moist wound environment, minimize the chance for damage to healthy periwound skin and provide extended wear for up to seven days.

90001 Hydrocolloid Dressing, 2 3/4" x 3 1/2" Oval, 5/bx, 20 bx/cs
 90002 Hydrocolloid Dressing, 4" x 4" Square, 5/bx, 20 bx/cs
 90003 Hydrocolloid Dressing, 4" x 4 3/4" Oval, 5/bx, 12 bx/cs
 90004 Hydrocolloid Dressing, 5 1/2" x 6 3/4" Oval, 3/bx, 20 bx/cs
 90005 Hydrocolloid Dressing, 6" x 6" Square, 3/bx, 20 bx/cs
 90007 Hydrocolloid Dressing, 4.9" x 5.5" Sacral Design, 6/bx, 4 bx/cs

3M™ TEGADERM™ THIN HYDROCOLLOID DRESSING

Sterile wound dressings which consist of a hydrocolloid adhesive with an outer clear adhesive cover film impermeable to liquids, bacteria and viruses. Available in oval and square shapes.

90021 Hydrocolloid Dressing, 2 3/4" x 3 1/2" Dressing, 4" x 4 3/4" Overall, Oval, 10/bx, 10 bx/cs
 90022 Hydrocolloid Dressing, 4" x 4" Square, 5/bx, 20 bx/cs
 90023 Hydrocolloid Dressing, 4" x 4 3/4" Dressing, 5" x 6" Overall, Oval, 10/bx, 6 bx/cs
 90024 Hydrocolloid Dressing, 5 1/2" x 6 3/4" Dressing, 6 3/4" x 8" Overall, Oval, 6/bx, 10 bx/cs

**HARTMANN-CONCO FLEXICOL® HYDROCOLLOID DRESSING**

Folds easily to conform to difficult-to-dress areas. Peel-off tabs and a thin beveled edge make application quick and prevent dressing roll-up. Handles wound drainage with a combination of absorption and evaporation - making it ideal for treating light to moderately draining chronic and acute wounds. Stays secure days after it's first applied so re-application or securing with tape is eliminated. Thin film exterior and moist core allow for higher and quicker fluid uptake - leaving virtually no residue in the wound. Latex free.

48600000	2" x 2" Hydrocolloid Dressing, 20/bx
48610000	4" x 4" Hydrocolloid Dressing, 10/bx
48620000	6" x 7" Sacral Hydrocolloid Dressing, 5/bx
48630000	3" x 5" Heel Hydrocolloid Dressing, 10/bx
48640000	4" x 4" Thin Hydrocolloid Dressing, 10/bx
48650000	6" x 6" thin, Used to protect the hip bone of bed-ridden patients, Latex Free, 5/bx
48660000	6" x 6" standard thickness, Used to protect the hip bone of bed-ridden patients if pressure ulcers form, Latex Free, 5/bx

J&J NU-DERM™ HYDROCOLLOID WOUND DRESSING

Sterile dressings designed to promote a moist wound environment encouraging natural debridement and enabling granulation under optimal conditions. The dressing material interacts with wound exudate to form a soft gel. The matrix formulation of the hydrocolloid allows it to be removed with the dressing, with little or no damage to newly formed tissue. Features include low-friction backing, patented tapered edges, low profile tray design, skin-friendly adhesive, semi-transparent, patented center-hinge design sacral dressing and patented heel/elbow dressing design.

HCB102	Border Hydrocolloid Dressing (film backing), 2" x 2", 20/tray, 5 tray/cs
HCB106	Border Hydrocolloid Dressing (film backing), 6" x 6", 5/tray, 4 tray/cs
HCF204	Standard Hydrocolloid Dressing (foam backing), 4" x 4", 5/tray, 10 tray/cs
HCF208	Standard Hydrocolloid Dressing (foam backing), 8" x 8", 5/tray, 4 tray/cs
HCS100	Border Sacrum Hydrocolloid Dressing (film backing), 6" x 7", 5/tray, 4 tray/cs
HCT101	Thin Hydrocolloid Dressing (film backing), 4" x 4", 10/tray, 10 tray/cs

**COVIDIEN/KENDALL ULTEC HYDROCOLLOID DRESSING**

High absorption level. Optimum peel strength. Minimal adhesive flow. Provide an optimal healing environment. For minimal to moderate exudate. Sterile.

8884475400	4" x 4" Hydrocolloid Dressing, 50/cs
8884475602	6" x 6" Hydrocolloid Dressing, 30/cs
8884475800	6" x 8" Hydrocolloid Dressing, 30/cs
8884475826	8" x 8" Hydrocolloid Dressing, 30/cs

SMITH & NEPHEW CUTINOVA® HYDRO DRESSING

CUTINOVA® Hydro is a unique, hydro-active polyurethane matrix with an embedded super-absorber providing maximum exudate absorption. The special structure of this dressing offers the unique action of absorbing water from the wound fluid. Features: Waterproof, Helps prevent bacterial contamination, Water vapor permeable, Highly absorbent even under compression, Does not break down in the wound, Conformable, Semi-transparent.


66047441	2" x 2 3/8" Hydro Dressing, 10/bx, 10 bx/cs (474410000004)
66047443	4" x 4" Hydro Dressing, 5/bx, 10 bx/cs
66047445	6" x 8" Hydro Dressing, 3/bx, 16 bx/cs (474450000004)

DERMA SCIENCES HONEYCOLLOID DRESSING

31222	Dressing, Honeycolloid, 2" x 2", 10/bx, 10 bx/cs
31245	Dressing, Honeycolloid, 4" x 5", 10/bx, 5 bx/cs



483000



483200

SMITH & NEPHEW REPLICARE® HYDROCOLLOID DRESSINGS

Combines absorbency with all the advantages of a film dressing. A balance of adhesion and elasticity promotes excellent adherence over joints, while allowing the dressing to move more freely with the patient. Superior conformability provides greater durability and patient comfort.

483000	1 1/2" x 2 1/2" Hydrocolloid Dressing, 30/pkg, 6 pkg/cs
483100	4" x 4" Hydrocolloid Dressing, 5/pkg, 8 pkg/cs
483200	6" x 6" Hydrocolloid Dressing, 5/pkg, 8 pkg/cs
483300	8" x 8" Hydrocolloid Dressing, 5/pkg, 6 pkg/cs

**SMITH & NEPHEW REPLICARE® THIN HYDROCOLLOID DRESSINGS**

Maintains a moist wound environment for lightly exuding wounds such as skin tears and abrasions, as well as partial thickness wounds. A conformable, comfortable and cost-effective alternative for end-stage healing.

59484000	2" x 2 3/4" Hydrocolloid Dressing, 10/pkg, 15 pkg/cs
59484100	3 1/2" x 5 1/2" Hydrocolloid Dressing, 10/pkg, 12 pkg/cs
59484200	6" x 8" Hydrocolloid Dressing, 5/pkg, 10 pkg/cs



SMITH & NEPHEW REPLICARE® ULTRA DRESSING

59484600	4" x 4" Ultra Dressing, 5/bx, 8 bx/cs
59484700	6" x 6" Ultra Dressing, 5/bx, 8 bx/cs
59484900	7" x 8" Ultra Dressing, 5/bx, 6 bx/cs

3M™ TEGADERM™ HYDROGEL WOUND FILLER

Amorphous hydrogel, formulated to help provide a moist wound healing environment that has been shown to enhance wound healing. Available in single-dose tubes.

91110	Wound Filler, Hydrogel, 0.528 oz, 10/bx
91111	Wound Filler, Hydrogel, 0.88 oz, 10/bx



HARTMANN-CONCO AQUACLEAR® GEL SHEET DRESSING

Clear gel dressing for moist wound treatment. Transparent, hydrocellular gel dressings made of absorbent polyurethane polymers covered with a semi-permeable polyurethane film to prevent water and micro-organism contamination. The three-dimensional gel composition has a water content of 60%. Aquaclear Adhesive has an adhesive film border to ensure that it stays in place and to prevent bacterial penetration. Latex free.

900795	2" x 3" Gel Sheet Dressing, 5/bx
900796	4" x 4" Gel Sheet Dressing, 5/bx
900797	3" x 4" Adhesive Gel Sheet Dressing, 5/bx
900798	Dressing Gel Sheet Dressing, 5" x 5", Adhesive, 5/bx



J&J NU-GEL™ WOUND DRESSING

An occlusive non-adherent hydrogel dressing for light to moderate drainage. Maintains a moist environment to support the wound healing process. Helps cool, soothe and relieve pain. Ideal for burn-related wound sites and for partial-thickness chronic and traumatic wounds. Sterile.

2498	6" x 8" Wound Dressing, 5/bx, 6 bx/cs
------	---------------------------------------



J&J TIELLE™ HYDROPOLYMER DRESSING

Provides an exudate handling system for low to moderate exuding wounds. During use, the absorbent island gently expands as it absorbs exudate, maintaining a moist environment. Water resistant, permeable polyurethane backing. Gentle adhesive border. Nonwoven wicking layer. May be left in place for 5-7 days.

MTL100EN	2 3/4" x 3 1/2" Dressing, 10/bx, 5 bx/cs (was #2439) or MTL100
MTL101EN	4 1/4" x 4 1/4" Dressing, 10/bx, 5 bx/cs
MTL102	5 7/8" x 5 7/8" Dressing, 5/bx, 5 bx/cs
MTL103	7" x 7" Dressing, 5/bx, 5 bx/cs (was #2442)
MTL104	7" x 7" Sacrum Dressing, 5/bx, 5 bx/cs
MTP501	4 1/4" x 4 1/4" Dressing, 10/bx, 5 bx/cs (was #5440)
MTP502	5 7/8" x 7 3/4" Tielles™ Plus Dressing, 5/bx, 5 bx/cs (was #5441)
MTP505	5 7/8" x 5 7/8" Tielles™ Plus Dressing, 10/bx, 5 bx/cs (was #5444)
MTP506	5 7/8" x 5 7/8" Tielles™ Plus Sacrum Dressing, 10/bx, 5 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - HYPERGEL®

A water-based, hypertonic saline gel that softens and debrides eschar. 20% sodium chloride solution in gel form.

360500	0.17 oz Tube Hypertonic Gel, 10/bx, 4 bx/cs
361500	0.50 oz Tube Hypertonic Gel, 10/bx, 4 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPIFORM®

Breathable, comfortable, and waterproof. No extra fixation is needed. Thin, flexible, and discreet. Can be worn during all daily activities. Can be reapplied. One dressing can be used up to seven days.

293190	1.6" x 10" Self-Adherent Silicone Dressing, 5/bx, 5 bx/cs
293290	2" x 3" Self-Adherent Silicone Dressing, 5/bx
293490	4" x 7" Self-Adherent Silicone Dressing, 5/bx

MOLNLYCKE TENDRA® WOUND MANAGEMENT - NORMIGEL®

370500	0.17 oz Tube Isotonic Saline Gel, 10/bx, 4 bx/cs
371500	0.50 oz Tube Isotonic Saline Gel, 10/bx, 4 bx/cs

HARTMANN-CONCO HYDROTUL® HYDROACTIVE IMPREGNATED DRESSING

499581	2" x 2", Hydroactive Dressing, Latex Free, 10/bx
499583	4" x 5", Hydroactive Dressing, Latex Free, 10/bx
499584	6" x 8", Hydroactive Dressing, Latex Free, 10/bx

SMITH & NEPHEW IODOSORB & IODOFLEX

Bed wound preparation with 0.9% Cadexomer Iodine, developed to build a platform for the treatment of chronic wounds.

6602124014	Iodosorb Wound Gel, 10gm tube (0.9% Cadexomer Iodine), 4/bx, 12 bx/cs
6602125040	Iodosorb Wound Gel, 40gm tube (0.9% Cadexomer Iodine), 12/cs
6602133005	Iodoflex Wound Pad, 5 x 5gm (0.9% Cadexomer iodine), 5/bx, 12 bx/cs
6602134010	Iodoflex Wound Pad, 3 x 10gm (0.9% Cadexomer iodine), 3/bx, 12bx/cs

HARTMANN-CONCO AQUAFIL™ HYDROGEL IMPREGNATED DRESSINGS

49610000	2" x 2", Sterile, Latex Free, 10/bx
49620000	4" x 4", Sterile, Latex Free, 10/bx

Dressings/Gauzes



SMITH & NEPHEW CICA-CARE™ ADHESIVE SILICONE GEL SHEETS

A flexible, adhesive, semi-occlusive silicone sheet indicated for the management of hypertrophic and keloid scars. Combines the benefits of silicone gel sheets utilized in scar therapy with the advantages of self-adhesiveness and durability. The results are improved patient compliance. Can be washed and re-used for up to four weeks.

66250707 5" x 6" Silicone Gel Sheet, 10/pkg



INTEGRITY INTEGRA-GEL™ IMPREGNATED HYDROGEL WOUND DRESSINGS

Sterile, 100% cotton, 12-ply gauze dressing impregnated with a clean, "cool to the touch", greaseless, SUPER CLEAR water-based aloe vera and glycerin hydrogel formulation. Used as a primary dressing, directly on the wound. Stays moist longer, which helps to reduce the frequency of dressing changes and protects against wound dehydration. Conforms well, non-adherent to the wound bed. Latex free.

40-0220-B 2" x 2", 1 dressing/pkg, 30 pkg/bx
40-0440-B 4" x 4", 1 dressing/pkg, 15 pkg/bx



SMITH & NEPHEW FLEXIGEL® HYDROGEL SHEET DRESSINGS

A moist wound dressing which is made of polyacrylamide in a matrix with embedded hydrophilic polysaccharide particles. Used on moderately exuding wounds because of its capability of absorbing up to 5 times its own weight. It is moisture vapor permeable and has a cooling, soothing effect on first and second degree burns and partial thickness wounds.

59500000 2" x 2" Hydrogel Sheet Dressing, 5/pkg, 5 pkg/cs
59500200 4" x 4" Hydrogel Sheet Dressing, 5/pkg, 5 pkg/cs



SMITH & NEPHEW FLEXIGEL® STRANDS® ABSORBENT WOUND DRESSING

Cut into strips to form an even more absorbent dressing. Naturally conforms more easily to wound contours and surface, which helps create a moist wound environment. Increases the surface area in contact with the dressing, dramatically increasing the absorption of exudate. Easy to use and remove. Won't shred or break down in the wound. Removal is in one piece.

59500400 6 gms Absorbent Wound Dressing, 10/pkg, 4 pkg/cs



66027308



66027311

SMITH & NEPHEW INTRASITE™ HYDROGEL WOUND GEL

A sterile amorphous hydrogel, gently rehydrates necrotic tissue, facilitating autolytic debridement in the initial management of chronic wounds. Absorbs exudate from the wound surface, creates a moist wound environment and reduces the frequency of dressing changes. Comes in an easy-to-use Applipak delivery system. Contains no preservatives, and can be used throughout the wound management process.

66027308 8g (.28 oz) Applipak HydrogelWound Gel, 10/pkg, 4 pkg/cs
66027311 15g (.52 oz) Applipak HydrogelWound Gel, 10/pkg, 4 pkg/cs
66027313 25g (.88 oz) Applipak HydrogelWound Gel, 10/pkg, 4 pkg/cs



SMITH & NEPHEW SOLOSITE® GEL CONFORMABLE WOUND DRESSING

A clean wound dressing manufactured using non-woven gauze thoroughly saturated. The gel stays on the dressing and maintains intimate contact with the wound surface to create a moist wound healing environment. With the water-swellaable formula it can absorb excess exudate, allowing for fewer dressing changes. Preservatives are added to the produce to assist in preventing bacterial growth in the packaged dressing.


59482300 2" x 2" Wound Dressing, 10/pkg, 10 pkg/cs
59482340 2" x 2" Wound Dressing, 100/pkg
59482400 4" x 4" Wound Dressing, 10/pkg, 10 pkg/cs
59482440 4" x 4" Wound Dressing, 100/pkg



SMITH & NEPHEW SOLOSITE® WOUND GEL

An amorphous hydrogel creates a moist wound environment and is non-cytotoxic. The polymer used is water swellaable allowing it to absorb exudate yet remain gel-like until saturated. This allows to remain on the wound longer than gels using water-miscible (soluble) polymers.

449600 3 oz. Tube Wound Gel, 12/cs
449700 7 oz. Pushbutton Applicator Wound Gel, 6/cs



SOUTHWEST ELASTO-GEL™ PLUS ISLAND WOUND DRESSING

A convenient, quick, one step method of applying to a wound to give a waterproof covering. Utilizes a polyurethane barrier film with good moisture vapor transmission, which is coated with a mildly aggressive adhesive. It is designed to adhere firmly in place. Bacteriostatic/fungistatic and highly absorbent. Helps prevent the formation of slough, debris, necrotic tissue and eschar and will not stick to or dissolve into the wound.

DR8075 4" x 4" Sterile Island Dressing w/barrier film, 5/bx, 20 bx/cs



DR8050



DR8250

SOUTHWEST ELASTO-GEL™ PLUS WOUND DRESSING

An occlusive hydrogel sheet dressing with Sta-Fix tape included in the dressing pouch. The tape is included for added convenience when a water resistant tape is required to attach the dressing to the wound site. The Sta-Fix tape is a non-woven, polyester base tape with an open structure and coated with adhesive by a special process to produce a highly porous, breathable tape. The use of Sta-Fix tape allows the dressing to be adhered to difficult and complex contours. Utilizes a high concentration of glycerine (65%) as a water soluble humectant in a mixture of water and polymer matrix.

DR8050 4" x 4" Wound Dressing w/tape, 5/bx

DR8250 2" x 3" Wound Dressing w/tape, 5/bx



SOUTHWEST ELASTO-GEL™ TOE-AID

A sterile, bacteriostatic, absorbent hydrogel dressing attached to a unique "T" shaped border adhesive of our Sta-Fix tape. The Elasto-Gel™ pad has been shown to prevent and control the formation of slough, debris, necrotic tissue and eschar at the wound site. This high glycerine pad will not dry out and will not stick to or dissolve into the wound. The pad absorbs excess exudate and helps preserve and maintain a clean, moist, infection-free wound bed. When used as a post-operative dressing after toenail removal, Toe-Aid helps reduce the frequency of dressing changes, negates daily soaks, nearly eliminates pain and helps assure a positive clinical outcome. The water resistant border adhesive helps assure a secure placement.

DR8450 Toe-Aid/T-Shaped Dressing w/1.25" x 1.25" Gel Area, 3/bx



SOUTHWEST GOLDDUST™ HYDROGEL WOUND FILLER

A high absorbent hydrophilic polymer capable of absorbing large amounts of exudate and retaining the exudate in the matrix even under high pressure. May be used as dry granules on highly exudating wounds, but must be monitored for the potential of over-drying the tissue.

DR9300 GoldDust Super Absorbent Wound Filler, 3 grams, 10/bx, 10 bx/cs



SOUTHWEST STIMULEN™ ENHANCED COLLAGEN WOUNDCARE-LOTION

A thixotropic liquid that will form a gel when kept undisturbed for several hours at temperatures below 72°F (21°C). It becomes liquid immediately with agitation, at slightly higher temperatures or when applied to the skin. It is designed to moisturize and condition the skin. It may improve skin problems, such as: skin rash or irritation caused by various diseases or exposure to various chemicals, soaps and detergents, microbes, or other common problems - dry hands, feet and elbows, athletes foot, jock itch, etc. It may be used to soften the skin that may be causing wrinkles, callous and corns. May be effective for improving scar remodeling.

ST9555 Lotion, 2 oz/btl



SOUTHWEST ELASTO-GEL™ FACE MASK DRESSING

FM7500 Face Mask Dressing w/stretch material securing device, 2 Gel Masks & 2 Covers/bx, 5 bx/cs



SOUTHWEST ELASTO-GEL™ HORSE SHOE SHAPED DRESSING

High glycerine content hydrogel dressing with a unique horse shoe shaped design to cover the buttocks area, including the coccyx and ischials, which are the most common areas for pressure and pressure ulcers. Used as a wound management dressing as well as a preventative padding to help prevent the formation of pressure ulcers. The dressing utilizes a water resistant, hypoallergenic adhesive border tape system to hold the dressing firmly in place. Bacteriostatic/fungistatic, highly absorbent, helps prevent the formation of slough, debris, necrotic tissue and eschar and will not stick to or dissolve into the wound.

DR9051 8" x 8" Horseshoe-Shaped Dressing, 5/bx, 5 bx/cs


Dressings/Gauzes



SOUTHWEST ELASTO-GEL™ WOUND CARE

A superior absorbent dressing. Bacteriostatic reduces the change for infection. Will not damage tissue thanks to the patented glycerine content. Designed with unique features which provide for greater conformability, cushioning and padding compared to other available dressings. Mildly adhesive which allows for true inspection of the wound any time. It does not liquefy into the wound and can be cut to conform to hard to dress areas.

DR3000	3" Circle Dressing, 5/bx, 40 bx/cs
DR8000	4" x 4" Wound Dressing w/out tape, 5/bx
DR8200	2" x 3" Wound Dressing w/out tape, 5/bx
DR8500	5" x 5" Wound Dressing w/out tape, 5/bx
DR8600	6" x 8" Wound Dressing w/out tape, 5/bx
DR8700	12" x 12" Wound Dressing w/out tape, 5/bx, 5 bx/cs
DR8800	8" x 16" Wound Dressing w/out tape, 5/bx, 5 bx/cs



PRO ADVANTAGE® TRANSPARENT FRAME-STYLE FILM DRESSING

This multipurpose, transparent film dressing provides a thin flexible covering to help protect the wound surface or secure catheters and IVs. The dressings are water-resistant to enable patient bathing and transparent to permit observation without removing the dressing. Their moisture vapor permeability allows excess moisture to vent while maintaining a moist wound environment conducive to healing. The latex-free adhesive is gentle on the skin. Sterile.

P157402	Transparent Frame-Style Film Dressing, 2.375" x 2.75", 100/bx
P157404	Transparent Frame-Style Film Dressing, 4" x 4.75", 50/bx

3M™ TEGADERM™ MATRIX DRESSING WITH PHI™ TECHNOLOGY

Helping stalled wounds progress toward healing. A unique, affordable wound care solution that aids in the management of hard-to-heal wounds. It is designed to help jump-start wounds, enabling the wound to progress through the normal healing cascade. Indicated for use with chronic, non-infected wounds including: diabetic ulcers, pressure ulcers (stages II-IV), and venous stasis ulcers.

90570	Dressing, 2" x 2", 10/bx, 4 bx/cs
-------	-----------------------------------



3M™ TEGADERM™ + PAD FILM DRESSING WITH NON-ADHERENT PAD

A sterile, waterproof bacterial barrier which consists of a non-adherent absorbent pad bonded to a larger thin film transparent dressing. Breathable film.

3582	2" x 2 3/4" Transparent Dressing Size, 1" x 1 1/2" Pad Size, 50/bx, 4 bx/cs
3584	2 3/8" x 4" Transparent Dressing Size, 1" x 2 3/8" Pad Size, 50/bx, 4 bx/cs
3586	3 1/2" x 4" Transparent Dressing Size, 1 3/4" x 2 3/8" Pad Size, 25/bx, 4 bx/cs
3587	3 1/2" x 4 1/8" Transparent Dressing Size, 1 3/4" x 2 3/8" Pad Size, 25/bx, 4 bx/cs
3588	6" x 6" Transparent Dressing Size, 4" x 4" Pad Size, 25/bx, 4 bx/cs
3589	3 1/2" x 6" Transparent Dressing Size, 1 3/4" x 4" Pad Size, 25/bx, 4 bx/cs
3590	3 1/2" x 8" Transparent Dressing Size, 1 3/4" x 6" Pad Size, 25/bx, 4 bx/cs
3591	3 1/2" x 10" Transparent Dressing Size, 1 3/4" x 8" Pad Size, 25/bx, 4 bx/cs
3593	3 1/2" x 13 3/4" Transparent Dressing Size, 1 3/4" x 11 3/4" Pad Size, 25/bx, 4 bx/cs

3M™ TEGADERM™ CHG CHLORHESIDINE GLUCONATE IV SECUREMENT DRESSING

Proven to be effective at reducing skin flora. Clear, allowing continuous visualization of insertion site. Integrated design of dressing and gel pad reduces application steps and minimizes application error. Absorbs fluid.

1657	CHG Dressing, 3.5" x 4.5", 25/bx, 4 bx/cs
1658	CHG Dressing, 4" x 4.75", 25/bx, 4 bx/cs
1659	CHG Dressing, 4" x 6 1/8", 25/bx, 4 bx/cs

3M™ TEGADERM™ NON-ADHERENT CONTACT LAYER

5642	Dressing 3" x 4", Tegaderm™ Non-Adherent Contact Layer, Non-Irritating & Hypoallergenic, 10/bx, 4 bx/cs
5643	Dressing 3" x 8", Tegaderm™ Non-Adherent Contact Layer, Non-Irritating & Hypoallergenic, 10/bx, 4 bx/cs
5644	Dressing 8" x 10", Tegaderm™ Non-Adherent Contact Layer, Non-Irritating & Hypoallergenic, 5/bx, 4 bx/cs



1622W



1626



1630

3M™ TEGADERM™ TRANSPARENT FILM DRESSING FRAME STYLE

A breathable dressing, impermeable to liquids, bacteria and viruses*. "Picture-frame" delivery — preferred by caregivers. Secure adhesion, et gentle to the skin, breathable film, allows oxygen in and moisture vapor out, barrier to external contaminants and comprehensive clinical documentation. *Laboratory testing has proven 3M™ Tegaderm™ and Tegaderm™ HP Transparent Dressings provide a viral (HIV-1 and HBV) barrier while dressings remain intact without leakage.

1614	Transparent Dressing w/border, 2 3/8" x 2 3/4", 100/bx, 4 bx/cs
1616	Transparent Dressing w/border, 4" x 4 3/4", 50/bx, 4 bx/cs
1622W	Easy Application Frame Style Transparent Dressing, 1 3/4" x 1 3/4", 100/bx, 4 bx/cs
1624W	Easy Application Frame Style Transparent Dressing w/label, 2 3/8" x 2 3/4", 100/bx, 4 bx/cs
1626	Original Frame Style Transparent Dressing, 4" x 4 3/4", 50/bx, 4 bx/cs
1626W	Easy Application Frame Style Transparent Dressing w/label, 4" x 4 3/4", 50/bx, 4 bx/cs
1627	Original Frame Style Transparent Dressing, 4" x 10", 20/bx, 4 bx/cs
1628	Original Frame Style Transparent Dressing, 6" x 8", 10/bx, 8 bx/cs
1629	Original Frame Style Transparent Dressing, 8" x 12", 10/bx, 8 bx/cs
1630	Picture Frame Style Transparent Dressing, 4" x 4 1/2" w/label, 50/bx, 4 bx/cs
1634	Original Frame Style Transparent Dressing, 2 3/8" x 2 3/4", 100/bx, 4 bx/cs
9505W	Easy Application Frame Style Transparent Dressing w/label, 2 3/8" x 2 3/4", 20/bx, 10 bx/cs
9506W	Easy Application Frame Style Transparent Dressing w/label, 4" x 4 3/4", 10/bx, 10 bx/cs

3M™ TEGADERM™ AG MESH DRESSING WITH SILVER

A fast acting, long lasting antimicrobial barrier in an affordable silver dressing. Clinically effective fast sustained microbial barrier control for up to 7 days. Patient Friendly non-cytotoxic and non-staining. Versatile easy to cut, multiple sizes. Clinician Friendly easy to use, familiar, gauze-like porous material.

90500	2" x 2" Dressing, 5/bx, 4 bx/cs
90501	4" x 5" Dressing, 5/bx, 4 bx/cs
90502	4" x 8" Dressing, 5/bx, 4 bx/cs
90503	8" x 8" Dressing, 3/bx, 4 bx/cs

3M™ TEGADERM™ TRANSPARENT FILM ROLL

16002	Dressing, Transparent Film Roll, 2" x 11 yd, 4/cs
16004	Dressing, Transparent Film Roll, 4" x 11 yds, 4 rl/cs
16006	Dressing, Transparent Film Roll, 6" x 11 yds, 4 rl/cs



3M™ TEGADERM™ TRANSPARENT FILM DRESSING FIRST AID STYLE

A sterile, waterproof, clean thin film dressing. Designed to cover and protect I.V. and wound sites. Sterile dressing provides a bacterial and viral barrier to outside contaminants. Breathable film allows for exchange of moisture vapor and oxygen. Notched design provides better seal around the catheter. Unique frame style makes application easy and allow for precise dressing placement. Transparent window allows easy site monitoring without the need to change dressings.

1620	First-Aid Style Transparent Dressing, 2 3/8" x 2 3/4", 100/bx, 6 bx/cs
1621	First-Aid Style Transparent Dressing, 4" x 4 3/4", 50/bx, 4 bx/cs



3M™ TEGADERM™ HP (HOLDING POWER) TRANSPARENT FILM DRESSING

Providing all of the benefits of Tegaderm™ dressings, Tegaderm™ HP dressings have a special adhesive and film that differentiate it from standard Tegaderm™ film dressings. The hydrophilic nature of Tegaderm™ HP Film adhesive makes it more adherent and useful in moist conditions and difficult-to-dress areas. It offers holding power in challenging clinical conditions, reducing unscheduled dressing changes.

9519HP	2 3/8" x 2 3/8" w/secure strip, 100/bx, 4 bx/cs
9534HP	2 3/8" x 2 3/4", 100/bx, 4 bx/cs
9536HP	4" x 4 3/4", 50/bx, 4 bx/cs
9543HP	4 1/2" x 4 3/4", Sacral, 12/bx, 4 bx/cs
9545HP	2 1/8" x 2 1/2", Oval, 50/bx, 6 bx/cs
9546HP	4" x 4 1/2", Oval, 50/bx, 4 bx/cs
9548HP	5 1/2" x 6 1/2", Oval, 10/bx, 8 bx/cs



3M™ TEGADERM™ IV TRANSPARENT FILM DRESSING WITH BORDER

Tegaderm™ Advanced IV Dressings are designed for IV applications. The dressings have reinforced borders and notches designed to create a better seal around the catheter, extend wear time, reduce edge lift, and offer convenience and versatility because of pre-cut sterile tape strips.


1610	2" x 2 1/4" IV Transparent Dressing, 100/bx, 4 bx/cs
1633	2 3/4" x 3 1/4" Transparent Dressing, 100/bx, 4 bx/cs
1635	3 1/2" x 4 1/4" Transparent Dressing, 50/bx, 4 bx/cs
1650	4" x 6 1/8" Transparent Dressing, 25/bx, 4 bx/cs
1655	3 1/2" x 4 1/2" IV Transparent Dressing, 50/bx, 4 bx/cs

Dressings/Gauzes

3M™ TEGADERM™ TRANSPARENT FILM DRESSING WITH BORDER

A sterile, waterproof, clean thin film dressing. Designed to cover and protect I.V. and wound sites. Sterile dressing provides a bacterial and viral barrier to outside contaminants. Breathable film allows for exchange of moisture vapor and oxygen. Notched design provides better seal around the catheter. Unique frame style makes application easy and allow for precise dressing placement. Transparent window allows easy site monitoring without the need to change dressings.

- 1614 Transparent Dressing w/border, 2 3/8" x 2 3/4", 100/bx, 4 bx/cs
1616 Transparent Dressing w/border, 4" x 4 3/4", 50/bx, 4 bx/cs



3M™ CAVILON™ NO-STING BARRIER FILM

A sting-free, alcohol-free liquid film that dries quickly to form a breathable, transparent coating on the skin. Designed to protect intact or damaged skin from urine, feces, other body fluids, tape trauma, and friction.

- 3343 Small Foam Applicator, 1.0ml, 25/bx, 4 bx/cs
3344 Wipes, 1.0ml, 25/bx, 4 bx/cs
3345 Large Foam Applicator, 3.0ml, 25/bx, 4 bx/cs
3346 Pump Spray Bottle, 28.0ml, 12/cs

3M™ STERI-DRAPE™ ABSORBENT PAD

Made of a plastic film with a highly fluid absorbent layer on top with an adhesive strip.

- 9095 29" x 20" Absorbent Pad, 50/bx, 4 bx/cs

COVIDIEN/KENDALL POLYSKIN II TRANSPARENT DRESSINGS

Ideal for I.V. sites, donor sites, burns, ulcers and sutures. Simple, effective delivery system allows one-handed application. Highly permeable to oxygen and moisture vapor. Outstanding barrier to bacteria and fluid contaminants.

- 6640 Transparent Adhesive Dressing, 2" x 2 3/4", 100/bx, 4 bx/cs
6640B Transparent Adhesive Dressing, Non-Sterile, 2" x 2.75", 100/bx, 12 bx/cs
6641 Transparent Adhesive Dressing, 4" x 4 3/4", 50/bx, 4 bx/cs
6641B Transparent Adhesive Dressing, Non-Sterile, 4" x 4.75", 50/bx, 12 bx/cs
6642 Transparent Adhesive Dressing, 6" x 8", 10/bx, 4 bx/cs
6647 Transparent Adhesive Dressing, 4" x 8", 20/bx, 6 bx/cs
6648 Transparent Adhesive Dressing, 8" x 10", 10/bx, 8 bx/cs
6651 Transparent Adhesive Dressing, 1 1/2" x 1 1/2", 100/bx, 4 bx/cs

DUKAL TRANSPARENT SEMI-PERMEABLE DRESSINGS

With an easy two step application, this is an excellent dressing for IV sites, minor burns, and wounds. They keep contamination out and promote healing while allowing for moisture vapor transfer. This product features a removable patient label for the dressing or chart. Similar to Tegaderm.

- 4045 Dressing, 1 1/2" x 2", 1/pk, 100 pk/bx
4067 Dressing, 2 3/8" x 2 3/4", 1 pk, 100 pk/bx
4112 Dressing, 4" x 4 3/4", 1/pk, 50 pk/bx
4152 Dressing, 6" x 8", 1/pk, 10 pk/bx



DYNAREX VIEWGUARD TRANSPARENT DRESSINGS

Easy to apply. Conforms to difficult body contours. Provides a breathable barrier that lets air in while keeping moisture vapors and fluid out.

- 3642 2 3/8" x 2 3/4", Sterile, 100/bx, 4bx/cs
3643 4" x 4 3/4", Sterile, 50/bx, 4bx/cs
3644 6" x 8", Sterile, 10/bx, 8bx/cs
3646 2-3/8" x 2-3/4", 100/bx, 4bx/cs
3647 4" x 4-3/4", 50/bx, 4bx/cs
3648 6" x 8", 50/bx, 4bx/cs

J&J BIOCLUSIVE™ SELECT TRANSPARENT DRESSING

This semi-occlusive, bacterial/viral barrier helps protect the site from external contamination and permits constant monitoring. Patented frame delivery allows a one-hand aseptic application. The MVP Select Dressing is a high moisture vapor permeable film dressing that is indicated for invasive access devices only.

- 2455 2 3/4" X 2 3/8" Transparent Dressing, 50/bx, 4 bx/cs
2457 4" x 5" Transparent Dressing, 50/bx, 4 bx/cs
2474 1 3/4" x 2 3/4" Transparent Dressing, 100/bx, 4 bx/cs
2475 3" x 4" Transparent Dressing, 50/bx, 4 bx/cs



J&J BIOCLUSIVE™ TRANSPARENT DRESSING

Moisture vapor permeable dressing creates a moist wound healing environment which has less scarring than dry conventional gauze. Permits air circulation while maintaining a barrier against viruses, bacteria and fluid. Transparent for continuous wound monitoring, which is useful for management of dermal ulcers, burns, donor sites and surgical incisions. Easy application and removal.

- 2460 1 1/2" x 1 1/2" Transparent Dressing, 100/bx, 4 bx/cs
2461 2" x 3" Transparent Dressing, 100/bx, 4 bx/cs
2463 4" x 5" Transparent Dressing, 50/bx, 4 bx/cs
2465 5" x 7" Transparent Dressing, 20/bx, 5 bx/cs
2467 4" x 10" Transparent Dressing, 20/bx, 6 bx/cs
2469 8" x 10" Transparent Dressing, 10/bx, 8 bx/cs



SMITH & NEPHEW OPSITE™ IV3000 TRANSPARENT ADHESIVE FILM DRESSINGS

The patented Reactive™ film is 3 to 8 times more permeable to water vapor than standard film dressings. Under most situations, the unique film allows the I.V. site to remain dry. The risk of catheter-related infections may be reduced, unscheduled dressing changes are minimized, and patient comfort is enhanced.

- | | |
|------|--|
| 4007 | 2 3/8" x 2 3/4" 1-Hand Peripheral Film Dressing, 100/pkg, 9 pkg/cs |
| 4008 | 4" x 4 3/4" 1-Hand Central Film Dressing, 50/pkg, 4 pkg/cs |
| 4649 | 4" x 8" Tape Handles Epidural Film Dressing, 50/pkg, 4 pkg/cs |
| 4654 | 2" x 1 1/2" 1-Hand Pediatric Film Dressing, 100/pkg, 10 pkg/cs |
| 4925 | 4" x 5" Tape Handles Central Film Dressing, 10/bx, 10 bx/cs |
| 4973 | 4" x 5" Tape Handles Central Film Dressing, 50/bx, 4 bx/cs |



SMITH & NEPHEW OPSITE™ TRANSPARENT ADHESIVE FILM DRESSINGS

Transparent polyurethane membrane which creates a moist environment by trapping the wound exudate. Because the dressing is waterproof, and aids in the prevention of bacterial contamination, it reduces the risk of infection from airborne bacteria.

- | | |
|------|---|
| 4542 | 11" x 4" Film Dressing, 10/pkg, 6 pkg/cs |
| 4963 | 5 1/2" x 4" Film Dressing, 10/pkg, 10 pkg/cs |
| 4967 | 5 1/2" x 10" Film Dressing, 20/pkg, 6 pkg/cs |
| 4975 | 5 1/2" x 4" Film Dressing, 50/pkg, 4 pkg/cs |
| 4986 | 11" x 6" Film Dressing, 10/pkg, 6 pkg/cs |
| 4987 | 11" x 11 3/4" Film Dressing, 10/pkg, 6 pkg/cs |
| 4988 | 11" x 17 3/4" Film Dressing, 10/pkg, 6 pkg/cs |
| 4989 | 17 3/4" x 21 5/8" Film Dressing, 10/pkg, 4 pkg/cs |



SMITH & NEPHEW OPSITE® FLEXIFIX TRANSPARENT FILM ROLLS

A waterproof, transparent film dressing which provides excellent adherence even while patients bathe and shower. It comes in roll form and aids in the prevention of bacterial contamination, thus reducing the risk of secondary infection. Flexifix provides high moisture vapor permeability and allows the skin to breathe and may reduce the risk of maceration.

- | | |
|----------|--|
| 66000040 | 2" x 11 yds Transparent Film Roll, 24/cs |
| 66000041 | 4" x 11 yds Transparent Film Roll, 12/cs |



SMITH & NEPHEW OPSITE™ FLEXIGRID™ TRANSPARENT ADHESIVE FILM DRESSINGS

Easy one-hand delivery allows this dressing to be used for both I.V. and wound care.

- | | |
|----------|---|
| 66024628 | 2 3/8" x 2 3/4" Peripheral Transparent Adhesive Film Dressing, 100/pkg |
| 66024629 | 4" x 4 3/4" Central Line Transparent Adhesive Film Dressing, 10/pkg, 8 pkg/cs |
| 66024630 | 4" x 4 3/4" Central Line Transparent Adhesive Film Dressing, 50/pkg |
| 66024631 | 6" x 8" Transparent Adhesive Film Dressing, 10/pkg, 10 pkg/cs |
| 66024632 | 4 3/4" x 10" Transparent Adhesive Film Dressing, 20/pkg |



3M™ STERI-DRAPE™ WOUND EDGE PROTECTOR

A wound edge protector that separates edges of the incision from contact with visceral contents, instruments and gloves during the procedure, while allowing for visualization of the wound edges. Four adhesive patches. Semi-rigid ring.

- | | |
|------|--|
| 1073 | Wound Edge Protector, 35" x 35", Ring Diameter: 4 3/4", 10/bx, 4 bx/cs |
| 1074 | Wound Edge Protector, 35" x 35", Ring Diameter: 6 5/8", 10/bx, 4 bx/cs |
| 1075 | Wound Edge Protector, 35" x 35", Ring Diameter: 8 5/8", 10/bx, 4 bx/cs |
| 1076 | Wound Edge Protector, 35" x 35", Ring Diameter: 10", 10/bx, 4 bx/cs |



SMITH & NEPHEW PRIMAPORE™ I.V. ADHESIVE DRESSINGS

A one-step island dressing that helps prevent site contamination and avoids the multiple steps of building a gauze and tape dressing. This water-resistant, non-woven polyester dressing helps reduce the potential of moisture contamination and the strike-through effect. Conformable and very gentle on the skin.

- | | |
|----------|--|
| 66000317 | 4" x 3 1/8" Central IV Adhesive Dressing, 20/pkg, 10 pkg/cs |
| 66000318 | 6" x 3 1/8" Central IV Adhesive Dressing, 20/pkg, 10 pkg/cs |
| 7133 | 2.75" x 2", Pad Size 40mm x 27mm, Peripheral IV Adhesive Dressing, 100/pkg |

Dressings/Gauzes

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEFILM®

Breathable and transparent. Easy application. Helps protect from friction injuries. Barrier for leakage and bacteria. Secure fixation. Gentle to skin and wound site. Bacteria proof.

270600	2.4" x 2.6" Adhesive Polyurethane Film Dressing, 100/bx, 4 bx/cs
271500	4" x 5" Adhesive Polyurethane Film Dressing, 70/bx, 3 bx/cs
272500	4" x 10" Adhesive Polyurethane Film Dressing, 20/bx, 5 bx/cs
273000	6" x 8.5" Adhesive Polyurethane Film Dressing, 10/bx, 8 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEFIX®

Skin-friendly, comfortable to wear. Pre-measured protection backing allows for convenient measurement/cutting. Conforms to body contours, allowing free, unrestricted movement. Air permeable.

310250	1" x 11 yds Fixation Fabric Dressing, 39/cs
310500	2" x 11 yds Fixation Fabric Dressing, 20/cs
311000	4" x 11 yds Fixation Fabric Dressing, 20/cs
311500	6" x 11 yds Fixation Fabric Dressing, 20/cs
312000	8" x 11 yds Fixation Fabric Dressing, 10/cs

HARTMANN-CONCO HYDROFILM® TRANSPARENT FILM DRESSING

685755	2.4" x 2.75", Transparent Dressing, Latex Free, 10/bx
685756	2.4" x 2.75", Transparent Dressing, Latex Free, 100/bx
685757	4" x 5", Transparent Dressing, Latex Free, 10/bx
685758	4" x 5", Transparent Dressing, Latex Free, 100/bx
685759	4" x 6", Transparent Dressing, Latex Free, 10/bx
685760	4" x 6", Transparent Dressing, Latex Free, 100/bx
685761	6" x 8", Transparent Dressing, Latex Free, 10/bx
685762	6" x 8", Transparent Dressing, Latex Free, 50/bx
685770	2" x 2.8", Transparent Dressing w/absorbent pad, Latex Free, 5/bx
685771	2" x 2.8", Transparent Dressing w/absorbent pad, Latex Free, 50/bx
685772	3.5" x 4", Transparent Dressing w/absorbent pad, Latex Free, 5/bx
685773	3.5" x 4", Transparent Dressing w/absorbent pad, Latex Free, 50/bx
685774	3.5" x 6", Transparent Dressing w/absorbent pad, Latex Free, 5/bx
685775	3.5" x 6", Transparent Dressing w/absorbent pad, Latex Free, 25/bx

SMITH & NEPHEW OPSITE™ IV3000 FRAME DELIVERY ADHESIVE FILM DRESSINGS

The Large frame border opening allows for easy molding around larger catheters and connectors. Large tabs make handling the dressing easy, even with gloves. NOW LATEX FREE.


59410082	2.38" x 2.75" IV Frame Dressing, 100/bx, 10 bx/cs
59410882	4" x 4.75" IV Frame Dressing, 50/bx, 4 bx/cs



SMITH & NEPHEW SKIN-PREP® PROTECTIVE DRESSINGS

This protective liquid dressing may reduce adhesive removal trauma during removal of tapes and films. Allows skin to breathe so tapes and films may adhere better. May also assist in increasing the interval between dressing changes. The No-Sting Skin-Prep® does not contain alcohol and will reduce the potential for stinging.

420200	Spray-On 4 oz. Liquid Dressing, 12/cs
420400	Wipes Liquid Dressing, 50/pkg, 20 pkg/cs
420500	Swabs Liquid Dressing, 50/pkg, 10 pkg/cs
59420600	No-Sting Wipes Liquid Dressing, 50/pkg, 20 pkg/cs
59420700	No-Sting Swabs Liquid Dressing, 50/pkg, 10 pkg/cs



PRO ADVANTAGE® CALCIUM ALGINATE DRESSING

Calcium alginate is a wound exudate absorber for moderately to heavily draining wounds. It may be appropriately used for Stage II and IV pressure ulcers, venous ulcers, arterial and diabetic ulcers, donor sites, superficial burn wounds, surgical incisions, abrasions, lacerations and traumatic wounds. Available in flat dressings or wound packing rope.

P157702	Calcium Alginate Dressings, 2" x 2", 10/bx
P157704	Calcium Alginate Dressings, 4.25" x 4.25", 10/bx
P157712	Calcium Alginate Dressings, 12" Ribbon, 5/bx

J&J FIBRACOL™ COLLAGEN-ALGINATE WOUND DRESSING

Collagen compound breaks down gradually and maintains a moist wound environment. Alginate creates a hydrogel. Nonadherent and maintains integrity when wet. Can be cut to size. Sterile.

2981	2" x 2" Wound Dressing, 12/bx, 6 bx/cs
2982	4" x 4 3/8" Wound Dressing, 12/bx, 6 bx/cs
2983	4" x 8 3/4" Wound Dressing, 6/bx, 6 bx/cs
2984	3/8" x 3/8" x 15 3/4" Wound Dressing, 6/bx, 6 bx/cs

J&J MEDICAL SILVERCEL ANTIMICROBIAL ALGINATE DRESSING


SILVERCEL Antimicrobial Alginate Dressing is a sterile, non-woven pad composed of high tensile strength alginate, carboxymethylcellulose (CMC) and silver coated nylon fibers. Contains elemental silver (8%) as a sustained release formulation. The sustained release of silver aids the dressing in creating an antimicrobial environment. Antimicrobial properties are built-in through the use of X-STATIC® Silver Fibers.

800112	1" x 12", Sterile Dressing, 5/bx, 5 bx/cs
800202	2" x 2", Sterile Dressing, 10/bx, 5 bx/cs
800404	4.25" x 4.25", Sterile Dressing, 10/bx, 5 bx/cs
800408	4" x 8", Sterile Dressing, 5/bx, 5 bx/cs

J&J NU-DERM™ ALGINATE WOUND DRESSING

A sterile, non-woven pad consisting of a co-spun high G (guluronic acid) alginate (85%) and carboxy methyl cellulose (CMC) fiber (15%) with spin finish. Designed to be fast gelling and to allow removal intact. Helps create a moist wound healing environment. Assists in autolytic debridement. Maintains integrity when wet. High tensile strength - dry. Product and packaging are latex free. Outstanding absorbent capacity allows for fewer dressing changes.

AWD112	1" x 12" Rope Alginate Wound Dressing, 5/bx, 5 bx/cs
AWD202	2" x 2" Alginate Wound Dressing, 10/bx, 5 bx/cs
AWD404	4" x 4" Alginate Wound Dressing, 10/bx, 5 bx/cs
AWD408	4" x 8" Alginate Wound Dressing, 10/bx, 5 bx/cs



3M™ TEGADERM™ HI & HG ALGINATE DRESSING

Sterile, highly conformable, and absorbent primary wound care dressings for use on moderately to heavily exuding wounds. Tegaderm HI alginate dressing (high integrity) and Tegaderm HG alginate dressing (high gelling) are available in a variety of sizes.

90112	High-Integrity Alginate Dressing, 4" x 4", 10/bx, 5 bx/cs
90120	High-Integrity Alginate Dressing, 12" Rope, 5/bx, 4 bx/cs
90212	High-Gelling Alginate Dressing, 4" x 4", 10/bx, 5 bx/cs
90220	High-Gelling Alginate Dressing, 12" Rope, 5/bx, 4 bx/cs



HARTMANN-CONCO SORBALGON® CALCIUM ALGINATE DRESSING

A soft, highly absorbent (absorbs more than 14 times its weight) dressing that quickly forms a hydrophilic gel to protect and cover the wound bed. While many alginates dissolve, this dressing keeps its form so it's easy to remove, generally on one piece. Latex free.

49200000	2" x 2" Calcium Alginate Dressing, 10/bx
49210000	4" x 4" Calcium Alginate Dressing, 10/bx
49220000	12" Ribbon Calcium Alginate Dressing, 10/bx
49230000	4" x 8", Used for very large wounds or several smaller wounds, Latex Free, 5/bx



COVIDIEN/KENDALL CURASORB CALCIUM ALGINATE DRESSINGS

Absorbs exudate from wound and creates a protective cushioning gel and optimal moist healing environment. Protective gel binds exudate and helps prevent skin maceration. Easy, painless dressing removal minimizes tissue trauma and maximizes patient comfort. Absorbs approximately 20 times dressing weight.

9231	12" Rope Calcium Alginate Dressing, Sterile 1s, 5/bx, 4 bx/cs
9232	2" x 2" Calcium Alginate Dressing, Sterile 1s, 10/bx, 100 bx/cs
9233	4" x 4" Calcium Alginate Dressing, Sterile 1s, 10/bx, 50 bx/cs
9236	4" x 4" Plus Calcium Alginate Dressing, Sterile 1s, 10/bx, 50 bx/cs
9238	8" x 4" Calcium Alginate Dressing, Sterile 1s, 5/bx, 50 bx/cs
9239	6" x 10" Calcium Alginate Dressing, Sterile 1s, 10/bx, 50 bx/cs
9240	4" x 5 1/2" Calcium Alginate Dressing, Sterile 1s, 10/bx, 5 bx/cs
9242	Calcium Alginate Dressing w/zinc, 12" x 24", Sterile 1's, 5/cs
9243	24" Rope Calcium Alginate Dressing, Sterile 1s, 5/bx, 20 bx/cs
9244	36" Rope Calcium Alginate Dressing, Sterile 1s, 5/bx, 20 bx/cs
9351	12" Rope Calcium Alginate Dressing w/zinc, Sterile 1s, 5/bx, 20 bx/cs
9354	2" x 2" Calcium Alginate Dressing w/zinc, Sterile 1s, 10/bx, 100 bx/cs
9355	4" x 4" Calcium Alginate Dressing w/zinc, Sterile 1s, 10/bx, 5 bx/cs
9356	4" x 8" Calcium Alginate Dressing w/zinc, Sterile 1s, 5/bx, 10 bx/cs

MOLNLYCKE TENDRA® MELGISORB® AG ALGINATE DRESSING

This product can be used under compression bandages. Because it contains alginate, it assists in supporting the control of minor bleeding in superficial wounds. Suitable for use, under medical supervision on wounds that are critically colonized. Indicated for external use only. Absorbs up to 45% more than market leading silver hydrofiber dressing. High fiber strength allows removal in one piece. Rapid and sustained antimicrobial effect. Sustained silver release up to 4 days.

255050	Dressing Silver Containing Alginate, 2" x 2", 10/bx, 10 bx/cs
255100	Dressing Silver Containing Alginate, 4" x 4", 10/bx, 10 bx/cs
255150	Dressing Silver Containing Alginate, 6" x 6", 10/bx, 10 bx/cs
255200	Dressing Silver Containing Alginate, 8" x 12", 5/bx, 10 bx/cs
255600	Dressing Silver Containing Alginate, 1.2" x 18", 10/bx, 10 bx/cs



SMITH & NEPHEW ALGISITE™ M CALCIUM ALGINATE DRESSING

Made of fast-gelling, high mannuronic acid fibers from calcium alginate, which is extracted and purified from a natural, renewable source - seaweed. When exudate comes into contact with the dressing, it forms a hydrophilic gel that remains in intimate contact with the wound surface, helping to create and maintain a moist wound environment.

59480100	2" x 2" Calcium Alginate Dressing, 10/pkg, 10 pkg/cs
59480200	4" x 4" Calcium Alginate Dressing, 10/pkg, 10 pkg/cs
59480300	6" x 8" Calcium Alginate Dressing, 10/pkg, 6 pkg/cs
59480400	3/4" x 12" Calcium Alginate Dressing, 10/pkg, 6 pkg/cs


MOLNLYCKE TENDRA® WOUND MANAGEMENT - MELGISORB®

Non-adherent to moist wound bed. High gel blocking. Helps maintain a moist wound environment. Reduces the risk of maceration. Absorbs moderate to high amounts of exudate.

250600	2" x 2" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
251100	4" x 4" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
251500	4" x 8" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
253000	12.5" Rope Absorbent Calcium Alginate Dressing, 5/bx, 10 bx/cs

**Brand names
you can trust.**

Dressings/Gauzes



PRO ADVANTAGE® WATERPROOF COMPOSITE ISLAND DRESSING

This dressing has four distinct layers designed to cover all stages of wound care. Layer 1, a non-adherent contact layer, protects the wound site and will not stick to granulation tissue. Layer 2, a soft pad, absorbs drainage. Layer 3, a non-woven backing with an adhesive border, securely holds the contact pad in place. Layer 4, a semi-occlusive, polyurethane film, helps deter external contaminants and maintains an optimal wound environment. The vapor permeable film coating allows showering and protects the wound from water and contamination. It also protects clothing and bed linens from blood or exudate strike-through. Sterile.

- P157304 Waterproof Composite Island Dressing, 4" x 4" w/2" x 2" pad, 10/bx
 P157306 Waterproof Composite Dressing, 6" x 6" w/4" x 4" pad, 10/bx



3M™ MEDIPORE™ +PAD SOFT CLOTH ADHESIVE WOUND DRESSING

Sterile all-in-one island dressing that features a breathable soft cloth tape backing with a unique non-adherent pad.

- 3562 2" x 2 3/4" Wound Dressing Size, 1" x 1 1/2" Pad Size, 50/bx, 4 bx/cs
 3564 2 3/8" x 4" Wound Dressing Size, 1" x 2 3/8" Pad Size, 50/bx, 4 bx/cs
 3566 3 1/2" x 4" Wound Dressing Size, 1 3/4" x 2 3/8" Pad Size, 25/bx, 4 bx/cs
 3568 Wound Dressing, 6" x 6", Pad Size 4" x 4 1/2", 25/bx, 4 bx/cs
 3569 3 1/2" x 6" Wound Dressing Size, 1 3/4" x 4" Pad Size, 25/bx, 4 bx/cs
 3570 3 1/2" x 8" Wound Dressing Size, 1 3/4" x 6" Pad Size, 25/bx, 4 bx/cs
 3571 3 1/2" x 10" Wound Dressing Size, 1 3/4" x 8" Pad Size, 25/bx, 4 bx/cs
 3573 3 1/2" x 13 3/4" Wound Dressing Size, 1 3/4" x 11 3/4" Pad Size, 25/bx, 4 bx/cs

DERMA SCIENCES COVERALL™ 5 N 1 ISLAND DRESSING

- 89066 Bordered Wound Dressing, 6" x 6", Non-Sterile, 30/bg, 4 bg/cs
 89144 4" x 4" Island Dressing, 25/bx, 4 bx/cs
 89166 6" x 6" Island Dressing, 25/bx, 4 bx/cs
 89168 6" x 8" Island Dressing, 25/bx, 4 bx/cs

DERMA SCIENCES XTRASORB™ SUPER ABSORBENT DRESSING

Flat edges conforms well to healthy skin and facilitate application. The elastic wound contact layer prevents adhesion to the wound. Absorbent core with super absorbent polymer absorbs large quantities of wound exudate, forming a gel that binds and locks the exudate away from the wound. A blue strike-through protective outer layer protects clothing and bed sheets against soiling.

- 89545 4" x 5" Sterile Dressing, Absorbent, 10/bx, 10 bx/cs
 89569 6" x 9" Sterile Dressing, Absorbent, 10/bx, 5 bx/cs



HARTMANN-CONCO COSMOPORE® ADHESIVE WOUND DRESSING

A non-adherent wound contact layer that absorbs exudate for a cleaner wound bed. 100% cotton, air and water permeable. Non-woven adhesive cover has rounded edges to prevent the dressing from becoming detached even after being worn for prolonged periods of time. Adhesive border protects the skin and wound area from harmful contamination and trauma. Latex free.

- 48000000 2" x 3.5", Perfect size for Peripheral IV Dressing or arthroscopic surgery sites, Latex Free, 50/bx, 8bx/cs
 48100000 3 1/8" x 4" Wound Dressing, 25/bx, 8 bx/cs
 48200000 3 1/8" x 6" Wound Dressing, 25/bx, 8 bx/cs
 48300000 4" x 8" Wound Dressing, 25/bx, 8 bx/cs
 48400000 4" x 10" Wound Dressing, 25/bx, 8 bx/cs
 48500000 4" x 13 3/4" Wound Dressing, 25/bx, 8 bx/cs
 49500000 4" x 4" Wound Dressing, 25/bx, 8 bx/cs
 49510000 6" x 6" Wound Dressing, 25/bx, 8 bx/cs



SMITH & NEPHEW CONFORMANT 2® WOUND VEILS

Designed to minimize adherence and save nursing time. Permeable and flexible, it can be used over most wound sites under an absorbent outer dressing. Made entirely from high-density polyethylene, it tears easily and uniformly in the direction of the grain for dressing customization. Can be used with a variety of topical agents.


- 5955044 4" x 4" Wound Veils, Sheets, 48/cs
 59551212 12" x 12" Wound Veils, Sheets, 48/cs
 59551224 12" x 24" Wound Veils, Sheets, 48/cs
 59552436 24" x 36" Wound Veils, Sheets, 50/cs
 5955305 3" x 5 yds Wound Veils, Rolls, 20/cs
 5955412 4" x 12" Wound Veils, Sheets, 48/cs
 5955602 6" x 2 yds Wound Veils, Rolls, 20/cs



SMITH & NEPHEW COVRSITE® COVER DRESSINGS

Indicated for the management of acute or chronic wounds. Designed as a secondary cover dressing for gels, gel sheets, alginates, wound fillers and non-adhesive foams. Can also be used as an alternative to tape and gauze over central line dressings. The gentle adhesive is designed so that the dressing can be changed daily. Replaces traditional tape and gauze with a single step dressing. Water-resistant, extensible, and conformable.

59714000	4" x 4" Cover Dressing, 10/pkg, 10 pkg/cs
59714100	4" x 4" Cover Dressing, 30/pkg, 10 pkg/cs
59714300	6" x 6" Cover Dressing, 10/pkg, 10 pkg/cs
59714400	6" x 6" Cover Dressing, 30/pkg, 10 pkg/cs



SMITH & NEPHEW COVRSITE® PLUS COMPOSITE DRESSING

A waterproof version of CovrSite® which is indicated as a primary or secondary dressing in the management of acute or chronic wounds. The waterproof film layer covering the dressing helps protect the wound from contamination by urine or feces.


59715000	4" x 4" Composite Dressing, 10/pkg, 10 pkg/cs
59715100	6" x 6" Composite Dressing, 10/pkg, 10 pkg/cs



SMITH & NEPHEW OPSITE® POST-OP COMPOSITE DRESSINGS

Optimal dressing for post-operative wound care. This sterile, island-style dressing provides a waterproof barrier, and aids in the prevention of bacterial contamination while allowing the skin to breathe. Wound drainage can be monitored through the clear film without removing the dressing, promoting maximum wear time.

66000708	2 1/2" x 2" Composite Dressing, 100/bx
66000709	3 3/4" x 3 3/8" Composite Dressing, 20/bx, 10 bx/cs
66000710	4 3/4" x 4" Composite Dressing, 10/bx, 10 bx/cs
66000712	6 1/8" x 3 3/8" Composite Dressing, 20/bx, 10 bx/cs
66000713	8" x 4" Composite Dressing, 20/bx, 10 bx/cs
66000714	10" x 4" Composite Dressing, 20/bx, 10 bx/cs
66000715	11 3/4" x 4" Composite Dressing, 20/bx, 8 bx/cs
66000716	13 3/4" x 4" Composite Dressing, 20/bx, 8 bx/cs



SMITH & NEPHEW PROFORE™ WOUND CONTACT LAYER

Allows exudate to pass freely to an absorptive dressing or bandage layer, thereby preventing exudate pooling and the possibility of maceration. Assists in preventing wound trauma during removal. The contact layer protects the wound from absorbent pad materials, while the natural viscose rayon monofilament produces no lint or breakdown residues to complicate wound cleansing.

66000701	Wound Contact Layer, 5 1/2" x 8", 50/pkg, 4 pkg/cs
----------	--

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPORE PRO®

670890	Absorbent Shower Proof Adhesive Dressing, 2.5" x 3", Latex Free, Sterile, 60/bx, 8 bx/cs
670990	Absorbent Shower Proof Adhesive Dressing, 3.6" x 4", Latex Free, Sterile, 60/bx, 8 bx/cs
671090	Absorbent Shower Proof Adhesive Dressing, 3.6" x 6", Latex Free, Sterile, 40/bx, 8 bx/cs
671190	Absorbent Shower Proof Adhesive Dressing, 3.6" x 8", Latex Free, Sterile, 30/bx, 6 bx/cs
671290	Absorbent Shower Proof Adhesive Dressing, 3.6" x 10", Latex Free, Sterile, 30/bx, 6 bx/cs
671390	Absorbent Shower Proof Adhesive Dressing, 3.6" x 12", Latex Free, Sterile, 30/bx, 9 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPORE®

Skin-friendly, water-based adhesive. Low-adherent wound contact layer minimizes risk of adherence to the wound. Fluid-repellent backing. Easy-to-remove protection backing allows for convenient application. Air permeable.


670800	2.5" x 3" Absorbent Island Dressing, 60/bx, 8 bx/cs
670900	3.6" x 4" Absorbent Island Dressing, 50/bx, 8 bx/cs
671000	3.6" x 6" Absorbent Island Dressing, 50/bx, 8 bx/cs
671100	3.6" x 8" Absorbent Island Dressing, 30/bx, 6 bx/cs
671200	3.6" x 10" Absorbent Island Dressing, 30/bx, 6 bx/cs
671300	3.6" x 12" Absorbent Island Dressing, 30/bx, 9 bx/cs
671400	3.6" x 14" Absorbent Island Dressing, 30/bx, 6 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MEPILEX® HEEL

Mepilex® Heel with Safetac soft silicone technology is shaped to fit any heel so there is no need for measuring or cutting. Its high conformability ensures maximum skin contact over the entire dressing area. The Safetac layer provides a good seal to reduce the risk of maceration and minimizes the trauma and pain at dressing changes.

288100	5" x 8", 5/bx
--------	---------------

Dressings/Gauzes



PRO ADVANTAGE® BORDERED ISLAND GAUZE DRESSINGS

An absorbent island dressing constructed with a heavy-duty absorbent pad and conformable, non-woven fabric adhesive backing. The porous structure of the Sontara® non-woven fabric ensures air and water vapor permeability. The fabric is stretchable, conforms well to body contours and allows body movements when in place. The adhesive has been designed to be gentle to the skin but at the same time provide secure fixation. The wound pad can absorb moderate amounts of wound fluids and the low-adherent surface reduces the risk of adherence to the wounds.

- P157204 Bordered Island Gauze Dressing, 4" x 4" w/ 2" x 2" Pad, 25/bx
 P157206 Bordered Island Gauze Dressing, 6" x 6" w/ 4" x 4" Pad, 25/bx

DUKAL STERILE BURN GAUZE PADS

Can be used in conjunction with medications such as silver nitrate to enhance healing and facilitate cleaner wound environment. Whether being used for debridement or to dress a grafted area or donor site, these dressings are the highest quality burn dressings available on the market today.

- 9170 Burn Gauze Pad, 18" x 18", 10-Ply, Sterile 1s, 2 x1 Mesh, Cut Edges, 100/cs

DERMA SCIENCES DUSOFT DRESSING

- 84122 Dressing, Non Woven, 2" x 2", Sterile 2's, 4-Ply, 25 pkgs of 2/bx, 60 bx/cs
 84133 Dressing, Non Woven, 3" x 3", Sterile 2's, 4-Ply, 25 pkgs of 2/bx, 48 bx/cs
 84144 Dressing, Non Woven, 4" x 4", Sterile 2's, 4-Ply, 25 pkgs of 2/bx, 24 bx/cs

DERMA SCIENCES DUSOFT SPONGE

- 94122 Sponge, Drain, Non Woven, Non Sterile, 2" x 2", 4-Ply, 200/slv, 40 slv/cs
 94133 Sponge, Drain, Non Woven, Non Sterile, 3" x 3", 4-Ply, 200/slv, 20 slv/cs
 94144 Sponge, Drain, Non Woven, Non Sterile, 4" x 4", 4-Ply, 200/slv, 10 slv/cs

SMITH & NEPHEW PRIMAPORE™ SPECIALTY ABSORBENT DRESSINGS

A soft, conformable, water-resistant, non-adherent island dressing that ensures patient comfort and immediate acceptability. It is an economical all-in-one tape and gauze product with a first-aid dressing application ready for aseptic application. The acrylic adhesive, evenly spread on the non-woven backing surface, simplifies safe and secure fixation over the wound site.

- 66000319 8" x 4" Absorbent Dressing, 20/pkg, 10 pkg/cs
 66000321 11 3/4" x 4" Absorbent Dressing, 20/pkg, 10 pkg/cs
 66007140 13 3/4" x 4" Absorbent Dressing, 20/bx, 10 bx/cs



SMITH & NEPHEW EXU-DRY® WOUND DRESSINGS & GARMENTS

Uniquely designed to minimize trauma, save nursing time and lower total dressing costs. Replaces non-adherent dressings, gauze and abdominal pads. Designed to minimize adherence and is more comfortable for patients during the between dressing changes. Its Anti-Shear Layer™ helps reduce friction and shearing. The one-piece, multi-layer sealed construction eliminates frayed edges, loose threads and linting. Will not disturb new granulation tissue. Non-occlusive, soft and pliable. Can be used as a wet or dry dressing on clean or contaminated wounds. Can be used with a variety of topical agents.

- 5999003 3" Wound Dressing Disc, 50/bx, 2 bx/cs
 5999003S 3" Wound Dressing Slit Disc, Medium Absorbency, 100/cs
 5999004120 4" x 6" Wound Dressing, Full Absorbency, 120/cs
 5999006 Leg/Arm 6" x 9" Dressing, Infant, 48/cs
 5999009 Leg/Arm 9" x 15" Dressing, Toddler, 30/cs
 5999018 15" x 18" Wound Dressing, 30/cs
 5999024 15" x 24" Wound Dressing, 30/cs
 5999028 20" x 28" Wound Dressing, 20/cs
 5999034 3" x 4" Wound Dressing, Full Absorbency, 100/cs
 5999034S 3" x 4" Wound Dressing Slit Tube Dressing, Full Absorbency, 50/cs
 5999101 4" x 6", Slit Tube Dressing, 100/cs
 5999BP1 Adult Buttocks Dressing, 10/cs
 5999DEL Elbow/Knee/Heel Dressing, 20/cs
 5999FDL Large Boot/Foot Dressing, 20/cs
 5999FDM Medium (Child) Boot/Foot Dressing, 20/cs
 5999FM1 Adult Face Dressing, 20/cs
 5999INC9 3" x 9" Incision Dressing, 50/cs
 5999L72 36" x 72" Sheet, Non-Permeable, 15/cs
 5999L73 36" x 72" Sheet, Non-Permeable w/straps, 15/cs
 5999L74 36" x 72" Sheet, Permeable w/straps quilted, 15/cs
 5999L75 20" x 28" Sheet, Crib Permeable w/straps, 20/cs
 5999LJ1 Large Burn Jacket, 10/cs
 5999LPA Arm Dressing, 20/cs
 5999LPAS Arm Dressing w/shoulder, Ties, 10/cs
 5999LPG Large Padded Hand Dressing, Velcro Closures, 20/cs
 5999LPL Leg Dressing, 20/cs
 5999LPS Adult Scalp Dressing, 20/cs
 5999LV1 Large Burn Vest, 20/cs
 5999M04 4" x 6" Wound Dressing, 100/cs
 5999M06 6" x 9" Wound Dressing, 50/cs
 5999M09 9" x 15" Wound Dressing, 30/cs
 5999M18 15" x 18" Wound Dressing, 30/cs
 5999M24 15" x 24" Wound Dressing, 30/cs
 5999M28 20" x 28" Wound Dressing, 30/cs
 5999M36 24" x 36" Pad, Non-Permeable Heavy Absorbency, 15/cs
 5999M37 24" x 36" Pad, Non-Permeable Medium Absorbency, 24/cs
 5999MJ1 Medium Burn Jacket, 10/cs
 5999MPG Medium Padded Hand Dressing, Velcro Closures, 20/cs
 5999NKB6 Adult Padded Neck Dressing w/bib, 6-Ply, 10/cs
 5999PTM 2" x 3" Slit Tube Dressing, 50/cs
 5999SJ1 Small Burn Jacket, 10/cs
 5999SPG Small (Child) Padded Hand Dressing, Velcro Closures, 20/cs
 5999SV1 Small (Child) Vest, 20/cs

**MOLNLYCKE TENDRA® WOUND DRESSING - MEPITEL®**

Non-adherent to moist wound beds yet adheres ultra-gently to dry skin. Minimizes trauma and pain during dressing changes. Can remain in place for several days, which ensures undisturbed wound healing. Enables less frequent dressing changes. Manages exudate, thereby minimizing the risk of peri-wound maceration. Helps prevent or minimize damage to newly formed delicate tissue. Transparent for easy wound inspection during application and during wear. Conforms well to body contours.

- 282000 Border, Scrum, 7.2" x 7.2", Self-Adherent Soft Silicone Border Foam Dressing, 5/bx, 8 bx/cs
 284000 Border, Sacrum, 9.2" x 9.2", Self-Adherent Soft Silicone Boarder Foam Dressing, 5/bx, 5 bx/cs
 290590 2" x 3" Non-Adherent Silicone Dressing, 10/bx, 6 bx/cs
 290790 3" x 4" Non-Adherent Silicone Dressing, 10/bx, 5 bx/cs
 291090 4" x 8" Non-Adherent Silicone Dressing, 10/bx
 292005 8" x 12" Non-Adherent Silicone Dressing, 5/bx, 6 bx/cs
 295800 Border, Post-Op, 4" x 8", Self-Adherent Soft Silicone Border Foam Dressing, 5/bx, 7 bx/cs
 295900 Border, Post-Op, 4" x 12", Self-Adherent Soft Silicone Border Foam Dressing, 5/bx, 5 bx/cs



360500



265324



270600



310250



250600



293190



294190



284090



285280



370500



288100

MOLNLYCKE TENDRA® WOUND MANAGEMENT - ALLDRESS®

Low adherence to the wound. Helps maintain a moist wound environment. Minimizes the risk of maceration. Waterproof/bacterial barrier. Provides protection from contamination.

- 360500 0.17 oz Tube Hypertonic Gel, 10/bx, 4 bx/cs
 361500 0.50 oz Tube Hypertonic Gel, 10/bx, 4 bx/cs
 265324 4" x 4" Composite Dressing, 10/bx, 10 bx/cs
 265329 Composite Dressing, 4" x 4", 10/bx, 7 bx/cs
 265346 6" x 6" Composite Dressing, 10/bx, 10 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - ALLDRESS® (Continued)

- 265349 Composite Dressing, 6" x 6", 10/bx, 10 bx/cs
 265368 6" x 8" Composite Dressing, 10/bx, 10 bx/cs
 265369 Composite Dressing, 6" x 8", 10/bx, 12 bx/cs
 270600 2.4" x 2.6" Adhesive Polyurethane Film Dressing, 100/bx, 4 bx/cs
 271500 4" x 5" Adhesive Polyurethane Film Dressing, 70/bx, 3 bx/cs
 272500 4" x 10" Adhesive Polyurethane Film Dressing, 20/bx, 5 bx/cs
 273000 6" x 8.5" Adhesive Polyurethane Film Dressing, 10/bx, 8 bx/cs
 310250 1" x 11 yds Fixation Fabric Dressing, 39/cs
 310500 2" x 11 yds Fixation Fabric Dressing, 20/cs
 311000 4" x 11 yds Fixation Fabric Dressing, 20/cs
 311500 6" x 11 yds Fixation Fabric Dressing, 20/cs
 312000 8" x 11 yds Fixation Fabric Dressing, 10/cs
 250600 2" x 2" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
 251100 4" x 4" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
 251500 4" x 8" Absorbent Calcium Alginate Dressing, 10/bx, 5 bx/cs
 253000 12.5" Rope Absorbent Calcium Alginate Dressing, 5/bx, 10 bx/cs
 293190 1.6" x 10" Self-Adherent Silicone Dressing, 5/bx, 5 bx/cs
 293290 2" x 3" Self-Adherent Silicone Dressing, 5/bx
 293490 4" x 7" Self-Adherent Silicone Dressing, 5/bx
 294190 4" x 4" Self-Adherent Absorbent Foam Dressing, 5/bx, 12 bx/cs
 294290 4" x 8" Self-Adherent Absorbent Foam Dressing, 5/bx
 294390 6" x 6" Self-Adherent Absorbent Foam Dressing, 5/bx, 6 bx/cs
 294490 8" x 8" Self-Adherent Absorbent Foam Dressing, 5/bx, 6 bx/cs
 294592 8" x 20" Silicone Soft Transfer Foam Dressing, 2/bx, 6 bx/cs
 294890 6" x 8" Silicone Soft Transfer Foam Dressing, 5/bx, 7 bx/cs
 295200 3" x 3" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 14bx/cs
 295300 4" x 4" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 10 bx/cs
 295400 6" x 6" Self-Adherent Absorbent Foam Dressing w/border, 5/bx
 295600 6" x 8" Self-Adherent Absorbent Foam Dressing w/border, 5/bx, 9 bx/cs
 284090 2.4" x 3.4" Soft Silicone Thin Foam Dressing, 5/bx, 14 bx/cs
 284190 4" x 4" Soft Silicone Thin Foam Dressing, 5/bx, 10 bx/cs
 284390 6" x 6" Soft Silicone Thin Foam Dressing, 5/bx, 10 bx/cs
 284590 8" x 20" Soft Silicone Thin Foam Dressing, 2/bx, 6 bx/cs
 670890 Absorbent Shower Proof Adhesive Dressing, 2.5" x 3", Latex Free, Sterile, 60/bx, 8 bx/cs
 670990 Absorbent Shower Proof Adhesive Dressing, 3.6" x 4", Latex Free, Sterile, 60/bx, 8 bx/cs
 671090 Absorbent Shower Proof Adhesive Dressing, 3.6" x 6", Latex Free, Sterile, 40/bx, 8 bx/cs
 671190 Absorbent Shower Proof Adhesive Dressing, 3.6" x 8", Latex Free, Sterile, 30/bx, 6 bx/cs
 671290 Absorbent Shower Proof Adhesive Dressing, 3.6" x 10", Latex Free, Sterile, 30/bx, 6 bx/cs
 671390 Absorbent Shower Proof Adhesive Dressing, 3.6" x 12", Latex Free, Sterile, 30/bx, 9 bx/cs
 670800 2.5" x 3" Absorbent Island Dressing, 60/bx, 8 bx/cs
 670900 3.6" x 4" Absorbent Island Dressing, 50/bx, 8 bx/cs
 671000 3.6" x 6" Absorbent Island Dressing, 50/bx, 8 bx/cs
 671100 3.6" x 8" Absorbent Island Dressing, 30/bx, 6 bx/cs
 671200 3.6" x 10" Absorbent Island Dressing, 30/bx, 6 bx/cs
 671300 3.6" x 12" Absorbent Island Dressing, 30/bx, 9 bx/cs
 671400 3.6" x 14" Absorbent Island Dressing, 30/bx, 6 bx/cs
 285280 3/4" x 39" Ribbon Impregnated Wound Dressing, 10/bx, 6 bx/cs
 285580 4" x 4" Impregnated Wound Dressing, 30/bx, 8 bx/cs
 285780 6" x 6" Impregnated Wound Dressing, 30/bx, 8 bx/cs
 286080 8" x 8" Impregnated Wound Dressing, 30/bx, 8 bx/cs
 370500 0.17 oz Tube Isotonic Saline Gel, 10/bx, 4 bx/cs
 371500 0.50 oz Tube Isotonic Saline Gel, 10/bx, 4 bx/cs
 288100 5" x 8", 5/bx

Dressings/Gauzes



DUKAL EMERGENCY PRESSURE DRESSING

The DUKAL Emergency Pressure Dressing is an all-in-one dressing. It is specifically designed with a highly absorbent pad that offers maximum compression when held in place with the attached bandage roll.

1050 Emergency Pressure Dressing, 1/pk, 50 pk/cs

EXEL ZORBAND™ PRESSURE BANDAGES

Whenever a large bore needle is used these bandages can help to reduce pressure time to stop bleeding. The absorbent pad soaks up blood. The swelling action of pad acts as compress.

26834 Large Bandage, 100/bx, 10 bx/cs
26836 X-Large Bandage, 100/bx, 10 bx/cs

COVIDIEN/KENDALL CURITY MULTI-TRAUMA DRESSING

Soft, nonwoven facing with highly absorbent filings. Superior padding and protection for major wounds. Available in specialized emergency room put-ups.

1967 Multi-Trauma Dressing, 10" x 30", Sterile, 1/bg, 50 bg/cs



AMD-RITMED ABDOMINAL PADS

A7059 5" x 9" Abdominal Pads, Sterile 1s, Sealed Ends, 20/pk, 20 pk/cs
A7081 8" x 10" Abdominal Pads, Sterile 1s, Sealed Ends, 24/pk, 15 pk/cs
A7087 8" x 7.5" Abdominal Pads, Sterile 1s, Sealed Ends, 12/pk, 20 pk/cs
A7159 5" x 9" Abdominal Pads, Non-Sterile, Sealed Ends, 576/cs
A7181 8" x 10" Abdominal Pads, Non-Sterile, Sealed Ends, 432/cs
A7187 8" x 7.5" Abdominal Pads, Non-Sterile, Sealed Ends, 576/cs

DERMA SCIENCES DUPAD® ABDOMINAL PADS

Filled with cellulose fibers, making them super-absorbent. The fibers quickly absorb and disperse fluids laterally, to prevent pooling and maceration. The backing is constructed of a hydrophobic material, preventing strike-through. All edges are sealed to prevent linting. Latex free.

87059 Abdominal Pad, 5" x 9", Sterile, 25/bx, 8 bx/cs (200 ea/cs)
87088 Abdominal Pad, 8" x 7.5", Sterile, 25/bx, 8 bx/cs (200 ea/cs)
87810 Abdominal Pad, 8" x 10", Sterile, 25/bx, 8 bx/cs (200 ea/cs)



5590



5940

DUKAL ABD PADS

Highly absorbent, large dressing with a moisture resistant barrier provides for superior care of large wounds and prevents exudate from wicking through.

5590 Pad, Sterile, 5" x 9", 1/pk, 25 pk/tray, 16 tray/cs
5810 Pad, Sterile, 8" x 10", 1/pk, 20 pk/tray, 16 tray/cs
5875 Pad, Sterile, 8" x 7 1/2", 1/pk, 20 pk/tray, 12 tray/cs
5940 Pad, Non-Sterile, 5" x 9", 25/bg, 16 bg/cs
5942 Pad, Nonsterile, ABD, 8" x 7.5" 25/bg, 16 bg/cs
5943 Pad, Nonsterile, ABD, 8" x 10", 25/bg, 16 bg/cs
5944 Pad, Nonsterile, ABD, 10" x 12", 25/bg, 10 bg/cs
5945 Pad, Nonsterile, ABD, 12" x 16", 25/bg, 6 bg/cs



DYNAREX COMBINE PADS

3501 Sterile Combine Pad, 5" x 9", 1/pch, 20/tray, 20 trays/cs
3502 Sterile Combine Pad, 8" x 7 1/2", 1/pch, 12/tray, 20 trays/cs
3503 Sterile Combine Pad, 8" x 10", 1/pch, 24/tray, 15 trays/cs
3511 Non-Sterile Combine Pad, 5" x 9", Bulk, 576/cs
3512 Non-Sterile Combine Pad, 8" x 9 1/2", Bulk, 576/cs
3513 Non-Sterile Combine Pad, 8" x 10", Bulk, 432/cs

COVIDIEN/KENDALL CURITY ABDOMINAL PADS

All four edges sealed to prevent lint residue and leaking. Soft nonwoven layer for patient comfort and fluff filler for absorbency.

30066 12" x 16" Abdominal Pads, Non-Sterile, Bulk, 144/cs
68659 24" x 8" Abdominal Pads, Non-Sterile, Bulk, 144/cs
69043 12" x 10" Abdominal Pads, Non-Sterile, Bulk, 288/cs

TIDI COMBINATION PADDING - ABD PADS

Quick wicking, high wet strength padding with non-woven cover. Sterile padding has absorbent lamino filler, non-sterile padding has non-absorbent back. Available in either an all absorbent fill (AAF) or with a non-absorbent backing (NAB) to impede fluid strike through.

979553 8" x 10" Combination Padding, Non-Sterile, NAB, 432/cs
979555 10" x 12" Combination Padding, Non-Sterile, NAB, 288/cs
979556 12" x 16" Combination Padding, Non-Sterile, NAB, 144/cs
979560 Rolls Combination Padding, Non-Sterile, NAB, 8" x 20 yds, 12/cs
979562 10" x 24" Combination Padding, Non-Sterile, NAB, 144/cs
979580 Rolls Combination Padding, Non-Sterile, NAB, 12" x 30 yds

J&J SURGIPAD™ COMBINE DRESSING

Multi-layer abdominal pad that helps provide protection and padding for major surgical wounds, while maximizing absorbency, wicking and dispersal of wound exudate. Sealed on all four sides. Green indicator line on back side for placement away from the wound. Sterile.

005511 5" x 9" Sterile Combine Dressings, 12/bx, 24 bx/cs



COVIDIEN/KENDALL TENDERSORB WET PRUF ABDOMINAL PADS

Consists of three key layers: soft outer nonwoven layer, fluff filler to absorb and disperse fluid and wet-proof moisture barrier to retard fluid strike-through. Unique wet-proof moisture barrier helps minimize exposure to blood-borne pathogens. Absorbs more fluid than ordinary ABD pads for fewer dressing changes. All four edges sealed to prevent lint residue and leaking.

8190A Abdominal Pads, Non-Sterile, Bulk, 5" x 9", 880/cs
8192A Abdominal Pads, Non-Sterile, Bulk, 7 1/2" x 8", 648/cs
8194A Abdominal Pad, Non-Sterile, Bulk, 8" x 10", 432/cs
9190A 5" x 9" Abdominal Pads, Sterile, 1s, 36/tray, 12 tray/cs
9192A 7 1/2" x 8" Abdominal Pads, Sterile, 1s, 18/tray, 12 tray/cs
9194A 8" x 10" Abdominal Pads, Sterile, 1s, 18/tray, 12 tray/cs



DUKAL TRAUMA DRESSING

The multiple-purpose Trauma Dressing offers greater absorbency, protection and padding to the site. Featuring a moisture resistant barrier this pad will wick quickly and prevent exudates from pooling.

1030TD Trauma Dressing, 10" x 30", 1/pk, 25 pk/cs
3050 Trauma Dressing, 12" x 30", 1/pk, 25 pk/cs

AMBRA LE ROY EMERGENCY BANDAGE

Sterile trauma wound dressing combines an elastic bandage, non-adhering gauze pad, pressure applicator and quick closure to control bleeding in emergency situations.

70403-B Emergency Bandage, 4"W, 10/bx
70606-B Emergency Bandage, 6"W, 10/bx
70810-B Emergency Bandage, Abdominal Pad, 8"W, 10/bx

DYNAREX MULTI-TRAUMA DRESSING

Sterile.

3531 10" x 30" Multi-Trauma Dressing, 50/cs
3532 12" x 30" Multi-Trauma Dressing, 50/cs

TIDI TRAUMA DRESSINGS

979564 All Absorbent Trauma Dressing, Sterile 1s, 10" x 30", 50/cs



ALBA STERILE PETROLATUM GAUZE DRESSINGS

Clings and conforms to all body contours without sticking to wound. Provides soothing non-adherent protection with more convenience and less mess. Non-toxic, non-sensitive and non-irritating. An oil-resistant paper wrap folded over the dressing protects it within the peel-open foil-lined pouch and greatly simplifies handling. The dressing is easily taken out of its package without sticking or pad distortion. The folded edges of the inner wrap pull neatly away from the dressing with minimum adherence. The dressings are useful as an initial layer in surgical dressings and on skin graft sites, for soothing treatment of first degree burns and abrasions, and as wound packing.

204 1/2" x 72" Gauze Dressing, 50/bx, 4 bx/cs
205 1" x 36" Gauze Dressing, 50/bx, 4 bx/cs
206 3" x 9" Gauze Dressing, 50/bx, 4 bx/cs
207 3" x 18" Gauze Dressing, 50/bx, 4 bx/cs
208 3" x 36" Gauze Dressing, 50/bx, 4 bx/cs
209 6" x 36" Gauze Dressing, 50/bx, 4 bx/cs
210 1" x 8" Gauze Dressing, 50/bx, 4 bx/cs



ALBA XEROFORM PETROLATUM GAUZE DRESSINGS

Clings and conforms to all body contours without sticking to wound. Provides soothing non-adherent protection with more convenience and less mess. Includes a special deodorizing agent - 3% bismuth tribomophenate. Non-toxic, non-sensitive and non-irritating. An oil-resistant paper wrap folded over the dressing protects it within the peel-open foil-lined pouch and greatly simplifies handling. Whether using forceps or gloved hands, the dressing is easily taken out of its package without sticking or pad distortion. The folded edges of the inner wrap pull neatly away from the dressing with minimum adherence. The dressings are useful as an initial layer in surgical dressings and on skin graft sites, for soothing treatment of first degree burns and abrasions, and as wound packing.

212 1" x 8" Gauze Dressing, 50/bx, 4 bx/cs
213 5" x 9" Gauze Dressing, 50/bx, 4 bx/cs

DERMA SCIENCES PETROLATUM IMPREGNATED GAUZE

20049 Petrolatum Impregnated Gauze, 1/2" x 72", 12/bx, 6 bx/cs
20053 Petrolatum Impregnated Gauze, 1" x 36", 12/bx, 6 bx/cs
20056 Petrolatum Impregnated Gauze, 3" x 9", 12/bx, 6 bx/cs
20059 Petrolatum Impregnated Gauze, 3" x 18", 12/bx, 6 bx/cs
20062 Petrolatum Impregnated Gauze, 3" x 36", 12/bx, 6 bx/cs
20065 Petrolatum Impregnated Gauze, 6" x 36", 12/bx, 6 bx/cs
77030 Petrolatum Wound Closure, 1" x 8", Sterile, 50/bx, 4 bx/cs

SMITH & NEPHEW COLACTIVE AG™ DRESSINGS

A sterile collagen dressing that contains 90% porcine-derived collagen and 10% sodium alginate. When applied to the wound bed, ColActive™ forms a moist wound environment in which fibroblasts and macrophages can easily move and function more efficiently. Absorbs more than thirty times its weight in wound drainage.

66003002 Collagen Dressings w/silver chloride, 2" x 2", Sterile, 10/bx
66003003 Collagen Dressings w/silver chloride, 4" x 4", Sterile, 10/bx

Dressings/Gauzes

DERMA SCIENCES SHUR-CONFORM® OIL EMULSION DRESSING

Non-adherent and sterile. Constructed of knitted cellulose acetate fabric, coated with a formulated petrolatum emulsion. Conforms to body contours while reducing the risk of adhering to wounds and disrupting the healing process. Packaged in clear poly pouches.

- 77041 3" x 3" Oil Emulsion Dressing, 50/bx, 12 bx/cs
- 77042 3" x 8" Oil Emulsion Dressing, 24/bx, 6 bx/cs
- 77043 3" x 16" Oil Emulsion Dressing, 36/bx, 6 bx/cs
- 77044 3" x 8" (3) Oil Emulsion Dressing, 36/bx, 6 bx/cs



DERMA SCIENCES XEROFORM PETROLATUM GAUZE

- 20068 Xeroform Petrolatum Gauze, 5" x 9", Sterile, 50/bx, 4 bx/cs
- 71164 Xeroform Petrolatum Gauze, 2" x 2", Sterile, 25/bx, 6 bx/cs
- 71165 Xeroform Petrolatum Gauze, 4" x 4", Sterile, 25/bx, 6 bx/cs
- 77034 Xeroform Petrolatum Gauze, 1" x 8", Sterile, 50/bx, 4 bx/cs

DERMA SCIENCES XEROFORM/PETROLATUM GAUZE DRESSINGS OVERWRAPS

- 99815 1/2" x 72" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99816 1" x 36" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99817 3" x 9" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99818 3" x 18" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99819 3" x 36" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99820 6" x 36" Petrolatum Gauze Overwrap, Sterile, 12/bx, 6 bx/cs
- 99821 1" x 8" Petrolatum Gauze Overwrap, Sterile, 50/bx, 4 bx/cs
- 99822 5" x 9" Xeroform Petrolatum Gauze Overwrap, Sterile, 25/bx, 6 bx/cs
- 99823 1" x 8" Xeroform Petrolatum Gauze Overwrap, Sterile, 50/bx, 4 bx/cs



20181



20241



20401

SMITH & NEPHEW ACTICOAT™ WITH NANOCRYSTALLINE SILVER

- 20101 4" x 4" Dressing, 12/bx
- 20111 Acticoat™ Moisture Control, 2" x 2" Dressing, 10/bx
- 20141 4" x 5" Dressing, 5/bx
- 20151 Acticoat™, 5" x 5" Dressing, 5/bx, 5 bx/cs
- 20181 3/4" x 12" Rope, 5/bx, 12 bx/cs
- 20201 4" x 8" Dressing, 12/bx
- 20241 6" x 6" Dressing, 5/bx, 4 bx/cs
- 20301 8" x 16" Dressing, 6/bx, 4 bx/cs
- 20311 Acticoat™ Moisture Control, 4" x 4" Dressing, 10/bx
- 20381 4" x 5" Absorbent Dressing, 5/bx
- 20401 16" x 16" Dressing, 6/bx, 2 bx/cs
- 20411 Acticoat™ Moisture Control, 4" x 8" Dressing, 10/bx
- 20501 4" x 48" Roll, 6/bx, 2 bx/cs
- 20601 Acticoat™, 2" x 2" Dressing 5/bx, 20 bx/cs



INTEGRITY OIL EMULSION DRESSINGS

Knitted, high porosity fabric with uniform weave evenly impregnated with U.S.P. White Petrolatum in an emulsion blend. Provides non-adherence without occlusion. An excellent dressing for draining wounds. Stretches and conforms to the wound site. Packaged sterile in a convenient tamperproof, peel-open package.

- 25-0316 Dressing, Sterile, 3" x 16", 36pkg/bx, 6bx/cs
- 25-0330 Dressing, Sterile, 3" x 3", 50pkg/bx, 12bx/cs
- 25-0381 Dressing, Sterile, 3" x 8", 24pkg/bx, 6bx/cs



INTEGRITY PETROLATUM GAUZE

U.S.P. absorbent fine mesh gauze impregnated with U.S.P. White Petrolatum. Occlusive and non-adhering. It helps seal around wound, stoma, and tube sites, and aids in the prevention of air leaks and unwanted fluid loss. Our special manufacturing process results in an even distribution of clean, smooth, white petrolatum. U.S.P. White Petrolatum Gauze is available sterile in a handy, single peel-open package.

- 16-1136 Gauze, Sterile, 1" x 36", 12pkg/bx, 12bx/cs
- 16-1180 Gauze, Sterile, 1" x 8", 12pkg/bx, 12bx/cs
- 16-1272 Gauze, Sterile, .5" x 72", 12pkg/bx, 12bx/cs
- 16-1318 Gauze, Sterile, 3" x 18", 12pkg/bx, 12bx/cs
- 16-1336 Gauze, Sterile, 3" x 36", 12pkg/bx, 12bx/cs
- 16-1636 Gauze, Sterile, 6" x 36", 12pkg/bx, 12bx/cs
- 25-1390 Gauze, Sterile, 3" x 9", 12pkg/bx, 12bx/cs



INTEGRITY XEROFORM PETROLATUM GAUZE

U.S.P. absorbent fine mesh gauze impregnated with U.S.P. White Petrolatum in a blend with 3% Bismuth Tribromophenate - a medicating and deodorizing, occlusive and non-adhering dressing. Helps prevent external contamination. Packaged in a sterile, peel-open package.

- 24-0220 Gauze, Xeroform Petrolatum, Sterile, 2" x 2", 25/bx, 6 bx/cs
- 24-0440 Gauze, Xeroform Petrolatum, Sterile, 4" x 4", 25/bx, 6 bx/cs
- 24-1590 Gauze, Xeroform Petrolatum, Sterile, 5" x 9", 12pkg/bx, 12bx/cs
- 24-3180 Gauze, Xeroform Petrolatum, Sterile, 1" x 8", 12pkg/bx, 12bx/cs

**COVIDIEN/KENDALL VASELINE PETROLATUM GAUZE DRESSING**

Fine mesh, absorbent gauze impregnated with white petrolatum. Non-adherent. Clings and conforms to all body contours. Petrolatum impregnant remains moist, nontoxic, nonsensitizing and nonirritating. Sterile.

8884411600	1/2" x 72" Strip in Tube, 72/cs
8884412600	1" x 36" Strip in Peelable Foil Packs, 72/cs
8884413605	3" x 9" Petrolatum Gauze, 50/bx, 4 bx/cs
8884414600	3" x 18" Strip in Peelable Foil Packs, 72/cs
8884415600	3" x 36" Strip in Peelable Foil Packs, 72/cs
8884416600	6" x 36" Strip in Peelable Foil Packs, 72/cs
8884417601	1" x 8" Cision Dressing Foil Packs, 200/cs
8884421600	Petrolatum Gauze, 1/2" x 72", Tube Foil Overwrap Peelable Foil Packs, 12/bx, 6 bx/cs
8884422600	1" x 36" Strip in Overwrap Peelable Foil Packs, 72/cs
8884423600	3" x 9" Strip in Overwrap Peelable Foil Packs, 72/cs
8884424600	1" x 18" Strip in Overwrap Peelable Foil Packs, 72/cs
8884425600	3" x 36" Strip in Overwrap Peelable Foil Packs, 72/cs
8884426600	6" x 6" Strip in Overwrap Peelable Foil Packs, 72/cs
8884427601	1" x 8" Cision Dressing in Overwrap, 200/cs

**COVIDIEN/KENDALL XEROFORM PETROLATUM GAUZE DRESSING**

3% Bismuth Tribromophenate in a special petrolatum blend on fine mesh gauze. Non-adherent. Clings and conforms to all body contours. Bacteriostatic action. Sterile.

8884431302	1" x 8" Strip in Overwrap Peelable Foil Packs, 200/cs
8884431605	5" x 9" Strip in Overwrap Peelable Foil Packs, 200/cs
8884432000	4" x 3 yds Petrolatum Gauze Dressing, 6/bx, 6 bx/cs
8884433301	Petrolatum Gauze, 1" x 8", Patch Peelable Foil Packs, 50/bx, 4 bx/cs
8884433400	2" x 2" Patch Peelable Foil Packs, 150/cs
8884433500	4" x 4" Patch Peelable Foil Packs, 150/cs
8884433605	5" x 9" Petrolatum Gauze Dressing, 50/bx, 4 bx/cs

MOLNLYCKE TENDRA® WOUND MANAGEMENT - MESALT®

Helps stimulate the cleansing of moist necrosis (slough) in draining and infected wounds. Absorbs exudate, bacteria and moist necrotic material.

285280	3/4" x 39" Ribbon Impregnated Wound Dressing, 10/bx, 6 bx/cs
285580	4" x 4" Impregnated Wound Dressing, 30/bx, 8 bx/cs
285780	6" x 6" Impregnated Wound Dressing, 30/bx, 8 bx/cs
286080	8" x 8" Impregnated Wound Dressing, 30/bx, 8 bx/cs

HARTMANN-CONCO ATRAUMAN® AG SILVER IMPREGNATED DRESSING

499571	2" x 2", Anti-Microbial Dressing, Latex Free, 10/bx
499573	4" x 4", Anti-Microbial Dressing, Latex Free, 10/bx
499575	4" x 8", Anti-Microbial Dressing, Latex Free, 10/bx

SMITH & NEPHEW CUTICERIN® GAUZE

66045502	8" x 16" Gauze, 25/bx, 6bx/cs
66045503	3" x 8" Gauze, 3/pk, 25pk/bx, 6bx/cs
66045560	3" x 3" Gauze, 50/bx, 6bx/cs
66045561	3" x 8" Gauze, 50/bx, 6bx/cs

INTEGRITY PROTODERM™ XEROFORM PETROLATUM PIN SITE DRESSING

Specifically designed to help manage external fixation device pin site wounds. Used as a primary layer to cover and protect wounds from environmental contamination with mild deodorizing. The unique die-cut individual dressings, allows for easy access and precise pin/wire placement. It is the only product packaged 10 ready-to-use dressings in one, sterile pouch resulting in significantly lower cost per dressing and ease of use. They are easily pulled apart from the single packaged strip of 10. Each sterile pouch contains ten 1" x 1" square dressings. Each dressing is slit and slotted to fit easily around external fixation pins and wires.

XSQ-10	Dressing, Sterile, 1" x 1", 10/bx, 5 bx/cs
--------	--

**HARTMANN-CONCO ECONO-PASTE® PLUS CALAMINE CONFORMING ZINC-OXIDE PASTE BANDAGE**

Latex-free, non-sterile, cotton zinc-oxide, impregnated gauze bandage with added soothing benefits of calamine. Non-hardening emulsion. Bandage stretches and conforms easily to body contours and has finished selvage edges. Individually rolled on plastic core and vacuum-sealed in foil to guarantee a long shelf life.

47310000	Zinc-Oxide Paste w/calamine on a flexible gauze bandage, 3" x 10 yds, 12rl/cs
47410000	Zinc-Oxide Paste w/calamine on a flexible gauze bandage, 4" x 10 yds, 12rl/cs

**HARTMANN-CONCO ECONO-PASTE® CONFORMING ZINC-OXIDE PASTE BANDAGE**

Latex free, non-sterile cotton zinc-oxide impregnated gauze bandage with finished selvage edges. Soothing, non-hardening, zinc-oxide-based emulsion. Stretches and conforms easily to body contours. Individually rolled on plastic cores and vacuum-sealed in foil to guarantee a long shelf life.

47300000	3" x 10 yds Zinc-Oxide Paste Bandage, 12/cs
47400000	4" x 10 yds Zinc-Oxide Paste Bandage, 12/cs

Dressings/Gauzes



GRAHAM FIELD MEDICOPASTE® BANDAGE

Zinc impregnated premium gauze bandage. Aids in healing of venous stasis ulcers, reduces and controls edema, and improves blood flow by providing support for subcutaneous tissues. Also acts as a "flexible cast" in the treatment of sprains, strains, and dislocations. Works as a protective covering for skin grafts of the lower extremities. Manufactured at our Graham-Field facility since 1947. The brand medical professionals rely on for superior performance and quality.

1565-3	3" x 10 yds Medicopaste Bandage, 12/bx
1565-4	4" x 10 yds Medicopaste Bandage, 12/bx
1565C-3	3" x 10 yds Medicopaste Bandage w/calamine, 12/bx
1565C-4	4" x 10 yds Medicopaste Bandage w/calamine, 12/bx

COVIDIEN/KENDALL TENDERWRAP UNNA BOOT BANDAGE

Formulated with nonhardening zinc oxide paste which allows the bandage to remain flexible encouraging patient ambulation and compliance. Gelatin free, eliminating the need for preservatives which in turn reduces the chance of adverse skin reactions. Available with the soothing relief of Calamine. Nonsterile. Soft, nonraveling, 100% cotton gauze. Conforms well to leg contours for a comfortable fit. Helps deliver firm, even compression therapy for a multitude of bandaging needs.

8033	Unna Boot Bandage, 3" x 10 yds, 1 r/bx, 12 bx/cs
8034	Unna Boot Bandage, 4" x 10 yds, 1 r/bx, 12 bx/cs
8035	Unna Boot Bandage, 3" x 10 yds w/calamine, 1 r/bx, 12 bx/cs
8036	Unna Boot Bandage, 4" x 10 yds w/calamine, 1 r/bx, 12 bx/cs

DERMA SCIENCES PRIMER FLEX™ UNNA BOOT

A gauze bandage that is stretchable and non-rigid. The bandage is evenly impregnated with a non-hardening paste of zinc oxide. Ideal for treating: stasis leg ulcers, lymphatic edema, thrombophlebitis, eczema of the leg, minor fractures, sports injuries, sprains and strains. Contains no gelatin, which prevents the hardening or caking of the bandage. Made without preservatives, no expiration date. Remains flexible. Packaged in a unique vacuum-sealed pouch.

GL-300-1F	3" x 10 yds, 12/cs
GL-400-1F	4" x 10 yds, 12/cs



300-1



300-1C

DERMA SCIENCES PRIMER® UNNA BOOT

A gauze bandage that has been evenly impregnated with a non-hardening paste of zinc-oxide. Available with or without calamine. Ideal for treating: venous stasis leg ulcers, lymphatic edema, minor sprains, strains, fractures, thrombophlebitis and skin grafts on burn wounds. Contains no gelatin which prevents the hardening or caking of the bandage. Made without preservatives. Made of non-raveling threadlocked gauze. Packaged in a unique vacuum-sealed pouch.

300-1	3" x 10 yd Boot, 12/cs
300-1C	3" x 10 yd Boot w/calamine, 12/cs
400-1	4" x 10 yd Boot, 12/cs
400-1C	4" x 10 yd Boot w/calamine, 12/cs



DERMA SCIENCES UNNA-PAK® TWO-LAYER SYSTEM

A convenient system containing two high performance bandages that together form a two-layer Unna Boot. All of the advantages and uses described for the Primer® Unna Boot with the additional benefit of a flesh-colored secondary dressing that is lightweight, strong and easy to apply. System available in two bandage widths, 3" and 4".

GL200-3	3" x 10 yds Unna Boot w/3" x 5 yds elastic bandage
GL200-4	4" x 10 yds Unna Boot w/4" x 5 yds elastic bandage



DYNAREX UNNA BOOT BANDAGE

3453	3" x 10 yds Unna Boot Bandage, 12/cs
3454	4" x 10 yds Unna Boot Bandage, 12/cs
3455	3" x 10 yds Unna Boot Bandage w/calamine, 12/cs
3456	4" x 10 yds Unna Boot Bandage w/calamine, 12/cs



SMITH & NEPHEW VISCOPASTE™ PB7 ZINC PASTE BANDAGE

Used in the application of unna boots for the management of venous leg ulcers. Also indicated for the management of skin tears and stasis eczema. Soothes dry skin and does not contain calamine. Made from a light open-weave fabric, is simple to apply and easy to remove.

4956	3" x 10 yds Zinc Paste Bandage, 12/pkg, 4 pkg/cs
------	--



ALBA PACKING STRIPS

Woven gauze fabric resists linting. Used for packing or as drainage conduits in nasal or sinus packing. The addition of iodoforn provides an odor-resisting packing strip ideal for draining wounds. Sterile in unbreakable poly bottles with barrier caps and easily identified tamper evident seal.

260	Plain Packing Strips, 1/4" x 5 yds, 1/btl, 12 btl/cs
261	Plain Packing Strips, 1/2" x 5 yds Packing Strips, 1/btl, 12 btl/cs
262	Plain Packing Strips, 1" x 5 yds, 1/btl, 12 btl/cs
263	Plain Packing Strips, 2" x 5 yds, 1/btl, 12 btl/bx
270	Iodoform 5% Packing Strips, 1/4" x 5 yds, 1/btl, 12 btl/cs
271	Iodoform 5% Packing Strips, 1/2" x 5 yds, 1/btl, 12 btl/cs
272	Iodoform 5% Packing Strips, 1" x 5 yds, 1/btl, 12 btl/bx
273	Iodoform 5% Packing Strips, 2" x 5 yds, 1/btl, 12 btl/bx

**COVIDIEN/KENDALL CURITY PACKING STRIPS**

All natural 100% cotton. Fine mesh gauze ideal for wet-to-dry packing. Available in plain or iodoform varieties. Guaranteed sterility.

7631	Packing Strips, Plain, 1/4 x 5 yds, 12/cs
7632	Packing Strips, Plain, 1/2 x 5 yds, 12/cs
7633	Packing Strips, Plain, 1 x 5 yds, 12/cs
7634	Packing Strips, Plain, 2 x 5 yds, 12/cs
7831	Packing Strips, Iodoform, 1/4 x 5 yds, 12/cs
7832	CURITY Packing Strips, Iodoform, 1/2 x 5 yds, 12/cs
7833	Packing Strips, Iodoform, 1" x 5 yds, 12/cs
7834	Packing Strips, Iodoform, 2 x 5 yds, 12/cs

DERMA SCIENCES PAK-ITS® GAUZE PACKING STRIPS - PLAIN & IODOFORM

Can be used dry to absorb drainage in exudating wounds or soaked with a medication of choice. Pre-cut sizes suit different wound types. Saves time and money. A secondary dressing is needed to hold the strips in place and to protect the wound from outside contaminants. 5 yards per bottle and available in multiple widths. Latex free.

59121	Plain Gauze Packing Strip, 1/4" x 5 yds, 12/cs
59146	Iodoform Packing Strip, 1/4" x 5 yds, 12/cs
59221	Plain Gauze Packing Strip, 1/2" x 5 yds, 12/cs

DERMA SCIENCES PAK-ITS® GAUZE PACKING STRIPS - PLAIN & IODOFORM W/STRIP SYSTEM

Can be used dry to absorb drainage in exudating wounds or soaked with a medication of choice. Pre-cut sizes suit different wound types. Saves time and money. A secondary dressing is needed to hold the strips in place and to protect the wound from outside contaminants. 5 yards per bottle and available in multiple widths. The bottle has special patent-pending strip delivery system that ensures packing material is kept clean, tangle free, with the end always easy to locate at the top of the bottle. Latex free.

59120	Plain Gauze Packing Strip, 1/4" x 5 yds, 12/cs
59145	Iodoform Packing Strip, 1/4" x 5 yds, 12/cs
59220	Plain Packing Strip, 1/2" x 5 yds, 12/cs
59245	Iodoform Packing Strip, 1/2" x 5 yds, 12/cs
59320	Plain Packing Strip, 1" x 5 yds, 12/cs
59345	Iodoform Packing Strip, 1" x 5 yds, 12/cs
59420	Plain Packing Strip, 2" x 5 yds, 12/cs
59445	Iodoform Packing Strip, 2" x 5 yds, 12/cs

DUKAL SURGICAL GAUZE & PACKING

100% cotton woven gauze that are designed to be of the highest quality for use in the OR.

1320	Section Sponges, Vaginal Packing, 4" x 36", Sterile 1s, 8-Ply, Prewashed, X-Ray Detectable, 2 x1 Mesh, 100/cs
1336X	Section Sponges, Vaginal Packing, 4" x 36", Non-Sterile, 8-Ply, X-Ray Detectable, 2 x1 Mesh, 200/cs
1339S	Section Sponges, Vaginal Packing, 2" x 36", Sterile 1s, 4-Ply, X-Ray Detectable, 28 x 24 Mesh, Edges Folded In, 100/cs
9902	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 1/2" x 100 yds, 4-Ply, 10/cs
9902X	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 1/2" x 100 yds, 4-Ply, X-Ray Detectable, 10/cs
9903	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 1" x 100 yds, 4-Ply, 10/cs
9903X	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 1" x 100 yds, 4-Ply, X-Ray Detectable, 10/cs
9904	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 2" x 100 yds, 4-Ply, 10/cs
9904X	Gauze Packing, Non-Sterile, 28 x 24 Mesh, 2" x 100 yds, 4-Ply, X-Ray Detectable, 10/cs

**DYNAREX PACKING STRIPS**

Sterile.

3401	1/4" x 5 yds Plain Packing Strips, 1 strip/btl, 12 btl/cs
3402	1/2" x 5 yds Plain Packing Strips, 1 strip/btl, 12 btl/cs
3403	1" x 5 yds Plain Packing Strips, 1 strip/btl, 12 btl/cs
3404	2" x 5 yds Plain Packing Strips, 1 strip/btl, 12 btl/cs
3411	1/4" x 5 yds Iodoform Packing Strips, 1 strip/btl, 12 btl/cs
3412	1/2" x 5 yds Iodoform Packing Strips, 1 strip/btl, 12 btl/cs
3413	1" x 5 yds Iodoform Packing Strips, 1 strip/btl, 12 btl/cs
3414	2" x 5 yds Iodoform Packing Strips, 1 strip/btl, 12 btl/cs

**INTEGRITY PACKING STRIPS**

Packing Strips are made of a woven gauze fabric which resists linting. They are used for packing or as drainage conduits in wound care management applications, such as nasal or sinus packing. The addition of Iodoform provides an odor-resisting packing strip that is ideal for draining wounds. Packaged in unbreakable plastic bottles with vapor barrier caps and easily identified tamper evident container.

60-1005	Packing Strips, 1"W x 5yds, 12bottle/cs
60-2005	Packing Strips, 2"W x 5yds, 12bottle/cs
60-2505	Packing Strips, 1/4"W x 5yds, 12bottle/cs
60-5005	Packing Strips, 1/2"W x 5yds, 12bottle/cs
61-1005	Packing Strips, Iodoform, 1"W x 5yds, 12bottle/cs
61-2005	Packing Strips, Iodoform, 2"W x 5yds, 12bottle/cs
61-2505	Packing Strips, Iodoform, 1/4"W x 5yds, 12bottle/cs
61-5005	Packing Strips, Iodoform, 1/2"W x 5yds, 12bottle/cs

**AMD-RITMED LAPAROTOMY SPONGES**

Premium quality mesh gauze. X-ray detectable with colored tape for better visibility. Cotton loop.

1505	12" x 12" Laparotomy Sponges, 4-Ply, Prewashed, 5/pk, 40 pk/cs
1515	18" x 18" Laparotomy Sponges, 4-Ply, Prewashed, 5/pk, 40 pk/cs
1515T	Sponge, Laparotomy, 18" x 18", 4-Ply, Prewashed, 5/plastic tub, 40 tub/cs
1525	4" x 18" Laparotomy Sponges, 4-Ply, Prewashed, 5/pk, 40 pk/cs

Dressings/Gauzes

COVIDIEN/KENDALL CURITY LAPAROTOMY SPONGES

Designed for maximum absorbency. Sterile sponges are pre-washed 100% cotton and disposable. X-ray detectable, multiple sizes; banded in 5's.

6022	Lap Sponges, 18" x 18", Sterile, Tray Pack, 100 trays/cs
6034	Lap Sponges, 12" x 12", Sterile, Tray Pack, 100 trays/cs
6045	Lap Sponges, 4" x 18", Sterile, Tray Pack, 100 trays/cs
6056	Lap Sponges, 8" x 36", Sterile, Tray Pack, 100 trays/cs
6522	Lap Sponges, 18" x 18", Sterile, Soft Pouch Pack, 200/cs
6534	Lap Sponges, 12" x 12", Sterile, Soft Pouch Pack, 200/cs
6545	Lap Sponges, 4" x 17", Sterile, Soft Pouch Pack, 200/cs
6556	Lap Sponges, 8" x 36", Sterile, Soft Pouch Pack, 200/cs
6822	Lap Sponges, 18" x 18", Non-Sterile, Bulk Pack, 400/cs
6834	Lap Sponges, 12" x 12", Non-Sterile, Bulk Pack, 750/cs

COVIDIEN/KENDALL LAPAROTOMY SPONGES

Designed for maximum absorbency. Sterile, pre-washed 100% cotton, disposable. X-ray detectable. Kerlix sponges have dual-colored radiopaque handles that remain visible when saturated with blood. Multiple sizes; banded in 5's.

4031	Lap Sponges, 17" x 26", Sterile, Tray Pack, 100 trays/cs
4034	Lap Sponges, 16" x 17", Sterile, Tray Pack, 100 trays/cs
4037	Lap Sponges, 8" x 32", Sterile, Tray Pack, 100 trays/cs
4048	Lap Sponges, 3" x 17", Sterile, Tray Pack, 100 trays/cs
4069	Lap Sponges, 8" x 11", Sterile, Tray Pack, 100 trays/cs

DUKAL LAPAROTOMY SPONGES

The 100% cotton lap sponges are cross stitched with super locked edges to minimize linting and loose threads. Manufactured to exceed expectations, the high quality cotton and tight weave ensures that sponges offer superior absorption. Disposable and pre-washed.

10-0004	Laparotomy Sponge, 4" x 18", Sterile, X-Ray Detectable, Prewashed, Softpack, 5/pk, 20 pk/cs
10-0012	Laparotomy Sponge, 12" x 12", Sterile, X-Ray Detectable, Prewashed, Softpack, 5/pk, 20 pk/cs
10-0018	Laparotomy Sponge, 18" x 18", Sterile, X-Ray Detectable, Prewashed, Softpack, 5/pk, 20 pk/cs
10-0036	Laparotomy Sponge, 8" x 36", Sterile, X-Ray Detectable, Prewashed, Softpack, 5/pk, 20 pk/cs
10-1836-3	Laparotomy Sponge, 18" x 36", Sterile, X-Ray Detectable, Prewashed, Softpack, 3/pk, 20 pk/cs
10-7004	Laparotomy Sponge, 4" x 18", Non-Sterile, X-Ray Detectable, Prewashed, 5/pk, 40 pk/cs
10-7018	Laparotomy Sponge, 18" x 18", Non-Sterile, X-Ray Detectable, Prewashed, 5/pk, 40 pk/cs
10-7036	Laparotomy Sponge, 8" x 36", Non-Sterile, X-Ray Detectable, Prewashed, 5/pk, 40 pk/cs
99-0004	Laparotomy Sponge, 4" x 18", Sterile, Prewashed, X-Ray Detectable, Softpack, 5/pk, 40 pk/cs
99-0012	Laparotomy Sponge, 12" x 12", Sterile, Prewashed, X-Ray Detectable, Softpack, 5/pk, 40 pk/cs
99-0018	Laparotomy Sponge, 18" x 18", Sterile, Prewashed, X-Ray Detectable, Softpack, 5/pk, 40 pk/cs
99-0036	Laparotomy Sponge, 8" x 36", Sterile, Prewashed, X-Ray Detectable, Softpack, 5/pk, 40 pk/cs
99-1220	Laparotomy Sponge, 12" x 12", Sterile, Prewashed, X-Ray Detectable, w/plastic rings, Softpack, 5/pk, 40 pk/cs
99-1221	Laparotomy Sponge, 18" x 18", Sterile, Prewashed, X-Ray Detectable w/plastic ring, Softpack, 5/pk, 40 pk/cs
99-9012	Laparotomy Sponge, 12" x 12", Sterile, X-Ray Detectable, No Loop, Prewashed, 5/pk, 40 pk/cs
99-9018	Laparotomy Sponge, 18" x 18", Sterile, X-Ray Detectable, No Loop, Prewashed, Softpack, 5/pk, 40 pk/cs
KP-7012	Laparotomy Sponge, 12" x 12", Non-Sterile, Banded in 5s - No Wrap, Prewashed, 5/bdl, 280 bdl/cs
KP-7018	Laparotomy Sponge, 18" x 18", Non-Sterile, Banded in 5s - No Wrap, Prewashed, 5/bdl, 160 bdl/cs

AMD-RITMED DERMAFILL BIOCELLULOSE DRESSING

A naturally occurring biosynthetic skin substitute. 100% natural cellulose that stimulates the growth of skin cells. Your answer to the closure of many chronic vascular, venous and pressure ulcers that never showed improvement. Especially formulated for non-complex battlefield wounds and long-term care patients. Amazing results on minor skin abrasions and 2nd degree burns. Eases pain and dramatically improved healing like no other product on the market.

DI-10-DIS	Dressing, 3" x 4", Sterile 1's, 12 sheet/bx
DI-16-DIS	Dressing, 4" x 6", Sterile 1's, 12 sheet/bx
DI-21-DIS	Dressing, 6" x 8", Sterile 1's, 12 sheet/bx



BIRCHWOOD A.R.D® ANOPERINEAL DRESSINGS

Unique butterfly shape keeps dressing in place without tape. Soft, smooth, lint-free material is highly absorbent to wound drainage and perspiration usually associated with anorectal surgery, pruritus and other perianal irritations. Excellent for mild anal incontinence. Also maybe used to apply topical medications. Designed for quick, easy dispensing. Packaging: 24/box, 12 boxes/case.

14-7031-224 A.R.D® Anoperineal Dressings, 24/bx, 12 bx/cs



BIRCHWOOD FULLER SHIELD®

Designed to hold anal, perianal and sacral dressings comfortably in place without tape, and provide protection against soiling of clothing or bed linens. Superior anatomical design and soft, pliable fabric provides a lightweight garment with minimal bulk. Laundry safe. Open front design facilitates male urination and front Velcro fasteners simplify dressing changes on bed-ridden patients. Dressings can be changed without completely removing shield. Latex free. Single unisex size 24-48.

14-7012-1 Fuller Shield® Perineal Binder, Universal Size 24-48, 12/cs

ETHICON SURGICEL NU-KNIT

Absorbable hemostat. Sterile.

1940	1" x 1" Surgicel Nu-Knit, 24/cs
1941	1" x 3 1/2" Surgicel Nu-Knit, 10/cs
1943	3" x 4" Surgicel Nu-Knit, 24/cs
1946	6" x 9" Surgicel Nu-Knit, 10/cs

ETHICON ABSORBABLE SURGICEL

Sterile.

1951	2" x 14" Absorbable Surgicel, 24/cs
1952	4" x 8" Absorbable Surgicel, 24/cs
1953	2" x 3" Absorbable Surgicel, 24/cs
1955	1/2" x 2" Absorbable Surgicel, 24/cs

ETHICON SURGICEL FIBRILLAR

Absorbable hemostat. Sterile.

1961	1" x 2" Surgicel Fibrillar, 10/cs
1962	2" x 4" Surgicel Fibrillar, 10/cs
1963	4" x 4" Surgicel Fibrillar, 10/cs

J&J PROMOGRAN™ MATRIX WOUND DRESSING

A unique advanced wound care dressing comprised of a sterile, freeze-dried matrix composite of 45% ORC and 55% collagen. Indicated for multiple types of chronic wounds including but not limited to, diabetic, venous and pressure ulcers.

PG004	Dressing, 4.34" Hexagon Sheet, 10/bx, 4 bx/cs
PG019	Dressing, 19.1" Hexagon Sheet, 10/bx, 4 bx/cs

J&J PROMOGRAN™ PRISMA MATRIX WOUND DRESSING

Designed to provide protection and growth by addressing the continuous changes in the wound microenvironment. The unique combination of 55% Collagen, 44% ORC and 1% silver-ORC absorbs destructive components of chronic wound fluid and creates the optimal environment for cellular growth. PRISMA Matrix is a simple-to-use primary dressing and should be considered early as part of good ulcer care.

MA028	Dressing, 4.34 sq in., 10/bx, 4 bx/cs
MA123	Dressing, 19.1 sq in., 10/bx, 4 bx/cs



AMD-RITMED TRACH & IV NON-WOVEN DRESSING SPONGES

A9210	4" x 4" Drain Sponge, 6-Ply, Sterile 2s, 50/pk, 12 pk/cs
A9210-NP	4" x 4" No-Split Sponge, 6-Ply, Sterile 2s, 50/pk, 12 pk/cs
A9220	2" x 2" IV Sponge, 6-Ply, Sterile 2s, 70/pk, 20 pk/cs
A9220-NP	2" x 2" No-Split Sponge, 6-Ply, Sterile 2s, 70/pk, 20 pk/cs



BUSSE TRACHEOSTOMY DRESSING

Sterile. Precut gauze of 4" x 4".

707	Sterile Tracheostomy Dressing, 50/bx, 6 bx/cs
-----	---

DUKAL I.V. & T-DRAIN SPONGES

These rayon/poly blend sponges are pre-cut to assist in providing a snug fit around catheters, tubes, I.V.'s, and tracheotomies. They assist in providing a lean environment at the site that is soft and gentle against the skin. As the precut slit wicks moisture away from the site. T-Drain sponges are also available in woven gauze.

7265	IV Sponge, 2" x 2" Sterile, 6-Ply, 2/pk, 35 pk/bx, 20 bx/cs
7465	Gauze Sponge, 4" x 4", Sterile, 6-Ply, Trach Drain, 2/pk, 25 pk/bx, 12 bx/cs
7468	T-Drain Sponge, 4" x 4" Sterile, 16 Ply, 1/pk, 25 pk/tray, 12 tray/cs



DYNAREX DRAIN SPONGES

Sterile.

4610	2" x 2" I.V. Sponge, 6 Ply, 2/pch, 35 pch/tray, 20 tray/cs
4611	4" x 4" Drain Sponge, 6 Ply, 2/pch, 25 pch/tray, 12 tray/cs



DYNAREX TRACHEOSTOMY SPONGES

Sterile.

4612	4" x 4" Tracheostomy Sponge, 12 Ply, 50/tray, 6 trays/cs
------	--



COVIDIEN/KENDALL EXCILON DRAIN & I.V. SPONGES

Superior absorbency and fast wicking action help reduce infections and decrease the frequency of dressing changes. Unique pre-cut T-split conforms snugly around drains, catheters, chest tubes, I.V.'s and tracheotomies. Sterile. No loose threads to fall into wound. Eliminates need for sterile scissors, conserves nursing time.

7086	4" x 4" Drain Sponge, Sterile 2s, 6-Ply, 50/tray, 12 tray/cs
7087	I.V. Sponge, 2" x 2", 6 Ply, 70/tray, 20 tray/cs



COVIDIEN/KENDALL SORB-IT DRAIN & I.V. SPONGES

Lower cost alternative to premium EXCILON Drain and I.V. Sponges. Pre-cut notch conforms well around drains, catheters, chest tubes, I.V.'s and tracheotomies. Sterile. Virtually lint-free. Eliminates need for sterile scissors, conserves nursing time.

6241	2" x 2" I.V. Sponge, Sterile 2s, 6 Ply, 70/tray, 20 tray/cs
6242P	4" x 4" Drain Sponge, Sterile 2s, 6-Ply, 50/tray, 12 tray/cs

DERMA SCIENCES TRACHEOTOMY/DRAIN SPONGES

An added fenestration makes these sponges ideal for dressing around a tracheotomy tube or any type of tube on the body. Available in a soft 6-ply non-woven or in standard 12-ply gauze. The non-woven provides superior absorption and are low-linting, fast-wicking, and non-adherent to wound beds. The woven version is for when natural breathability is desired. Latex free.

84916	Sponge, Tracheotomy Drainage, Sterile, 4" x 4", 6-Ply, 2's, 25 pk/bx, 12 bx/cs (600 ea/cs)
84918	Sponge, Tracheotomy Drainage, Pediatric, Sterile, 2" x 2", 6-Ply, 2's, 35 pk/bx, 20 bx/cs (1400 ea/cs)



DYNAREX TRACHEOSTOMY CARE KIT

Sterile. Contains: 3-compartment plastic tray, (1) pair sterile gloves, (1) trachea tube brush, (1) waterproof drape, (4) 4"x4" gauze sponges, (2) pipe cleaners, (2) 6" cotton tip applicators, (1) pre-cut gauze tracheostomy dressing, (1) 30" twill tape.

4601	Tracheostomy Care Kit, 20/cs
------	------------------------------

Dressings/Gauzes

AMD-RITMED FLUFF SPONGES

A highly absorbent primary dressing with softness and bulkiness that surrounds, protects and cushions the wound site. Krimped weave construction with a diamond fold shape that offers cushioning and additional coverage of the wound. Sterile - available in formed plastic trays or pouches. Made of the highest quality fluff dried absorbent cotton.

A66750	6" x 6.75", Sterile, 10/hard tray, 48 tray/cs
A66752	6" x 6.75", Sterile, 2/pouch, 20 pouches/bx, 12bx/cs
A66755	6" x 6.75", Sterile, 5/hard tray, 120 tray/cs

DUKAL SUPER SPONGES

This multi-layer 100% fluff cotton sponge is highly absorbent and durable which allows it to be used in a number of applications. The thick, diagonally folded construction adds loft, which enhances the overall performance.

2334	Sponge, Diagonal Measurement: 6" x 6.75" (Medium), Non-Sterile, 100/bg, 6 bg/cs
2463	Sponge, Sterile 2s, Medium, Soft Pouch, 2/pk, 20 pk/tray, 12 tray/cs
2487	Sponge, Sterile 5s, Medium, Hard Tray, 5/tray, 120 tray/cs
2495	Sponge, Sterile 10s, Medium, 10/pk, 48 pk/cs



NUTRAMAX GAUZE FLUFFS

Used when greater loft and absorbency is needed in surgical prepping, padding, wound cleansing, dressing or to apply pressure. Widely used as postoperative dressing in mastectomy procedures.

11124	Large Gauze Fluff, Non-Sterile, 30" x 18", 500/cs
11126	Super Gauze Fluff, Non-Sterile, 36" x 36", 500/cs
11127	Super Gauze Fluff, Sterile, 36" x 36", 1 ply, 250/cs
11130	Super/Super Gauze Fluff, Sterile, 36" x 36", 2 ply, 100/cs

TIDI POST-OP SPONGES

Multi-layered for excellent absorption, dispersion and capacity. The center layer of absorbent cotton holds many times its weight in body fluids. Plus, each sponge contains a full-size inner sheet of natural cellulose. The outer layer of soft, absorbent cotton gauze has deep infolded edges to protect against loose threads.

908242	3" x 3" Post-Op Sponge, Non-Sterile, 100/bg, 4000/cs
908244	4" x 4" Post-Op Sponge, Non-Sterile, 100/bg, 2000/cs
908253	4" x 3" Post-Op Sponge, Sterile 2s, 25 pkg/tray, 1200/cs
908254	4" x 4" Post-Op Sponge, Sterile 2s, 25 pkg/tray, 1200/cs



5340



5360

DUKAL TOP SPONGES - NON-WOVEN FACING COVER SPONGE

The Top Sponge is an excellent example of DUKAL quality. This high performance multi-layer sponge features a soft rayon/poly outer facing. It may add protection to the site while assisting in the healing process by offering fast wicking and enhanced absorbency.

5340	Top Sponge, 4" x 4", Non-Woven, Non-Sterile, 100/bg, 20 bg/cs
5350	Top Sponge, 4" x 3", Non-Woven, Sterile, 2/pk, 25 pk/bx, 24 bx/cs
5360	Top Sponge, 4" x 4", Non-Woven, Sterile, 2/pk, 25 pk/bx, 24 bx/cs



COVIDIEN/KENDALL KERLIX SUPER SPONGES

Prewashed, fluff-dried woven gauze, with crinkle-weave pattern and diamond-fold construction offers superior loft, bulk and high absorbency. Superior primary dressing designed to transfer exudate away from wound into secondary dressing. Absorbency, strength and virtual lint-free performance ideal for pre-op prepping and cleansing, debriding and packing. Measured diagonally, mediums are 6" x 6 3/4"; large are 7 3/4" x 8 3/4"; x-large are 9 3/4" x 10 3/4".

1167	Kerlix Super Sponge, Med, 6" x 6 3/4", Bulk, 600/cs
1272	Kerlix Super Sponge, Large, 7 3/4" x 8 3/4", Bulk, 400/cs
2585	Kerlix Super Sponge, Med, 6" x 6 3/4", Sterile 2's, 480/cs
3085	Super Sponges, Sterile 5's in Soft Pouch Package, Medium 50/tray, 12 tray/cs
4308	Kerlix Super Sponge, Lg, 7 3/4" x 8 3/4", 400/cs
6035	Kerlix Super Sponge, X-Lg, 9 3/4" x 10 3/4", Bulk, 200/cs
7310	Kerlix Super Sponge, Medium, 6" x 6 3/4", Sterile 10's, 480/cs



DYNAREX KRINKLE SUPER SPONGES

Sterile.

3163	6" x 6 3/4" Sponge, 2/pk, 20 pk/tray, 12 trays/cs
------	---

J&J TOPPER™ DRESSING SPONGES

Fabric facing permits drainage to disperse into, rather than across pad. Absorbent.

005501	4" x 4" Sterile Dressing Sponges, 24/bx, 24 bx/cs
--------	---

HARTMANN-CONCO STERILUX® SUPER SPONGES

56700000	6" x 6.75", Non-Sterile Sponge, 100/slv, 6 slv/cs
56710000	7.75" x 8.75", Non-Sterile Sponge, 100/slv, 4 slv/cs
56720000	6" x 6.75", Sterile Sponge, 2/pk, 80/bx, 6 bx/cs
56730000	7.75" x 8.75", Sterile Sponge, 2/pk, 40/bx, 10 bx/cs

**DYNAREX POST-OP SPONGES**

3346	4" x 3" Post-Op Sponge, 8 Ply, 2/pch, 25 pch/tray, 24 trays/cs
3347	4" x 4" Post-Op Sponge, 8 Ply, 2/pch, 25 pch/tray, 24 trays/cs

DUKAL CHERRY SPONGES

100% cotton woven gauze that are X-ray detectable and designed to be of the highest quality for use in the O.R.

63401-P	Cherry Sponge, C-5 Holder, Sterile 5s, 5/pk, 40 pk/cs
---------	---

DUKAL KITTNER SPONGES

100% woven gauze that are X-ray detectable for use in the O.R. Packaging ensures accuracy and traceability of count.

60401-P	Kittner Dissectors, 0.25" x 0.56", C-5 Foam Holder, Sterile 5s, 5/pk, 40 pk/cs
60404	Kittner Dissectors, Sterile 5s, 0.25" x 0.56", No Holder, 5/pk, 40 pk/cs

DUKAL O.R. SPECIALTY SPONGES

These high quality sponges are used in applications ranging from packing and prepping to light sponging and dissecting.

76408	Flatstick Sponges, X-Ray Detectable, Non-Sterile, 1 1/2", 100/pk, 5 pk/cs
76410	Flatstick Sponges, X-Ray Detectable, Non-Sterile, 2" Sq, 100/pk, 5 pk/cs

DUKAL PEANUT SPONGES

Made of 100% cotton and are X-ray detectable. Available in a variety of packaging choices which ensure accurate counts. Perfect for minor absorption during surgery.

61401-P	Peanut Sponges, C-5 Holder, Sterile 5s, 5/pk, 40 pk/cs
61402	Peanut Sponges, Pinned, Sterile 5s, 5/pk, 40 pk/cs
61404	Peanut Sponges, Sterile 5s, No Holder, 5/pk, 40 pk/cs
61405-5	Peanut Sponges, Non-Sterile in C-5 Holders, 5/pk, 50 pk/cs

DUKAL POST-OP SPONGES

This gauze facing multi-layer sponge can be used in place of a number of traditional sponges by offering faster wicking and more absorbency than a single sponge.

5434	Sponge, 4" x 3" Sterile Post-Op, Gauze Facing, 2/pk, 25 pk/bx, 24 bx/cs
5444	Sponge, 4" x 4" Sterile Post-Op, Gauze Facing, 2/pk, 25 pk/bx, 24 bx/cs

DUKAL PREP SPONGES

100% cotton woven gauze and are designed to be of the highest quality for use in the O.R.

66406	Prep Sponges, Non-Sterile, 1 1/4", 100/pk, 5 pk/cs
66407	Prep Sponges, Non-Sterile, 1 1/2", 100/pk, 5 pk/cs

DUKAL ROUND STICK SPONGES

100% cotton woven gauze and are designed to be of the highest quality for use in the O.R.

75405	Roundstick Sponges, Cotton Filled, X-Ray Detectable, Large 1", Non-Sterile, 100/pk, 5 pk/cs
75406	Roundstick Sponges, Cotton-Filled, X-Ray Detectable, Super 1 1/2", Non-Sterile, 100/pk, 5 pk/cs
75426	Roundstick Sponges, Cotton-Filled, X-Ray Detectable, Super 1 1/2, Sterile 10s, 10/pk, 20 pk/cs

DUKAL TONSIL SPONGES

Designed for high absorption. The 100% cotton woven gauze are double strung for easy retrieval and are X-ray detectable to ensure visibility. Packaging ensures count accuracy.

74436	Double Strung Tonsil Sponges, 15" Blue Thread, Sterile 5s, Medium, 7/8", 5/pk, 100 pk/cs
74437	Double Strung Tonsil Sponges, 15" Blue Thread, Sterile 5s, Large, 1", 5/pk, 100 pk/cs



10244



10604

NUTRAMAX TONSIL & ADENOID SPONGES

Made from select surgical cotton fiber to maximize absorbency. Prevent dripping and reduce risk of contamination from lint or loose threads.

10244	Single Strung Tonsil Sponge, Non-Sterile Bulk 500's, Medium, 7/8", 500/cs
10245	Single Strung Tonsil Sponge, Non-Sterile Bulk 500's, Large, 1", 500/cs
10247	Single Strung Tonsil Sponge, Non-Sterile Bulk 500's, X-Large, 1 1/4", 500/cs
10604	Double Strung Tonsil Sponge w/15" blue thread, Sterile 5's, Medium, 7/8", 500/cs
10606	Double Strung Tonsil Sponge w/15" blue thread, Sterile 5's, Large, 1", 500/cs
10608	Double Strung Tonsil Sponge w/15" blue thread, Sterile 5's, X-Large, 1 1/4", 500/cs

Dressings/Gauzes



19-2000



20-3000




20-3002



20-3012



20-3022



20-3042

MILTEX ENT PACKING

Nose packings are made in the USA. High quality with guaranteed performance for one time use.

19-2000	Otitis External Wick 9mm diameter x 15mm long, 50/bx
19-2002	Otitis External Wick 9mm diameter x 24mm long, 10/bx
19-2004	Ear Packing w/string, 12mm diameter x 15cm long, 10/bx
20-3000	Small Standard Nasal Dressing w/string, 4.5cm long x 1.5cm thick x 2cm height, 10/bx
20-3002	Small Standard Nasal Dressing w/airway & string, 4.5cm long x 1.5cm thick x 2cm height, 10/bx
20-3004	Large Standard Nasal Dressing w/string, 8.0cm long x 1.5cm thick x 2cm height, 10/bx
20-3006	Large Standard Nasal Dressing w/airway & string, 8.0cm long x 1.5cm thick x 2cm height, 10/bx
20-3010	Thin Anatomical Nasal Dressing w/string 8.0cm long x .75cm thick x 3cm height, 10/bx
20-3012	Anatomical Nasal Dressing w/string 8.0cm long x 1.5cm thick x 3cm height, 10/bx
20-3014	Anatomical Nasal Dressing w/airway & string 8.0cm long x 1.5cm thick x 3cm height, 10/bx
20-3020	Small Epistaxis Dressing w/string, 5.5cm long x 1.5cm thick x 2.5cm height, 20/bx
20-3022	Large Epistaxis Dressing w/string, 10cm long x 1.5cm thick x 2.5cm height, 10/bx
20-3030	Sinus Dressing B-Type w/string, 2.5cm long x .9cm thick x 2cm height, 10/bx
20-3032	Sinus Dressing K-Type w/string, 3.5cm long x .9cm thick x 1.2cm height, 10/bx
20-3034	Sinus Dressing T-Type w/string, 4.5cm long x .9cm thick x 2cm height, 10/bx
20-3036	Slim Sinus Dressing B-Type w/string, 2.5cm long x .6cm thick x 2cm height, 10/bx
20-3038	Slim Sinus Dressing K-Type w/string, 3.5cm long x .6cm thick x 1.2cm height, 10/bx
20-3040	Slim Sinus Dressing T-Type w/string, 4.5cm long x .6cm thick x 2cm height, 10/bx
20-3042	Sinus Dressing Kit, 10/bx



10533



10544



10545



10546

NUTRAMAX DISSECTING SPONGES

Small, firm gauze sponges for dissection or absorption in delicate surgical procedures. Locked-in X-ray detectable element.

10533	Kittner Blunt Dissectors, 9/16" x 1/4", C-5® Holder, Sterile 5s, 500/cs
10534	Kittner Blunt Dissectors, 9/16" x 1/4", C-5® Holder, Non-Sterile 5s, 500/cs
10544	Cherry Sponges, 1/2", C-5® Holder, Non-Sterile 5s, 500/cs
10545	Peanut Sponges, 3/8", C-5® Holder, Non-Sterile 5s, 500/cs
10546	Round Dissectors, 3/8" x 1/4", C-5® Holder, Non-Sterile 5s, 500/cs
10560	Round Dissectors, 3/8" x 1/4", Sterile 5s, 500/cs
10563	Round Dissectors, 3/8" x 1/4", C-5® Holder, Sterile 5s, 500/cs
10660	Peanut Sponges, 3/8", Sterile 5s, 500/cs
10662	Peanut Sponges, 3/8", Pinned, Sterile 10s, 500/cs
10663	Peanut Sponges, 3/8", C-5® Holder, Sterile 5s, 500/cs
10667	Peanut Sponges, 3/8", Non-Sterile, Bulk 100s, 500/cs
10670	Cherry Sponges, 1/2", Sterile 5s, 500/cs
10673	Cherry Sponges, 1/2", C-5® Holder, Sterile 5s, 500/cs

**NUTRAMAX STICK SPONGES**

All gauze or gauze covered, X-ray detectable stick sponges for prepping or wound cleaning.

10628	Super Gauze-Covered Stick Sponge, 1 1/2", Sterile 10s, 500/cs
10644	Large Gauze-Covered Stick Sponge, 1", Non-Sterile, 500/cs
10648	Super Gauze-Covered Stick Sponge, 1 1/2", Non-Sterile, 500/cs



25-1500-50AS



25-2500-50AS

SHIPPET SPEEDI-WICK® EAR PACK


PVA Expandacell ear canal pack for treatment of otitis external. Keeps antibiotics in contact with ear lining and prevents stenosis of the external ear canal. Absorbs postoperative fluids. Sterile.

25-1500-50AS	Speedi-Wick Long, 15mm x 10mm, 50/bx
25-2500-50AS	Speedi-Wick X-Long, 25mm x 10mm, 50/bx

**SHIPPERT SUPER PAK® POSTERIOR NASAL PACK**


PVA Expandacell nasal pack used in the control and management of posterior bleeds. Less painful and expensive than balloons for posterior bleeds. Tan anterior flange prevents posterior creep. Large foam posterior tip applies pressure to Sphenopalatine Artery.

11W-SU800-08AS Super Pak Posterior Pack, 8 x 1.5 x 4cm, 8/bx

**SHIPPERT TURBO PAK® NASAL PACK**


PVA Expandacell nasal pack used to control bleeding after inferior turbinate surgery. Longitudinal stitching that limits expansion and applies uniform pressure the full length of the pack along the turbinate. Sterile.

11W-TU800-08AS Turbo Pak, Large, 8 x 1 x 2cm, 8/bx


**SHIPPERT VAC WICK® NASAL PACK**

PVA Expandacell nasal suction pack with Y connector for aspiration of blood and mucus. Saline or medication can be administered via the indwelling tube. Sterile.

11-VACWICK-05 Vac Wick Nasal Suction Pack, 5/bx



11-A0500-08AS




11-LT800-08AS

SHIPPERT VENTI-PAK® ANATOMICAL NASAL PACK


PVA Expandacell nasal foam pack with encased angled airway tube for control and management of post-operative epistaxis. Anatomically contoured to conform to the nasal cavity. Sterile.

11-A0500-08AS	Anatomical Nasal w/out airway, Pediatric, 5.6 x 1 x 2.2cm, 8/bx
11-A0800-08AS	Anatomical Nasal w/out airway, Large, 8 x 1.5 x 3.2cm, 8/bx
11-AT800-08AS	Anatomical Nasal w/airway, Large, 8 x 1.5 x 3.2cm, 8/bx
11-LT800-08AS	Venti Pak Wideline Nasal w/airway, Large, 8 x 1.5 x 2cm, 8/bx
11-MT500-08AS	Venti Pak Wideline Nasal w/airway, Medium, 5 x 1.5 x 2cm, 8/bx

**SHIPPERT PLUS PAK® SINUS PACK W/SIDEOUT APPLICATOR**


PVA Expandacell foam with Sideout applicator for ease of insertion. Management of post-op bleeding of ethmoid or maxillary sinus. Prevents lateralization of the turbinate. Gentle but firm pressure to the sinus mucosa. Foam sewn full length for extra strength and string attached for easy removal. Sterile.

30-PLUS2400-10PAS	Plus Pak Sinus, Small, 2.4 x .8 x 1.0cm, 1pair/pouch, 10pouch/bx
30-PLUS2800-10PAS	Plus Pak Sinus, Medium, 2.8 x .8 x 1.2cm, 1pair/pouch, 10pouch/bx
30-PLUS3200-10PAS	Plus Pak Sinus, Large, 3.2 x 1 x 1.2cm, 1pair/pouch, 10pouch/bx

**SHIPPERT POST-OP EAR PACK**

Expands, absorbs and keeps antibiotics in contact with tissue. Small pore size allows less tissue ingrowth and minimizes scar tissue stricture of canal. String attached for removal. Sterile.

25-P12-10	Post-Op Ear Pack, Small, 12mm dia. x 24mm length, 10/bx
25-P15-10	Post-Op Ear Pack, Large, 15mm dia. x 24mm length, 10/bx

**SHIPPERT PURE PAK® NASAL PACK**

For tamponade of anterior or middle segment epistaxis. Expands and gives gentle but firm pressure to the nasal mucosa after trauma or surgery. String attached for easy removal. X-large size for posterior bleeds and Speedi Pedi for pediatric patients. Sterile PVA Expandacell foam.

11S-NS1000-08AS	Pure Pak Nasal Slimline Singles, X-Large, 10 x 1 x 2cm, 8/bx
11S-NS300-08AS	Pure Pak Nasal Slimline Singles, Small, 3 x 1 x 2cm, 8/bx
11S-NS500-08AS	Pure Pak Nasal Slimline Singles, Medium, 5 x 1 x 2cm, 8/bx
11S-NS800-08AS	Pure Pak Nasal Slimline Singles, Large, 8 x 1 x 2cm, 8/bx
11W-NS1000-08AS	Pure Pak Nasal Wideline Singles, X-Large, 10 x 1.5 x 2cm, 8/bx
11W-NS300-08AS	Pure Pak Nasal Wideline Singles, Small, 3 x 1.5 x 2cm, 8/bx
11W-NS500-08AS	Pure Pak Nasal Wideline Singles, Medium, 5 x 1.5 x 2cm, 8/bx
11W-NS800-08AS	Pure Pak Nasal Wideline Singles, Large, 8 x 1.5 x 2cm, 8/bx

**SHIPPert SLIK PAK® LAMINATED NASAL PACK**


Encased in double polymer sleeve for epistaxis management of post-operative bleeding. Polymer sleeve simplifies insertion and removal and reduces ingrowth. Pack absorbs fluid, expands and applies gentle but firm pressure to the nasal mucosa. String attached for removal. Sterile.

11-ANATSLK-08AS	Slik Pak Anatomical, Large, 8 x 1.5 x 3.2cm, 8/bx
11-ATSLK800-08AS	Slik Pak Anatomical w/tube, 8 x 1.5 x 3.2cm, 8/bx
11-SKS500-08AS	Slik Pak Slimline, Medium, 5 x 1 x 2cm, 8/bx
11-SKS800-08AS	Slik Pak Slimline, Large, 8 x 1 x 2cm, 8/bx
11-SKW800-08AS	Slik Pak Wideline, Large, 8 x 1.5 x 2cm, 8/bx
11-SKWT800-08AS	Slik Pak Wideline, Large, w/tube, 8 x 1.5 x 2cm, 8/bx

**SHIPPert SLIK PAK® LAMINATED SINUS PACK**


Encased in double polymer sleeve for control of hemostasis after F.E.S.S. or maxillary sinus surgery. Polymer sleeve simplifies insertion and removal and reduces ingrowth. Pack absorbs fluid, expands and applies gentle but firm pressure to the sinus mucosa. String attached for removal. Sterile.

30-SLK3200-10AS	Slik Pak Non-Stick Sinus Pack, Large, 3.2 x 1 x 1.2cm, 1/pouch, 10pouch/bx
-----------------	--


**SHIPPert SPEEDI PAK® SINUS PACK**

Management of post-op bleeding of ethmoid or maxillary sinus. Prevents lateralization of the turbinate. Gentle but firm pressure to the sinus mucosa. Foam sewn full length for extra strength and string attached for easy removal. Sterile PVA Expandacell foam.

30-2400-10PAS	Speedi Pak, Small, 2.4 x 8 x 1cm, 1pair/pouch, 10pouch/bx
30-2800-10PAS	Speedi Pak, Medium, 2.8 x 8 x 1.2cm, 1pair/pouch, 10pouch/bx
30-3200-10PAS	Speedi Pak, Large, 3.2 x 1 x 1.2cm, 1pair/pouch, 10pouch/bx
30-APP2800-10AS	Speedi Pak Sinus with Sideout Applicator, Medium, 2.8 x 8 x 1.2cm, 1/pouch, 10pouch/bx
30-APP3200-10AS	Speedi Pak Sinus with Sideout Applicator, Large, 3.2 x 1 x 1.2cm, 1/pouch, 10pouch/bx



30-B120-10P



30-K120-10P



30-T60-10P

SHIPPert B, K, AND T-TYPE SINUS PACKS

Sinus packs are comparable or identical to Merocel products. Used to control post-operative bleeding of the ethmoid sinuses. Expand and give gentle but firm pressure to the sinus mucosa. Can be used in the maxillary, sphenoid or frontal sinus. Prevents lateralization of the turbinate. String attached for easy removal. Sterile PVA Expandacell foam.

30-B120-10P	B-Type Regular 2.5 x 1.2 x 2cm, 1pair/pouch, 10pouch/bx
30-B60-10P	B-Type Ultra Slim, 2.5 x 6 x 2cm, 1pair/pouch, 10pouch/bx
30-B90-10P	B-Type Slim 2.5 x .9 x 2cm, 1pair/pouch, 10pouch/bx
30-K120-10P	K-Type Regular, 3.5 x 1.2 x 1.2cm, 1pair/pouch, 10pouch/bx
30-K60-10P	K-Type Ultra Slim, 3.5 x 6 x 1.2cm, 1pair/pouch, 10pouch/bx
30-K90-10P	K-Type Slim, 3.5 x 9 x 1.2cm, 1pair/pouch, 10pouch/bx
30-T60-10P	T-Type Ultra Slim, 4.5 x 6 x 2cm, 1pair/pouch, 10pouch/bx
30-T90-10P	T-Type Slim, 4.5 x 9 x 2cm, 1pair/pouch, 10pouch/bx

**SHIPPert MEDWICK® SINUS IRRIGATION SYSTEM**


PVA Expandacell nasal pack with indwelling tube and luer lock system for delivery of medication or irrigation to the sinus mucosa while pack is in place. Control of hemostasis of sinus hemostasis where irrigation or medication is needed topically. Sterile.

30-MED3200-08AS	Medwick Sinus, Large, 3.2 x 1 x 1.2cm, 1/pouch, 8pouch/bx
-----------------	---

**SHIPPert PLUS PAK® NASAL PACK W/ADVANTAGE APPLICATOR**


PVA Expandacell foam with AdVantage applicator for ease of insertion and blood containment concerns. Epistaxis management that expands and gives gentle but firm pressure to the nasal mucosa after trauma or surgery. String attached for easy removal. Sterile.

11S-PLUS1000-08AS	Plus Pak with Advantage Applicator Slimline Singles, X-Large, 10 x 1 x 2cm, 8/bx
11S-PLUS300-08AS	Plus Pak with Advantage Applicator Slimline Singles, Small, 3 x 1 x 2cm, 8/bx
11S-PLUS500-08AS	Plus Pak with Advantage Applicator Slimline Singles, Medium, 5 x 1 x 2cm, 8/bx
11S-PLUS800-08AS	Plus Pak with Advantage Applicator Slimline Singles, Large, 8 x 1 x 2cm, 8/bx
11W-PLUS1000-08AS	Plus Pak with Advantage Applicator Wideline Singles, X-Large, 10 x 1.5 x 2cm, 8/bx
11W-PLUS300-08AS	Plus Pak with Advantage Applicator Wideline Singles, Small, 3 x 1.5 x 2cm, 8/bx
11W-PLUS500-08AS	Plus Pak with Advantage Applicator Wideline Singles, Medium, 5 x 1.5 x 2cm, 8/bx
11W-PLUS800-08AS	Plus Pak with Advantage Applicator Wideline Singles, Large, 8 x 1.5 x 2cm, 8/bx


**SHIPPERT 30 SECOND RHINO ROCKET® NASAL PACK**

Used post trauma or surgery to control epistaxis. The patented applicator aids in insertion and minimizes blood concerns. When in contact with moisture, the PVA foam absorbs blood, swells to 6 times its compressed size and gives gentle but firm pressure to the nasal mucosa. 30 seconds to insert and string attached for easy removal. Sterile PVA Expandacell foam.

11-OB-OTCMED-20	Over the counter Rhino Rockets Slim, Medium, 20/bx
11OB-OTCSM20	Over the counter Rhino Rockets Slim, Small, 20/bx
11S-S0300-08AS	Rhino Rockets Slim Singles, Small, 3 x 1 x 2cm, 8/bx
11S-S0500-08AS	Rhino Rockets Slim Singles, Medium, 5 x 1 x 2cm, 8/bx
11S-S0800-08AS	Rhino Rockets Slim Singles, Large, 8 x 1 x 2cm, 8/bx
11S-S0800-4PAS	Rhino Rocket Slimline Doubles Large, 8 x 1 x 2cm, 4pairs/bx
11S-S1000-08AS	Rhino Rockets Slim Singles, X-Large, 10 x 1 x 2cm, 8/bx
11W-W0300-08AS	Rhino Rockets Wide Singles, Small, 3 x 1.5 x 2cm, 8/bx
11W-W0500-08AS	Rhino Rockets Wide Singles, Medium, 5 x 1.5 x 2cm, 8/bx
11W-W0800-08AS	Rhino Rockets Wide Singles, Large, 8 x 1.5 x 2cm, 8/bx
11W-W1000-08AS	Rhino Rockets Wide Singles, X-Large, 10 x 1.5 x 2cm, 8/bx



11-M0450-10




11S-D800-10

SHIPPERT M, D, AND P TYPE NASAL PACKS

Comparable or identical to Merocel products. Used for tamponade of anterior or middle segment epistaxis. Expand and give gentle but firm pressure to the nasal mucosa after trauma or surgery. String attached for easy removal. Sterile PVA Expandacell foam.

11-M0450-10	M-Type, Small, 4.5 x 1.5 x 2cm, 10/bx
11-M0800-10	M-Type, Large, 8 x 1.5 x 2cm, 10/bx
11-MSC800-10	M-Type, Large, Single Compressed, 8 x 1.5 x 2cm, 10/bx
11-MT450-10	M-Type, Small, w/airway, 4.5 x 1.5 x 2cm, 10/bx
11-MT800-10	M-Type Large w/airway, 8 x 1.5 x 2cm, 10/bx
11-MTSC800-10	M-Type, Large, w/air Single Compressed, 8 x 1.5 x 2cm, 10/bx
11-P0550-10	P-Type, Small, Epistaxis, 5.5 x 1.5 x 2.5cm, 10/bx
11-P1000-10	P-Type, X-Large, Epistaxis, 10 x 1.5 x 2.5cm, 10/bx
11S-D800-10	D-Type Thin Anatomical 8 x 1 x 3cm, 10/bx
11S-DT550-10	D-Type Thin Pediatric Anat w/air, 5.5 x 1 x 2.5cm, 10/bx
11W-D800-10	D-Type Wideline Anatomical, 8 x 1.5 x 3cm, 10/bx
11W-DT800-10	D-Type Wideline Anatomical w/air, 8 x 1.5 x 3cm, 10/bx

**SHIPPERT MEDWICK® NASAL IRRIGATION PACK**

PVA Expandacell nasal pack with indwelling tube and luer lock system for delivery of medication or irrigation to the nasal mucosa while pack is in place. Control of hemostasis in nasal surgery where irrigation or medication is needed topically. Sterile.

11W-M800-08AS	Medwick Nasal, Large, 8 x 1.5 x 2cm, 8/bx
---------------	---

COVIDIEN/KENDALL CURITY EYE PADS

Constructed of long fiber, absorbent cotton filler and covered with fine mesh gauze on both sides. Sterile and ready to use. Oval, contoured shape. One per package in convenient dispensing carton.

2841	Oval Eye Pad, 1 5/8" x 2 5/8", Sterile, 1/pk, 50/dispensing box, 12 bx/cs
------	---

DERMA SCIENCES SURECARE® EYE PADS

Constructed of absorbent cotton filler covered with fine mesh gauze on both sides. Sterile and ready to use, this pad is ideal for applying medication or to cover and protect they eye. Latex free.

82911	Eye Pad, Sterile, 2 5/8" x 1 5/8" Oval, 1's, 25/tray, 24 tray/cs
-------	--

**DUKAL EYE PADS**

Soft and gentle to the touch and offer protection to the site. The sealed edges provide minimal linting and the easy place tab on the back ensures clean placement. They are individually wrapped and packaged in a dispenser box for your convenience.

841	Eye Pad, Oval, 1 5/8" x 2 5/8", Sterile, 50/bx, 12 bx/cs
841B	Eye Pad, Oval, 1 5/8" x 2 5/8", Sterile, 1000 pads/cs

**AMD-RITMED OVAL EYE PADS**

A1110	1 5/8" x 2 5/8 Oval Eye Pad, 50/pk, 12 pk/cs
A1111	2 1/8" x 2 5/8 Oval Eye Pad, 50/pk, 12 pk/cs

**DYNAREX OVAL EYE PADS**

3359	2 5/8" x 1 5/8" Eye Pad, Sterile Bulk Pack, 1000/cs
3360	2 5/8" x 1 5/8" Eye Pad, Sterile, 50/bx, 12 bx/cs

BD VISISORB™ ABSORBENT STICKS

581060	6.6cm x 5mm Absorbent Stick, 10/pk, 25 pk/bx
--------	--

Dressings/Gauzes - Tapes

BD VISISPEAR™ EYE SPONGE

581089 7cm Eye Sponge, 5/pkg, 20 pk/bx

LOOK EASY SHIELD™ SYSTEM

Includes an Eye Garter® cover, a durable aluminum eye shield and a foam adhesive for application. System limits light exposure to the eye, but still allows the patient to see. Ideal for patients after cataract surgery.


G101-100RB	Eye Garter Only, Adult & Youth, Royal Blue, 100/pkg
G101-50	Eye Garter Only, Assorted Colors, 50/pkg
G102-50	Eye Garter Only, White, 50/pkg
G102-500	Eye Garter Only, White, 500/bx
G251-25	Adult w/colored cover w/adhesive, 25/pkg
G252-25	Adult w/white cover w/adhesive, 25/pkg
G257-25	Easy Shield Adhesives, Adult/Youth, 25/pkg



DERMA SCIENCES OVAL EYE PADS

Sterile pads readily fit contours of eye cavity. Cover of fine mesh seamless gauze, with thick cotton fiber filler.

17570	2 1/8" x 2 5/8" Oval Eye Pads, 50/bx, 24 bx/cs
17571	2 1/8" x 2 5/8" Oval Eye Pads, Bulk, 500/cs
17576	1 5/8" x 2 5/8" Oval Eye Pads, 50/bx, 24 bx/cs
9315000	Eye Strips, Non-Sterile, 3/8" x 5", 2/card, 1250 cards/cs



PRO ADVANTAGE® RETENTION TAPE

These dressings are the ideal cover dressings made from Sontara® non-woven fabric. The permeable structure of the fabric ensures oxygen permeability. The fabric is stretchable, conforms well to body contours and allows body movements when in place. The fabric is non-absorbent and provides a barrier for blood and fluid strike-through. The adhesive is gentle to the skin but is very secure when in place. Because it is air permeable, it can be used to provide additional fixation for other adhesive dressings. Retention tape sizes provide added versatility.

P157102	Retention Tape, Perf
P157104	
P157106	Retention Tape, Perf



4209



4210

SMITH & NEPHEW HYPAFIX™ DRESSING RETENTION ROLLS

Secures dressings, tubes, catheters and instruments quickly and easily even in hard-to-tape areas. A comfortable, non-woven dressing retention tape that stretches for easy application and patient comfort. Water-resistant with a non-irritating adhesive that acts as an excellent secondary dressing for gauze, foams, gels, etc.

4209	2" x 10 yds Dressing Retention Roll, 24/cs
4210	4" x 10 yds Dressing Retention Roll, 24/cs
4211	6" x 10 yds Dressing Retention Roll, 24/cs
4215	2" x 2 yds Dressing Retention Roll, 36/cs
4216	4" x 2 yds Dressing Retention Roll, 18/cs
4217	6" x 2 yds Dressing Retention Roll, 12/cs



3M™ BLENDERM™ SURGICAL TAPE

Occlusive, transparent plastic tape. Transparent to allow monitoring. Waterproof. Flexible for hard-to-tape areas. Hypoallergenic. Occlusive test patches.

1525-0	1/2" x 5 yd Surgical Tape, 24 rl/bx, 10 bx/cs
1525-1	1" x 5 yd Surgical Tape, 12 rl/bx, 10 bx/cs
1525-2	2" x 5 yd Surgical Tape, 6 rl/bx, 10 bx/cs



3M™ CLOTH ADHESIVE TAPE

Combines high backing strength with good adhesion. Ideal for securing medical devices and for light immobilization support.


2950-0	1/2" x 10 yd Cloth Adhesive Tape, 24/bx, 10 bx/cs
2950-1	1" x 10 yd Cloth Adhesive Tape, 12/bx, 10 bx/cs
2950-2	2" x 10 yd Cloth Adhesive Tape, 6/bx, 10 bx/cs
2950-3	3" x 10 yd Cloth Adhesive Tape, 4/bx, 10 bx/cs



3M™ DURAPORE™ SURGICAL TAPE

"Silk-like" cloth tape with strong adhesion for securing dressings or devices. Adheres well with good initial stick. Comfortable, versatile, high strength. Allows skin to breathe. Hypoallergenic. Convenient bi-directional tear.

1538-0	1/2" x 10 yd Surgical Tape, 24 r/bx
1538-1	1" x 10 yd Surgical Tape, 12 r/bx, 10 bx/cs
1538-2	2" x 10 yd Surgical Tape, 6 r/bx, 10 bx/cs
1538-3	3" x 10 yd Surgical Tape, 4 r/bx, 10 bx/cs
1538S-1	Single Use, 1" x 1 1/2 yd Surgical Tape, 100 r/bx, 5 bx/cs
1538S-2	Single Use, 2" x 1 1/2 yd Surgical Tape, 50 r/bx, 5 bx/cs



3M™ MEDIPORE™ H SOFT CLOTH SURGICAL TAPE

A stronger holding power for more challenging applications. Porous and breathable. Soft, stretchy tape is hypoallergenic. Perforated rolls are easy to tear.

2861	1" x 10 yd Cloth Surgical Tape, 2 rolls/pk, 12 pk/cs
2862	2" x 10 yd Cloth Surgical Tape, 12/cs
2862S	2" x 2 yd Cloth Surgical Tape (short roll), 48/cs
2863	3" x 10 yd Cloth Surgical Tape, 12/cs
2864	4" x 10 yd Cloth Surgical Tape, 12/cs
2864S	4" x 2 yd Cloth Surgical Tape (short roll), 24/cs
2866	6" x 10 yd Cloth Surgical Tape, 12/cs
2866S	6" x 2 yd Cloth Surgical Tape (short roll), 16/cs
2868	8" x 10 yd Cloth Surgical Tape, 6/cs



3M™ MEDIPORE™ SOFT CLOTH PRE-CUT DRESSING COVERS

Soft, stretchy and conformable. Easy to use and easy on skin. Fast, easy pad delivery saves dressing time and materials. Excellent adhesion, yet gentle to the skin. Highly conformable, with moderate three-way stretch. Hypoallergenic.

2954	Dressing Covers, 3 7/8" x 4 5/8", 3 sheets/pad, 25 pads/bx, 4 bx/cs
2955	Dressing Covers, 3 7/8" x 7 7/8", 3 sheets/pad, 25 pads/bx, 4 bx/cs
2956	Dressing Covers, 5 7/8" x 5 7/8", 3 sheets/pad, 25 pads/bx, 4 bx/cs
2957	Dressing Covers, 5 7/8" x 11", 3 sheets/pad, 25 pads/bx, 4 bx/cs
2958	Dressing Covers, 7 7/8" x 11", 3 sheets/pad, 25 pads/bx, 4 bx/cs

3M™ MEDIPORE™ SOFT CLOTH SURGICAL TAPE

Soft, stretchy and conformable. Easy to use and easy on skin. Easy-tear perforated rolls. No fussing with a paper liner. Easy to customize to the length you need. Handles easily without curling or sticking to itself; less waste from misapplication.

2961	1" x 10 yd Cloth Surgical Tape, 2/pk, 12 pk/cs
2962	2" x 10 yd Cloth Surgical Tape, 12/cs
2962S	2" x 2 yd Cloth Surgical Tape (short roll), 48/cs
2963	3" x 10 yd Cloth Surgical Tape, 12/cs
2964	4" x 10 yd Cloth Surgical Tape, 12/cs
2964S	4" x 2 yd Cloth Surgical Tape (short roll), 24/cs
2966	6" x 10 yd Cloth Surgical Tape, 12/cs
2966S	6" x 2 yd Cloth Surgical Tape (short roll), 16/cs
2968	8" x 10 yd Cloth Surgical Tape, 6/cs



3M™ MICROFOAM™ SURGICAL TAPES & STERILE TAPE PATCH

Highly conformable elastic foam tape for compression applications, or securing dressings on challenging areas. Secure adhesion conforms to irregularly contoured sites. Gives with swelling. Can be used to provide compression. Water resistant and hypoallergenic.

1528-1	1" x 5 1/2 yd (stretched) Surgical Tape, 12 r/bx, 6 bx/cs
1528-2	2" x 5 1/2 yd (stretched) Surgical Tape, 6 r/bx, 6 bx/cs
1528-3	3" x 5 1/2 yd (stretched) Surgical Tape, 4 r/bx, 6 bx/cs
1528-4	4" x 5 1/2 yd (stretched) Surgical Tape, 3 r/bx, 6 bx/cs
1562	4" x 7" Surgical Tape, 1 Individually Packaged, Sterile Tape Patch, 25/bx, 4 bx/cs



3M™ MICROPOR™ SURGICAL TAPES

A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. An economical, general purpose, breathable surgical tape. Available in white or tan.

1530-0	1/2" x 10 yd Paper Surgical Tape, 24 r/bx, 10 bx/cs
1530-1	1" x 10 yd Paper Surgical Tape, 12 r/bx, 10 bx/cs
1530-2	2" x 10 yd Paper Surgical Tape, 6 r/bx, 10 bx/cs
1530-3	3" x 10 yd Paper Surgical Tape, 4 r/bx, 10 bx/cs
1530S-1	Single Use, 1" x 1 1/2 yd Paper Surgical Tape, 100 r/bx, 5 bx/cs
1530S-2	Single Use, 2" x 1 1/2 yd Paper Surgical Tape, 50 r/bx, 5 bx/cs
1533-0	Tan, 1/2" x 10 yd Paper Surgical Tape, 24 r/bx, 10 bx/cs
1533-1	Tan, 1" x 10 yd Paper Surgical Tape, 12 r/bx, 10 bx/cs
1533-2	Tan, 2" x 10 yd Paper Surgical Tape, 6 r/bx, 10 bx/cs
1535-0	Dispenser Pack, 1/2" x 10 yd Paper Surgical Tape, 24 r/bx, 10 bx/cs
1535-1	Dispenser Pack, 1" x 10 yd Paper Surgical Tape, 12 r/bx, 10 bx/cs
1535-2	Dispenser Pack, 2" x 10 yd Paper Surgical Tape, 6 r/bx, 10 bx/cs
1535-3	Dispenser Pack, 3" x 10 yd Paper Surgical Tape, 4 r/bx, 10 bx/cs

**3M™ TRANSPORE™ SURGICAL TAPE**

Transparent, easy to tear, perforated plastic tape for securing dressings or devices. Hypoallergenic. Easy to handle with gloves. Perforated for easy bi-directional tear.

1527-0	1/2" x 10 yd Surgical Tape, 24 rl/bx, 10 bx/cs
1527-1	1" x 10 yd Surgical Tape, 12 rl/bx, 10 bx/cs
1527-2	2" x 10 yd Surgical Tape, 6 rl/bx, 10 bx/cs
1527-3	3" x 10 yd Surgical Tape, 4 rl/bx, 10 bx/cs
1527S-1	Single Use, 1" x 1 1/2 yd Surgical Tape, 100 rl/bx, 5 bx/cs
1527S-2	Single Use, 2" x 1 1/2 yd Surgical Tape, 50 rl/bx, 5 bx/cs

3M™ TRANSPORE™ WHITE DRESSING TAPE

1534-0	White Dressing Tape, 1/2" x 10 yd, 24/bx, 10 bx/cs
1534-1	White Dressing Tape, 1" x 10 yd, 12/bx, 10 bx/cs
1534-2	White Dressing Tape, 2" x 10 yd, 6/bx, 10 bx/cs
1534-3	White Dressing Tape, 3" x 10 yd, 4/bx, 10 bx/cs

**AMD-RITMED HYPOALLERGENIC CLEAR SURGICAL TAPE**

A5005	1/2" x 10 yds Clear Surgical Tape, 24/bx, 12 bx/cs
A5010	1" x 10 yds Clear Surgical Tape, 12/bx, 12 bx/cs
A5020	2" x 10 yds Clear Surgical Tape, 6/bx, 12 bx/cs
A5030	3" x 10 yds Clear Surgical Tape, 4/bx, 12 bx/cs

**AMD-RITMED HYPOALLERGENIC PAPER SURGICAL TAPE**

A5105	1/2" x 10 yds Paper Surgical Tape, 24/bx, 12 bx/cs
A5110	1" x 10 yds Paper Surgical Tape, 12/bx, 12 bx/cs
A5120	2" x 10 yds Paper Surgical Tape, 6/bx, 12 bx/cs
A5130	3" x 10 yds Paper Surgical Tape, 4/bx, 12 bx/cs

**AMD-RITMED HYPOALLERGENIC SILK SURGICAL TAPE**

A5205	1/2" x 10 yds Silk Surgical Tape, 24/bx, 12 bx/cs
A5210	1" x 10 yds Silk Surgical Tape, 12/bx, 12 bx/cs
A5220	2" x 10 yds Silk Surgical Tape, 6/bx, 12 bx/cs
A5230	3" x 10 yds Silk Surgical Tape, 4/bx, 12 bx/cs

**DUKAL SURGICAL TAPE - CLOTH**

A good quality durable tape that can be ripped in both directions. This strong, highly adhesive tape was designed for bilateral tearing without the use of scissors and is non-irritating to provide maximum comfort. This durable tape will leave minimal residue on the skin. These tapes are hypoallergenic and latex free.

C110	Surgical Tape, 1" X 10 yds, 12 rl/bx, 12 bx/cs
C210	Surgical Tape, 2" x 10 yds, 6 rl/bx, 12 bx/cs
C310	Surgical Tape, 3" x 10 yds, 4 rl/bx, 12 bx/cs
C50	Surgical Tape, 1/2" x 10 yds, 24 rl/bx, 12 bx/cs

DUKAL WATERPROOF TAPE

Conveniently offered in plastic spools and shells and it can be cleaned by soap and water. It is ideal for conventional strapping, physical therapy, dressing adhesion, and for the EMS market.

AT2316	Tape, 1/2" x 5 yds, Waterproof, 24 rl/bx, 18 bx/cs
AT2323	Tape, 1/2" x 2.5yds, Waterproof, 36 rl/bx, 16 bx/cs

**DUKAL SURGICAL TAPE - PAPER**

P110	Surgical Tape, 1" x 10 yds, 12 rl/bx, 12 bx/cs
P210	Surgical Tape, 2" x 10 yds, 6 rl/bx, 12 bx/cs
P310	Surgical Tape, 3" x 10 yds, 4 rl/bx, 12 bx/cs
P50	Surgical Tape, 1/2" x 10 yds, 24 rl/bx, 12 bx/cs

**DUKAL SURGICAL TAPE - TRANSPARENT**

T110	Surgical Tape, 1" x 10 yds, 12 rl/bx, 12 bx/cs
T210	Surgical Tape, 2" x 10 yds, 6 rl/bx, 12 bx/cs
T310	Surgical Tape, 3" x 10 yds, 4 rl/bx, 12 bx/cs
T50	Surgical Tape, 1/2" x 10 yds, 24 rl/bx, 12 bx/cs

**DYNAREX SURGICAL TAPE**

Hypoallergenic. 10 yards length.

3551	Paper Surgical Tape, 1/2", 24/bx, 12 bx/cs
3552	Paper Surgical Tape, 1", 12/bx, 12 bx/cs
3553	Paper Surgical Tape, 2", 6/bx, 12 bx/cs
3554	Paper Surgical Tape, 3", 4/bx, 12 bx/cs
3561	Cloth Surgical Tape, 1/2", 24/bx, 12 bx/cs
3562	Cloth Surgical Tape, 1", 12/bx, 12 bx/cs
3563	Cloth Surgical Tape, 2", 6/bx, 12 bx/cs
3564	Cloth Surgical Tape, 3", 4/bx, 12 bx/cs
3571	Transparent Surgical Tape, 1/2", 24/bx, 12 bx/cs
3572	Transparent Surgical Tape, 1", 12/bx, 12 bx/cs
3573	Transparent Surgical Tape, 2", 6/bx, 12 bx/cs
3574	Transparent Surgical Tape, 3", 4/bx, 12 bx/cs



HARTMANN-CONCO AC-TAPE PLUS® LATEX FREE ADHESIVE TAPE

All-cotton, latex free tape that provides superior support, compression and fixation. Perfect for a variety of applications especially athletic conditions. All-cotton elasticity allows it to conform to every contour yet stay securely in place. Wide core and feathered edge.

64210000	2" x 5 yds Adhesive Tape, 6/bx, 12 bx/cs
64310000	3" x 5 yds Adhesive Tape, 4/bx, 12 bx/cs
64410000	4" x 5 yds Adhesive Tape, 6/bx, 6 bx/cs
65210000	2" x 5 yds Adhesive Tape, Team Pack, 24/cs
65310000	3" x 5 yds Adhesive Tape, Team Pack, 16/cs



HARTMANN-CONCO PRO'S CHOICE® ADHESIVE & COHESIVE ATHLETIC STRETCH TAPE

The first choice of athletic trainers, physicians and therapists who need a durable elastic, hand-tear tape for multiple applications. Constructed with superior tensile strength and support. Stays securely in place, even during the most rigorous activity, its lightweight fabric blend makes it cool and comfortable to wear. Ideal for a range of taping and strapping applications. Also tears easily, leaving a clean, crisp edge. Easy unwind allows for consistent taping all the way to the core. Pro's Choice Adhesive is latex free. NOTE: Pro's Choice Cohesive contains latex.

40195001	1.5" x 6 yds, White, 32rl/cs
40195002	2" x 6 yds, White, 24rl/cs
40195003	3" x 6 yds, White, 16rl/cs
40195012	2" x 6 yds, Black, 24rl/cs
40195013	3" x 6 yds, Black, 16rl/cs
64710000	1 1/2" x 7.5 yds Stretch Tape, 32/cs
64720000	2" x 7.5 yds Stretch Tape, 24/cs
64730000	3" x 7.5 yds Stretch Tape, 16/cs
64740000	1 1/2" x 5 yds Stretch Tape, 32/cs
64750000	2" x 5 yds Stretch Tape, 24/cs
64760000	3" x 5 yds Stretch Tape, 16/cs
64770000	1" x 5 yds Stretch Tape, 48/cs



HARTMANN-CONCO OMNIFIX® NON-WOVEN DRESSING RETENTION TAPE

Designed to make application and removal easy. Soft, secure, flexible and latex free. Permeable to air and water vapor. Ideal for patients with hypersensitive skin and a tendency toward skin allergies. Perfect to use for securing wound dressings, catheters and drainage lines. Won't restrict movement and does not absorb X-rays.

36200000	2" x 2 yds Retention Tape, 24/cs
36210000	2" x 10 yds Retention Tape, 24/cs
36400000	4" x 2 yds Retention Tape, 24/cs
36410000	4" x 10 yds Retention Tape, 12/cs
36610000	6" x 10 yds Retention Tape, 12/cs

J&J COACH® SPORTS CARE BANDAGES

004975	Bandage, Elastic, 2" x 2.2 yds, Self-Adhering, 24/cs
004976	Bandage, Elastic, 3" x 2.2 yds, Self-Adhering, 24/cs
007902	1 1/2" x 10 yds Sports Tape, 6/bx, 4 bx/cs
007949	1 1/2" x 10 yds Value Pack Sports Tape, 4/bx, 12 bx/cs



J&J ELASTIKON™ ELASTIC TAPE

Helps maintain constant controlled pressure on large draining wounds. Convenient custom fit around hard to tape areas. Ideal for securing dynamic pressure dressing. 5 yard lengths.

005170	2" x 5 yds Elastic Tape, 24 rolls/cs
005172	1" x 5 yds Elastic Tape, 12 rl/bx, 12 bx/cs
005174	2" x 5 yds Elastic Tape, 6 rl/bx, 12 bx/cs
005175	3" x 5 yds Elastic Tape, 4 rl/bx, 12 bx/cs
005177	4" x 5 yds Elastic Tape, 6 rl/bx, 6 bx/cs

J&J FIRST AID TAPES DISPENSER

004811	1" x 5 yds Paper First Aid Tape, 6/bx, 4 bx/cs
004831	1/2" x 5 yds Cloth First Aid Tape, 6/bx, 8 bx/cs

J&J FIRST AID TAPES NON-DISPENSER

004815	1" x 10 yds Clear First Aid Tape, Blister Card, 6/bx, 8 bx/cs
004900	1" Hurt-Free™ First Aid Tape, Blister Card, 6/bx, 4 bx/cs
004902	2" Hurt-Free™ First Aid Tape, Blister Card, 4/bx, 6 bx/cs
004979	Tape, Cloth, First Aid, All Purpose, 1" x 10 yds, No Dispenser, Blister Card, 48/cs
004980	Tape, Cloth, First Aid, All Purpose, 2" x 10 yd, No Dispenser, Blister Card, 48/cs
004981	Tape, Paper, First Aid, Non-Irritating, 1" x 10yds, No Dispenser, Blister Card, 48/cs



J&J SECURE-COMFORT™ HYPOALLERGENIC CLOTH TAPE

Flexible and conforming for contouring applications, as well as x-ray transparent to allow easy viewing of the wound site. Perforated to provide custom tearability for individual taping techniques in IV therapy, tube and general dressing securement. Latex free adhesive. Plastic backing. Slight opaque/matte finish.

005162	1" x 5 yds Cloth Tape, 6/bx, 4 bx/cs
--------	--------------------------------------

**J&J WATERPROOF TAPE**

Holds securely in presence of moisture. Prevents moisture from reaching patient's skin. Unwinds uniformly and adheres. Can be cleaned with soap and water. Easy to apply. Can be written on. 10 yard lengths.

005131	1/2" x 10 yds Waterproof Tape, 24/bx, 12 bx/cs
005134	2" x 10 yds Waterproof Tape, 6/bx, 12 bx/cs
005135	3" x 10 yds Waterproof Tape, 4/bx, 12 bx/cs

**J&J ZONAS™ POROUS TAPE**

Strong, easy to remove from cloth tape for immobilization of sprains, support strapping, patient positioning and restraining. 10 yard lengths.

005103	1/2" x 10 yds Porous Tape, 24 rl/bx, 12 bx/cs
005104	1" x 10 yds Porous Tape, 12 rl/bx, 12 bx/cs
005106	2" x 10 yds Porous Tape, 6 rl/bx, 12 bx/cs
005107	3" x 10 yds Porous Tape, 4 rl/bx, 12 bx/cs

**COVIDIEN/KENDALL SILK TAPE**

The most economical, high-strength and high-adhesion tape. Excellent for securing regular and bulky dressings, tubing, and immobilization.

7137	Silk Tape, Hypoallergenic, 1/2" x 10 yds, 24/bx, 12 bx/cs
7138	Silk Tape, Hypoallergenic, 1" x 10 yds, 12/bx, 12 bx/cs
7138S	Silk Tape, Hypoallergenic, 1" x 1.5 yds, Single Roll, 100/bx, 5 bx/cs
7139	Silk Tape, Hypoallergenic, 2" x 10 yds, 6/bx, 12 bx/cs
7139S	Silk Tape, Hypoallergenic, 2" x 1.5 yds, Single Roll, 50/bx, 5 bx/cs
7140	Silk Tape, Hypoallergenic, 3" x 10 yds, 4/bx, 12 bx/cs

**COVIDIEN/KENDALL STANDARD POROUS TAPE**

Adhesive tape with linear porosity for maximum breathability. Designed to provide high tensile strength and aggressive adhesion. Pinked edges allow easy tearing.

2304	Standard Porous Tape, 1/2", 24/bx, 12 bx/cs
2531	Standard Porous Tape, 1", 12/bx, 12 bx/cs
3027	Standard Porous Tape, 1 1/2", 8/bx, 12 bx/cs
3615	Standard Porous Tape, 4", 8/bx, 12 bx/cs
5806	Standard Porous Tape, Tan, 1", 12/bx, 12 bx/cs
5807	Standard Porous Tape, Tan, 1.5", 8/bx, 12 bx/cs
5808	Standard Porous Tape, Tan, 2", 6/bx, 12 bx/cs
6613	Standard Porous Tape, 2", 6/bx, 12 bx/cs
7046	Standard Porous Tape, 3", 4/bx, 12 bx/cs

COVIDIEN/KENDALL PAPER TAPE

Excellent for repeated applications on sensitive, fragile skin and elderly patients. Hypoallergenic, breathable and gentle to the skin. Minimal adhesive residue upon removal.

1596C	Hypoallergenic Paper Tape, 1/2" x 10yd, 24 rl/bx, 12 bx/cs
1596T	Paper Tape, Hypoallergenic, 0.5" x 10 yds, Tan, 24/bx, 10 bx/cs
1914	Hypoallergenic Paper Tape, 1" x 10 yds, 12 rl/bx, 12 bx/cs
1914S1	Paper Tape, Hypoallergenic, 1" x 1.5 yds, 100/bx, 5 bx/cs
1914T	Paper Tape, Hypoallergenic, 1" x 10 yds, Tan, 12/bx, 10 bx/cs
2419	Hypoallergenic Paper Tape, 2" x 10 yds, 6 rl/bx, 12 bx/cs
2419S2	Paper Tape, Hypoallergenic, 2" x 1.5 yds, Tan, 50/bx, 5 bx/cs
2419T	Paper Tape, Hypoallergenic, 2" x 10 yds, Tan, 6/bx, 10 bx/cs
3394C	Hypoallergenic Paper Tape, 3" x 10 yds, 4 rl/bx, 12 bx/cs

**COVIDIEN/KENDALL WATERPROOF TAPE**

Recognized leader in the field. Ideal for general taping needs and outpatient use. Waterproof and oil resistant tape with aggressive adhesion and high tensile strength.

3063	Waterproof Tape 1/2", 24/bx, 12 bx/cs
3142	Waterproof Tape 1", 12/bx, 12 bx/cs
3267	Waterproof Tape 2", 6/bx, 12 bx/cs
3354	Waterproof Tape 3", 4/bx, 12 bx/cs

COVIDIEN/KENDALL BABY TAPE™ AND HYDRO TAPE™ ADHESIVE TAPES

Clear adhesive hydrogel formula gently yet firmly secures catheters, I.V. lines, and lead wires to infant's skin. Tapes may be repositioned and reapplied without irritating or damaging skin. Available in individually packaged tape strips or on a roll.

MI00677	Adhesive Hydrogel Tape, Strips, 50strips/pk, 10 pk/cs
---------	---

**MOLNLYCKE MEPITAC®**

Mepitac is highly conformable, easy to apply and remove. The Safetac layer assures secure fixation and atraumatic removal from skin. Mepitac can be lifted and repositioned if inspection and repositioning are necessary. Mepitac is an ideal fixation for patients with fragile skin. Mepitac is suitable for patients requiring repeated application and removal of tape over the same area e.g. neonates and dialysis patients.

298300	.75" x 118" (2cm x 3cm) Soft Silicone Tape, 1rl/bx
298400	1.5" x 59" (4cm x 1.5cm) Soft Silicone Tape, 1rl/bx



NUTRAMAX WATERPROOF ADHESIVE TAPE

2025033	1/2" x 2 1/2 yds Adhesive Tape, 36 dz rolls/cs
2035033	1/2" x 5 yds Adhesive Tape, 48 rl/cs
2040033	1" x 5 yds Adhesive Tape, 48 rl/cs
2045033	2" x 5 yds Adhesive Tape, 72 rl/cs
2050033	2" x 5 yds Adhesive Tape, Tri-Cut, 72 rl/cs
2060033	1/2" x 10 yds Adhesive Tape, 36 rl/cs



SOUTHWEST STA-FIX TAPE™

Very thin, non-woven polyester base with an open structure and coated with adhesive by a special process to produce a highly porous, breathable tape. The hypo-allergenic adhesive of moderate strength gives a tape which will firmly adhere without damage to tissue upon removal. The non-woven polyester material has a slight amount of stretch which allows the tape to conform to complex contour and difficult areas to tape. The stretch reduces irritation and pain, especially when used with Elasto-Gel™ dressings. Although Sta-Fix is breathable, it is water resistant.

TP100	1" x 400" Tape, 10 rolls/cs
TP200	2" x 400" Tape, 10/cs



HARTMANN-CONCO AC-TAPE® ELASTIC ADHESIVE BANDAGES

Woven with a high twist cotton yarn and combined with a latex free adhesive base. Provides excellent support, compression and fixation. Its all-cotton elasticity allows it to conform to every contour, and its non-latex adhesive helps to reduce skin irritation and damage. The elasticity doesn't restrict movement - making it ideal to use during athletic activities. Finished edge.

64100000	1" x 5 yd Adhesive Tape, 12 rl/bx, 12 bx/cs
64200000	2" x 5 yd Adhesive Tape, 6 rl/bx, 12 bx/cs
64300000	3" x 5 yd Adhesive Tape, 4 rl/bx, 12 bx/cs
64400000	4" x 5 yd Adhesive Tape, 6 rl/bx, 6 bx/cs
65200000	2" x 5 yds Adhesive Tape, Team Pack, 24/cs
65300000	3" x 5 yds Adhesive Tape, Team Pack, 16/cs



COVIDIEN/KENDALL CURITY DRESSING KIT

Components sequentially packaged in accordance with aseptic wound care techniques. All-in-one convenience improves nursing efficiency. All kits contain a waterproof drape, waste bag with tie, and one roll of medical tape. One of each component and sterile unless otherwise indicated.

80159	Includes: Iodophor Triple Swabs, (4) 4" x 4" Gauze Sponges, (2) 5" x 9" ABD Pad, (2) Forceps, Sharp/Blunt Scissors, Safety Pin, (1) Pair Gloves, 20/cs
-------	--

MEDICAL ACTION DRESSING CHANGE TRAYS

For cleansing and dressing central venous catheter insertion sites. Each tray contains: (2) latex gloves, (3) alcohol swab sticks, (3) PVP swab sticks, PVP ointment and gauze. All trays are CSR wrapped.

2614	Dressing Change Tray w/Tegaderm dressing, mask, 2" x 2" gauze, 4" x 4" gauze, 20/cs
2615C	Dressing Change Tray w/Tegaderm dressing, mask, 2" x 2" gauze, 4" x 4" gauze, 20/cs

BD PERSIST™ DRESSING CHANGE KITS

Basic kit includes: (4) 2" x 2" gauze sponge, (1) 2" x 2" non-woven split sponge, (1) alcohol swabstick, (1) cotton-tipped applicator, (1) dressing change label, (1) CSR wrap, (1) pr gloves-exam style, (1) pr gloves-wallet style, (1) mask, (1) Persist™ swabstick (triple), (1) waste bag, (1) plastic tie, (3) reinforce tape strips, (1) roll Transpore® tape, (1) tape measure.

386508	Basic Kit w/(1) Large OpSite™ IV3000 Dressing, 30/cs
--------	--

BD PERSIST™ DRESSING CHANGE MINI KITS

Basic kit includes: (4) 2" x 2" gauze sponge, (1) alcohol swabstick, (1) cotton-tipped applicator, (1) dressing change label, (1) pr gloves-wallet style, (1) mask, (1) Persist™ swabstick (triple), (1) reinforce tape strip (3), (1) tape measure.

386510	Basic Kit w/(1) Oval Tegaderm™ Dressing, 30/cs
--------	--

BUSSE CENTRAL LINE DRESSING CHANGE KIT WITH TEGADERM™ DRESSING

Kit contains: (1) polylined drape, (1) mask, (1) pr gloves, (3) isopropyl alcohol swabsticks, (3) iodophor PVP swabsticks, (1) skin protectant swabstick, (1) 2" x 2" non-woven dressing, (1) 2" x 2" gauze dressing, (1) 3" x 3" gauze dressing, (1) 4" x 5" Tegaderm™ transplant dressing, (1) roll medical tape.

834	Central Line Dressing Change Kit, 50/cs
-----	---

Wound Care Trays



BUSSE CENTRAL LINE DRESSING CHANGE TRAY WITH TEGADERM™ DRESSING

Tray with transparent dressing. Each Tray Contains: 1 CSR Overwrap, 1 Mask, 1 Pair Gloves, 1 Non-Adhering Dressing, 2" x 2", 1 Non-Adhering Split Dressing, 2" x 2", 1 Gauze Dressing, 2" x 2", 1 Gauze Dressing, 4" x 4", 1 Iodophor PVP Swabstick (triples), 1 Povidone Iodine Ointment, 1 Transparent Dressing, 4" x 5", 1 Roll Medical Tape, 1 Dressing Change Label, 1 Isopropyl Alcohol Swabsticks (triples), 1 Tray

829 Tray w/Tegaderm™ Dressing & Isopropyl Alcohol, 30/cs

BUSSE CENTRAL LINE KIT WITH CHLORAPREP® ONE-STEP APPLICATOR & TEGADERM™ DRESSING

835 Central Line Kit w/Tegaderm™ Dressing & ChloraPrep® One-Step Applicator, 30/cs

BUSSE CENTRAL VENOUS CATHETER DRESSING CHANGE KIT

Kit contains: (1) polylined drape, (1) pr gloves, (1) iodophor PVP ointment, (3) isopropyl alcohol swabsticks, (3) iodophor PVP swabsticks, (1) skin protectant swabstick, (1) 2" x 2" non-woven dressing, (1) 2" x 2" non-woven split dressing, (1) 3" x 3" gauze dressing, (1) 4" x 5" Tegaderm™ transplant dressing, (1) roll medical tape, (1) dressing change label.

832 Central Venous Catheter Dressing Change Kit, 50/cs



BUSSE DRESSING CHANGE TRAY WITH INSTRUMENTS

Each set contains: 1 Polylined Sterile Field, 1 Poly Bag w/twist tie, 1 Safety Pin, 1 Iodophor PVP Swabstick (triples), 4 Gauze Sponges, 1 Sharp/Blunt Scissors, 1 Roll Medical Tape, 1 Transfer Forceps, 1 Posi-Grip Forceps, 1 Abdominal Pad, 1 Tray

759 Dressing Change Tray, 20/cs



BUSSE PRIMARY DRESSING CHANGE TRAY

Each set includes: 1 three compartment tray; 1 pair gloves; 1 polylined sterile field; 1 abdominal pad; 8 4"x4" gauze sponges; 1 100ml container sterile saline solution; 1 polybag w/twist tie; 1 roll plastic tape; 1 sharp/blunt scissors; 1 stretch gauze dressing; 1 dressing change label.

1443 Dressing Change Tray, 20/cs



AMSINO AMSURE® TRACHEOSTOMY CARE TRAY

Sterile. Latex free.

AS860 Tracheostomy Care Tray Contains: 500mL Tray, Cleaning Brush, (4) 4x4 Gauze Sponges, (4) Pipe Cleaners, Vinyl Gloves, (2) Cotton-Tip Applicators, Twill Tape, Trach Dressing, Waterproof Drape, Popup Basin, 20/cs

AS861 Tracheostomy Care Tray Contains: 2-Compartment Tray, Removable Inner Tray, Cleaning Brush, (4) 4x4 Gauze Sponges, (4) Pipe Cleaners, Vinyl Gloves, (2) Cotton-Tip Applicators, Twill Tape, Trach Dressing, Waterproof Drape, 20/cs



713



715

BUSSE TRACHEOSTOMY CARE KIT

Sterile. Includes all necessary components (except cleansing solutions) for care and cleaning of a tracheostomy cannula, sequentially packed in order of use. Each kit contains: 1 3-Compartment Plastic Tray, 1 Pair Gloves, 1 Cleaning Brush, 4 Gauze Sponges, 4" x 4", 4 Large Pipe Cleaners, 2 Cotton-Tipped Applicators, 1 Cotton Twill Tape, 1 tracheostomy Dressing, 1 Large Polylined Drape, 16" x 18".

713 Tracheostomy Kit w/removable basin, 20/cs

715 Tracheostomy Kit, 20/cs



BUSSE TRACHEOSTOMY CARE SET


Sterile. Set includes: tray w/2 compartments; 1 pair gloves; 1 pre-cut gauze tracheostomy dressing; 1 polylined towel/drape; Four 4"x4" gauze sponges; 2 pipe cleaners; 3 cotton tipped applicators; flexible nylon bristle brush; 30" twill tape; 1 foam neck pad; 2 100ml packets hydrogen peroxide; 1 100ml container sterile saline solution.

800 Tracheostomy Care Set w/hydrogen peroxide, 1 Pair Gloves, 1 Polylined Drape, 20/cs

BUSSE TRACHEOSTOMY CARE SET WITH DRESSING

Sterile. Set includes: 3-compartment tray; 1 pre-cut gauze tracheostomy dressing; 2 polylined towel/draper; 4 4"x4" gauze sponges; 6 pipe cleaners; 2 cotton tipped applicators; flexible nylon bristle brush; 30" twill tape.

710 Tracheostomy Care Set w/pre-cut gauze dressing, 24/cs



BUSSE TRACHEOSTOMY CARE SET WITH FORCEPS

Sterile. Set includes: 3-compartment tray; 1 pair gloves; 1 pre-cut gauze tracheostomy dressing; 2 polylined towel/drapes; 1 pair forceps; 3 4"x4" gauze sponges; 3 pipe cleaners; 3 cotton tipped applicators; flexible nylon bristle brush; 30" twill tape.

- 703 Tracheostomy Care Set w/forceps & gauze dressing, 20/cs
799 Tracheostomy Care Set w/out gloves, 20/cs



BUSSE TRACHEOSTOMY CARE SET WITH GLOVES

Sterile. Set includes: 3 compartment tray; 1 pair gloves; 1 pre-cut gauze tracheostomy dressing; 2 polylined towel/drapes; 4 4"x4" gauze sponges, 6 pipe cleaners, 2 cotton tipped applicators; flexible nylon bristle brush; 30" twill tape.

- 711 Tracheostomy Care Set w/extra towel & 2 extra pipe cleaners, 24/cs



BUSSE TRACHEOSTOMY CARE SET WITH HYDROGEN PEROXIDE

Each kit contains: 1 3-Compartment Plastic Tray, 1 Pair of Gloves, 1 Tracheostomy Dressing, 1 Polylined Sterile Field, 4 Gauze Dressings, 4" x 4", 2 Pipe Cleaners, 2 Cotton-Tipped Applicators, 1 Cleaning Brush, 1 Twill Tape, 2 Packets of Hydrogen Peroxide, 1 oz.

- 714 Tracheostomy Care Set w/hydrogen peroxide, 20/cs



COVIDIEN/KENDALL TRACHEOSTOMY CARE TRAYS

Sterile, single-use. Compact and disposable. All necessary components for care and cleaning of trach site.

- 42201 Tray Includes: (2) Blue Nitrile Gloves, CSR Wrap, (3) Solution Basins, 6" Nylon Bristle Trach Tube Brush, Drape, (4) 4"x4" 8-Ply Gauze Sponges, (3) Pipe Cleaners, (3) 4"x4" Dressings, 30" Twill Tape, (3) 6" Cotton Tipped Applicators, 20/cs
47800 Standard Trach Care Tray w/removable basin, Includes: (2) Blue Nitrile Gloves, Trach Brush, Drape, 36" Twill Tape, (2) Cotton Tipped Applicators, (2) Pipe Cleaners, (4) 4" x 4" Gauze Sponges, 20/cs
47802 Standard Trach Care Tray w/bottle of sterile saline, Includes: (2) Blue Nitrile Glove, (1) Trach Sponge, (2) Packs Hydrogen Peroxide, (1) Drape, No Removable Basin, 24/cs
47815 Standard Trach Care Tray w/14 FR Suction Catheter, 20 trays/cs
47835 Standard Trach Care Tray w/plastic forceps, 20 trays/cs
47885 Standard Trach Care Tray w/hydrogen peroxide, 20 trays/cs
47890 Economy Trach Care Kit, 20 trays/cs
47892 Soft Pack Trach Care Kit, 50/cs
49000 Velcro Tracheostomy Tube Holder, Non-Sterile, 100/cs

RUSCH TRACH CARE KIT

- 47100 Trach Care Kit, 20/bx (on contract)

MEDICAL ACTION GENT-L-KARE® TRACHEOSTOMY CARE TRAYS

Designed for cleansing and dressing tracheostomies. All supplies are packaged in order of use to promote sterile technique. All trays contain a waterproof drape, (2) cotton-tipped applicators, tracheostomy brush, pipe cleaners, (4) 4" x 4" gauze, tracheostomy dressing, 1/2" x 36" twill tape.

- 4118 Deluxe Tracheostomy Care Tray w/gloves, 24/cs
4120 Deluxe Tracheostomy Care Tray w/gloves, (2) Waterproof Drapes & Extra Pipe Cleaners, 24/cs
4128 Tracheostomy Care Tray w/gloves & Hydrogen Peroxide, 20/cs
4129 Tracheostomy Care Tray w/gloves, Hydrogen Peroxide, Saline & Tyvek Lid, 20/cs



BUSSE SUTURING KIT WITH SATIN INSTRUMENTS

Each tray contains: 1 Halsey Needle Holder, Fine Point, 5"; 1 Iris Scissors, Straight; 1 Adson Forceps 1" x 2", Fine Point, Satin; 2 2 oz Medicine Cups; 5 2" x 2" Gauze Sponges, 8 ply; 8 4" x 4" Gauze Sponges, 8 ply; 1 18" x 26" Fenestrated Drape; 1 Absorbent Towel, 19 1/2" x 23", Blue. Latex free.

- 747 Suturing Kit w/satin instruments, 20/cs

Quality supplies at a reasonable cost.

Wound Care Trays



DYNAREX SUTURE REMOVAL KIT

Sterile in plastic tray. Contains: Littauer Scissors, 4" Metal Forceps, Gauze Sponge.

4521 Suture Removal Kit, 50/cs

MEDICAL ACTION GENT-L-KARE® LACERATION TRAYS

Assembled with quality components, perfect for wound closures. Arranged in order of use. Select from satin finished or mirror finished stainless steel instruments. All trays are CSR wrapped. Trays contain: 4" x 4" gauze, 2" x 2" gauze, 3 different size needles, 12cc syringe, 2 oz medicine cup, Webster needle holder, Iris scissors, Mosquito forceps, Adson forceps, fenestrated drape, towel/drape.

- 2681 Laceration Tray w/mirror finish instruments & Paper Towel/Drape, 20/cs
- 2683 Laceration Tray w/satin finish instruments & Cloth Towel/Drape, 20/cs
- 2687 Laceration Tray w/mirror finish instruments & Cloth Towel/Drape, 20/cs

MEDICAL ACTION GENT-L-KARE® STERILE SUTURE REMOVAL KITS

Available with wire formed or stainless steel instruments. Conveniently packed in a peel back tray.

- 2632 Suture Removal Kit Tray Includes: (1) SS Iris Scissor, (1) SS Thumb Forceps, (1) Alcohol Prep Pad, (1) 2" x 2" Gauze, 50/cs
- 4131 Suture Removal Kit Tray Includes: (1) WF Littauer Scissors, (1) Metal Forceps, (1) 2" x 2" Gauze, 50/cs

COVIDIEN/KENDALL IRRIGATION TRAYS

Tray includes: 500cc graduated container, waterproof drape, alcohol prep, protector cap.

- 3685 60cc Piston Syringe Irrigation Tray, (Mini-Tray), 20/cs
- 67800 60cc Bulb Syringe 1200cc Irrigation Tray, Lidded, 20/cs
- 67803 60cc Bulb Syringe 1200cc Irrigation Tray, w/Iodophor PVP prep pad, CSR Wrapped, 20/cs
- 68800 60cc Piston Syringe 1200cc Irrigation Tray w/lid, 20/cs
- 68803 60cc Piston Syringe 1200cc Tray w/Iodophor PVP Prep Pad, CSR Wrapped, 20/cs
- 68820 60cc Piston Syringe 1200cc Tray w/100mL Sterile Saline, Lidded, 20/cs



AMSINO AMSURE® ENTERAL FEEDING/IRRIGATION KITS & TRAYS

Latex free.

- AS121 Bulb Irrigation Kit Includes: 500cc Graduated Container, 60cc Bulb Irrigation Syringe, Patient ID Label, Small Tube Adapter, Packaged in a Resealable IV Pole Bag, 30/cs
- AS126 Piston Irrigation Kit Includes: 500cc Graduated Container, 60cc Thumb Control Ring Syringe, Patient ID Label, Small Tube Adapter, Packaged in a Dust Cover, 30/cs
- AS127 Piston Irrigation Kit Includes: 500cc Graduated Container, 60cc Thumb Control Ring Syringe, Patient ID Label, Small Tube Adapter, Packaged in a Resealable IV Pole Bag, 30/cs



AS130



AS136

AMSINO AMSURE® STERILE IRRIGATION TRAY

Latex free.

- AS130 Bulb Irrigation Tray Includes: 1000cc Outer Tray, 500cc Graduated Container, 60cc Bulb Irrigation Syringe, Alcohol Prep Pad, Large Moisture-Proof Underpad, 20/cs
- AS136 Piston Irrigation Tray Includes 1000cc Outer Tray, 500cc Graduated Container, 60cc Thumb Control Ring Syringe, Alcohol Prep Pad, Large Moisture-Proof Underpad, 20/cs

MEDICAL ACTION GENT-L-KARE® IRRIGATION TRAYS

Ideal for irrigating wounds, catheters, etc. All trays contain a 60 cc syringe, graduated container, tip protector, waterproof drape and alcohol prep pad. Graduated container is marked in 50 cc (ml) and 1 oz increments up to 500 cc's.

- 4093 Irrigation Tray w/piston syringe, 20/cs



1581



2081

3M™ COBAN™ SELF-ADHERENT WRAP

Self-adherent elastic wrap functions like a tape, but sticks only to itself. Can be used to secure, protect, compress and support. Sticks to itself; no pins or clips required. Also available with or without latex.


1581	1" x 5 yd Self-Adherent Wrap, Tan, 5/pk, 6 pk/cs
1581B	1" x 5 yd Self-Adherent Wrap, Blue, 30/cs
1582	2" x 5 yd Self-Adherent Wrap, Tan, 36/cs
1582B	2" x 5 yd Self-Adherent Wrap, Blue, 36/cs
1583	3" x 5 yd Self-Adherent Wrap, Tan, 24/cs
1583A	3" x 5 yd Self-Adherent Wrap, Assortment Pack, 2 rolls of Red, Blue, Green, Yellow, White, Tan, 12/cs
1583B	3" x 5 yd Self-Adherent Wrap, Blue, 24/cs
1583G	3" x 5 yd Self-Adherent Wrap, Green, 24/cs
1583N	3" x 5 yd Self-Adherent Wrap, Neon Assortment Pack, 2 rolls of Neon Pink, Neon Green, Neon Orange, Purple, Black, White, 12/cs
1583R	3" x 5 yd Self-Adherent Wrap, Red, 24/cs
1583S	3" x 5 yd Self-Adherent Wrap, Tan, Sterile, 24/cs
1583V	3" x 5 yd Self-Adherent Wrap, White, 24/cs
1583Y	3" x 5 yd Self-Adherent Wrap, Yellow, 24/cs
1584	4" x 5 yd Self-Adherent Wrap, Tan, 18/cs
1584B	4" x 5 yd Self-Adherent Wrap, Blue, 18/cs
1584L	4" x 6.5 yd Self-Adherent Wrap (longer length), Tan, 18/cs
1584S	4" x 5 yd Self-Adherent Wrap, Tan, Sterile, 18/cs
1586	6" x 5 yd Self-Adherent Wrap, Tan, 12/cs
1586S	6" x 5 yd Self-Adherent Wrap, Tan, Sterile, 12/cs
2081	1" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Non-Sterile, 5/pk, 6 pk/cs
2082	2" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Non-Sterile, 36/cs
2083	3" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Non-Sterile, 24/cs
2083S	3" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Sterile, 24/cs
2084	4" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Non-Sterile, 18/cs
2084L	4" x 6.5 yd Self-Adherent Wrap (longer length), Latex Free, Tan, Non-Sterile, 18/cs
2084S	4" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Sterile, 18/cs
2086	6" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Non-Sterile, 12/cs
2086S	6" x 5 yd Self-Adherent Wrap, Latex Free, Tan, Sterile, 12/cs

Fast, easy access to the products you need.



HARTMANN-CONCO CO-LASTIC® COHESIVE ELASTIC BANDAGES - LATEX FREE

45110000	Bandage, Cohesive, Elastic, 1.5" x 5 yd., 8 Assorted Colors, Latex Free, 48/cs
45120000	Bandage, Cohesive, Elastic, 1" x 5 yd., Tan, Latex Free, 30/cs
45150000	Bandage, Cohesive Elastic, 1.5" x 5 yd., Tan, Latex Free, 48/cs
45200000	Bandage, Cohesive, Elastic, 2" x 5 yd., Tan, Latex Free, 36/cs
45300000	Bandage, Cohesive, Elastic, 3" x 5 yd., Tan, Latex Free, 24/cs
45400000	Bandage, Cohesive, Elastic, 4" x 5 yd., Tan, Latex Free, 18/cs
45600000	Bandage, Cohesive, Elastic, 6" x 5 yd., Tan, Latex Free, 12/cs




DUKAL COHESIVE BANDAGES

Can be easily torn by hand and are an excellent product when additional support is needed. This durable self adhesive wrap is well suited for use in a variety of applications and does not require clips or pins.

8015AS	Bandage, Cohesive, 1", Non-Sterile, Assorted Colors, 5 yds/rl, 30 rl/cs
8015DB	Bandage, Cohesive, 1", Non-Sterile, Dark Blue, 5 yds/rl, 30 rl/cs
8015T	Bandage, Cohesive, 1", Non-Sterile, Tan, 5 yds/rl, 30 rl/cs
8015W	Bandage, Cohesive, 1", Non-Sterile, White, 5 yds/rl, 30 rl/cs
8025AS	Bandage, Cohesive, 2", Non-Sterile, Assorted Colors, 5 yds/rl, 36 rl/cs
8025DB	Bandage, Cohesive, 2", Non-Sterile, Dark Blue, 5 yds/rl, 36 rl/cs
8025T	Bandage, Cohesive, 2", Non-Sterile, Tan, 5 yds/rl, 36 rl/cs
8035AS	Bandage, Cohesive, 3", Non-Sterile, Assorted Colors, 5 yds/rl, 24 rl/cs
8035DB	Bandage, Cohesive, 3", Non-Sterile, Dark Blue, 5 yds/rl, 24 rl/cs
8035T	Bandage, Cohesive, 3", Non-Sterile, Tan, 5 yds/rl, 24 rl/cs
8045AS	Bandage, Cohesive, 4", Non-Sterile, Assorted, 5 yds/rl, 18 rl/cs
8045T	Bandage, Cohesive, 4", Non-Sterile, Tan, 5 yds/rl, 18 rl/cs
8155AS	Bandage, Cohesive, 1.5", Non-Sterile, Assorted Colors, 5 yds/rl, 48 rl/cs
8155DB	Bandage, Cohesive, 1.5", Non-Sterile, Dark Blue, 5 yds/rl, 48 rl/cs
8155R	Bandage, Cohesive, 1.5", Non-Sterile, Red, 5 yds/rl, 48 rl/cs
8155T	Bandage, Cohesive, 1.5", Non-Sterile, Tan, 5 yds/rl, 48 rl/cs
8155W	Bandage, Cohesive, 1.5", Non-Sterile, White, 5 yds/rl, 48 rl/cs

Elastic Bandages



AMD-RITMED COHESIVE BANDAGES

A highly adaptable "no-slip" bandage that sticks to itself but not to skin or other materials. Assorted colors to meet all your patient's needs. Easy one-hand tear ability. No scissors needed. Convenient dispenser box for easy storage and distribution. Medium, but adjustable, compression. Also available in blue, red, white and assorted (rainbow).

A1061-T	1" x 5 yds, Non-Sterile, Tan, Individually Wrapped, 72/bx
A1061-TLF	1" x 5 yds, Non-Sterile, Tan, Individually Wrapped, Latex Free, 72/bx
A2061-LF	2" x 5 yds, Non-Sterile, Tan, Individually Wrapped, Latex Free, 36/bx
A2061-T	2" x 5 yds, Non-Sterile, Tan, Individually Wrapped, 36/bx
A3061-LF	3" x 5 yds, Non-Sterile, Tan, Individually Wrapped, Latex Free, 24/bx
A3061-T	3" x 5 yds, Non-Sterile, Tan, Individually Wrapped, 24/bx
A4061-LF	4" x 5 yds, Non-Sterile, Tan, Individually Wrapped, Latex Free, 18/bx
A4061-T	4" x 5 yds, Non-Sterile, Tan, Individually Wrapped, 18/bx
A6061-LF	6" x 5 yds, Non-Sterile, Tan, Individually Wrapped, Latex Free, 12/bx
A6061-T	6" x 5 yds, Non-Sterile, Tan, Individually Wrapped, 12/bx

**Your sales representative
can assist you with
product information and
current pricing.**



COVIDIEN/KENDALL COHESIVE BANDAGES TAPE

Self-adherent wraps made with a non-woven material and elastic fibers placed lengthwise provide elasticity. Elastic wrap contains a cohesive material allowing it to stick to itself but not to other materials or skin. Many options available latex and latex-free, sterile and non-sterile, and a variety of widths and colors.

4581BC	Cohesive Bandages Tape, 1" x 5 yds, Bulk, Blue, 30 rl/cs
4581C	Cohesive Bandages Tape, 1" x 5 yds, Bulk, Tan, 30 rl/cs
4581LF	Cohesive Bandages Tape, 1" x 5 yds, Latex Free, Non-Sterile, 5/pk 6 pk/cs
4582BC	Cohesive Bandages Tape, 2" x 5 yds, Bulk, Blue, 30 rl/cs
4582BLF	Cohesive Bandages Tape, 2" x 5 yds, Latex Free, Non-Sterile, Black, 36/cs
4582C	Cohesive Bandages Tape, 2" x 5 yds, Bulk, Tan, 30 rl/cs
4582LF	Cohesive Bandages Tape, 1" x 5 yds, Latex Free, Non-Sterile, 36/cs
4582WLF	Cohesive Bandages Tape, 2" x 5 yds, Latex Free, Non-Sterile, White, 36/cs
4583A	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Color Assortment, 12/cs
4583BC	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Blue, 24 rl/cs
4583BLF	Cohesive Bandages Tape, 3" x 5 yds, Latex Free, Non-Sterile, Black, 24/cs
4583C	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Tan, 24 rl/cs
4583G	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Green, 24/cs
4583LF	Cohesive Bandages Tape, 3" x 5 yds, Latex Free, Non-Sterile, 24/cs
4583N	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Neon, 12/cs
4583R	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Red, 24/cs
4583SC	Cohesive Bandages Tape, 3" x 5 yds, Tan, 24 rl/cs
4583SLF	Cohesive Bandages Tape, 3" x 5 yds, Latex Free, Sterile, 24/cs
4583W	Cohesive Bandages Tape, 3" x 5 yds, Bulk, White, 24/cs
4583WLF	Cohesive Bandages Tape, 3" x 5 yds, Latex Free, Non-Sterile, White, 24/cs
4583Y	Cohesive Bandages Tape, 3" x 5 yds, Bulk, Yellow, 24/cs
4584BC	Cohesive Bandages Tape, 4" x 5 yds, Bulk, Blue, 18 rl/cs
4584C	Cohesive Bandages Tape, 4" x 5 yds, Bulk, Tan, 18 rl/cs
4584LC	Cohesive Bandages, 4" x 6.5 yds, Bulk, Tan, 18/cs
4584LF	Cohesive Bandages Tape, 4" x 5 yds, Latex Free, Non-Sterile, 18/cs
4584LLF	Cohesive Bandages Tape, 4" x 6.55 yds, Latex Free, Non-Sterile, 18/cs
4584SC	Cohesive Bandages Tape, 4" x 5 yds, Tan, 18 rl/cs
4584SLF	Cohesive Bandages Tape, 4" x 5 yds, Latex Free, Sterile, 18/cs
4586C	Cohesive Bandages Tape, 6" x 5 yds, Tan, 12 rl/cs
4586LF	Cohesive Bandages Tape, 6" x 5 yds, Latex Free, Non-Sterile, 12/cs
4586SC	Cohesive Bandages Tape, 6" x 5 yds, Tan, 12 rl/cs
4586SLF	Cohesive Bandages Tape, 6" x 5 yds, Latex Free, Sterile, 12/cs



DYNAREX SENSI WRAP, SELF ADHERENT BANDAGE ROLLS

A lightweight compression bandage which sticks to itself, but not to other materials or skin. Easily torn without scissors. Will not slip. Individually poly bagged. Available in tan or in an assorted rainbow pack.

- 3171 Sensi-Wrap, Self-Adherent, 1" x 5yds, Tan, Contains Latex, 2/bg, 30bg/cs
- 3172 Sensi-Wrap, Self-Adherent, 2" x 5yds, Tan, Contains Latex, 1/bg, 36bg/cs
- 3173 Sensi-Wrap, Self-Adherent, 3" x 5yds, Tan, Contains Latex, 1/bg, 24bg/cs
- 3174 Sensi-Wrap, Self-Adherent, 4" x 5yds, Tan, Contains Latex, 1/bg, 18bg/cs
- 3176 Sensi-Wrap, Self-Adherent, 6" x 5yds, Tan, Contains Latex, 12/cs
- 3181 Rainbow Packs, 1" x 5yds, Rainbow, 5/color, Contains Latex, 30/cs
- 3182 Rainbow Packs, 2" x 5yds, Rainbow, 5/color, Contains Latex, 36/cs
- 3183 Rainbow Packs, 3" x 5yds, Rainbow, 5/color, Contains Latex, 24/cs
- 3184 Rainbow Packs, 4" x 5yds, Rainbow, 5/color, Contains Latex, 18/cs
- 3186 Rainbow Packs, 6" x 5yds, Rainbow, 5/color, Contains Latex, 12/cs
- 3187 Sensi-Wrap, Self-Adherent, 1" x 5yds, Tan, Latex-Free, 2/bg, 30 bg/cs
- 3188 Sensi-Wrap, Self-Adherent, 2" x 5yds, Tan, Latex-Free, 36/cs
- 3189 Sensi-Wrap, Self-Adherent, 3" x 5yds, Tan, Latex-Free, 24/cs
- 3190 Sensi-Wrap, Self-Adherent, 4" x 5yds, Tan, Latex-Free, 18/cs
- 3191 Sensi-Wrap, Self-Adherent, 6" x 5yds, Tan, Latex-Free, 12/cs



44110000



44120000

HARTMANN-CONCO CO-LASTIC® LF LATEX FREE COHESIVE ELASTIC BANDAGES

High quality, cohesive elastic wrap, versatile enough to use as a securing, protective compression and support bandage. Designed with both you and your patients in mind. It sticks securely to itself, but never to your patients' hair or skin. Provides even, controlled compression and readily conforms to every anatomical structure. It is water resistant, lightweight and comfortable. Hand-tears to desired length and has a folded starting edge.

- 45210000 Bandage, Cohesive, Elastic, 2" x 5 yd., 8 Assorted Colors, Latex Free, 36/cs
- 45310000 Bandage, Cohesive, Elastic, 3" x 5 yd., 8 Assorted Colors, Latex Free, 24/cs



HARTMANN-CONCO CO-WRAP® SELF-ADHERENT SECURING WRAP

A unique blend of cotton and polymer fibers that adheres securely to itself - not skin, hair or clothing. Its special construction also makes it porous so it's cool and comfortable to wear no matter how long it's worn. An easy-to-use, versatile and is effective even where sensitive tissue or wet environments are a consideration. Ideal as a secondary or securing bandage and for overwrapping primary dressings. Unlike bandages that require clips, tapes or pins, it won't shift and slide. Ideal for delicate or traumatized skin when tape cannot be used. CAUTION: This product contains natural rubber latex which may cause allergic reactions.

- 29100000 1" x 5 yds Bandage, 24/bx, 4 bx/cs
- 29150000 1 1/2" x 5 yds Bandage, Bulk, 125/cs
- 29200000 2" x 5 yds Bandage, 12/bx, 8 bx/cs
- 29300000 3" x 5 yds Bandage, 12/bx, 8 bx/cs
- 29400000 4" x 5 yds Bandage, 12/bx, 8 bx/cs



HARTMANN-CONCO MEDI-RIP® SELF-ADHERENT BANDAGE

Self-adherent elastic support and compression bandage. Woven of 99% cotton yarn, it is porous, lightweight and cool. Sticks to itself, but not to skin or hair. Conforms readily to body contours for ease of application and maximum comfort. Can be torn easily to exact length needed. CAUTION: This product contains natural rubber latex which may cause allergic reactions.

- 25100000 1" x 5 yds Bandage, Non-Sterile, Tan, 24 rl/bx, 4 bx/cs
- 25200000 2" x 5 yds Bandage, Non-Sterile, Tan, 12 rl/bx, 8 bx/cs
- 25300000 3" x 5 yds Bandage, Non-Sterile, Tan, 12 rl/bx, 8 bx/cs
- 25320000 2" x 5 yds Bandage, Non-Sterile, Blue, 12/bx, 8 bx/cs
- 25330000 3" x 5 yds Bandage, Non-Sterile, Blue, 12/bx, 8 bx/cs
- 25400000 4" x 5 yds Bandage, Non-Sterile, Tan, 12 rl/bx, 8 bx/cs
- 25600000 6" x 5 yds Bandage, Non-Sterile, Tan, 12 rl/bx, 4 bx/cs
- 25810000 1 1/2" x 5 yds Bandage, Non-Sterile, Team Pack, Tan, 32/cs
- 25820000 2" x 5 yds Bandage, Non-Sterile, Team Pack, Tan, 24/cs
- 25830000 3" x 5 yds Bandage, Non-Sterile, Team Pack, Tan, 16/cs
- 70005599 1 1/2" x 5 yds Bandage, Non-Sterile, Tan, 250/cs



SOUTHWEST FINGER BOBS

An elastic cotton cloth covering comprised of two rolled bandage rings. It is ready to apply so there is no need for tape, scissors or applicator. Finger Bobs stay firmly in place, are comfortable to wear and are ideal for holding a dressing in place. Ideal for protecting fingers and toes! Available in 6 assorted colors per package.

- FB7605 Medium Finger Bobs, 50/bg
- FB7610 Large Finger Bobs, 25/bg

Elastic Bandages

AMBRA LE ROY ORTHOPEDIC ELASTIC BANDAGE

A superior quality, non-sterile woven compression bandage. It is manufactured with cotton, polyester and latex-free elastic yarns. The densely woven construction maximizes comfort and durability, making the elastic bandage resistant to harsh medications and ointments, and allows it to stand up to frequent laundering and sterilization. Ideal for orthopedic applications and also useful for the support of prosthetic devices, reduction fractures, and support of limbs post-injury or surgery. It is recommended for strains and sprains, and for the treatment of tendonitis, tenosynovitis, varicosities, edema and poor peripheral circulation. It comes with elastic stretch clips.

- 71250 Premium Elastic Bandage, Orthopedic, 2" x 5 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs
- 71350 Premium Elastic Bandage, Orthopedic, 3" x 5 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs
- 71410 Premium Elastic Bandage, Orthopedic, 4" x 10 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs
- 71450 Premium Elastic Bandage, Orthopedic, 4" x 5 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs
- 71610 Premium Elastic Bandage, Orthopedic, 6" x 10 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs
- 71650 Premium Elastic Bandage, Orthopedic, 6" x 5 yds (Stretched) w/elastic stretch clips, Beige, Latex-Free, 10/bx, 5 bx/cs



AMBRA LE ROY MIDLASTIC ELASTIC BANDAGE

A high quality, non-sterile woven compression bandage. It is manufactured with cotton, polyester and latex-free elastic yarns. The controlled stretch in this reinforced elastic bandage makes it an excellent choice for applications requiring firm, consistent tension. It is a good choice for securing dressings, splints and ice packs and also provides good compression and support in the treatment of edema, sprains, and injured ligament. It comes with elastic stretch clips.

- 72250 Premium Elastic Bandage, 2" x 5 yds (Stretched), w/elastic stretch clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 72350 Premium Elastic Bandage, 3" x 5 yds (Stretched), w/elastic stretch clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 72450 Premium Elastic Bandage, 4" x 5 yds (Stretched), w/elastic stretch clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 72650 Premium Elastic Bandage, 6" x 5 yds (Stretched), w/elastic stretch clips, Tan, Latex-Free, 10/bx, 5 bx/cs



AMBRA LE ROY SUPRA-GRIP ELASTIC BANDAGE

A premium quality, non-sterile knitted compression bandage. The controlled stretch provides support that is consistent and firm. No clips are necessary to fasten this latex-free elastic bandage, which is secured with a 1" Velcro brand start strip and 2" Velcro closure strip. Use this bandage where the convenience and efficiency of Velcro closure is required along with the benefit of adjustable compression levels.

- 90250 Premium Elastic Bandage, 2" x 5 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs
- 90350 Premium Elastic Bandage, 3" x 5 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs
- 90411 Premium Elastic Bandage, 4" x 11 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs
- 90455 Premium Elastic Bandage, 4" x 5 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs
- 90611 Premium Elastic Bandage, 6" x 11 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs
- 90655 Premium Elastic Bandage, 6" x 5 yds (Stretched), w/double Velcro closure, Tan/White Honeycomb, Latex-Free, 10/bx, 5 bx/cs



AMBRA LE ROY VALUELASTIC ELASTIC BANDAGE

A good, general-use non-sterile woven compression bandage. Manufactured with cotton, polyester and latex-free elastic yarns, the elastic bandage is lightweight and economical because it can be laundered and sterilized. It secures with aluminum clips and is durable and comfortable, with good stretch, compression and balance.

- 73250 Economy Elastic Bandage, 2" x 5 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 73350 Economy Elastic Bandage, 3" x 5 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 73410 Economy Elastic Bandage, 4" x 10 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 73450 Economy Elastic Bandage, 4" x 5 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 73610 Economy Elastic Bandage, 6" x 10 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs
- 73650 Economy Elastic Bandage, 6" x 5 yds (Stretched), w/standard clips, Tan, Latex-Free, 10/bx, 5 bx/cs

COVIDIEN/KENDALL TENSOR ELASTIC BANDAGES W/REMOVABLE CLIPS

Nonsterile, individually wrapped. Ideal for pressure dressings, vascular and muscular support. Unique crocheted construction allows extra air circulation and superior breathability. Provides firm compression and comfortable uniform support. Detachable clip provides convenient, exact positioning.

4202	Elastic Bandage w/removable clips, 2" x 4.1 yds., 12/bg, 12 bg/cs
4203	Elastic Bandage w/removable clips, 3" x 4.1 yds., 12/bg, 12 bg/cs
4204	Elastic Bandage w/removable clips, 4" x 4.5 yds., 12/bg, 12 bg/cs
4206	Elastic Bandage w/removable clips, 6" x 4.5 yds., 12/bg, 6 bg/cs

BD ACE® BANDAGE CLIPS

207406	Bandage Clips, 12/bx, 12 bx/cs
--------	--------------------------------

BD ACE® BRAND ELASTIC BANDAGES

207336	2" Rubber Elastic Bandages, 10/bx, 5 bx/cs
207337	3" Rubber Elastic Bandages, 10/bx, 5 bx/cs
207338	4" Rubber Elastic Bandages, 10/bx, 5 bx/cs
207339	6" Rubber Elastic Bandages, 10/bx, 5 bx/cs
207430	2" Elastic Bandages, 10/bx, 5 bx/cs
207432	3" Elastic Bandages, 10/bx, 5 bx/cs
207433	4" Elastic Bandages, 10/bx, 5 bx/cs
207435	6" Elastic Bandages, 10/bx, 5 bx/cs

BD PEG® SELF-ADHERING ELASTIC BANDAGES

207419	1" Elastic Bandages, 20/bx, 5 bx/cs
207420	2" Elastic Bandages, 10/bx, 5 bx/cs
207422	3" Elastic Bandages, 10/bx, 5 bx/cs
207423	4" Elastic Bandages, 10/bx, 5 bx/cs
207425	6" Elastic Bandages, 10/bx, 5 bx/cs

DERMA SCIENCES FLEXOPLAST® ELASTIC ADHESIVE BANDAGE

An all elastic adhesive bandage that provides compression. Adapts easily to contours of the body. Used in combination with Primer® which effectively manages venous leg ulcers. Also provides support for fractures and muscle injuries. Sticks to itself and will not slip, bunch up or loosen.

GL-100-3	3" x 5 1/2 yds Elastic Bandage, 12/cs
GL100-1	1" x 5 yds (approx) Stretched Adhesive Bandage, 12/bx
GL100-2	2" x 5 yds (approx) Stretched Adhesive Bandage, 6/bx
GL100-3	3" x 5 yds (approx) Stretched Adhesive Bandage, 4/bx
GL100-3R	3" x 5 yds (approx) Stretched Adhesive Bandage (reversed adhesion), 12/bx
GL100-4	4" x 5 yds (approx) Stretched Adhesive Bandage, 6/bx



650-6062-0043



650-6062-0044



650-6066-0046

DMI® RUBBER REINFORCED ELASTIC BANDAGES

Maximum support with constant even pressure. Bound edges prevents fraying. Washable.

650-6062-0043	3" x 5 yds Bandage, 12/cs
650-6062-0044	4" x 5 yds Bandage, 12/cs
650-6066-0046	6" x 5 yds Bandage, 12/cs

DUKAL ELASTIC BANDAGE

Provides consistent support and are a good product to use when even compression is needed. All rolls are individually wrapped for clean storage.

502	Elastic Bandage, 2", 10/bx, 5 bx/cs
502LF	Elastic Bandage, 2", Latex Free, 10/bx, 5 bx/cs
503	Elastic Bandage, 3", 10/bx, 5 bx/cs
503LF	Elastic Bandage, 3", Latex Free, 10/bx, 5 bx/cs
504	Elastic Bandage, 4", 10/bx, 5 bx/cs
504LF	Elastic Bandage, 4", Latex Free, 10/bx, 5 bx/cs
506	Elastic Bandage, 6", 10/bx, 5 bx/cs
506LF	Elastic Bandage, 6", Latex Free, 10/bx, 5 bx/cs
512LF	Elastic Bandage, 2", Latex Free, Double Velcro, 10/bx, 5 bx/cs
513LF	Elastic Bandage, 3", Latex Free, Double Velcro, 10/bx, 5 bx/cs
514LF	Elastic Bandage, 4", Latex Free, Double Velcro, 10/bx, 5 bx/cs
516LF	Elastic Bandage, 6", Latex Free, Double Velcro, 10/bx, 5 bx/cs



DYNAREX ELASTIC BANDAGES, LATEX FREE

A 100% latex free elastic bandage, providing outstanding compression and stretch. Individually poly-wrapped with 2 clips. 4.5 yards, stretched.

3662	2", 10/bx, 5bx/cs
3663	3", 10/bx, 5bx/cs
3664	4", 10/bx, 5bx/cs
3666	6", 10/bx, 5bx/cs

GRAHAM FIELD ELASTIC BANDAGE

Manufactured from cotton/poly blend and heat-resistant latex rubber. Non-sterile, individually wrapped 5 yard (approx.) stretched length in convenient dispenser box.

1971E-2	2"
1971E-3	3"
1971E-4	4"
1971E-6	6"



HARTMANN-CONCO CEB® LATEX FREE COTTON ELASTIC BANDAGES

The premium choice when 100% cotton elastic bandages are required. High-quality cotton yarns are woven to form a bandage which is porous, hypo-allergenic elastic with a non-raveling feathered edge. Rubberless weave reduces risk of constriction. Sterilizable and reusable. Withstands washing without affecting elasticity. 5 yard length. Latex free.

11200000	CEB, 2" x 5 yd Bandage, 10/pk, 6 pk/cs
11300000	CEB, 3" x 5 yd Bandage, 10/pk, 6 pk/cs
11400000	CEB, 4" x 5 yd Bandage, 10/pk, 6 pk/cs
11600000	CEB, 6" x 5 yd Bandage, 10/pk, 6 pk/cs

Elastic Bandages



HARTMANN-CONCO COMPERM® LF TUBULAR ELASTIC BANDAGES

Provides 360° of uniform compression and support with minimal risk of constriction. The 100% cotton-knit-tubular construction reinforced with latex-free elastomers creates a bandage which is easy to apply and comfortable to wear. Stays in place without clips or tape. Compression and support remain consistent - even after reapplication. Offers superior support for edema reduction, post-cast support and as coverage for a wound dressing.

83010000	Size A, 1 5/8" x 11 yds Latex Free Tubular Bandage
83020000	Size B, 2 1/2" x 11 yds Latex Free Tubular Bandage
83030000	Size C, 2 3/4" x 11 yds Latex Free Tubular Bandage
83040000	Size D, 3" x 11 yds Latex Free Tubular Bandage
83050000	Size E, 3 1/2" x 11 yds Latex Free Tubular Bandage
83060000	Size F, 4" x 11 yds Latex Free Tubular Bandage
83070000	Size G, 5" x 11 yds Latex Free Tubular Bandage
83080000	Size J, 7" x 11 yds Latex Free Tubular Bandage
83090000	Size K, 8 1/2" x 11 yds Latex Free Tubular Bandage
83100000	Size L, 13" x 11 yds, Latex Free, Tubular Bandage

HARTMANN-CONCO COMPERM® -PD LF LATEX FREE PADDED TUBULAR BANDAGE

83120000	Small, 3" x 13', 1l/cs
83130000	Large, 4" x 13 ft Padded Tubular Bandage



HARTMANN-CONCO CONTEX® LF LATEX FREE REINFORCED ELASTIC BANDAGE

A latex free, premium grade, reinforced elastic bandage offering consistent and dependable tension and firm support. Its high cotton content and latex-free elastomer make it cool and comfortable to wear. For treating sprains, strains, tendonitis, torn or injured ligaments and edema reduction.

50200000	2" x 5 yds Bandage, Sterile, 10/cs
50300000	3" x 5 yds Bandage, Sterile, 10/cs
50400000	4" x 5 yds Bandage, Sterile, 10/cs
50600000	6" x 5 yds Bandage, Sterile, 10/cs
54200000	2" x 5 yds Bandage, 10/pk, 6 pk/cs
54300000	3" x 5 yds Bandage, 10/pk, 6 pk/cs
54400000	4" x 5 yds Bandage, 10/pk, 6 pk/cs
54410000	4" x 10 yds Bandage, 6/cs
54600000	6" x 5 yds Bandage, 10/pk, 6 pk/cs
54610000	6" x 10 yds Bandage, 6/cs



HARTMANN-CONCO DELUXE® LF LATEX FREE ELASTIC BANDAGES

Orthopedic grade, latex free bandage that provides consistent elasticity, support and compression. Made of the finest quality yarns that resist deterioration from ointments and medications and can withstand numerous washings and sterilization. For use treating varicosities, edema, poor circulation, tendonitis, tenosynovitis, sprains, strains and fracture reduction.

27200000	2" x 5.5 yds Elastic Bandage, Sterile, 20/cs
27300000	3" x 5.5 yds Elastic Bandage, Sterile, 20/cs
27400000	4" x 5.5 yds Elastic Bandage, Sterile, 20/cs
27600000	6" x 5.5 yds Elastic Bandage, Sterile, 20/cs
28400000	4" x 11 yds Elastic Bandage, Sterile, 20/cs
28600000	6" x 11 yds Elastic Bandage, Sterile, 20/cs
38200000	2" x 5 1/2 yds Bandage, Deluxe, 10/pk, 6 pk/cs
38300000	3" x 5 1/2 yds Bandage, Deluxe, 10/pk, 6 pk/cs
38400000	4" x 5 1/2 yds Bandage, Deluxe, 10/pk, 6 pk/cs
38410000	4" x 11 yds Bandage, 6/cs
38600000	6" x 5 1/2 yds Bandage, Deluxe, 10/pk, 6 pk/cs
38610000	6" x 11 yds Bandage, 6/cs
39200000	2" x 5 yds Bandage, Deluxe "480", 10/pk, 6 pk/cs
39300000	3" x 5 yds Bandage, Deluxe "480", 10/pk, 6 pk/cs
39400000	4" x 5 yds Bandage, Deluxe "480", 10/pk, 6 pk/cs
39600000	6" x 5 yds Bandage, Deluxe "480", 10/pk, 6 pk/cs



HARTMANN-CONCO ECONO-WRAP® LF REINFORCED ELASTIC BANDAGE

Offers balanced stretch and compression. Lightweight and ideal for general use. Perfect for securing splints, dressings, ice packs and in the treatment of varicosities, sprains and strains. Latex Free.

33200000	2" x 4.5 yd Bandage, 10/pk, 6 pk/cs
33300000	3" x 4.5 yd Bandage, 10/pk, 6 pk/cs
33400000	4" x 4.5 yd Bandage, 10/pk, 6 pk/cs
33600000	6" x 4.5 yd Bandage, 10/pk, 6 pk/cs



HARTMANN-CONCO LATEX FREE CONFORMING STRETCH BANDAGE

Unique weave clings to itself and stretches to conform. Polyester yarns create a strong bandage that provides light compression and conforms to body contours. Woven finished edges minimize lint and loose threads. In multi-layers, weave provides padding and protection. Latex free.

80100000	1" x 4.1 yds Bandage, Non-Sterile, 24/bg, 4 bg/cs
80200000	2" x 4.1 yds Bandage, Non-Sterile, 12/bg, 8 bg/cs
80300000	3" x 4.1 yds Bandage, Non-Sterile, 12/bg, 8 bg/cs
80400000	4" x 4.1 yds Bandage, Non-Sterile, 12/bg, 8 bg/cs
80600000	6" x 4.1 yds Bandage, Non-Sterile, 6/bg, 8 bg/cs
81200000	2" x 4.1 yds Bandage, Sterile, 12/bg, 8 bg/cs
81300000	3" x 4.1 yds Bandage, Sterile, 12/bg, 8 bg/cs
81400000	4" x 4.1 yds Bandage, Sterile, 12/bg, 8 bg/cs
81600000	6" x 4.1 yds Bandage, Sterile, 12/bg, 4 bg/cs



SMITH & NEPHEW PROFORE™ FOUR LAYER BANDAGING SYSTEM

Comprised of a wound contact layer and four bandages, has been specifically designed for the management of venous leg ulcers. Designed to maintain effective levels of compression after application. May be left in place for up to one week. More frequent changes may be necessary especially in the first few weeks of management when the ulcers may be draining heavily. Contains everything necessary for the application of graduated compression therapy. May reduce preparation time, simplifies stock control and ensures the correct components are always available.

66020016 Four Layer Bandaging System, 8/cs

DERMA SCIENCES BANDNET STRETCH BANDAGE

BA2502 Bandage, Stretch, Size 2, 25 yds/bx

DERMA SCIENCES BANDNET ELASTIC NET DRESSING

- BA2501 Small Elastic Dressing for Fingers & Toes, 25yds
- BA2503 Small Elastic Dressing for Knee, Foot, Elbow, and Hand, 25yds
- BA2504 Large Elastic Dressing for Knee, Foot, Elbow, and Hand, 25yds, ea
- BA2505 Small Elastic Dressing for Calf and Thigh, 25yds
- BA2505.5 Pediatric Head, Thigh, Shoulder, Elastic Net Dressing, 25yds, ea
- BA2506 Elastic Dressing for Axilla, Abdomen, Adult Head, and Chest, 25yds, ea
- BA2507 Small Elastic Dressing for Axilla, Abdomen, and Chest, 25yds, ea
- BA2508 Medium Elastic Dressing for Axilla, Abdomen, and Chest, 25yds, ea
- BA2509 Large Elastic Dressing for Axilla, Abdomen, and Chest, 25yds, ea
- BA2510 X-Large Elastic Dressing for Axilla, Abdomen, and Chest, 25yds, ea
- BA2511 Axilla - 250# Plus, Chest, Abdomen, Elastic Net Dressing, 25yds, ea



DERMA SCIENCES SURGIGRIP® TUBULAR ELASTIC SUPPORT BANDAGE - LATEX FREE

A latex-free tubular elastic bandage designed to provide support and compression. Made of a soft cotton knit interlaced with elastic threads. Perfect for sprains and strains, support following cast removal, soft-tissue injuries, joints and muscles, post-graft pressure for burn patients and edema management/control. Convenient alternative to flat wrap bandages. Available in 11 yard (10 meter) lengths and has 8-12mm Hg compression at 50% stretch. Can be applied to sensitive areas using an applicator. Stays in place with no clips or tape. Machine washable. Available in a wide range of sizes.

- GLB10 2 1/2" Wide Tubular Bandage, Small Hands & Limbs
- GLC10 2 3/4" Wide Tubular Bandage, Adult Hands, Arms & Legs
- GLD10 3" Wide Tubular Bandage, Large Arms & Legs
- GLE10 3 1/2" Wide Tubular Bandage, Legs or Small Thighs
- GLF10 4" Wide Tubular Bandage, Large Knees or Thighs
- GLG10 4 1/2" Wide Tubular Bandage, Large Thighs
- GLJ10 6 3/4" Wide Tubular Bandage, Small Trunk
- GLK10 8" Wide Tubular Bandage, Medium Trunk
- GLL10 12 3/4" Wide Tubular Bandage, Large Trunk



GL501



GL600



GL620



GLLF2501



GLLF2507

DERMA SCIENCES SURGILAST® TUBULAR ELASTIC BANDAGE RETAINER

A tubular elastic stretch net designed to serve as a secondary dressing, applying gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings. Easy to customize. Allows maximum air-flow. Available in a wide range of sizes.

- GL10-1 Size 1 Tubular Bandage Retainer, Fingers, Toes, Wrist, 10 yd/ea
- GL10-10 Size 10 Tubular Bandage Retainer, Large Chest, Back, Perineum, Axilla, 10 yd/ea
- GL10-11 Size 11 Tubular Bandage Retainer, X-Large: Chest, Back, Perineum, Axilla, 10 yd/ea
- GL10-12 Size 12 Tubular Bandage Retainer, XX-Large: Chest, Back, Perineum, Axilla, 10 yd/ea
- GL10-2 Size 2 Tubular Bandage Retainer, Small: Hand, Arm, Leg, Foot, 10 yd/bx
- GL10-3 Size 3 Tubular Bandage Retainer, Medium, Hand, Arm, Leg, Foot, 10 yd/ea
- GL10-4 Size 4 Tubular Bandage Retainer, Large: Hand, Arm, Leg, Foot, 10 yd/bx
- GL10-5 Size 5 Tubular Bandage Retainer, Small: Head, Shoulder, Thigh, 10 yd/bx
- GL10-6 Size 6 Tubular Bandage Retainer, Medium: Head, Shoulder, Thigh, 10 yd/bx
- GL10-7 Size 7 Tubular Bandage Retainer, Large: Head, Shoulder, Thigh, 10 yd/bx
- GL10-8 Size 8 Tubular Bandage Retainer, Small: Chest, Back, Perineum, Axilla, 10 yd/bx
- GL10-9 Size 9 Tubular Bandage Retainer, Medium: Chest, Back, Perineum, Axilla, 10 yd/bx
- GL2522 Size 22 Tubular Bandage Retainer, Special Sizing, 25 yd/bx
- GL501 Size 1 Tubular Bandage Retainer, Fingers, Toes, Wrist, 50 yd/bx
- GL502 Size 2 Tubular Bandage Retainer, Small: Hand, Arm, Leg, Foot, 50 yd/bx
- GL503 Size 3 Tubular Bandage Retainer, Medium: Hand, Arm, Leg, Foot, 50 yd/bx
- GL504 Size 4 Tubular Bandage Retainer, Large: Hand, Arm, Leg, Foot, 50 yd/bx
- GL505 Size 5 Tubular Bandage Retainer, Small: Head, Shoulder, Thigh, 50 yd/bx
- GL506 Size 5.5 Tubular Bandage Retainer, Medium: Head, Shoulder, Thigh, 50 yd/bx
- GL507 Size 6 Tubular Bandage Retainer, Large: Head, Shoulder, Thigh, 50 yd/bx
- GL508 Size 7 Tubular Bandage Retainer, Small: Chest, Back, Perineum, Axilla, 50 yd/bx
- GL509 Size 8 Tubular Bandage Retainer, Medium: Chest, Back, Perineum, Axilla, 50 yd/bx
- GL510 Size 9 Tubular Bandage Retainer, Large: Chest, Back, Perineum, Axilla, 50 yd/bx
- GL511 Size 10 Tubular Bandage Retainer, X-Large: Chest, Back, Perineum, Axilla, 50 yd/ea
- GL512 Special Sizing 7 1/2" Tubular Bandage Retainer, 50 yd/ea

Elastic Bandages

DERMA SCIENCES SURGILAST® TUBULAR ELASTIC BANDAGE RETAINER (Continued)

GL600	Head Tubular Bandage Retainer, 20/bx
GL620	Perineum Panty, Small/Medium, 20/bx
GL622	Perineum Panty, Large/X-Large, 20/bx
GL641	Knee, Foot, Elbow, Hand Tubular Bandage Retainer, 20/bx
GL701	Size 1 Tubular Bandage Retainer, Fingers, Toes, Wrists (WLM401), 25 yd/ea
GL702	Size 2 Tubular Bandage Retainer, Small: Hand, Arm, Leg, Foot (WLM402), 25 yd/ea
GL703	Size 3 Tubular Bandage Retainer, Medium: Hand, Arm, Leg, Foot (WLM403), 25 yd/ea
GL704	Size 4 Tubular Bandage Retainer, Large: Hand, Arm, Leg, Foot (WLM404), 25 yd/ea
GL705	Size 5 Tubular Bandage Retainer, Small: Head, Shoulder, Thigh (WLM405), 25 yd/ea
GL706	Size 5.5 Tubular Bandage Retainer, Medium: Head, Shoulder, Thigh (WLM406), 25 yd/ea
GL707	Size 6 Tubular Bandage Retainer, Large: Head, Shoulder, Thigh (WLM407), 25 yd/ea
GL708	Size 7 Tubular Bandage Retainer, Small: Chest, Back, Perineum, Axilla (WLM408), 25 yd/ea
GL709	Size 8 Tubular Bandage Retainer, Medium: Chest, Back, Perineum, Axilla (WLM409), 25 yd/ea
GL710	Size 9 Tubular Bandage Retainer, Large: Chest, Back, Perineum, Axilla (WLM410), 25 yd/ea
GL711	Size 10 Tubular Bandage Retainer, X-Large: Chest, Back, Perineum, Axilla, 25 yd/ea
GL712	Size 11 Tubular Bandage Retainer, XX-Large: Chest, Back, Perineum, Axilla, 25 yd/ea
GL713	Size 12 Tubular Bandage Retainer, Obese: (cir stretch to 42"), 25 yd/ea
GL714	Size 13 Tubular Bandage Retainer, X-Obese: (cir stretch to 46"), 25 yd/ea
GL715	Size 14 Tubular Bandage Retainer, XX-Obese (cir stretch to 90"), 25 yd/ea
GL720	Size 3 Tubular Bandage Retainer, Medium: Hand, Arm, Leg, Foot, 10 yd/ea
GL722	Size 5 Tubular Bandage Retainer, Small: Head, Shoulder, Thigh, 10 yd/ea
GL724	Size 6 Tubular Bandage Retainer, Large: Head, Shoulder, Thigh, 10 yd/ea
GL726	Size 10 Tubular Bandage Retainer, X-Large: Chest, Back, Perineum, Axilla, 10 yd/ea
GLLF2501	Size 1 Tubular Bandage Retainer, Fingers, Toes & Wrists (WLM501), Latex Free, 25 yd/ea
GLLF2502	Size 2 Tubular Bandage Retainer, Small: Hand, Arm, Leg, Foot (WLM502), Latex Free, 25 yd/ea
GLLF2503	Size 3 Tubular Bandage Retainer, Medium: Hand, Arm, Leg, Foot (WLM503), Latex Free, 25 yd/ea
GLLF2504	Size 4 Tubular Bandage Retainer, Large: Hand, Arm, Leg, Foot (WLM504), Latex Free, 25 yd/ea
GLLF2505	Size 5 Tubular Bandage Retainer, Small: Head, Shoulder, Thigh (WLM505), Latex Free, 25 yd/ea
GLLF2505.5	Size 5.5 Tubular Bandage Retainer, Medium: Head, Shoulder, Thigh, Latex Free, 25 yd/ea
GLLF2506	Size 6 Tubular Bandage Retainer, Large: Head, Shoulder, Thigh (WLM506), Latex Free, 25 yd/bx
GLLF2507	Size 7 Tubular Bandage Retainer, Small: Chest, Back, Perineum, Axilla (WLM507), Latex Free, 25 yd/bx
GLLF2508	Size 8 Tubular Bandage Retainer, Medium: Chest, Back, Perineum, Axilla (WLM508), Latex Free, 25 yd/bx
GLLF2509	Size 9 Tubular Bandage Retainer, Large: Chest, Back, Perineum, Axilla (WLM509), Latex Free, 25 yd/bx
GLLF2510	Size 10 Tubular Bandage Retainer, X-Large: Chest, Back, Perineum, Axilla (WLM510), Latex Free, 25 yd/bx
GLLF2511	Size 11 Tubular Bandage Retainer, XX-Large: Chest, Back, Perineum, Axilla, Latex Free, 25 yd/bx
GLLF2512	Size 12 Tubular Bandage Retainer, Obese (cir. stretch to 42"), Latex Free, 25 yd/bx
GLLF2514	Size 14 Tubular Bandage Retainer, XX-Obese: (cir. stretch to 90"), Latex Free, 25 yd/bx



DERMA SCIENCES SURGITUBE® FOR USE WITH APPLICATORS

Cotton seamless tubular bandage system for neat, comfortable, professional dressings. Easy to use, efficient and functional design. Just wrap, anchor, and wrap again. Use as many layers as needed.

GL105F	Size 1, 5/8", Flesh, Small Fingers & Toes, Includes Aluminum Splint Type Applicator, 5 yd/bx
GL105W	Size 1, 5/8", White, Small Fingers & Toes, Includes Aluminum Splint Type Applicator, 5 yd/bx
GL110W	Size 1, 5/8", White, Small Fingers & Toes, Includes Aluminum Splint Type Applicator, 10 yd/bx
GL205F	Size 2, 7/8" x 5 yds, Flesh, Large Fingers & Toes, Includes Aluminum Splint Type Applicator, 12 rl/bx
GL205W	Size 2, 7/8", White, Large Fingers & Toes, Includes Aluminum Splint Type Applicator, 5 yd/bx, 12/cs
GL209	Size 1, Small Fingers and Toes, 5/8", Flesh, 50 yd/bx
GL210	Size 2, Large Fingers & Toes, 1", Flesh, 50 yd/bx
GL210F	Size 2, 7/8", Flesh, Large Fingers & Toes, Includes Aluminum Splint Type Applicator, 10 yd/bx
GL210W	Size 2, 7/8", White, Large Fingers & Toes, Includes Aluminum Splint Type Applicator, 10 yd/bx
GL211	Size 3, Hands, Wrist & Feet, 1 1/2", Flesh, 50 yd/bx
GL212	Size 4, Arms & Lower Leg, 2 5/8", Flesh, 50 yd/bx
GL219	Size 1, Small Fingers & Toes, 5/8" (WLM701), White, 50 yd/bx
GL220	Size 2, Large Fingers & Toes, 1" (WLM703), White, 50 yd/bx
GL221	Size 3, Hands, Wrists & Feet, 1 1/2" (WLM704), White, 50 yd/bx
GL222	Size 4, Arms & Lower Leg, 2 5/8" (WLM706), White, 50 yd/bx
GL223	Size 5, Legs, Thighs & Head, 3 5/8" (WLM707), White, 50 yd/bx
GL224	Size T-1, Torso, 5", White, 50 yd/bx
GL225	Size T-2, Torso, 7", White, 50 yd/bx
GL251	Size 1, Children's Fingers & Toes, Adult's Small Fingers & Toes, 5/8", White, For Use w/o Applicators, 15 yd/bx
GL252	Size 2, Children's Hands & Feet, Adult's Large Fingers & Toes, 7/8", White, For Use w/o Applicators, 15 yd/bx
GL252A	Size 2A, Children's Legs, Eye Compresses, Adult's Bulky Finger Dressings, Removable Bunion Pads, 1 1/8", White, For Use w/o Applicators, 15 yd/bx
GL305F	Size 3, Hands, Feet, Lower Legs & Shoulders, 1 1/2", Flesh, 5 yd/bx
GL305W	Size 3, Hands, Feet, Lower Legs & Shoulders, 1 1/2", White, 5 yd/bx
GL310W	Size 3, Hands, Feet, Lower Legs & Shoulders, 1 1/2", White, 10 yd/bx
GL405W	Size 4, Wrists & Elbows, 1 1/2", White, 5 yd/bx, Tight Weave

DERMA SCIENCES SURGITUBE® FOR USE WITHOUT APPLICATORS

Features a heavier knit construction than the gauze for use with an applicator. Simply pull onto the area to be bandaged.

GL241	Size 1P, Small Finger & Toes, 5/8", White, 50 yd/bx
GL242	Size 2P, Large Fingers & Toes, 7/8" (VLM702), White, 50 yd/bx
GL242A	Size 2AP, Infant's Legs, Bulky Fingers, 1 1/8", White, 50 yd/bx
GL243	Size 3P, Hands, Feet, Lower Legs, Shoulders, 1 1/2", White, 50 yd/bx
GL244	Size 4P, Wrists & Elbows, 1 1/2", Tight Weave, White, 50 yd/bx
GL245	Size 5P, Head, Shoulders, Breasts, Legs, 2" (WLM705), White, 50 yd/bx
GL246	Size 1P, Small Fingers & Toes, 5/8", Flesh, 50 yd/bx
GL247	Size 2P, Large Fingers & Toes, 7/8", Flesh, 50 yd/bx
GL248	Size 3P, Hands, Feet, Lower Legs, Shoulders, 1 1/2", Flesh, 50 yd/bx
GL260	Size 5P, Head, Shoulders, Breasts, Legs, 2", Flesh, 50 yd/bx

DERMA SCIENCES TUBE BANDAGE™ STARTER KIT

Surgitube - one 50 yd roll each of size 1 (GL-219), size 2 (GL-220), and size 3 (GL-221). Lucite tube applicators - one each of size 1 (GL-231P), size 2 (GL-232P), size 3 (GL-234P), and an aluminum applicator.

GL235 Tube Bandage Starter Kit



DERMA SCIENCES TUBULAR BANDAGE APPLICATORS

ALUAPP	Aluminum Splint Applicator For Use w/GL209 & GL219
GL227	Metal, Cylindrical Cage Applicator Set Containing Sizes 1,2,3
GL230	Metal, Cylindrical Cage Applicator For Surgitube Size 0, Pediatric
GL231	Metal, Cylindrical Cage Applicator For Surgitube Size 1
GL231P	Plastic, Cylindrical Cage Applicator For Surgitube Size 1
GL232	Metal, Cylindrical Cage Applicator For Surgitube Size 2
GL232P	Plastic, Cylindrical Cage Applicator For Surgitube Size 2
GL234	Metal, Cylindrical Cage Applicator For Surgitube Size 3
GL234P	Plastic, Cylindrical Cage Applicator For Surgitube Size 3
GL236	Metal, Cylindrical Cage Applicator For Surgitube Size 4



SMITH & NEPHEW PROFORE™ LATEX FREE MULTI-LAYER BANDAGING SYSTEM

66020626 Latex Free Multi-Layer Bandaging System, 8/cs

SMITH & NEPHEW PROFORE™ LITE MULTI-LAYER COMPRESSION BANDAGE SYSTEM

Comprised of a wound contact layer and three bandages. Specially developed for the management of mixed aetiology leg ulcers where assessment has identified that the patient has some degree of arterial impairment which prevents the use of full compression. Designed to maintain effective levels of compression after application. The system may be left in place for up to one week.

66000771 Multi-Layer Compression Bandage System, Latex Free, 8/cs

SMITH & NEPHEW PROGUIDE™ MULTI-LAYER HIGH-COMPRESSION BANDAGE SYSTEM

A multi-layer, high-compression bandage system, incorporating the patented, revolutionary "Vari-Stretch" Technology. May be used to effectively manage venous leg ulcers on mobile or immobile patients as its elastic nature provides sustained, graduated levels of compression even if patient is at rest. Available in 3 sizes to accommodate different ankle circumferences. Each system pack contains the following: 1 ProGuide wound contact layer, 1 ProGuide #1 natural padding bandage and 1 ProGuide #2 outer compression bandage.

66000780	Red Kit, Ankle Circumference 7 1/8" x 8.5", 18-22cm, 1kit/bx, 8 bx/cs
66000781	Yellow Kit, Ankle Circumference 8 5/8" x 11", 22-28cm, 1kit/bx, 8 bx/cs
66000782	Green Kit, Ankle Circumference 11" x 12 5/8", 28-32cm, 1kit/bx, 8 bx/cs



HARTMANN-CONCO ESMARK LF LATEX FREE BANDAGE

Provides consistent 400% elongation for ideal compression and application. An ideal choice for use in surgical procedures. Composed of special elastomeric compounds. Contains no rubber proteins, resins or fatty acids, to deliver maximum protection for latex-sensitive patients and health professionals. Individually wrapped and sterile.

18310000	3" x 9 ft Bandage, 20/cs
18320000	3" x 12 ft Bandage, 20/cs
18410000	4" x 9 ft Bandage, 20/cs
18420000	4" x 12 ft Bandage, 20/cs
18610000	6" x 9 ft Bandage, 20/cs
18620000	6" x 12 ft Bandage, 20/cs
V18310073	3" x 9' Esmark Sterile Bandage, 20 rl/cs
V18320073	3" x 12' Esmark Sterile Bandage, 20 rl/cs
V18410073	4" x 9' Esmark Sterile Bandage, 20 rl/cs
V18420073	4" x 12' Esmark Sterile Bandage, 20 rl/cs
V18610073	NovaPlus Bandage, 6" x 9', Unstretched, Latex Free, 20 rl/cs
V18620073	6" x 12' Esmark Sterile Bandage, 20 rl/cs

Elastic Bandages



HARTMANN-CONCO EZE-BAND® LATEX FREE ELASTIC BANDAGE WITH SELF CLOSURE

EZe-Band is an orthopedic grade knitted bandage that uses no pins, clips, or tape. Front and end self-closures make EZe-Band fast and easy to use again and again. For controlled compression and support. Full range of sizes in regular and double length sterile and non-sterile. Latex free.

59120000	2" x 5 yd, Bandage, Non-Sterile, 10/pk, 6 pk/cs
59130000	3" x 5 yd, Bandage, Non-Sterile, 10/pk, 6 pk/cs
59140000	4" x 5.5 yd, Bandage, Non-Sterile, 10/pk, 6 pk/cs
59160000	6" x 5.5 yd, Bandage, Non-Sterile, 10/pk, 6 pk/cs
59180000	4" x 11 yd, Bandage, Non-Sterile, 12/bx, 3 bx/cs
59190000	6" x 11 yd, Bandage, Non-Sterile, 12/bx, 2 bx/cs
59720000	2" x 5 yd, Bandage, Sterile, 36/cs
59730000	3" x 5 yd, Bandage, Sterile, 36/cs
59740000	4" x 5.5 yd, Bandage, Sterile, 36/cs
59760000	6" x 5.5 yd, Bandage, Sterile, 36/cs
59790000	6" x 11 yd, Bandage, Sterile, 24/cs



HARTMANN-CONCO FLEXICON® LATEX FREE CONFORMING STRETCH BANDAGE

Conforming stretch bandage that holds dressings securely in place without slipping. Unique knit structure clings to itself. Uniform knitting of polyester yarns creates a gauze with stretchability that conforms even to the most difficult areas. Finished edges minimize lint and loose threads. Durable with great strength, yet it is soft and comfortable to wear. Blue application guideline for consistent wrapping. Latex free.

18200000	2" x 4.1 yds Bandage, Non-Sterile, Individually Wrapped, 20/bx, 5 bx/cs
18300000	3" x 4.1 yds Bandage, Non-Sterile, Individually Wrapped, 20/bx, 5 bx/cs
18400000	4" x 4.1 yds Bandage, Non-Sterile, Individually Wrapped, 20/bx, 5 bx/cs
18600000	6" x 4.1 yds Bandage, Non-Sterile, Individually Wrapped, 20/bx, 3 bx/cs
19100000	1" x 4.1 yds Bandage, Sterile, 12/bx, 8 bx/cs
19200000	2" x 4.1 yds Bandage, Sterile, 12/bx, 8 bx/cs
19300000	3" x 4.1 yds Bandage, Sterile, 12/bx, 8 bx/cs
19400000	4" x 4.1 yds Bandage, Sterile, 12/bx, 8 bx/cs
19600000	6" x 4.1 yds Bandage, Sterile, 12/bx, 4 bx/cs
22100000	1" x 4.1 yds Bandage, Non-Sterile, 24/bx, 4 bx/cs
22200000	2" x 4.1 yds Bandage, Non-Sterile, 12/bx, 8 bx/cs
22300000	3" x 4.1 yds Bandage, Non-Sterile, 12/bx, 8 bx/cs
22400000	4" x 4.1 yds Bandage, Non-Sterile, 12/bx, 8 bx/cs
22600000	6" x 4.1 yds Bandage, Non-Sterile, 6/bx, 8 bx/cs



HARTMANN-CONCO FOURPRESS® COMPRESSION BANDAGING SYSTEM

Used to effectively manage venous leg ulcers and related conditions. Composed of four bandages that work together. Provides effective sustained compression for up to seven days when applied as directed. Layer One is a non-woven padding composed of a blend of synthetic materials that absorbs wound exudate. This highly comfortable layer primarily provides protection to bony prominences. Layer Two is a crepe bandage which anchors and smoothes the padding layer while absorbing exudate. Layer Three is a compression bandage that conforms to any contour while providing about one half the compression of the system. Layer Four is a latex free cohesive bandage that sticks to itself and secures the bandages without tape for week-long compression.

43400000 4" Compression Bandaging System, 8 kits/cs



HARTMANN-CONCO LOPRESS® LATEX FREE COMPRESSION BANDAGE

A premium grade 100% cotton bandage that provides maximum resistance while reducing the risk of constriction. Designed for repeated use, thus it withstands numerous washings and resists deterioration from ointments and medications. Management of lymphedema and chronic venous insufficiencies, including post-operative treatment after vein stripping. Treatment of primary and secondary varicosis, edema, post-thrombotic syndrome, phlebotrombosis and thrombophlebitis. Can be sterilized. Latex free.

42200000	2.3" x 5.4 yds Compression Bandage, 12 rl/bx
42300000	3.1" x 5.4 yds Compression Bandage, 12 rl/bx
42400000	3.9" x 5.4 yds Compression Bandage, 12 rl/bx
42410000	3.9" x 10.9 yds Compression Bandage, 6 rl/bx
42500000	4.7" x 5.4 yds Compression Bandage, 12 rl/bx



HARTMANN-CONCO REB® LF REINFORCED ELASTIC BANDAGES

Offers balanced stretch and compression. Lightweight and ideal for general use. Perfect for securing splints, dressings, ice packs and in the treatment of varicosities, sprains and strains. The Sterile REB® is sealed in an easy-open Tyvek® pouch. Latex Free.

16200000	2" x 5 yd Bandage, 10/pk, 6 pk/cs
16300000	3" x 5 yd Bandage, 10/pk, 6 pk/cs
16400000	4" x 5 yd Bandage, 10/pk, 6 pk/cs
16410000	4" x 10 yd Bandage, 6/cs
16600000	6" x 5 yd Bandage, 10/pk, 6 pk/cs
16610000	6" x 10 yd Bandage, 6/cs
34200000	2" x 5 yd Bandage, Sterile, 10/cs
34300000	3" x 5 yd Bandage, Sterile, 10/cs
34400000	4" x 5 yd Bandage, Sterile, 10/cs
34600000	6" x 5 yd Bandage, Sterile, 10/cs
35410000	4" x 10 yd Bandage, Sterile, 10/cs
35610000	6" x 10 yd Bandage, Sterile, 10/cs



HARTMANN-CONCO SHUR-BAND® LF LATEX FREE SELF-CLOSURE ELASTIC BANDAGE

Latex free, premium grade support and compression bandage knitted with high quality yarns to provide superior performance. The non-latex component offers latex sensitive patients and healthcare professionals an alternative to the "standard" elastic bandage, and can be worn comfortably without risk of allergic reactions. A self-closure flap fastens to the bandage and holds it securely without clips, pins, or tape. Ice wrap is designed with a sewn pouch to hold standard size hot and cold therapy packs.

59520000	2" x 5 yds Bandage, 10/pk, 6 pk/cs
59530000	3" x 5 yds Bandage, 10/pk, 6 pk/cs
59540000	4" x 5 yds Bandage, 10/pk, 6 pk/cs
59560000	6" x 5 yds Bandage, 10/pk, 6 pk/cs
59570000	4" x 10 yds Bandage, 6/cs
59580000	6" x 10 yds Bandage, 6/cs
59880000	4" x 30" Bandage, Ice Wrap w/pouch for hot or cold pack, 12/cs
59890000	6" x 30" Bandage, Ice Wrap w/pouch for hot or cold pack, 12/cs
59920000	2" x 5 yds Bandage, Sterile, 10/cs
59930000	3" x 5 yds Bandage, Sterile, 10/cs
59940000	4" x 5 yds Bandage, Sterile, 10/cs
59960000	6" x 5 yds Bandage, Sterile, 10/cs



HARTMANN-CONCO THREEPRESS® COMPRESSION BANDAGING SYSTEM

Latex free multi-layer bandage system that is used in the management of venous leg ulcers. The system consists of three bandages that when applied provide 35-40 mmHg of sub-bandage pressure.

43420000	Compression Bandaging System Contains: (3) Bandages & (3) Tape Strips, 8/cs
----------	---



NUTRAMAX ESMARK BANDAGES

For exsanguination of a limb prior to tourniquet application. High quality latex free rubber provides uniform response along the length and width of the bandage. Durable and tear-resistant. Latex free.

12301	3" x 9 yds Esmark Bandage, 20/cs
12303	4" x 9 yds Esmark Bandage, 20/cs
12304	4" x 12 yds Esmark Bandage, 20/cs
12305	6" x 9 yds Esmark Bandage, 20/cs
12306	6" x 12 yds Esmark Bandage, 20/cs

PROCARE ELASTIC BANDAGES - SELF CLOSURE

Constructed of breathable elastic. Ideal for minor sprains and strains. Each bandage is individually wrapped and may be secured by hook and loop closure or with two clip fasteners. Latex free.

79-98233	2" x 5 yds Bandage, 10/pkg
79-98235	3" x 5 yds Bandage, 10/pkg
79-98237	4" x 5 yds Bandage, 10/pkg
79-98238	6" x 5 yds Bandage, 10/pkg

PROCARE FLEX-MASTER® NON-STERILE CLIP CLOSURE BANDAGE

Latex free, 100% cotton/rubber minimizes moisture retention and allergic reactions. Ideal for groin wraps, ice packs, knees, and ankles.

79-98843	2" x 5.5 yds Non-Sterile Bandage, 10/pkg
79-98845	3" x 5.5 yds Non-Sterile Bandage, 10/pkg
79-98847	4" x 5.5 yds Non-Sterile Bandage, 10/pkg
79-98848	6" x 5.5 yds Non-Sterile Bandage, 10/pkg
79-98857	4" x 11 yds Non-Sterile Bandage, 12/pkg
79-98859	6" x 11 yds Non-Sterile Bandage, 12/pkg

PROCARE FLEX-MASTER® STERILE CLIP CLOSURE BANDAGE

For post-surgical applications such as thoracic, knee, amputations, mastectomy, vascular, shoulder and burns. Latex free, cotton/elastic construction with clip closure.

79-98717	4" x 11 yds Sterile Bandage, 20/pkg
79-98718	6" x 11 yds Sterile Bandage, 20/pkg
79-98723	2" x 5.5 yds Sterile Bandage, 20/pkg
79-98725	3" x 5.5 yds Sterile Bandage, 20/pkg
79-98727	4" x 5.5 yds Sterile Bandage, 20/pkg
79-98728	6" x 5.5 yds Sterile Bandage, 20/pkg

PROCARE LATEX RUBBER BANDAGES

Constructed of natural latex and designed to force blood out of an extremity by compressing the superficial vessels.

79-98900	3" x 144" Ltex Bandage
79-98910	4" x 144" Latex Bandage
79-98920	6" x 144" Latex Bandage

PROCARE RUBBER ELASTIC BANDAGES

Cotton/elastic wrap stretches to a full 4.5 yds and secures with metal clip. Washable and autoclavable. Latex free.

79-98583	2" x 4.5 yds Bandage, 10/pkg
79-98585	3" x 4.5 yds Bandage, 10/pkg
79-98587	4" x 4.5 yds Bandage, 10/pkg
79-98588	6" x 4.5 yds Bandage, 10/pkg



SMITH & NEPHEW UNI-SOLVE® ADHESIVE REMOVER

Formulated to reduce adhesive trauma to the skin by thoroughly dissolving dressing, tape and appliance adhesives. Ideal for removing ostomy faceplate adhesive bandages as well as acrylic and hydrocolloid based adhesives. Does not contain CFCs (chlorinated fluoro carbons).

402300	Adhesive Remover Wipes, 50/pkg, 20 pkg/cs
59402500	8 oz. Bottle Adhesive Remover, 12/cs

Hand Sanitizers



BURNISHINE ADHESIVE REMOVER

Quickly and painlessly removes tape from skin. Safely removes tape from equipment.

R-6 16 oz Adhesive Remover, 12/cs



SMITH & NEPHEW I.V. PREP ANTISEPTIC

Provides a clear, one-step antiseptic for I.V. site preparation as well as a film which may help reduce skin trauma associated with adhesive removal when dressings are changed.

59421200 IV Prep Antiseptic Wipes, 50/pkg, 20 pkg/cs

GRAHAM PROFESSIONAL NAIL & COSMETIC PADS

Premium nail and cosmetic pad used by professional Nail and Spa Technicians. Features clean grip handle and protective barrier to protect the technicians nails. One pad capable of removing polish from all ten nails.

42950 Handsdown® Ultra Nail & Cosmetic Pads, 1.75" Round, White, 240/bg, 6 bg/cs



DYNAREX ADHESIVE TAPE REMOVER

1505 Adhesive Tape Remover, Mint Scented, 100/bx 10 bx/cs

GRAHAM PROFESSIONAL SANEK™ HAND SANITIZING DISPOSABLE WIPES

Will kill 99.9% of common germs on contact. For use in physician offices, hospitals, emergency medical settings, labs, clinics, dental offices, health clubs, schools, day care and veterinary facilities.

55727 Wipe, Hand Sanitizing, 80 count, 12/cs

55757 Wipe, Hand Sanitizing, Disposable, 80/canister, 12 can/cs



D43600



P13472



P15984

PDI SANI-HANDS® ALC INSTANT HAND SANITIZING WIPES

Pre-moistened wipe delivery system. 65.9% ethyl alcohol. Proven more effective than alcohol gels alone. The right amount of active at the right time. Friction of wipe removes soil and bacteria. Kills 99.99% of germs. Helps reduce the risk of cross-contamination and spread of infection. Cleans and sanitizes hands when soap and water are not available. Contains moisturizing aloe and Vitamin E which helps retain moisture in skin with no gummy residue. Dries quickly so no need to "wave" hands to air dry. Fragrance-free and dye-free which minimized allergic reactions and sensitivity. FDA food compliant. Convenient quick-pull canister and individual packet options. Canister features a unique deep-well lid. Sani-Bracket® wall and mobile equipment attachment system and floor stands dispenses wipes from wall, mobile equipment mounting or floor stand.

D43600 Handwipes, Pocket Packets (69.9% Alcohol), 100/bx, 10 bx/cs

P13472 Handwipes, Medium Canister (69.9% Alcohol), 135/can, 12 can/cs

P15984 Wipes, 6" x 7.5", Large Canister, 220/can, 6 can/cs

CROSSTEX SANI-DEX® ANTIMICROBIAL HAND WIPES

Demonstrated 99.99% effective against 28 bacteria and 5 viruses as well as gram positive and gram negative bacteria, fungi and yeast. 65.9% ethyl alcohol as recommended by CDC for maximum skin protection. Clinically proven to remove bacteria and soil from hands. Aloe Vera & Vitamin E - moisturizer/condition/soften skin & skin & combat dryness. Fragrance-free and dye-free.

SAND 6" x 7.5" Wipes, 135/tub, 12 tub/cs

GOJO PURELL® SANITIZING HAND WIPES

Enhance guest experience and satisfaction with economical towelettes that clean and sanitize. Individually wrapped, on-the-go amenity item for a wide variety of markets, including foodservice, travel and hospitality. Made with 62% ethyl alcohol the same active ingredient as PURELL Instant Hand Sanitizer.

9010-06 Non Linting, Durable Textured Wipe for Superior Cleaning, 175ct Popup Canister, Wiper Size 6" x 8", Tested & Approved for Use w/hands, Not Approved for Use as a Surface Disinfectant, 6 Canisters/cs

9011-12 Non Linting, Durable Textured Wipe for Superior Cleaning, 35ct Popup Canister, Wiper Size 7" x 8", Tested & Approved for Use w/hands, Not Approved for Use as a Surface Disinfectant, 12 Canisters/cs

9020-4M Individual Wrapped Wipes, 4000 Count Bulk Packed Shipper

9021-1M Individual Wrapped Wipes, 1000 Count Bulk Packed Shipper

9022-10 Individually Wrapped Wipes, 100 count box, 10 bx/cs

Accessories

9014-01 Wall Bracket For Use w/9010-06



MOLNLYCKE HIBISTAT® TOWELETTE

Provides rapid bacterial action and has a persistent antimicrobial effect against a wide range of microorganisms.

58717 Towelette, 5mL Packet, 50/bx, 8 bx/cs



P779114

P779118

PRO ADVANTAGE ALCOHOL BASED HAND SANITIZER

Kills germs with 62% alcohol. Rubs into skin quickly for clean, sanitized hands and skin. Skin conditioners prevent skin from drying out. Available in extra thick gel or foam.

P779115 1.5 oz foam, 24/cs
P779114 4oz oval bottle, gel, 24/cs
P779118 8 oz Bottle w/pump, gel, 12/cs



GOJO PURELL® SANITIZING STATION STAND

This free-standing unit (dispenser and refills sold separately) offers the combination of a PURELL TFX™ TouchFree Dispenser and a sturdy, moveable stand for use in high-traffic locations. Highly visible signage sends the message that you care about reducing the spread of germs. Easy to assemble, with an attractive design and built-in shield beneath the dispenser to keep surrounding areas clean. Ideal for lobbies, cafeterias, halls reception areas, restaurants, and other high-traffic areas. For use with Dispenser: 2720-12. Suggested refills: 5456-04, 5457-04, and 5392-02.

- 2424-DS Floor Stand w/"Healthy Hands Start Here" Sign, Purell® TFX™
- 2426-DS Table Top Stand w/"Healthy Hands Start Here" Sign, Purell® TFX™
- 2427-DS Universal Mount Stand w/"Healthy Hands Start Here" Sign, Purell® TFX™
- 2428-MB Bracket, Mask, For Use w/Purell TFX™ Stands, (NOTE: MASK BOXES ARE NOT AVAILABLE FROM GOJO)
- 2429-TB Bracket, Tissue, For Use w/Purell TFX™ Stands, (NOTE: TISSUE BOXES ARE NOT AVAILABLE FROM GOJO)
- 5398-02 Hand Sanitizer, Instant, Purell Skin Nourishing Foam, 2/cs



GOJO PURELL® FOAMING HAND SANITIZER

An alcohol-based, aerosol hand sanitizer. Formulated with 62% ethyl alcohol for broad spectrum germ kill, including MRSA & VRE. Requires only 3 grams of product per use. Latex & CHG compatible, dye free.

9698-12 15 oz. Can Purell® Foaming Hand Sanitizer, 12/cs

YOUR FIRST LINE OF DEFENSE

INFECTION PREVENTION

Sani-Cloth®. Sani-Hands®. Preventing the spread of infections is a team effort encompassing good surface disinfection and hand hygiene. PDI® provides trusted, cost efficient products to help you prevent the spread of infections in every situation.

PDI®. Our passion is prevention™.

Learn more at www.pdipdi.com



OUR PASSION IS PREVENTION™

© 2009 Professional Disposables International, Inc. ISO 9001:2000 Certified
Sani-Cloth®, Sani-Hands® ALC and Sani-Hands® for Kids are registered trademarks of Professional Disposables International, Inc.

Hand Sanitizers



2007-5C

2137-08

2149-08

2151-08



2156-08



2166-04



2237-04



2256-04



3005-60



3005-WRP



5192-03



5392-02



5451-04



5692-24



5792-04



7187-06



9605-24



9606-24



9625-04



9631-24



9637-12



9639-12



9646-12



9647-24



9648-24



9649-12



9651-24



9652-12



9654-12



9656-06



9657-12



9659-12



9669-12



9674-12



9682-24

GOJO PURELL® INSTANT HAND SANITIZER

- 2007-5C Hand Sanitizer, 0.5 oz. Bottle, Clear Original Formula, 500/cs
- 2137-08 1000mL Hand Sanitizer w/aloe refill, Green, 8/cs
- 2149-08 1000mL Refill Instant Hand Sanitizer, White, 8/cs
- 2151-08 Instant Hand Sanitizer w/Dermaglycerin System, 1000mL, NXT™ Refill, 8/cs

GOJO PURELL® INSTANT HAND SANITIZER (Continued)

- 2156-08 NXT™ Original Hand Sanitizer Refill, Clear, 1000mL, 8/cs
- 2166-04 Food Code Compliant, 1000mL, NXT™ Refill, 4/cs
- 2237-04 2000mL Hand Sanitizer w/aloe refill, Green, 4/cs
- 2256-04 2000mL NXT Formula Hand Sanitizer Refill, Clear, 4/cs
- 3005-60 0.5 oz. Bottle Hand Sanitizer, Clear Original Formula, 60/cs
- 3005-WRP Hand Sanitizer, Instant, Original, 0.5 oz., 36/cs
- 3023-12 Original Formula, 20 fl oz. Pump Bottle, Clear, 12/cs
- 5190-03 Instant Hand Sanitizer Foam (Non-Aersol), FMX12™, Clear, 1000mL, 3/cs
- 5192-03 Instant Hand Sanitizer Foam (Non-Aersol), FMX12™, Clear, 1200mL, 3/cs
- 5392-02 Instant Hand Sanitizer Foam (Non-Aersol), FMX™, Clear, 1200mL, 2/cs
- 5451-04 Instant Hand Sanitizer w/Dermaglycerin System, 1200mL, TFX™ Refill, 4/cs
- 5456-04 Instant Hand Sanitizer, Original Formula, 1200mL TFX Refill, Clear, 4/cs
- 5457-04 Instant Hand Sanitizer w/aloe, 1200mL TFX Refill, Clear, 4/cs
- 5466-04 Instant Hand Sanitizer Food Code Compliant, 1200mL TFX Refill, Clear, 4/cs
- 5592-02 Instant Hand Sanitizer Foam (Non-Aersol, FMX™, Clear, 1000mL, 2/cs
- 5692-24 Instant Hand Sanitizer Foam (Non-Aersol), 45mL, Pump Bottle, 24/cs
- 5792-04 Instant Hand Sanitizer Foam (Non-Aerosol), 535mL, Counter Top Pump Bottle, 4/cs
- 7187-06 1000mL Original Hand Sanitizer Refill, Clear, 6/cs
- 9600-DC1-BLK 8 fl oz. Pump Bottle w/holder (Purell Desk Caddy), Includes Educational Material, 12/cs
- 9600-EA-PL1 Hand Sanitizer 8 oz. Pump Bottle w/holder
- 9600-PL1 8 fl oz. Pump Bottle w/holder (Purell Pal), Includes Educational Material, 12/cs
- 9605-24 2 oz Squeeze Bottle w/flip-cap, 24/cs
- 9606-24 2 oz Pump Bottle, 24/cs
- 9624-24 3 fl oz. Bottle, Foliage Green w/disc-cap, 24/cs
- 9625-04 2 Liter Pump Bottle, Original Hand Sanitizer, Clear, 4/cs
- 9631-24 4.25 oz. Bottle Hand Sanitizer w/aloe, Green, 24/cs
- 9637-12 800mL Hand Sanitizer w/aloe refill, Green, 12/cs
- 9639-12 12 oz. Pump Bottle Hand Sanitizer w/aloe, Green, 12/cs
- 9646-12 Instant Hand Sanitizer w/Dermaglycerin System, 12 oz. Pump Bottle, 12/cs
- 9647-24 2 oz Squeeze Bottle Instant Hand Sanitizer, White, 24/cs
- 9648-24 Instant Hand Sanitizer w/Dermaglycerin System, 2 oz. Pump Bottle, 24/cs
- 9649-12 12 oz Pump Bottle Instant Hand Sanitizer, White, 12/cs
- 9651-24 4 oz. Bottle Original Hand Sanitizer, Clear, 24/cs
- 9652-12 8 oz. Pump Bottle, Original Hand Sanitizer, Clear, 12/cs
- 9654-12 250mL Original Hand Sanitizer Refill, Clear, 12/cs
- 9656-06 800mL Original Hand Sanitizer Refill, Clear, 6/cs
- 9657-12 800mL Original Hand Sanitizer Refill, Clear, 12/cs
- 9659-12 12 oz. Pump Bottle Original Hand Sanitizer, Clear, 12/cs
- 9669-12 12 fl oz, Food Code Compliant, 12/cs
- 9674-12 8 oz Bottle Instant Hand Sanitizer w/aloe, Yellow, 12/cs
- 9682-24 2 fl oz Bottle w/flip-cap, contains Aloe, 24/cs

GOJO PREMIUM FOAM ANTIBACTERIAL HANDWASH

Kills germs on hands while providing a rich, luxurious lather. Gentle on hands. Contains Vitamin E and aloe. Fresh fruit fragrance. Translucent apricot color.

- 5362-02 Handwash, Antibacterial Foam, Apricot Color, 1200mL Refill For Use w/TFX Touch Free Dispenser, 2/cs



3M™ AVAGARD™ D INSTANT HAND ANTISEPTIC

Provides rapid bactericidal action against a broad spectrum of microorganisms while helping to prevent skin breakdown and maintaining the skin's integrity. Contains 61% w/w/ethyl alcohol in a moisturizing base. Kills more than 99% of harmful bacteria within 15 seconds. Unique rich emollient base feels smooth on hands. Doesn't leave hands sticky or impact glove performance.

- 9221 Instant Hand Sanitizer Antiseptic, 88ml, 48/cs
- 9222 Instant Hand Sanitizer Antiseptic Pump Bottle, 500ml, 12/cs

3M™ AVAGARD™ SURGICAL & HEALTHCARE PERSONNEL HAND ANTISEPTIC

Chlorhexidine gluconate 1% solution and ethyl alcohol 61% w/w

- 9200 Hand Antiseptic, 16.9 fl oz, 8/cs
- 9210 Hand Antiseptic w/moisturizers, Personal Size 0.34 fl oz, 40/cs

Accessories

- 9201 Wall Bracket & Foot Pump
- 9202 Wall Bracket & Hand Pump
- 9204 Nail Cleaners, 150/bx, 6 bx/cs

3M™ AVAGARD™ FOAM INSTANT HAND ANTISEPTIC WITH MOISTURIZERS

A non-aerosol, foam hand antiseptic (62% w/w ethyl alcohol) that kills 99% of harmful bacteria in 15 seconds in vitro. Quickly converts to a liquid when rubbed into hands for even, thorough coverage. Designed to spread easily and dry fast. No special fees for aerosol can disposal. Gel-free formula does not contain thickening agents that can clog dispensers or create gel sludge on floors.

- 9251 Hand Antiseptic, Instant Foam w/moisturizers, 1.6fl oz (50mL), Pump Bottle, 25/cs
- 9252 Hand Antiseptic, Instant Foam w/moisturizers, 13.5 fl oz (400mL), Pump Bottle, 12/cs
- 9253 Hand Antiseptic, Instant Foam w/moisturizers, 33.8 fl oz (1000mL), Cartridge (For Use w/9254 wall dispensers) 6/cs

Accessories

- 9254 Wall Dispenser (for 9253 - 1000mL Cartridge), 12/cs (FREE if purchased w/item #9253)

KIMBERLY-CLARK HAND SANITIZER - KIMCARE

Kill germs on the go with Instant Hand Sanitizer gel. Provides instant degerming when water is not available. Its ethyl alcohol formula evaporates quickly, making it useful for just about anywhere, especially for health care and foodservice environments.

- 93060 8 oz. Pump Bottle, Sweet Citrus Fragrance, Clear, 12/cs
- 95150 Hand Sanitizer, 1.2 Liter, Clear, Sweet Citrus Fragrance, (NOTE: Can Be Used in the PPE Dispensing System) 4/cs



KIMBERLY-CLARK HAND SANITIZER

A 62% alcohol formula that contains aloe and skin conditioners to a moisturize and prevent dry skin. Dermatologist tested.

- 93056 8 oz. Instant Hand Sanitizer, White, 12/cs

CROSSTEX SANITYZE™ WATERLESS MOISTURIZING ANTIMICROBIAL GEL

Profound antimicrobial action in 15 seconds - as it dries. Kills up to 99% of bacteria without water. 61% Ethanol: recommended by CDC for maximum skin protection and best activity against viruses*. Long lasting: inhibits bacteria re-growth. Aloe Vera & Vitamin E - Latex compatible emollients prevents skin dryness and glove damage. One-Touch dispensing system to promote easy compliance.

- JWG Sanitizer, 4 oz., 24/cs
- JWGEL Sanitizer, 800 mL Press Pack, 12/cs
- JWGP Sanitizer, 18 oz., 16/cs
- PGEL One Touch Dispenser, 12/cs

TRIAD RICHMOND STAT-FOAM

- 11125R-RM4 Antiseptic Hand Foam, Waterless, 1000mL, 4/cs

TRIAD RICHMOND HEALTH-STAT FOAM

Waterless antiseptic hand foam. Kills 99.99% pathogenic organisms.

- 11105R-FM1000 Antiseptic Hand Foam, 1000mL Bottle, White, Fresh & Clean Fragrance, 8/cs



TRIAD ANTISEPTIC HAND GEL

No-rinse formulation is ideal for use when patient contact is frequent but soap and water are not readily available. Leaves no sticky or tacky residue. Flip-top dispenser. Conditions your skin with aloe vera and will leave your hands with a clean, fresh scent.

- 10-8604 4 oz Hand Gel, 48/cs
- 10-8608 8 oz Hand Gel, 24/cs
- 10-8616 16 oz Hand Gel w/pump, 12/cs
- 10-8628 Gallon Hand Gel, 4/cs
- 10-8672 5 gram Hand Gel, 12/bx, 6 bx/cs

Accessories

- 10-8699 Wall Mount Bracket For 16 oz Hand Gel



10-1624

10-1628

METREX VIONEX® NO RINSE GEL ANTISEPTIC HANDWASH

A waterless gel for hand cleaning when soap and water are not available. Kills germs and cleans skin without running water. Non-flammable and non-irritating. Contains skin conditioners that leave hands soft. Active ingredient: chloroxylenol (PCMX)

- 10-1624 VioNex No Rinse Gel 4 oz w/flip-top, 24/cs
- 10-1628 VioNex No Rinse Gel 18 oz w/pump, 12/cs

Hand Sanitizers - Soaps/Cleansers



METREX VIONEXUS™ NO-RINSE SPRAY ANTISEPTIC HANDWASH

An alcohol-based hand sanitizer for protection against bacteria when hands are not visibly soiled or contaminated. Meets CDC Guideline for Hand Hygiene in Healthcare Settings. Kills germs and cleans skin without running water. 65.87% ethanol and 0.10% benzalkonium chloride provide long-lasting protection against bacteria that can cause disease. Emollients help prevent dry, irritated skin and leave hands feeling soft. No water or towels needed: Eliminates cross contamination. Fast drying for quick donning of gloves.

- 10-1800 VioNexus No Rinse Spray Handwash, 1 Liter, 6/cs
10-1802 Vionex No Rinse Spray Handwash, 2 oz, 48/cs



NEW WORLD FRESHSCENT™ HAND SANITIZER

- HS4 4 oz Hand Sanitizer, 60/cs



632712303 632712942 632712949

BEAUMONT CITRUS II INSTANT HAND SANITIZING LOTION

Citrus II Instant Hand Sanitizing Lotion is a moisture-rich gel that kills 99.99% of most common disease-causing germs on the hands in seconds. Citrus II is ideal to use when soap and water hand washing is not available. And for extra safety, it should be used after hand-washing to kill any remaining germs. Works in 15 seconds. No water or towels needed. Non-toxic - leaves no harmful residue. Leaves hands soft and refreshed. So gentle, repeated use will not dry skin. Fresh, clean citrus fragrance.

- 632712303 8 oz. Pump Bottle, 12/cs
632712942 4 oz. Bottle, 24/cs
632712949 32 oz. Pump Bottle, 6/cs
632771121 0.5 oz. (Sample) Size, 72/cs

MAJOR HAND SANITIZER

- 241903 Instant Hand Sanitizer, 240mL, Compare to Purell®

MEDICOM SOFT-T CLEANSE

Increase your comfort with new Soft-T-Cleanse foaming hand sanitizer, providing effective, yet gentle all-day antiseptic cleansing of your skin! Our exclusive formula is Alcohol-Free, Fragrance-Free and Non-Irritating. Soft-T-Cleanse foam contains soothing all-natural moisturizing extract of Aloe Vera, Lavender and Chamomile. These extracts are widely recognized as healing agents for a variety of skin disorders, often associated with gloves allergies. Our patented anti-bacterial formula is clinically proven to kill 99.9% of bacteria instantly, protecting you from cross-contamination. Save Time and Money with our instant water-based foam that requires no water, soap or paper towels. Provides up to 3 times more applications per oz. compared to alcohol-based hand sanitizing gels. Soft-T-Cleanse foam is quick, safe and convenient to use between patients - simply apply, rub hand until dry & go!

- 7702C Hand Sanitizer, Foam, 1.69 oz, 48/cs
7706C Hand Sanitizer, Foam, 16 oz., 12/cs



17350



17351



17352



17354



17355



17357



17360



17374



18350



18352



18354



18355



18357



18360

SAFETEC ANTISEPTIC BIO-HAND CLEANER™

Waterless hand sanitizers are available in fresh or citrus scent and kill 99.9% of germs in 15 seconds without soap and water. It evaporates completely with no rinsing required. a.b.h.c.™ contains aloe vera to add moisture to the skin with each use, thus encouraging frequent hand-washing to aid in reducing the risk of cross-contamination. a.b.h.c.™ is available in a convenient pump dispenser, single use pouch and flip top dispenser bottle. a .b.h.c.™ helps you to meet OSHA and APIC hand-washing standards in the absence of soap and water.

- 17350 Cleaner, Fresh Scent, 4oz, 24/cs
17351 Cleaner, Fresh Scent, 1/2oz Pouch, 150/cs
17352 Cleaner, Fresh Scent, 3g Pouch, 100/bx, 10bx/cs
17354 Cleaner, Fresh Scent, 16oz Bottle, 12/cs
17355 Cleaner, Fresh Scent, 2oz Bottle, 24/cs
17357 Cleaner, Fresh Scent, 64oz Bottle, 8/cs
17360 Cleaner, Fresh Scent, 27.05oz Bag-In-Box, 12/cs
17374 Cleaner, Fresh Scent, 1.5 oz Tottle Bottle, 50/cs
17376 Cleaner, Fresh Scent, 0.9g, 25/bx, 36 bx/cs
17377 Cleaner, Fresh Scent, 0.9g, 144/bx, 12 bx/cs
18350 Cleaner, Citrus Scent, 4oz, 24/cs
18351 Cleaner, Citrus Scent, 1/2oz Pouch, 150/cs
18352 Cleaner, Citrus Scent, 3g Pouch, 100/bx, 10bx/cs
18354 Cleaner, Citrus Scent, 16oz Bottle, 12/cs
18355 Cleaner, Citrus Scent, 2oz Bottle, 24/cs
18357 Cleaner, Citrus Scent, 64oz, 8/cs
18360 Cleaner, Citrus Scent, 27.05oz Bag-In-Box, 12/cs

Accessories

- 27500 Bag-In-Box Dispenser, 12/cs

BEUTLICH HURRISEPT® ANTISEPTIC HAND GEL

HurriSept is an antiseptic hand gel that reduces the risk of bacterial contamination and infection. Contains 70% ethyl alcohol and comes in a 16fl oz. pump bottle. Also contains aloe vera gel to keep hands soft and moisturized and has a "berry fresh" scent. Product NDC# 0283-0304-16.

- 0283-0304-16 Antiseptic Hand Gel, 16oz Pump Bottle



01585



95862



01585



95862

DIAL® ANTIBACTERIAL HAND SANITIZER

Provides convenient, effective germ-kill that works to kill harmful disease-causing germs, yet is gentle on the hands. Kills 99.999% of pathogenic germs in as little as 15 seconds. Effective against bacteria, yeast & molds. Certified by NSF international for an E-2 rating. Latex and Non-latex compatible. Dermatologist tested and hypoallergenic. Fragrance and dye-free. Contains added moisturizers and skin conditioners to soften and smooth the hands.

- 00685 Hand Sanitizer, 4 oz., Fragrance Free, Flip Top Bottle, 24/cs
- 00688 Instant Hand Sanitizer Fragrance Free, Flip-Top Cap, 2 oz., 24/cs
- 01585 Instant Hand Sanitizer w/moisturizers, Pump, 7.5 oz., 12/cs
- 01L00463 Lanyard/Belt Clip For 2 oz Hand Sanitizer, 24/cs
- 95862 Hand Sanitizer Refill, Fragrance Free, 800 mL Flex Pak, 12/cs



THERMOFISHER SCIENTIFIC HANDCARE ALCOSCRUB® INSTANT ANTISEPTIC HAND CLEANSER

Contains emollients that makes it mild on hands. Highly effective against a broad spectrum of microorganisms. No sticky residue. Fast acting. Absorbs completely. No strong or lingering fragrance. 99.9% effective within 15 seconds or less.

- 42004-24-001 Instant Antiseptic Hand Cleanser, 4 oz Bottle w/flip-top cap, 24/cs
- 42016-12-001 Instant Antiseptic Hand Cleanser, 16 oz Bottle w/pump-top, 12/cs
- 42027-12-001 Instant Antiseptic Hand Cleanser, 27 oz (800ml) Bag-In-A-Box, 12/cs
- 42032-12-001 Instant Antiseptic Hand Cleanser, 32 oz Bottle w/pump-top, 12/cs

DERMA SCIENCES MYSOTROL® HAND SANITIZER

A convenient and highly practical healthcare personnel hand sanitizer that meets the OSHA protocol for a healthcare personnel handwash, reducing the risk of nosocomial infections. As a waterless, no-rinse hand sanitizer, it also evaporates quickly to save time. Lastly, it is formulated without isopropyl alcohol, eliminating the strong alcohol scent common to other hand sanitizers.

- S505-4Z Hand Sanitizer, 4 oz., 24/cs

SULTAN INSTAFRESH HAND SANITIZER

99.9% effective in 30 seconds against harmful pathogens such as MRSA, Staphylococcus aureus, Herpes Simplex Type 1, and Influenza A2. Moisturizes hands leaving them feeling soft and clean. No sticky residue or overpowering fragrance remains on hands.

- 95709 4 oz. Flip Top Bottle, 48/cs
- 95710 18 oz. Pump Bottle, 12/cs
- 95711 32 oz. Refill Bottle, 12/cs



PRO ADVANTAGE ANTIBACTERIAL SOAP

Antibacterial soap contains the active ingredient Triclosan, which is an effective germ killer. It is strong on bacteria, yet still mild on hands. Choice of liquid or foam soap. Foam contains a fresh, lavender fragrance.

- P778108 8 oz., liquid soap w/ pump, 12/cs
- P778128 Gallon liquid soap, 4/cs
- P778208 7.5 oz Foam dispenser (lavender fragrance), 12/cs



GOJO MICRELL® ANTIBACTERIAL LOTION SOAP

- 2157-08 1000mL Lotion Soap, 8/cs
- 7188-10 1000mL Lotion Soap Refill, Amber, 10/cs
- 9517-18 Traditional Bag-in-box 500mL, Amber, 18/cs
- 9752-12 8 oz Pump Bottle Lotion Soap, Amber, 12/cs
- 9755-04 Pour Gallon Lotion Soap, Amber, 4/cs
- 9756-06 800mL Lotion Soap Refill, Amber, 6/cs
- 9757-12 800mL Lotion Soap Refill, Amber, 12/cs
- 9759-12 12 oz Pump Bottle Lotion Soap, Amber, 12/cs

Soaps/Cleansers



2118-08

2218-04

4212-10


4213-12

4216-04

4218-12

GOJO PROVON® ANTIMICROBIAL LOTION SOAP

- 2118-08 NXT™ Lotion Soap Refill, Amber, 1000mL, 8/cs
- 2218-04 NXT™ Lotion Soap Refill, Amber, 2000mL, 4/cs
- 4212-10 1000mL Lotion Soap Refill, Amber, 10/cs
- 4213-12 12 oz Lotion Soap Pump Bottle, Amber, 12/cs
- 4216-04 Pour Gallon Lotion Soap, Amber, 4/cs
- 4218-12 800mL Lotion Soap Refill, Amber, 12/cs



2113-08

2213-04

4013-12

4017-10

GOJO PROVON® ENRICHED LOTION CLEANSER

- 2113-08 NXT® 1000mL Refill, Pink, 8/cs
- 2213-04 2000mL Lotion Cleanser Refill, Pink, 4/cs
- 4013-12 800mL Lotion Cleanser Refill, Pink, 12/cs
- 4017-10 1000mL Lotion Cleanser Refill, Pink, 10/cs



5186-03

5286-02

GOJO PROVON® FOAMING ANTIMICROBIAL HANDWASH

- 5186-03 PROVON® Foaming Antimicrobial Handwash w/moisturizers, Orange, 1.25L, 3/cs
- 5286-02 PROVON® Foaming Antimicrobial Handwash w/moisturizers, Orange, 2L, 2/cs



5185-03

5285-02

5385-02

5388-02

5785-04

GOJO PROVON® FOAMING HANDWASH

A general purpose handwash with a pleasant berry fragrance with moisturizers. Dermatologist tested.

- 5185-03 PROVON® Foaming Handwash w/moisturizers, Berry, 1.25L, 3/cs
- 5285-02 GOJO® Luxury Foam Antibacterial Handwash, Orange, 2000mL, 2/cs
- 5285-02 PROVON® Foaming Handwash w/moisturizers, Berry, 2L, 2/cs
- 5385-02 Provon® Foaming Handwash w/advanced moisturizers, 1200mL TFX, Refill, Clear, 2/cs
- 5388-02 Provon® Medicated Foaming Handwash w/advanced Moisturizers, 1200mL TFX, Refill, Clear, 2/cs
- 5785-04 535mL Counter Top Pump Bottle Foaming Handwash, Pink Formula, 4/cs



5188-03

5288-02

5788-04

GOJO PROVON® FOAMING MEDICATED HANDWASH

A general purpose antimicrobial handwash with moisturizers and triclosan. Kills most common germs that may cause disease. Dermatologist tested.

- 5188-03 PROVON® Foaming Medicated Handwash w/moisturizers & Triclosan, Blue, 1.25L, 3/cs
- 5288-02 PROVON® Foaming Medicated Handwash w/moisturizers & Triclosan, Blue, 2L, 2/cs
- 5788-04 535mL Counter Top Pump Bottle Foaming Handwash, Blue Formula, 4/cs



2119-08

4222-10

4223-12



4226-04

4228-12

GOJO PROVON® GENTLE LOTION SOAP

- 2119-08 NXT™ Lotion Soap Refill, Yellow, 1000mL, 8/cs
- 4222-10 1000mL Lotion Soap Refill, Yellow, 10/cs
- 4223-12 12 oz Pump Bottle Lotion Soap, Yellow, 12/cs
- 4226-04 Pour Gallon Lotion Soap, Yellow, 4/cs
- 4228-12 800mL Lotion Soap Refill, Yellow, 12/cs



2102-08

2202-04

4201-24

4202-10



4203-12

4208-12

GOJO PROVON® MEDICATED LOTION SOAP W/PCMX

- 2102-08 NXT™ Lotion Soap Refill, Gold, 1000mL, 8/cs
- 2202-04 NXT™ Lotion Soap Refill, Gold, 2000mL, 4/cs
- 4201-24 4 oz. Squeeze Bottle Lotion Soap, Gold, 24/cs
- 4202-10 1000mL Lotion Soap Refill, Gold, 10/cs
- 4203-12 12 oz. Pump Bottle Lotion Soap, Gold, 12/cs
- 4208-12 800mL Lotion Soap Refill, Gold, 12/cs



2158-08

4152-10

4153-12

GOJO PROVON® MEDICATED LOTION SOAP W/TRICLOSAN

- 2158-08 1000mL Lotion Soap, Gold, 8/cs
- 4152-10 1000mL Lotion Soap Refill, Gold, 10/cs
- 4153-12 12 oz. Pump Bottle Lotion Soap, Gold, 12/cs



2108-08

2208-04

4019-10

GOJO PROVON® MILD LOTION SOAP

- 2108-08 1000mL Mild Lotion Soap, Blue, 8/cs
- 2208-04 2000mL Mild Lotion Soap, Blue, 4/cs
- 4019-10 1000mL Lotion Soap Refill, Blue, 10/cs



5485-04

6185-03

GOJO PURELL® SURGICAL SCRUB

- 5485-04 Purell® Surgical Scrub 1200mL TFX™ Refill, Clear, 4/cs
- 6185-03 Purell® Surgical Scrub 1200mL Foot Pump Refill, Clear, 3/cs



9034-12

9102-12

GOJO® 800ML BAG-IN-BOX SYSTEM

- 9034-12 White, Bag-In-Box System Dispenser, 12/cs
- 9102-12 800mL Enriched Lotion Soap, Gold, 12/cs
- 9112-12 800mL Lotion Skin Cleanser, Pink, 12/cs
- 9142-12 800mL Moisturizing Hand Cream, White, 12/cs
- 9152-12 800mL SPA BATH® For Body & Hair, Pink, 12/cs
- 9212-12 800mL Ultra Mild Antimicrobial Lotion Soap w/PCMX, Amber, 12/cs

GOJO® BULK FLAT TOP GALLON SOAP PRODUCTS

- 1275-01 Flat Top Gallon Dispenser
- 1805-04 All-Purpose Flat Top Gallon Lotion Soap, Pink, 4/cs
- 1845-04 Thick Pink Antiseptic Flat Top Gallon Lotion Soap, Pink, 4/cs



0955-04

0957-12

GOJO® NATURAL ORANGE™ PUMICE HAND CLEANER

- 0955-04 One Gallon Hand Cleaner w/pump dispenser, White, 4/cs
- 0957-12 14 oz. Hand Cleaner GRIPIT™ Bottle w/nail brush, White, 12/cs

GOJO® NXT™ 1000ML LOTION SOAP PRODUCTS

Antimicrobial handwash designed for general purpose and medium-risk use. Kills most common germs that may cause disease. Contains the broad-spectrum antimicrobial agent chloroxylenol (PCMX) (0.3%), for excellent germ-killing efficacy of gram-positive and gram-negative microorganisms, including Methicillin-resistant Staphylococcus aureus (MRSA) and Vancomycin-resistant Enterococcus (VRE). Helps reduce the risk of acquiring infections by reducing disease-causing microbes. Produces a rich, clean-rinsing lather. Developed with mild skin cleaners, emollients, and conditioners including aloe and vitamins A&E. Dermatologist tested; low irritation and allergy potential.

- 2112-08 Ultra Mild Antimicrobial Lotion Soap w/chloroxylenol, 1000mL, Amber, 8/cs



GOJO® PRO™ 2000 BAG-IN-BOX SYSTEM

- 7220-04 Rich Pink Antibacterial Lotion Soap, 4/cs

GOJO® PRO™ 5000 BAG-IN-BOX SYSTEM

- 7530-02 Shower Up® Soap & Shampoo, Pink, 2/cs
- 7556-02 Natural Orange™ w/pumice hand cleaner, White, 2/cs

Soaps/Cleansers

GOJO 1250 ML FOAMING PRODUCTS

- 5161-03 GOJO® Luxury Foam Handwash, Pink, 1250mL, 3/cs
 5162-03 FMX-12™ Luxury Foam Antibacterial Handwash, Color Orange, Orange Blossom Fragrance, For use w/dispenser FMX-12™, 3/cs



91220

91241

91298



91320

91325

KIMBERLY-CLARK LOTION SOAP

- 91220 KimCare General Lotion Soap, Pink, 800mL, 12/cs
 91241 Triangle Lotion Soap, Green, 800mL, 12/cs
 91298 KimCare Antibacterial Skin Cleanser, 800mL, Strawberry Color, 12/cs
 91320 KimCare Cleanser Hair and Body Wash, 800mL, Golden Color, 12/cs
 91325 EUROBATH® Hair and Body Shampoo, Golden, 1 gal (pour), 4/cs
 91748 Kimcare Super Duty Hand Cleanser w/Grit 8 Liter, Herbal, 2/cs

KIMBERLY-CLARK LUXURY FOAM SOAP

The KIMCARE® Luxury Foam Soap system provides a pleasing hand washing experience anywhere, anytime, for everyone. It has everything you want in a soap system: An appealing, luxurious experience to encourage more frequent hand washing, easy maintenance and operational efficiencies. This general use formula contains mild cleaning ingredients and has a skin neutral pH. In addition to being pH balanced and containing mild cleaning ingredients, Item #91176 also contains emollients to help moisturize the skin. The formula in Item #91177 contains 0.5% triclosan that kills microorganisms to help reduce the spread of germs. Available in US only.

- 91175 Light Pink, Herbal Fragrance, 1000mL, For Use w/92044 (white) & 92045 (black) Dispensers, 6/cs
 91176 Light Blue, Citrus Fragrance w/moisturizers, 1000mL, For Use in the 92044 (White), & 92045 (Black) Dispensers, 6/cs
 91177 Light Green, Fragrance Free, Antibacterial, 1000mL, For Use in the 92044 or 92045 Soap Dispensers, 6/cs



91546

92519

KIMBERLY-CLARK ULTRA-KIND® HAND WASH

ULTRA-KIND® Clear Gentle Hand Wash is designed for people who frequently wash their hands. It's a dye-free formula, with aloe vera.

- 91546 800mL box, Clear, 12/cs
 92519 Kimcare Gentle Hand Wash, Clear (Dye-Free), w/aloe vera, 500mL, 18/cs

PURELL® Instant Hand Sanitizer Skin Nourishing Foam



**PROVEN.
PREFERRED.
EFFECTIVE.**

NEW PURELL® Instant Hand Sanitizer
Skin Nourishing Foam clinically proven
 to improve skin condition in 14 days.¹

The power of PURELL is also available in original gel formula, or gel with Aloe. **For more information, visit www.healthcare.GOJO.com.**

1 Skin moisture content. Clinical field evaluation #1, (February 11, 2008 to March 21, 2008), GOJO Industries, Inc. • Study ID 2008-01-F10176. Clinical evaluation #2, (October 9, 2007 to October 11, 2007), RCTS Study ID 2301 ©2009, GOJO Industries, Inc. All rights reserved. PURELL is a trademark of Johnson & Johnson and is used under license

KIMBERLY-CLARK KIMCARE® GENERAL PINK LOTION SOAP

Provides gentle cleaning and easy to rinse lather.

91438 Soap, Pink Lotion, 1000mL (1 Liter), 3/cs



10-1500



10-1504



10-1518



10-1580

METREX VIONEX® ANTIMICROBIAL LIQUID SOAP

Creamy antimicrobial liquid soap, ideal for frequent use. Effective germ reduction in 30 seconds. Lotion-based formula for extra softness. Mild and nondrying to skin, contains conditioners, emollients and softeners. Developed for the healthcare professional. Active ingredient: chloroxylenol (PCMX)

10-1500	Vionex Gallon Liquid Soap Refill, 4/cs
10-1504	Vionex 4 oz Bottle Liquid Soap w/flip top, 24/cs (special order)
10-1518	Vionex 18 oz Bottle Liquid Soap w/pump, 12/cs
10-1580	Vionex 800 ml Liquid Soap Bag in Box Refill, 12/cs



10-1810



10-1820

METREX VIONEXUS™ NO TOUCH DISPENSER

Pre-measured application increases coverage and minimizes waste. No water or towels needed which eliminates cross-contamination. For use with VioNexus No Rinse Spray

10-1810	No Touch Dispenser
---------	--------------------

Accessories

10-1811	Replacement Pump For Dispenser
10-1812	Replacement Key
10-1820	Tray

**METREX VIONEXUS™ ANTIMICROBIAL FOAMING SOAP**

An economical foaming soap gentle enough for repeated use. Helps reduce bacteria that can cause disease. Provides strong protection against germs without being harsh to skin. Will not dry out hands. Active ingredient: chloroxylenol (PCMX)

10-1900	VioNexus 1 Liter Antimicrobial Foaming Soap, 6/cs
---------	---

**METREX VIONEXUS™ FOAMING SOAP WITH VITAMIN E**

A gentle foaming soap with a refreshing plumeria-apple fragrance. Leaves hands feeling smooth and soft. pH-balanced and hypoallergenic. Non-irritating and mild enough for repeated use. Will not dry out hands.

10-2001	VioNexus 1 Liter Foaming Soap w/Vitamin E, 6/cs
10-2002	Vionex 2 oz Foaming Soap w/Vitamin E, 48/cs

TRIAD CASTILE SOAP

A clear, odorless, pure vegetable oil-based soap at 7.7% concentrate in deionized water. Liquid castile soap is used full strength for gentle cleansing or dilute in 1 or 2 quarts of warm water for enema use.

10-8050	Castile Soap, 1/3 oz pouch, 50/bx, 12 bx/cs
---------	---

**TRIAD SOOTHE & BATHE OATMEAL TREATMENT**

Provides soothing temporary relief of itching due to dry skin, rashes, eczema, psoriasis, poison ivy, hives, sunburn, prickly heat, chicken pox, etc. Each premeasured 1.5 oz packet contains 100% oatmeal preparation.

10-7100	Regular Oatmeal Treatment, 8/bx, 24/cs
10-7200	Oilated Oatmeal Treatment, 8/bx, 24/cs

TRIAD RICHMOND GENTLE BLUE LOTION SOAP

Lotion-based skin cleaners. Packaged in #2 polyethylene plastic bottles. Packaging: disposable/closed-systems: 540ml bottle, 80ml bottle, 1000ml bottle, 1000ml pouch, 2000ml pouch.

11101R-12800	800mL Lotion Soap, Hard Bottle, Floral Scent, Blue, 12/cs
11101R-540PMP	540mL Lotion Soap, Bottle w/pump, Floral Scent, Blue, 12/cs
11101R-800PMP	1000mL Lotion Soap, Bottle w/pump, Floral Scent, Blue, 12/cs
11101R-A1000	1000mL Lotion Soap, Pouch, Floral Scent, Blue, 10/cs
11101R-A1000X8	1000mL Lotion Soap, Pouch, Circular Tip, Floral Scent, Blue, 8/cs
11101R-A2000X4	2000mL Lotion Soap, Pouch, Circular Tip, Floral Scent, Blue, 4/cs
11101R-G1000	1000mL Lotion Soap, Pouch, Wide Tip, Floral Scent, Blue, 10/cs

Soaps/Cleaners

TRIAD RICHMOND HEALTH-STAT HANDWASH

Contains 0.3% triclosan, aloe vera and lanolin. Packaged in #2 polyethylene plastic bottles. Packaging: disposable/closed-systems: 4 oz bottle, 540ml bottle, 800ml bottle, 1000ml bottle, 1000ml pouch, 2000ml pouch.

11105R-12800	Antiseptic Handwash, 800mL Hard Bottle, Green, Fruit Fragrance, 12/cs
11105R-3640Z	Antiseptic Handwash, 4mL Bottle w/flip cap, Green, Fruit Fragrance, 36/cs
11105R-540PMP	Antiseptic Handwash, 540mL Bottle w/pump, Green, Fruit Fragrance, 12/cs
11105R-800PMP	Antiseptic Handwash, 1000mL Bottle w/pump, Green, Fruit Fragrance, 12/cs
11105R-A1000	Antiseptic Handwash, 1000mL Pouch, Green, Fruit Fragrance, 10/cs
11105R-A1000E	Antiseptic Handwash, 1000mL Pouch, White, Sandalwood Fragrance, 10/cs
11105R-A1000X8	Antiseptic Handwash, 1000mL Pouch w/circular tip, Green, Fruit Fragrance, 8/cs
11105R-A2000X4	Antiseptic Handwash, 2000mL Pouch w/circular tip, Green, Fruit Fragrance, 4/cs
11105R-G1000	Antiseptic Handwash, 1000mL Pouch w/wide tip, Green, Fruit Fragrance, 10/cs

TRIAD RICHMOND MICRO-GUARD HANDWASH

Healthcare personnel antiseptic handwash with enriched emollients. Contains 0.3% triclosan. Packaging: 1000ml pouch and 800ml bottle.

11106R-12800	Antiseptic Handwash, 800mL Hard Bottle, Pink Pearlescent, Floral Fragrance, 12/cs
11106R-540PMP	Antiseptic Handwash, 540mL Bottle w/pump, Pink Pearlescent, Floral Fragrance, 12/cs
11106R-800PMP	Antiseptic Handwash, 1000mL Bottle w/pump, Pink Pearlescent, Floral Fragrance, 12/cs
11106R-A1000	Antiseptic Handwash, 1000mL Pouch w/wide tip, Pink Pearlescent, Floral Fragrance, 10/cs
11106R-A1000X8	Antiseptic Handwash, 1000mL Pouch w/circular tip, Pink Pearlescent, Floral Fragrance, 8/cs
11106R-A2000X4	Antiseptic Handwash, 2000mL Pouch w/circular tip, Pink Pearlescent, Floral Fragrance, 4/cs
11106R-G1000	Antiseptic Handwash, 1000mL Pouch, Pink Pearlescent, Floral Fragrance, 10/cs

TRIAD RICHMOND NATURA-FOAM HOMEOPATHIC HANDWASH

Formulated with tea tree oil and menthol for natural antiseptic qualities.

11120R-FM1000	Handwash, Foaming, 1000mL, 8/cs
---------------	---------------------------------

**SULTAN LYSOL® I.C.™ ANTIMICROBIAL SOAP**

Contains 0.3% Triclosan, an antibacterial active ingredient proven effective against gram-positive and most gram-negative organisms. Mild formula rinses cleanly, leaving hands soft.

95701	Antimicrobial Soap, Gallon, 4/cs
95717	Antimicrobial Soap, 17.5 oz, 12/cs
SUL 95701	Antimicrobial Soap, Gallon, 4/cs
SUL 95717	Antimicrobial Soap, 17.5 oz, 12/cs



84022



84024



84029



95503

DIAL® ANTIMICROBIAL LIQUID SOAP WITH MOISTURIZERS & VITAMIN E

Broad spectrum antimicrobial effectiveness with extra moisturizers and conditioners to protect against dry skin. Ideal for frequent hand washing.

84022	Liquid Soap, 1 Gallon, 4/cs
84024	Liquid Soap, 7.5 oz D'Ecor Pump, 12/cs
84029	Liquid Soap Refill Cartridge, 1 Liter, 8/cs
95503	Liquid Soap Refill, 800 ml Flex Pak, 12/cs

**DIAL® BORAXO HEAVY DUTY LIQUID LOTION SOAP**

A deep-cleaning lotion soap that washes thoroughly while remaining gentle on the skin. pH balanced formula has a light, pleasant fragrance and produces a thick, rich lather in seconds.

02709	Boraxo, Liquid Lotion Soap, Heavy Duty, 1 Gallon, 4/cs
-------	--

Can't find what
you're looking for?
Call your
representative.



00034

00041

00395

00396



00397



00400



81035



81064



81075

DIAL® COMPLETE™ FOAMING HANDWASH

Provides a high level of broad-spectrum germ-killing action to help stop the spread of infection. Its patented Activated Triclosan formulation provides 25X better germ-kill than other antimicrobial soaps. Effective against MRSA, VRE, Streptococcus, E. coli and Salmonella. Its hypoallergenic foam is excellent for very frequent hand washing and contains soothing skin conditioners like aloe and glycerin to leave your hands feeling soft, smooth, and refreshed. Use in any commercial setting, including healthcare facilities, schools, offices, food service, daycare, and more.

- 00034 Handwash, Complete™ Antimicrobial Foaming, Foodservice, 1 Liter, 4/cs
- 00041 Dispenser w/placard, 800mL, 6/cs
- 00162 Dial Complete, 1 Liter Refill, 6/cs
- 00395 Dispenser, Black, 1 Liter, 6/cs
- 00396 Dispenser, Smoke, 1 Liter, 6/cs
- 00397 Dispenser, Blue, 1 Liter, 6/cs
- 00400 Dispenser Complete™ Foam, White, 1 Liter, 6/cs
- 81035 Handwash Refill Cartridge, 800mL, 8/cs
- 81064 Handwash Healthcare Pump Refill, 64 oz, 4/cs
- 81075 Handwash, 7.5 oz Table Top Pump, 12/cs



00095



00910



02401

DIAL® DEODORANT BAR SOAPS - RETAIL PACKAGING

- 00095 Dial Deodorant Bar, Size 1 1/2, Unwrapped, 500/cs
- 00098 Dial Deodorant Bar, Size 2 1/2, Unwrapped, 200/cs
- 00910 Dial Antibacterial Deodorant Bar, 3.5 oz, Gold Wrapped, 72/cs
- 02401 Dial Antibacterial Deodorant Bar, 4.0 oz, Gold Wrapped, 72/cs



80733



80790



84014



84019



88047



97501

DIAL® GOLD ANTIMICROBIAL LIQUID SOAP

Broad spectrum antimicrobial effectiveness with added skin conditioners for mildness.

- 80733 Liquid Soap, 1.5 oz, 72/cs
- 80790 Liquid Soap, 16 oz Pump, 12/cs
- 84014 Liquid Soap, 7.5 oz D'Écor Pump, 12/cs
- 84019 Liquid Soap Refill Cartridge, 1 Liter, 8/cs
- 88047 Liquid Soap, Gallon, 4/cs
- 97501 Liquid Soap Refill, 800 ml Flex Pak, 12/cs



DIAL® LURON® EMERALD LOTION SOAP

Formulated with added skin conditioners to help keep hands feeling soft and smooth after use. With its rich, viscous formula and creamy lather, Luron is the perfect lotion soap for any commercial setting.

- 84050 Lotion Soap Refill Cartridge, 1 Liter, 8/cs



02814



02851



02859

DIAL® LUROSEP® ANTIMICROBIAL LIQUID SOAP WITH PCMX®

Effective in removing gram-positive and gram-negative bacteria, yeast and molds. It is nonirritating and contains the active ingredient Chloroxylenol (PCMX).

- 02814 Liquid Soap, 1 Gallon, 4/cs
- 02851 Liquid Soap, 16 oz Pump, 12/cs
- 02859 Liquid Soap Refill Cartridge, 1 Liter, 8/cs

Soaps/Cleansers



00402



00403

DIAL® PURE & NATURAL FOAMING HAND SOAP

Foaming Lotion Soap with Aloe is a mild, pleasing soap that is gentle on the skin. Dermatologist-tested and hypoallergenic, rinses clean and is excellent for frequent handwashing. Perfect foaming lotion soap for any commercial setting, including healthcare facilities, schools, offices, restaurants, daycare and more.

00402 Dispenser, Smoke, 1 Liter, 6/cs

00403 Dispenser, Blue, 1 Liter, 6/cs



00629



00628



80626



98505

DIAL® PURE & NATURAL® HYPOALLERGENIC LIQUID SOAP

A mild liquid lotion soap that is gentle on the skin. With its fresh floral fragrance and rich lotion look and feel, Pure & Natural soap is a premium product with an emphasis on value. Perfect for any commercial setting including healthcare facilities, schools, offices, restaurants, daycare and more.

00164 Pure & Natural Hypoallergenic Foaming Lotion Soap, 1 Liter Refill, 6/cs

00629 Liquid Soap Refill Cartridge, 1 Liter, 8/cs

00628 Liquid Soap, 1 Gallon, 4/cs

80626 Liquid Soap, 16 oz Pump, 12/cs

98505 Liquid Soap Refill, 800ml Flex Pak, 12/cs



80784



82834



82838



82839



91502

DIAL® SENSITIVE SKIN ANTIMICROBIAL SOAP

Broad spectrum antimicrobial effectiveness. Extra mild, and formulated for very frequent hand washing. Hypoallergenic and dermatologist tested for even the most sensitive skin.

80784 Liquid Soap, 16 oz Pump, 12/cs

82834 Liquid Soap, 7.5 oz D&Ecor Pump, 12/cs

82838 Liquid Soap, 1 Gallon, 4/cs

82839 Liquid Soap Refill Cartridge, 1 Liter, 8/cs

91502 Liquid Soap Refill, 800 ml Flex Pak, 12/cs



DIAL® SURGICAL SCRUB

Developed specifically for hospital and healthcare professional use only. Dial Surgical Scrub contains 4% Chlorhexidine Gluconate, an active antimicrobial agent, in addition to 4% Isopropyl Alcohol. Use as a surgical scrub, healthcare personnel hand wash, pre-op and skin wound cleanser to help prevent cross-contamination. The mild base provides for low irritation and contains emollients for moisturizing the skin.

02841 1 Gal Surgical Scrub, 4/cs



DIAL® SWEETHEART® ANTIBACTERIAL SOAP

Recommended for hand washing to decrease bacteria on the skin. Use in any institutional environment. Contains Triclosan.

96507 Antibacterial Soap Refill, 800 ml Flex Pak, 12/cs



DIAL® SWEETHEART® PINK LOTION SOAP

An excellent cost-effective choice for use in almost any facility.

99506 Lotion Soap Refill, 800ml Flex Pak, 12/cs

CROSSTEX SANISEPT® ANTIMICROBIAL SOAP

Effective against gram-negative and gram positive bacteria. Effective against resident microbes hidden in skin layers as well as transient microbes - gentle, non-drying. 0.3% Triclosan gives persistent antimicrobial activity. Aloe Vera, Vitamin E, Keratin and Liposome's.

JSOD Soap, 18 oz., 16/cs

JSOG Soap, Gallon, 4/cs

MEDICAL ACTION CASTILE SOAP

24-10 1/3 oz. Packets Castile Soap, 50/bx, 12 bx/cs

MEDICAL CHEMICAL TINCTURE GREEN SOAP, USP

Tincture Green Soap, USP is a combination of vegetable soap, lavender oil, glycerin and ethyl alcohol. The soap is biodegradable and mildly antiseptic. Because of the ethyl alcohol, it is a good skin degreaser. Its antiseptic properties make it useful for preoperative cleansing particularly if followed by providone iodine or a similar germicide.

140B-1GL Tincture Green Soap, USP, 128 oz. Bottle

MEDICAL CHEMICAL TINCTURE SURGICAL SOAP

A green soap with similar properties to tincture green soap, USP, but manufactured with isopropyl alcohol instead of ethyl alcohol.

141B-1GL Tincture Surgical Soap, 128 oz. Bottle



BRBG



FBG2



S1



S15C



US1

NEW WORLD FRESHSCENT™ SOAPS

Bars contain the active ingredient triclosan to help fight bacteria and odor. Significantly reduces levels of bacteria with a special formula.

BRBG	#2 Bath Gel/Liquid Soap, Boston Round Bottles, 144/cs
FBG2	2 oz Bath Gel/Liquid Soap, 96/cs
FBG4	4 oz Bath Gel/Liquid Soap, 60/cs
FBS1	#1 Soap Face & Body, Vegetable Oil, 100/bx, 5 bx/cs
FBS12	#1/2 Soap, Vegetable Oil, 100/bx, 10 bx/cs
FBS15	#1.5 Soap Face & Body, Vegetable Oil, 50/bx, 10 bx/cs
FBS3	#3 Soap Face & Body, Vegetable Oil, 72/cs
FBS34	#3/4 Soap Face & Body, Vegetable Oil, 100/bx, 10 bx/cs
S1	#1 Deodorant Soap, 50/bx, 10 bx/cs
S1/2	#1/2 Deodorant Soap, 100/bx, 10 bx/cs
S15	1.5 oz Deodorant Soap, 50/bx, 10 bx/cs
S15C	1.5 oz Clear Soap, Clear Wrap, 100/bx, 5 bx/cs
S3	3 oz Deodorant Soap, 72/cs
S3/4	#3/4 Deodorant Soap, 100/bx, 10 bx/cs
S3C	3 oz Clear Soap, Clear Wrap, 144/cs
US1	#1 Unwrapped Soap, Vegetable Oil, 100/bx, 5 bx/cs
US1/2	#1/2 Unwrapped Soap, Vegetable Oil, 100/bx, 10 bx/cs
US15	#1.5 Unwrapped Soap, Vegetable Oil, 100/bx, 5 bx/cs
US3	3 oz Unwrapped Soap, Vegetable Oil, 144/cs
US34	#3/4 Unwrapped Antibacterial Soap, Vegetable Oil, 100/bx, 10 bx/cs

LAGASSE COLGATE-PALMOLIVE LIQUID SOFTSOAP

Premium liquid soap formulated for mild, effective skin cleaning. Conditioners and moisturizing agents smooth and protect skin while maintaining its natural pH balance.

CPC26012	7.5 oz Pump Bottle w/aloe, 12/cs
CPC26017	Soft-soap Antibacterial 7.5 oz Pump Bottle, 12/cs

LAGASSE JOHNSON DIVERSEY FULL SIZE BAR SOAP

Recommended by dermatologists. Ideal for infant and long-term care. Nonirritating, unique moisture cream formula is milder to skin than any soap.

DRKCB614243 Dove® 3.5 oz Beauty Bar, 48/cs

LAGASSE RECKITT BENCKISER ANTIMICROBIAL SOAP

Mild and non-drying handwash contains 0.3% Triclosan™. Specifically formulated for health care personnel. Broad-spectrum anti-microbial effectiveness helps reduce the spread of germs. Light, pleasant fragrance. NSF/USDA Classification E4.

REC95717 Lysol® Brand I. C.™ Antimicrobial Soap, 15.5 oz Pump Bottle, 12/cs

BD E-Z SCRUB™ ANTIMICROBIAL FOAM SOLUTION DISPENSER SYSTEM

372402	32 oz. Antimicrobial Foam Solution w/4% CHG, Foot Pump Foamer, 6/cs
372404	32 oz Antimicrobial Foam Solution w/3% PCMX, Foot Pump Foamer, 6/cs
372405	32 oz Antimicrobial Foam Solution w/0.5% povidone iodine, Foot Pump Foamer, 6/cs
372412	32 oz. Antimicrobial Foam Solution w/4% CHG, Hand Pump Foamer, 6/cs
372416	32 oz. Antimicrobial Foam Solution w/2% CHG, Hand Pump Foamer, 6/cs

Accessories

370051	Foot Pump Assembly For Use w/CHG, Red, 4/cs
370052	Foot Pump Assembly For Use w/PCMX, Blue, 4/cs
370053	Foot Pump Assembly For Use w/povidone iodine, Brown, 4/cs
370055	In-Line Air Filter, 12/bx, 12 bx/cs
370056	Foot Pump Wall Bracket w/adaptor, 6/cs
370057	Foot Pump Wall Bracket Adapter For AMSCO Sink, 6/cs



632670973



632670985



632671423

BEAUMONT CITRUS II ANTIBACTERIAL SOAP

Citrus II Antibacterial Hand Soap kills a broad spectrum of germs with a proven germ fighting ingredient. Like other Citrus II products, it has a refreshing, natural citrus scent. Citrus II's unique formula helps to keep you, your staff and customers healthy. Citrus II Antibacterial Soap provides a strong protection against germs, yet gentle enough to use repeatedly. Enriched with luxurious moisturizers.

632670973	8 oz. Pump Bottle, 12/cs
632670985	1 Gallon Bottle, 4/cs
632671423	32 oz. Pump Bottle, 6/cs

**BURNISHINE EPI-WASH ANTIMICROBIAL HAND SOAP**

Lotion soap protects against germs and leaves hands non-irritated even with repeated washing.

AT-1	16 oz Hand Soap, 12/cs
AT-2	Gallon Hand Soap, 4/cs

Soaps/Cleansers

**THERMOFISHER SCIENTIFIC HANDCARE SOFTCIDE® HAND SOAPS**

Unique formulations replenish the natural oils of the skin. Effective against a broad spectrum of microorganisms. pH balanced with the skin. Very mild ingredients.

21004-12-001	Extra Mild Antimicrobial Hand Soap (1% PCMX), 4 oz Bottle w/flip-top cap, 12/cs
21016-06-001	Extra Mild Antimicrobial Hand Soap (1% PCMX), 16 oz Bottle w/pump-top, 6/cs
21027-12-001	Extra Mild Antimicrobial Hand Soap (1% PCMX), 27 oz (800ml) Bag-In-A-Box, 12/cs
21032-06-001	Extra Mild Antimicrobial Hand Soap (1% PCMX), 32 oz Bottle w/pump-top, 6/cs
21128-04-001	Extra Mild Antimicrobial Hand Soap (1% PCMX), 128 oz/1 Gallon Refill Jug, 4/cs
24027-12-001	EC Antimicrobial Hand Soap (0.3% Triclosan), 27 oz (800ml) Bag-In-A-Box, 12/cs
24032-12-001	EC Antimicrobial Hand Soap (0.3% Triclosan), 32 oz Bottle w/pump top, 12/cs
25027-12-001	NA Non-Antimicrobial Hand Soap, 27 oz (800ml) Bag-In-A-Box, 12/cs
25032-12-001	NA Non-Antimicrobial Hand Soap, 32 oz Bottle w/pump top, 12/cs

THERMOFISHER SCIENTIFIC HANDCARE ACCESSORIES

11901-12-001	Pump-Top For 128 oz/1 gallon refill jug, 12/cs
21902-12-001	16 oz & 32 oz Wall Bracket, 12/cs
DISP027-12-001	27 oz (800ml) Bag-In-Box, Wall Mounted Dispenser, 12/cs

**PRO ADVANTAGE HAND & BODY MOISTURIZING LOTION**

Rich, soothing, pH balanced formula absorbs quickly into the skin with no greasy residue. Seals in the skin's natural oils. Perfect for sensitive skin. Orange blossom scent.

P771022	2 oz Bottle w/ flip top cap, 96/cs
P771024	4 oz Round Bottle w/ flip top cap, 60/cs
P771028	8 oz Bottle w/ flip top cap, 48/cs
P771128	Gallon Bottle, 4/cs

**3M™ CAVILON™ MOISTURIZING LOTION**

9205	Moisturizing Lotion, 16 oz Bottle, 12/cs
9215	Moisturizing Lotion, 2 oz Bottle, 48/cs

**3M™ CAVILON™ ONE-STEP SKIN CARE LOTION**

The lotion cleans, replenishes the natural lipids in the skin and leaves behind a skin protectant. May be used for general patient bathing or for cleaning urine, stool and other body fluids from the skin.

3383	8 oz Bottle Skin Care Lotion, 12/cs
------	-------------------------------------

FERNDALDE PRAX® LOTION

A topical lotion with proven anti-pruritic action. It is specially formulated as a gentle and safe way to treat pruritus. Applied locally in a 1% concentration, for the relief of discomfort and pain associated with itching dermatoses, hemorrhoids and rectal surgery, episiotomies, anogenital pruritus, moderate burns and sunburns.

0748-03	8 oz Lotion (minimum order 12 ea)
0748-04	4 oz Lotion (minimum order 12 ea)
0748-15	15ml Lotion



3826-12



3838-06

GOJO AQUELL® MOISTURIZING LOTION

3826-12	12 oz Lotion Pump Bottle, White, 12/cs
3838-06	500mL Lotion Refill (for use in 3800 dispenser only), White, 6/cs



2533-08



4231-24



4236-04

GOJO PROVON® MOISTURIZING HAND & BODY LOTION

- 2533-08 1000mL Moisturizing Hand & Body Lotion, White, 8/cs
- 4231-24 4 oz Lotion Squeeze Bottle, White, 24/cs
- 4234-12 12 fl oz. Pump Bottle, White, 12/cs
- 4236-04 Pour Gallon Lotion, White, 4/cs



4251-24



4252-24

GOJO PROVON® SKIN MOISTURIZER

- 4251-24 4 oz Moisturizer Squeeze Bottle, White, 24/cs
- 4252-24 8 oz Moisturizer Squeeze Bottle, White, 24/cs

J&J AVEENO® BODY MOISTURE LOTION

- 001081 2.5 oz. Tube Daily Moisturizing Lotion, 6/bx, 2 bx/cs



10-1304



10-1318



10-1399

METREX VIONEX® SKIN LOTION

A fresh-scented skin lotion for extra moisturizing protection. Non-greasy formulation. Soothes hands with aloe. Does not contain mineral oil or petroleum. Refreshing fragrance. Latex safe.

- 10-1304 Vionex 4 oz Skin Lotion w/flip top, 24/cs (special order)
- 10-1318 Vionex 18 oz Skin Lotion w/pump, 12/cs

Accessories

- 10-1398 Wall Mount Bracket, 12/cs
- 10-1399 Wall Mount Dispenser, 12/cs



BRL



FL2



PKL

NEW WORLD FRESHSCENT™ HAND & BODY LOTION

- BRL #2 Hand & Body Lotion, Boston Round Bottles, 144/cs
- FL2 2 oz Hand & Body Lotion, 96/cs
- FL4 4 oz Hand & Body Lotion, 60/cs
- FL8 8 oz Hand & Body Lotion, 36/cs
- FLCB4 4 oz Cocoa Butter Lotion, 60/cs
- PKL .25 oz Hand & Body Lotion, 500/bx, 2 bx/cs

SMITH & NEPHEW SECURA™ MOISTURIZING LOTION

Use daily to replenish moisture. Contains humectants and emollients to moisturize and condition. Fragrance free; gentle to sensitive skin. Maintains healthy skin.

- 59433400 8 oz. Bottle Secura™ Moisturizing Lotion, 24/cs



SMITH & NEPHEW UNICARE® MOISTURIZING LOTION

Soothes, conditions and moisturizes dry skin. pH-balanced and safe; non-irritating, non-sensitizing and non-toxic. A rich lotion, it can be applied easily and is quickly absorbed into the skin.

- 445000 8 oz. Bottle Moisturizing Lotion, 24/cs
- 445400 2 oz. Bottle Moisturizing Lotion, 144/cs

Please ask before
ordering a product based
on the pictured item.

Lotions/Moisturizer - Ointments/Creams

J&J LUBRIDERM

Packaging: 6/box, 2 boxes/case unless noted.

48816	Scented 6 oz Lubriderm
48826	Unscented 6 oz Lubriderm, 6/bx
48846	Scented 16 oz Lubriderm, 6/bx
48856	Unscented 16 oz. Lubriderm, 6/bx
48863	Unscented 3.3 oz Tube Lubriderm, 6/bx, 2 bx/cs
48881	Fragrance Free 1oz Lubriderm Daily Moisture, 72/cs

CROSSTEX ALOECARE PLUS 3® SKIN CARE LOTION

Aloe Vera, Vitamin E, keratin and Liposome's. Moisturizing, Non-greasy. Latex compatible. Nourishes & strengthens skin and nails.

JLD	Lotion, 18 oz., 16/cs
-----	-----------------------



3M™ CAVILON™ EMOLLIENT CREAM

An all-body emollient cream for moisturizing or preventing dry skin. Contains vitamins A, D, and E. pH balanced to preserve the skin's natural acid mantle.

3385	4 oz Tube Emollient Cream, 12/cs
------	----------------------------------



3M™ CAVILON™ FOOT EMOLLIENT CREAM

3386	4 oz Tube, Foot Cream, 12/cs
------	------------------------------



THERMOFISHER SCIENTIFIC HANDCARE SOFTGUARD® EXTRA-STRENGTH BARRIER HAND CREAM

Helps prevent skin irritation caused by latex gloves and frequent washing. Latex, PCMX and CHG compatible. Extra-strength antimicrobial preservative system to prevent cross-contamination. Petroleum, mineral oil and fragrance-free. Absorbs immediately. Non-greasy.

11002-24-001	Barrier Hand Cream, 2 oz Bottle w/flip-top cap, 24/cs
11016-06-001	Barrier Hand Cream, 16 oz Bottle w/pump-top, 6/cs
11027-12-001	Barrier Hand Cream, 27 oz (800ml) Bag-In-A-Box, 12/cs
11128-04-001	Barrier Hand Cream, 128 oz/1 Gallon Refill Jug, 4/cs
12003-12-001	Barrier Hand Cream, 3 oz Tube w/flip-top cap, 12/cs
12008-06-001	Barrier Hand Cream, 8 oz Tube w/flip-top cap, 6/cs

FERNDALE NOURIVA REPAIR™ MOISTURIZING CREAM

The first skin care cream for dry skin with both a rapid and a lasting effect because it adds the natural oils that the skin itself normally takes days to produce.

0865-30	Nouriva Repair™ Moisturizing Cream, 30gm
---------	--



COVIDIEN/KENDALL VASELINE CONSTANT CARE CONDITIONING CREAM

Non-greasy cream formulation that penetrates skin quickly. Leaves a light moisture barrier. Non-irritating.

8884737300	Conditioning Cream, 2 oz. Bottle, 12/cs
8884737400	Conditioning Cream, 8 oz. Bottle, 12/cs



PRO ADVANTAGE PERINEAL CREAM

Designed to prevent moisture, urine & fecal matter from contacting skin. Provides a non-greasy barrier and allows the skin to breathe without being occlusive. Rubs in clear to permit easy monitoring. Unique formula of dimethicone and zinc oxide helps protect the skin. Dimethicone is a silicone-based skin protectant that retains moisture in the skin, while zinc oxide protects against skin irritants. Vitamins A, D and E nourish the skin. Specially formulated without dyes or perfumes.

P770024	Single use packets .27 fl oz (8ml), 144/cs
P770035	4 oz. tube, 12/CS



SMITH & NEPHEW UNI DERM® MOISTURIZING CREAM

For perineal and all-over body use, the cream's emollient conditions and moisturizes dry skin. Non-irritating, non-sensitizing and non-toxic. The water-based pH-balanced formula smooths on easily and quickly vanishes into the skin. Non-greasy and non-staining.

443500 3 oz. Tube Moisturizing Cream, 24/cs

DERMA SCIENCES MEDIHONEY OINTMENT

31515 Medihoney Tube, 1.5 fl. oz. w/applicator, 12/cs
31535 Medihoney Tube, 3.5 fl. oz. w/applicator, 12/cs



3M™ CAVILON™ DURABLE BARRIER CREAM

A concentrated, therapeutic cream that provides durable, long-lasting protection from body fluids while it moisturizes the skin. Hypoallergenic and pH balanced.

3391 1 oz Tube Barrier Cream, 48/cs
3392 3.25 oz Tube Barrier Cream, 12/cs

BURNISHINE EPI-SAFE BARRIER CREAM

Bands with upper layer of skin to protect against latex glove irritation and germs. Each application lasts 3-4 hours.

AC-1 16 oz Barrier Cream, 12/cs

SMITH & NEPHEW SECURA™ DIMETHICONE PROTECTANT

Contains 5% dimethicone for maximum protection. Protects minor skin irritation associated with diaper rash and helps protect from wetness. Helps prevent and temporarily protects chafed, chapped, cracked, or windburned skin and lips. Absorbs quickly and is not greasy; will not wash away. Does not clog briefs. Propylene glycol moisturizes skin. Transparent application allows visual skin inspection. Moisturizes and conditions dry skin.

59432200 4 oz. Tube Secura™ Dimethicone Protectant, 12/cs

SMITH & NEPHEW SECURA™ EXTRA PROTECTIVE CREAM

Contains 30% zinc. For treatment and prevention of rash associated with diaper use or exposure to urine and feces. Helps seal out wetness. Provides a high level of protection for patients with severe diarrhea. Will not wash away after repeated exposure to urine and feces. Contains karaya to absorb moisturize and help the product to adhere to weepy macerated skin. Soothing, nonsensitizing and contains vitamin E. Maintains healthy skin.

59432400 3.25 oz. Tube Secura™ Extra Protective Cream, 24/cs
59432500 7.75 oz. Tube Secura™ Extra Protective Cream, 12/cs

SMITH & NEPHEW SECURA™ MOISTURIZING CREAM

Restores moisture to very dry skin. Conditions and soothes. Contains petrolatum and vitamin E. Non-staining. Use daily. Absorbs quickly and non-greasy.

59431900 3 oz. Tube Secura™ Moisturizing Cream, 24/cs
59432000 6.5 oz. Tube Secura™ Moisturizing Cream, 12/cs

SMITH & NEPHEW SECURA™ PROTECTIVE CREAM

Contains 10% zinc. For treatment and prevention of rash associated with diaper use or exposure to urine and feces. Helps seal out wetness. Clove oil helps mask odor. Pediatrician tested for safety. Aloe, vitamin E and allantoin moisturize and condition the skin.

59431100 1.75 oz. Tube Secura™ Protective Cream, 24/cs
59431200 2.75 oz. Tube Secura™ Protective Cream, 24/cs

SMITH & NEPHEW SECURA™ PROTECTIVE OINTMENT

Contains 98% petrolatum. For treatment and prevention of rash associated with diaper use or exposure to urine and feces. Helps seal out wetness. Clove oil helps mask odor. Vitamins A and E soothe and condition sensitive skin.

59431500 2.47 oz. Tube Secura™ Protective Ointment, 24/cs
59431600 5.6 oz. Tube Secura™ Protective Ointment, 12/cs
59434800 Unit Dose Secura™ Protective Ointment, 150/bx, 4 bx/cs

SMITH & NEPHEW SECURA™ SKIN CARE KITS

59434100 Secura™ EPC Skin Care Kit w/Secura™ 4 oz. Moisturizing Cleanser & Secura™ 3.25 oz. Extra Protective Cream, 24/cs
59434200 Secura™ Starter Kit w/Secura™ 8 oz. Moisturizing Cleanser, Secura™ 3 oz. Secura™ Moisturizing Cream & Secura™ 2.47 oz. Protective Ointment, 24/cs
59434300 Secura™ Personal Skin Care Kit w/Secura™ 8 oz. Personal Cleanser & Secura™ 2.47 oz. Protective Ointment, 24/cs
59434400 Secura™ Two-Step Kit w/Secura™ 8 oz. Moisturizing Cleanser & Secura™ 2.75 oz. Protective Cream, 24/cs

CROSSTEX SANICARE® BARRIER CREAM

Non-sticky-forms a molecular bond with skin. 0.3% Triclosan inhibits microbial growth. Absorbed immediately. Aloe Vera, Keratin, Liposome's & Vitamins A & E to nourish and heal skin. Protects normally porous skin 3-4 hours, even with repeated washing.

JBC Cream, 18 oz., 16/cs

J&J TUCKS OINTMENT

Promptly relieves itching and burning.

10531 Tucks Ointment (formerly Anusol) 1 oz., 24/cs (14 Digit Case UPC/GTIN# 00312547150316)
15037 Tucks HC Ointment (formerly Anusol) .7 oz., 24/cs (14 Digit Case UPC/GTIN# 00312547150378)

PRO ADVANTAGE PERINEUM CLEANSER

No-rinse spray wash for incontinence and stoma care. Pleasantly scented with water lily papaya. Its gentle, soothing formula quickly cleans, refreshes and eliminates odors. Rich in lanolin and aloe, it replenishes natural skin oils.

P772008 8 oz Bottle w/pump spray, 48/cs
P772128 Gallon w/empty spray bottle, 4/cs

SMITH & NEPHEW SECURA™ MOISTURIZING CLEANSER

No-rinse, one-step cleansing for perineum and body saves time. Gentle surfactants aid in removal of urine, feces, or other foreign material. Non-irritating. Contains an antimicrobial agent to reduce odor. Conditions and moisturizes skin. pH balanced to protect the skin's acid mantle. Pediatrician tested for safety. Contains aloe and glycerin to condition and moisturize skin.

59430800 4 oz. Bottle Secura™ Moisturizing Cleanser, 24/cs
59430900 8 oz. Bottle Secura™ Moisturizing Cleanser, 24/cs

SMITH & NEPHEW SECURA™ PERSONAL CLEANSER

No-rinse, one-step cleansing for perineum and body saves time. Gentle surfactants aid in removal of urine, feces, or other foreign material. Non-irritating. Contains an antimicrobial agent to reduce odor. Conditions and moisturizes skin. pH balanced to protect the skin's acid mantle.

59430400 8 oz. Bottle Secura™ Personal Cleanser, 24/cs
59430500 1 gal Bottle Secura™ Personal Cleanser, 4/cs

Perineal Wash



DIAL® PERINEAL CLEANSER

Specially formulated to help clean, refresh and wash away odors. Convenient, no rinse formula may be used for both perineal and total body cleansing. Dial Perineal Cleanser is mild, nonmedicated formula is enriched with natural skin conditioners, including aloe, to help soothe skin.

- 99600 Perineal Cleanser, 10 oz Bottles, Spray Pump, 20/cs
99650 Perineal Cleanser, Gallon, 4/cs



GOJO PROVON® ANTIBACTERIAL PERINEAL WASH

No rinse formula.

- 4432-12 12 oz Perineal Wash Spray Bottle, Clear, 12/cs
4432-12-MT Empty Bottle 10 oz, 12/cs
4434-24 4 oz Perineal Wash Spray Bottle, Clear, 24/cs
4436-04 Pour Gallon Perineal Wash, Clear, 4/cs



GOJO PROVON® PERINEAL WASH

No rinse formula.

- 4422-12 12 oz. Perineal Wash Spray Bottle, Green, 12/cs
4424-24 4 oz. Perineal Wash Spray Bottle, Green, 24/cs
4426-04 Pour Gallon Perineal Wash, Green, 4/cs

Accessories

- 4422-12-MT Empty Bottle 10 oz, 12/cs



NEW WORLD FRESHSCENT™ RINSE FREE PERINEAL WASH

- PNR128 Perineal Wash 1 Gallon, 4/cs
PNR8 8 oz Perineal Wash w/spray pump, 36/cs



BIRCHWOOD B-SURE® PERINEAL WASH

Use for cleaning the perineal area of fecal matter and urine due to incontinence. Safe for use as an all-over skin cleanser. Enriched with Aloe and Glycerin to soothe and moisturize irritated skin. Gentle pH-balanced formula relieves and protects sensitive skin.

- 92101 Perineal Wash, B-SURE®, 8oz Bottle, 12/cs
92103 Perineal Wash, B-SURE®, Gallon Bottle, 4/cs

DERMA SCIENCES IN BETWEEN® PERI-WASH

A versatile product for the quick and effective clean up of all body exudates, particularly those from incontinent episodes. It eliminates odor from urine and feces while gently and thoroughly cleansing the entire perineal skin area. It can be sprayed directly on the skin, soil mass, clothing, bed linens, floors, walls, furniture, carpets, and any other surface that needs to be clean and free of odors. Cleans, soothes, and freshens the skin. It leaves no soap film and does not require rinsing. Contains a special antimicrobial preserving system that aids in eliminating cross infection concerns.

- S201-1G Peri-Wash, Gallon Bottles, 4/cs
S201-8Z Peri-Wash, 8 oz Bottle, 24/cs



GOJO PROVON® 3N1 PERINEAL WASH CREAM

- 2110-08 Perineal Wash Cream, NXT® 1000mL, 8/cs



NEW WORLD PERI BOTTLE

- PERI 8 oz Peri Bottle, Individually Bagged, 50/cs